SevenGo (Duo) pro™/OptiOx™

English

Operating Instructions **SevenGo (Duo) pro™/OptiOx™** SG9/SG98

Français

Mode d'emploi **SevenGo (Duo) pro™/OptiOx™** SG9/SG98

Español

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DOMINIQUE DUTSCHER SAS

Table of Contents

1	Introduction	3
2	Safety Measures	4
2.1	Definitions of signal words and warning symbols	4
2.2	Product specific safety notes	4
3	Installation	6
3.1	Installing the batteries	6
3.2	Sensor preparation	6
3.3	OptiOx™ BOD adapter and protective guard	7
3.4	Fitting the wrist strap	7
3.5	SevenGo™ clip	7
4	Operating the pH/ORP/ion/DO meter SG98 and DO meter SG9	9
4.1	Meter layout	9
4.2	The display	10
4.3	Key controls	12
4.4	Using the softkeys	13
4.5	Navigating between menus	13
4.6	Navigating within a menu	13
4.7	Using the alphanumeric keypad	14
4.7.1	Alphanumeric input	14
4.7.2	Entering IDs/PIN	14
4.7.3	Editing values in a table	14
4.8	Calibration	15
4.8.1	Running a one-point pH/ion calibration (only SG98)	15
4.8.2	Running a multi-point pH/ion calibration (only SG98)	15
4.8.3	Automatic buffer recognition (only SG98)	15
4.8.4	Dissolved oxygen (DO) calibration with the InLab® OptiOx	16
4.9	Temperature compensation	17
5	Setup	18
5.1	Menu structure of setup	18
5.2	Sample ID	18
5.3	User ID	18
5.4	Data logging	18
5.5	System settings	19
5.6	Instrument self-test	20
6	Menus and settings	21
6.1	Menu structure of pH/ion (only SG98)	21
6.2	Menu structure of DO	21
6.3	Temperature settings	21
6.4	pH/ion calibration settings (only SG98)	21
6.4.1	Buffer groups / Standards	21
6.4.1.1	Predefined pH buffer groups	21
6.4.1.2	Customized pH buffer group	21
6.4.2	Calibration mode	22
6.4.3	Calibration reminder	22
6.5	pH/ion measurement settings (only SG98)	22
6.6	DO measurement settings	23
6.7	DO calibration reminder	24
6.8	Endpoint formats	24
6.9	Measurement limits	25

6.10	Sensor ID/SN	25
7	Data management	26
7.1	Menu structure of data menu	26
7.1.1	SG98	26
7.1.2	SG9	26
7.2	Measurement data	26
7.3	Calibration data	27
7.4	ISM data	27
8	Maintenance	29
8.1	Meter maintenance	29
8.2	pH electrode maintenance	29
8.3	InLab® OptiOx sensor maintenance	29
8.4	Troubleshooting InLab® OptiOx	30
8.5	Replacement of the OptiOx sensor cap	31
8.6	Interfering substances of the InLab® OptiOx sensor	31
9	Error messages	33
9.1	Error limits	35
10	Disposal	36
11	Sensors, solutions and accessories	37
12	Specifications	39
13	Appendix	42
13.1	Buffer tables	42
13.2	Solubility of oxygen in water as a function of temperature and salinity	44
14	Declaration of conformity	46

1 Introduction

Thank you for purchasing this METTLER TOLEDO meter. SevenGo Duo pro™ SG98 and SevenGo pro™ SG9 combined with the InLab® OptiOx optical dissolved oxygen sensor are not only easy-to-operate portable meters for precise measurements, they also contain many exciting features:

- **New ISM®** (Intelligent Sensor Management) **technology**: the meter automatically recognizes the sensor and transfers the last set of calibration data from the sensor chip to the meter. The last five calibrations as well as the initial calibration certificate are also stored on the sensor chip. These can be reviewed, transferred and printed. ISM® provides additional security and helps eliminate mistakes.
- **Multi-language graphical user interface** on a backlit display with intuitive menu guidance, making the operating instructions primarily a source of reference.
- **Easy switching** between the various parameters before or after the measurement.
- **IP67 rating – fully water proof**. The rating refers to the meter, the sensor and the connections. The meter is perfectly suited for indoor as well as outdoor use.

In addition to the new features, SevenGo Duo pro™ SG98 and SevenGo pro™ SG9 meters provide the same high quality standards as in all other SevenGo™ models:

- **Excellent ergonomics** – as if the meter is part of you.
- **Great flexibility** in the mode of operation and transport – the ultimate aid for all measurements in the plant as well as in the field.
- **RDO®** (Rugged Dissolved Oxygen) **technology**: The InLab® OptiOx optical dissolved oxygen sensor is based on the reliable RDO technology. Thanks to RDO, measuring dissolved oxygen is easier than ever before:
 - Stable results and rapid response time
 - Immediately ready for use – no polarization needed
 - Extremely easy handling, zero, maintenance: save time!
 - Suitable for an extremely wide range of applications

Conventions and symbols

Refers to an external document.

Note

For useful information about the product.

Elements of instructions

Instructions always contain action steps and can contain prerequisites, intermediate results and results. If an instruction contains more than one action step, the action steps are numbered.

- Prerequisites that must be fulfilled before the individual action steps can be executed.
- 1 Action step 1
 - Intermediate result
 - 2 Action step 2
 - Result

2 Safety Measures

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER A hazardous situation with high risk, resulting in death or severe injury if not avoided.

WARNING A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

CAUTION A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Notice

2.2 Product specific safety notes

Intended use

This instrument is designed for a wide range of applications in various areas and is suitable for measuring pH and dissolved oxygen.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

CAUTION

Environmental influences

Avoid the following environmental influences:

- Powerful vibrations
- Direct sunlight
- Atmospheric humidity greater than 80%
- Corrosive gas atmosphere
- Temperatures below 5 °C and above 40 °C
- Powerful electric or magnetic fields

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

WARNING

Explosion hazard due to spark formation, corrosion caused by the ingress of gases

The housing of the instrument is not gas tight. Never work in an environment subject to explosion hazards!

WARNING

Serious injury due to chemicals and solvents

When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

FCC Rules

This device complies with Part 15 of the FCC Rules and Radio Interference Requirements of the Canadian Department of Communications. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

3 Installation

Carefully unpack the meter. Keep the calibration certificate in a safe place.


3.1 Installing the batteries

NOTICE

Damage to the instrument due to the unsealed battery cover

The IP67 rating requires the battery compartment to be perfectly sealed. The sealing ring around the battery cover must be replaced if it is damaged in any way.

- 1 Slide the release button on the battery cover in the direction of the arrow.
- 2 Hold the lid with two fingers and remove it.
- 3 Insert the batteries in the battery compartment, as indicated by the arrows on the inside of the compartment.
- 4 Replace the battery cover and push back the button to fix the lid in place.

3.2 Sensor preparation


Preparation of a pH sensor

Follow the instructions given in the pH sensor manual.

Preparation of an InLab® OptiOx sensor

Note

The OptiOx™ sensor contains an internal clock which counts down the 365 day life span of a new sensor cap. The countdown begins once the OptiOx™ cap has been attached, the sensor connected to the device and the first measurement carried out. This process cannot be undone once the first measurement has been carried out.

- Remove the protective shipping cap from the sensor. Keep the protective shipping cap for later use. See **A**.
- Ensure that the two O-rings on the sensor are correctly positioned. See **B**.
- Line up the arrow on the OptiOx cap with the arrow on the OptiOx sensor. See **C**.
- Push the OptiOx cap onto the OptiOx sensor until the cap is firmly connected to the sensor. Do not turn the OptiOx cap. See **D**.

Note

Do not remove the OptiOx cap after installation until a cap replacement is necessary.

Connecting an IP67 sensor

To connect the IP67 sensor, make sure that the plugs are properly inserted. Twist the RCA (Cinch) / mini LTW plug to ease the attachment of the sensor.

Connecting an ISM® sensor

ISM® sensor

When connecting an ISM® sensor to the meter, one of the following conditions have to be met for the calibration data to be transferred automatically from the chip of the sensor into the meter and is used for further measurements. After attaching the ISM® sensor,


- the meter must be switched on.
- If the meter is already switched on, the **READ** key is pressed.
- If the meter is already switched on, the **CAL** key is pressed.

We strongly recommend you to switch off the meter when disconnecting an ISM sensor. In doing so, you make sure that the sensor is not removed while the instrument is reading data from or writing data to the ISM-chip of the sensor.

The **ISM** icon **ISM** appears on the display and the sensor ID of the sensor chip is registered and appears on the display.


3.3 OptiOx™ BOD adapter and protective guard

Installing the BOD adapter or the protective guard

The special BOD (biochemical oxygen demand) adapter enables quick and easy measurement in all current types of BOD bottles with the InLab® OptiOx.

Thanks to the adapter the sensor reaches into the bottle only as far as needed so that less water is being displaced during measurement. According to the EPA (Environmental Protection Agency, USA), stirring is not necessary thanks to RDO technology.


The robust protective guard made of stainless steel provides optimum protection even in harsh environments.

Thanks to the additional weight of the protective guard it also acts as a sinker taking the InLab® OptiOx to deeper measurement locations.

- Unscrew and remove the thread ring from the InLab® OptiOx. Keep it for later use.
- Slide the BOD adapter or the protective guard over the front of the InLab® OptiOx and screw it onto the sensor.

3.4 Fitting the wrist strap

- Fit the wrist strap as shown in the diagram.

3.5 SevenGo™ clip

The SevenGo™ clip is an electrode holder that can be placed next to the display on either side of the housing.

Note

The SevenGo™ clip cannot be used for InLab® OptiOx sensors.


- To mount the clip, remove the cover over the clip's fixing point using your thumbnail.
- Attach the clip by pressing it into the recess.
- Slide the shaft of the sensor into the clip from the top.
- Rotate the sensor around the clip's axis to switch between the storage and working positions.

DOMINIQUE DUTSCHER SAS



4 Operating the pH/ORP/Ion/DO meter SG98 and DO meter SG9

4.1 Meter layout

- 1 **BNC socket** for mV/pH signal input (only SG98)
- 2 **RCA (Cinch) socket** for pH temperature signal input (only SG98)
- 3 **Mini LTW socket** for DO and DO temperature signal input
- 4 **Slots** for attaching the wrist strap
- 5 **Fixing points** for SevenGo™ clip (both sides)
- 6 **Display**
- 7 **Battery cover**
- 8 **Rubber key pad**
- 9 **Bottom cap (blue)** over the field assistant's fixing point
- 10 **Rubber feet** fixing points

4.2 The display

- 1 **Battery status** icon
- 2 **Auto-off override** icon
- 3 **Data logging** icon (timed interval reading)
- 4 **Backlight** icon
- 5 **Routine mode** icon (user access rights are restricted)
- 6 Date and time
- 7 Measurement temperature
- 8 Endpoint format
- 9 Atmospheric pressure
- 10 Temperature compensation
 - **ATC**: Temperature sensor connected
 - **MTC**: no temperature sensor connected or detected
- 11 Number of data sets in memory
- 12 User ID
- 13 Softkey
- 14 Softkey

- 15 Softkey
- 16 Sample ID
- 17 Sensor ID

18 DO OptiOx sensor cap lifetime icon

New sensor cap

Sensor cap needs to be replaced in less than 6 months

Sensor cap needs to be replaced in less than 3 months

Sensor cap needs to be replaced in less than 1 month

Sensor cap needs to be replaced in less than 2 weeks

Sensor cap needs to be replaced in less than 2 days

19 pH electrode condition criteria (only SG98)

Slope: 95-105%
Offset: $\pm(0-15)$ mV
Electrode is in good condition

Slope: 94-90%
Offset: $\pm(15-35)$ mV
Electrode needs cleaning

Slope: 89-85%
Offset: $\pm(>35)$ mV
Electrode is faulty

20 ISM® sensor connected

Stability criterion (only SG98)

Strict

Medium

Fast


22 Warning messages

23 Buffer groups or standards

DOMINIQUE DUISCHER SAS

4.3 Key controls

Key	Press and release	Press and hold for 3 seconds
	Switch meter on or off	Switch meter on or off
	Start or endpoint measurement (measurement screen) Confirm input or start editing a table Exit setting and go back to measurement screen	Turn backlight on or off
	Start calibration	Review the last calibration data
	Switch mode in single channel (measurement screen) Discard setting and go back to previous menu (setting screens)	Switch between single and dual channel display (measurement screen) (only SG98)

Measurement modes

A single channel has to be selected first in order to switch the measurement mode (only SG98).

- Press and hold the **MODE** key to switch between the dual and single channel measurement screen (only SG98).
- Press and release the **MODE** key in the single channel display to change between the different measurement modes.

The sequence of the alternating measurement modes for pH/ion measurement (only SG98) is:

1. pH
2. mV
3. rel. mV
4. ion

For the DO measurement the sequence is:

1. saturation (%)

2. ppm
3. mg/L

4.4 Using the softkeys

The meter has three softkeys. The functions assigned to them change during operation depending on the application. The assignment is shown on the bottom line of the screen.

In the measurement screen, the three softkeys are assigned as follows:

Menu	Store	Data
Access meter settings	Save an endpointed measurement	Access data menu

The other softkey functions are:

	Move one position to the right	Edit	Edit table or value
	Move one position to the left	End	End calibration
	Scroll up in the menu	Yes	Confirm
	Scroll down in the menu	No	Reject
	Increase value	Review	Review selected data
	Decrease value	Save	Save data, setting or value
	Scroll to next data set in memory	Select	Select the highlighted function or setting
	Delete letters or numbers on alphanumeric keypad	Start	Begin the reference measurement
Delete	Delete selected data	Trans	Transfer selected data

4.5 Navigating between menus

The meter display consists of a measurement frame, softkeys, areas for status icons and underlying menu areas.

To access the menu areas and to navigate between them, use various softkeys (see "Using the softkeys").

- 1 Press **Menu**.

⇒ The **Setup** menu appears and **Sample ID** is highlighted.

- 2 Press  to highlight the **Setup** tab.
- 3 Press  to highlight the **pH/Ion** tab (only SG98).
- 4 Press  to highlight the **DO** tab.
- 5 Press **MODE/EXIT** to return to the measurement screen.

4.6 Navigating within a menu

This example is based on the **Setup** menu, but the procedure applies to the other menus as well.

- Press **Menu**.

⇒ The **Setup** menu appears and **Sample ID** is highlighted.

- Press  as often as needed to navigate to a menu item.
- Press **Select** to move deeper in the menu for the chosen operation.
- Continue navigating with ,  or **Select** until the final destination is reached within the menu.
- Press **MODE/EXIT** to go back to the previous menu.
— or —
- Press **READ** to return to the measurement screen directly.


4.7 Using the alphanumeric keypad

4.7.1 Alphanumeric input

The meter has a screen keypad for entering IDs, SNs and PINs. Both numbers and letters are allowed for these entries.

Note

- When entering a PIN, each character entered will be displayed as (*).

- Press **<** to move left to highlight number or letter, use **>** to move right, and **↓** to move downwards.
- Press **READ** to confirm the entry.
 - The line where the position of alphanumeric character is being entered blinks.
- To end and confirm entry, use softkeys to highlight screen key **OK**, and press **READ** to save the ID.
— or —
- To delete information, use softkeys to highlight **** and press **READ** to delete the previously entered character.
— or —
- Press **MODE/EXIT** to return to the upper level of the menu.
 - The entries are rejected.

4.7.2 Entering IDs/PIN

The three softkeys and **READ** key are used for navigating on the keypad and entering the ID/PIN.

Example: WATER

- If **A** is highlighted, press **↓** three times.
 - V** is highlighted.
- Press **>** once.
 - W** is highlighted.
- Press **READ** to enter **W**.
- Repositioning the highlighted bar to **A**, **T**, **E** and **R**, and press **READ** to enter each letter of sample ID in sequence as described in steps a - c.
- Reposition the highlighted bar to **OK**, and press **READ** to save the sample ID.

4.7.3 Editing values in a table

The meter has a feature, which allows the user to enter, edit or remove values in tables. (for example, temperature and buffer values for a customized buffer group). This is accomplished by using the softkeys on the display to navigate from cell to cell.

- Press **READ** to start editing the cell in the table.
 - The softkeys on the display change.
- Press **+** and **-** to enter the value and press **READ** to confirm.

- ➔ The softkeys change back to **↑** and **↓**.
- 3 Navigate to a cell and press **Delete** to remove a value.
- 4 To finish editing the table, navigate with the **↑** and **↓** to highlight **Save**.
- 5 Press **READ** to confirm the action and exit the menu.

4.8 Calibration

The meter allows you to perform pH and ion calibrations for up to 5 points (only SG98) and DO calibrations for up to 2 points.

4.8.1 Running a one-point pH/ion calibration (only SG98)

- 1 Press and hold **MODE** for 3 seconds to switch to single channel pH or ion measurement screen when in dual channel measurement.
- 2 Place the electrode in a calibration buffer and press **CAL**.
 - ➔ **Cal 1** appears on the display.
- 3 The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **READ**.
 - ➔ The relevant buffer value is shown on the display.
- 4 Press **End** to accept the calibration and return to sample measurement.
 - ➔ The calibration result (offset and slope for pH) is shown on the display.
— or —
- 5 Press **Save** to keep the calibration.
- 6 Press **EXIT** to reject the calibration.

Note

With the one-point calibration only the offset is adjusted. If the sensor was previously calibrated with a multipoint calibration the previously stored slope will remain. Otherwise, the theoretical slope (-59.16 mV/pH) will be used.

4.8.2 Running a multi-point pH/ion calibration (only SG98)

pH and ion calibrations can be run with this meter for up to 5 points.

- 1 Run the calibration as described in "Running a one-point pH/ion or one-point conductivity calibration" (steps a - c).
- 2 Rinse the electrode with deionized water.
- 3 Place the electrode in the next calibration buffer.
- 4 Press **CAL**.
 - ➔ **Cal 2** appears on the display. The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **READ**. The relevant buffer value is shown on the display.
- 5 Repeat the steps b - d for all calibration buffers.
- 6 Press **End** to end the calibration procedure.
 - ➔ Alternatively, the meter ends the calibration automatically when 5 calibrations are performed. The offset value and slope are shown on the display.
- 7 Press **Save** to keep the calibration.
- 8 Press **EXIT** to reject the calibration.

Note

Up to 5 calibrations can be saved for a dedicated sensor ID. The oldest calibration data is automatically overwritten with the current calibration data.

4.8.3 Automatic buffer recognition (only SG98)

The meter features automatic pH buffer recognition for the predefined buffer groups (see "Appendix"). The buffers within a buffer group are automatically recognized by the meter and displayed during calibration. This feature allows the calibration in any order within a predefined pH buffer group.

There is no automatic pH buffer recognition for customized buffer groups; in this case, the defined order of the buffers needs to be followed.

4.8.4 Dissolved oxygen (DO) calibration with the InLab® OptiOx

Under equilibrium conditions, the partial pressure of oxygen in air-saturated water is equal to the partial pressure of oxygen in water-saturated air. This means that an OptiOx sensor calibrated in water-saturated air will correctly read the partial pressure of oxygen in water samples. When measuring low concentration samples (less than 1 mg/L), a second calibration with a zero oxygen standard may be done.

Performing a one-point calibration

The first point of a DO calibration is always done in water-saturated air (100% O₂).

- 1 Remove the OptiOx calibration tube cap and remove the sponge from the cap.
- 2 Saturate the sponge with distilled water and squeeze the excess water out of the sponge.
- 3 Reassemble the OptiOx calibration tube.
- 4 Make sure that no water droplets are on the surface of the OptiOx sensor cap.
- 5 Slide the calibration tube over the front of the sensor until the calibration tube is firmly connected to the sensor.
- 6 Allow at least five minutes for the temperature to stabilize prior to calibration.
- 7 Press and hold **MODE** for 3 seconds to switch to single channel measurement screen when in dual channel measurement (only SG98).
- 8 Press **CAL**.
 - ➔ **Cal 1** appears on the display.
 - ➔ The meter endpoints according to the preselected endpoint mode, automatically after the signal has stabilized or after pressing **READ**. The standard value is shown on the display.
- 9 Press **End** to accept the calibration and return to sample measurement.
 - ➔ The calibration result is shown on the display.
- 10 Press **Exit** to reject the calibration.

Note

Under equilibrium conditions, the partial pressure of oxygen in air-saturated water is equal to the partial pressure of oxygen in water-saturated air.

Performing a two-point calibration

The second point of a DO calibration is done with a zero oxygen solution.

- 1 For the first calibration point follow the steps (a-h) as described above in "**Performing a one-point calibration**".
- 2 Remove the calibration tube.
- 3 Rinse the sensor with deionized water.
- 4 Prepare a zero oxygen solution and place the InLab® OptiOx into the bottle.
- 5 Allow at least five minutes for the sensor to equilibrate prior to calibration.
- 6 Press **CAL**.
 - ➔ **Cal 2** appears on the display.
 - ➔ The meter endpoints according to the preselected endpoint mode, automatically after the signal has stabilized or after pressing **READ**. The relevant buffer/standard value is shown on the display.
- 7 Press **End** to accept the calibration and return to sample measurement.
 - ➔ The calibration result is shown on the display.
- 8 Press **Exit** to reject the calibration.
- 9 Thoroughly rinse the sensor under running water and blot it dry with a lint-free tissue

Note

- Zero-point calibrations are frequently a source of error. Due to the very low zero current of METTLER TOLEDO sensors, a zero-point calibration is unnecessary even for measurement at low oxygen concentrations.

- If the sensor is sluggish or inaccurate after a zero point calibration, not all of the zero oxygen solution was removed from the sensor. A very thorough soaking and rinsing of the sensor in distilled water is required to remove all of the zero oxygen solution and restore the sensor performance.

4.9 Temperature compensation

We recommend the use of either a built-in or a separate temperature probe. If a temperature probe is used, ATC and the sample temperature are displayed. If no temperature sensor is used, MTC is displayed and the sample temperature should be entered manually. The meter accepts only an NTC 30 k Ω temperature sensor.

The meter calculates the temperature-adjusted electrode slope using this temperature and shows the temperature-compensated pH/ion value in the measurement display (only SG98).

DOMINIQUE DUTSCHER SAS

5 Setup

5.1 Menu structure of setup

The individual items of the menu setup are described on the pages following the list below.

1. Sample ID		4. System settings	
1. Enter Sample ID		1. Language	
2. Select Sample ID		2. Time and Date	
3. Delete Sample ID		3. Access Control	
2. User ID		4. Acoustic Signal	
1. Enter User ID		5. Routine/Expert Mode	
2. Select User ID		6. Screen Settings	
3. Delete User ID			1. Screen Contrast
3. Data Logging			2. Auto-Shutoff
1. Automatic Storage			3. Backlight Off
2. Manual Storage		5. Instrument Self-test	
3. Timed Interval Readings			

5.2 Sample ID

An alphanumeric sample ID with up to 12 characters can be **entered**. Alternatively, a previously entered sample ID can be **selected** from the list. If a sample ID has been entered, which is either purely numeric (for example, 123) or ends with a number (for example, WATER123), the following options are available:

1. <Auto Sequential> On
Using this setting will automatically increment the sample ID by 1 for each reading.
2. <Auto Sequential> Off
The sample ID is not incremented automatically.

A maximum of 5 sample IDs are stored in memory and listed for selection. If the maximum of 5 has already been entered, a sample ID can either be deleted manually or the oldest ID will be automatically overwritten by the new ID.

5.3 User ID

A user ID with up to 8 characters can be **entered**. Alternatively, a previously entered user ID can be **selected** from the list.

A maximum of 5 user IDs are stored in memory and listed for selection. If the maximum of 5 has already been entered, a user ID can either be deleted manually or the oldest ID will be automatically overwritten by the new ID.

5.4 Data logging

The meter stores up to 500 sets of measurement data in the memory. The number of data sets already stored in the memory is indicated by MXXX on the display. A message appears on the display when the memory is full. To save further measurements if the memory is full, data has to be deleted first. When measuring in dual channel mode (only SG98), both results will be stored separately. Therefore, the memory number in this case will increase by 2. You can select between automatic and manual storage or you can log your data into the memory in a user-defined interval:

1. **Automatic storage**
Stores every endpoint reading to the memory automatically.
2. **Manual storage**
If "Manual Storage" is set, **Store** appears on the display. Press **Store** to save endpoint readings.

The endpoint reading can only be stored once. When the data is stored, **Store** disappears from the measurement screen.

3. Timed interval readings

A reading is stored to memory every time after a certain interval (3 – 9999 s) defined in the menu has elapsed. When working in the timed-interval reading mode, the interval by can be defined by entering the seconds. The measurement series stops according to the selected endpoint format or manually by pressing **READ**. When timed-interval reading is "on", the **DL** icon  appears.

For readings lasting longer than 15 minutes, switch off the auto-shutoff function. The **Auto-off override** icon  appears on the display.

5.5 System settings

Note

The system settings menu is protected by a PIN. Upon delivery, the PIN is set to 000000 and is activated. Please change the PIN to prevent unauthorized access.

Language

The following languages are available for the system: English, German, French, Spanish, Italian, Portuguese, Chinese, Japanese, Korean and Russian.

Time and date

• Time

Two time display formats are available:

24-hour format (for example, 06:56 and 18:56)

12-hour format (for example, 06:56 AM and 06:56 PM)

• Date

Four date display formats are available:

28-11-2010 (day-month-year)

28-Nov-2010 (day-month-year)

28/11/2010 (day-month-year)

11-28-2010 (month-day-year)

Access control

System settings

PIN settings are available for:

- System settings
- Deleting data
- Instrument login

To enter the PIN, proceed as follows:

- 1 Switch PIN protection for the required access control ON. The window for entering an alphanumeric PIN appears.
- 2 Enter an alphanumeric PIN (max. 6 characters).
 - ➔ The input window for PIN verification appears.
- 3 Confirm PIN.

A maximum of 6 characters can be entered as PIN. In the factory default settings, the PIN for system settings and deleting data is set to 000000 and is activated, no instrument login password is set.

Acoustic signal

An acoustic signal can be switched on in the following three cases:

- Key is pressed
- Alarm/warning message appears
- Measurement is stable and has endpointed (stability signal appears)

Expert/Routine modes

The meter has two working modes:

- **Expert mode:** the factory default setting enables all functions of the meter.
- **Routine mode:** some of the menu settings are blocked.

The concept of the two working modes is a GLP feature that ensures that important settings and stored data cannot be deleted cannot be unintentionally changed under routine working conditions.

The meter only allows the following functions in the routine mode:

- Calibrating and measuring
- Editing user, sample and sensor IDs
- Editing the MTC temperature
- Editing system-settings (PIN-protected)
- Storing and viewing
- Running the instrument self-test

Screen settings

Screen contrast

The screen contrast can be set from levels 1 to 6.

Auto-shutoff

The meter will auto-shutoff when no key is pressed in a preset time to save the battery life. The time can be set (5 min, 10 min, 30 min, 1 hour, 2 hour) for the meter auto-shutoff or set to "Never" to disable this feature. If "never" is selected, the **Auto-off override** icon  appears on the display and you need to manually switch off the meter by pressing **ON/OFF**.

Backlight off

If the backlight feature is activated (**Backlight** icon  on the display), the backlight switches on with a key press and switches off again when no key is pressed for a preset time to save the battery life. The time can be set (10 s, 15 s, 30 s, 1 minute) after which backlighting automatically switches off, or set to "Never" to leave the backlighting always switched on.

- Press and hold the **Backlight** key to deactivate backlighting.
 - ➔ The **Backlight** icon  disappears from the display.

5.6 Instrument self-test

The instrument self-test requires user interaction.

- 1 In the **Setup** menu, select "6. Instrument Self-test".
 - ➔ Selecting the menu item starts the self-test routine.
- 2 Press the function keys on the keypad one by one in any order.
 - ➔ The self-test result is displayed after a few seconds.
 - ➔ The meter returns to the system settings menu automatically.

Note

- The user needs to finish pressing all seven keys within two minutes, otherwise "Self-test failed!" appears and the procedure has to be repeated.
- If error messages repeatedly appear, contact METTLER TOLEDO Service.

6 Menus and settings

6.1 Menu structure of pH/ion (only SG98)

1.	Temperature Settings		3.	Measurement Settings
	1. Set MTC Temperature			1. Measurement Resolution
	2. Temperature Unit			2. Stability Criterion
2.	Calibration Settings			3. Ion Measurement Unit
	1. Buffer Group/Standards			4. Rel. mV Offset
	2. Calibration Mode		4.	Endpoint Formats
	3. Calibration Reminder		5.	Measurement Limits
			6.	Sensor ID/SN

6.2 Menu structure of DO

1.	Temperature Settings			2. Barometric Compensation
	1. Set MTC Temperature			3. Barometric Unit
	2. Temperature Unit		4.	Endpoint Formats
2.	Calibration Reminder		5.	Measurement Limits
3.	Measurement Settings		6.	Sensor ID/SN
	1. Salinity Compensation			

6.3 Temperature settings

- **Set MTC temperature**

If the meter does not detect a temperature probe, **MTC** appears on the display. In this case the sample temperature should be entered manually. An **MTC** value between -30 °C and 130 °C can be entered.

- **Temperature unit**

Select the temperature unit: °C or °F. The temperature value is automatically converted between the two units.

6.4 pH/ion calibration settings (only SG98)

6.4.1 Buffer groups / Standards

6.4.1.1 Predefined pH buffer groups

One of seven predefined buffer groups can be selected:

B1	1.68	4.01	7.00	10.01		(at 25°C)	Mettler US
B2	2.00	4.01	9.00	9.21	11.00	(at 25°C)	Mettler Europe
B3	2.00	4.00	7.00	9.00	12.00	(at 20°C)	Standard Merck buffer
B4	1.679	4.008	6.865	9.180		(at 25°C)	JIS Z 8802
B5	1.680	4.008	6.865	9.184	12.454	(at 25°C)	DIN19266
B6	1.09	4.65	6.79	9.23	12.75	(at 25°C)	DIN19267
B7	1.680	4.003	6.864	9.182	12.460	(at 25°C)	Chinese

Temperature tables for these buffers are programmed in the meter and can be found in the "Appendix".

6.4.1.2 Customized pH buffer group

A set of user-defined pH buffers with up to 5 different temperatures for each buffer can be created. The temperature difference between pH buffers must be at least 5 °C and the difference between the pH values must be at least 1.

When switching from predefined buffer group to customized buffer group, press **Save** in the table even if no values have changed.

Ion standards

Concentrations for up to 5 standards with one standard temperature can be defined (see "pH/ion measurement settings"). Five concentration units are available:

- mmol/L
- mol/L
- ppm
- mg/L
- %

6.4.2 Calibration mode

Two calibration modes are offered:

- **Segmented:** the calibration curve is made up of linear segments joining the individual calibration points. If high accuracy is required, the segment method is recommended.
- **Linear:** the calibration curve is determined using linear regression. This method is recommended for samples with widely varying values.

Note

These settings applies to both pH and Ion calibration.

6.4.3 Calibration reminder

When the calibration reminder is "On", the user is reminded to perform a new calibration after a certain user-defined interval (maximum 9999 h) has elapsed.

- Press **READ** to save the interval and another screen appears to select calibration expiration date.

Four different time spans can be programmed. In all four cases, a warning message appears that the electrode should be calibrated.

- **Immediately**
The meter is immediately blocked for measurement when the predefined interval has elapsed.
- **Reminder + 1h**
The meter is blocked for measurement 1 hour after the predefined interval has elapsed.
- **Reminder + 2h**
The meter is blocked for measurement 2 hours after the predefined interval has elapsed.
- **Continue Reading**
The user can continue measuring when the predefined interval has elapsed.

6.5 pH/ion measurement settings (only SG98)

Measurement Resolution

The resolution for pH and mV needs to be set for the display. Up to 3 decimal places can be chosen depending on the unit of measurement (see table below).

On display	Description	Option
X.XXX	three decimal places	pH
X.XX	two decimal places	pH
X.X	one decimal place	pH, mV
X	no decimal places	mV

In the ion mode, the measurement resolution depends on the concentration and the unit of the measured ion.

Stability Criterion

The **Stability** icon appears according to the following stability criteria:

Stability criteria for pH and mV measurement

Strict

The measured signal should not change by more than 0.03 mV in 8 seconds or by more than 0.1 mV in 30 seconds.

Medium

The measured signal should not change by more than 0.1 mV in 6 seconds.

Fast

The measured signal should not change by more than 0.6 mV in 4 seconds.

Stability criterion for ion measurement

Strict

The measured signal should not change by more than 0.03 mV in 12 seconds or by more than 0.08 mV in 26 seconds.

Medium

The measured signal should not change by more than 0.08 mV in 8 seconds.

Fast

The measured signal should not change by more than 0.3 mV in 4 seconds.

Ion Measurement Units

The unit (mmol/L, mol/L, ppm, mg/L or %) for measurements and calibration can be defined.

Note

In some cases, changing units requires the user to first recalibrate before starting a measurement, otherwise an error message will appear.

The units of measurement are divided into two groups: **1.** mmol/L, mol/L and **2.** ppm, mg/L, %. Changing within a group doesn't require recalibration, but changing between the two groups does.

Rel. mV Offset

In the rel. mV mode the offset value is subtracted from the measured value. Either an offset value can be entered or it can be determined by measuring the mV of a reference sample.

1) Enter offset value

Enter an offset value in mV between -1999.9 and +1999.9 mV.

2) Test a reference sample

- 1 Place electrode into the reference sample.
- 2 Press **Start** to begin the reference measurement and wait until the measurement display freezes.
— or —
- 3 Press **READ** to manually end the measurement.
- 4 Press **Save** to enter the measured mV value as offset into the meter.

6.6 DO measurement settings

Salinity Compensation

A salinity correction is required in the oxygen concentration measurement of samples with a salt content of more than 1 ppt. The meter corrects automatically after you have entered the ion concentration in this menu (salinity range from 0.0 to 42.0 ppt = 42 g/L).

Barometric Compensation

A barometric compensation is required in the oxygen concentration measurement and calibration.

- **Auto:** the meter can automatically determine the current pressure by means of an integrated pressure sensor.
- **Manual:** the user can enter the absolute ambient atmospheric pressure.

Barometric Unit

There are four barometric units available for barometric reading:

mbar hPa Torr atm

Every unit can be calculated from the standard unit Pa using the following formulas: mbar = hPa = 100 Pa
Torr = 133.322 Pa atm = 101325 Pa

6.7 DO calibration reminder

For details on the **DO calibration reminder**, see "pH/ion calibration settings".

6.8 Endpoint formats

Auto

With the automatic endpoint the selected stability criterion determines the end of an individual reading depending on the behavior of the sensor used. This ensures an easy, quick and precise measurement.

- 1 Place sensor in the sample.
- 2 Press **READ**.
 - ➔ **A** appears on the display.
 - ➔ The measurement ends automatically when the measured value is stable. \overline{A} appears.
 - ➔ If **READ** is pressed before the signal is stable, the endpoint format changes to manual \overline{M} .

Manual

Unlike **Auto**, user interaction is required to stop the measurement reading in manual mode.

- 1 Place sensor in the sample.
- 2 Press **READ**.
 - ➔ **M** appears on the display.
 - ➔ \overline{M} appears on the display to signalize measurement stability.
- 3 Press **READ** to end the measurement. \overline{M} appears.

Timed

The measurement stops after the set time, which can be set between 5 s and 3600 s.

- 1 Place sensor in the sample.
- 2 Press **READ**.
 - ➔ **T** appears on the display.
 - ➔ \overline{T} appears on the display to signalize measurement stability.
 - ➔ The measurement ends automatically when the set time period expires. \overline{T} appears.
 - ➔ If **READ** is pressed before the signal is stable, the endpoint format changes to manual \overline{M} .

Information on the display

The following symbols appear in the display, depending on the endpoint setting.

Preselected format	Start of measurement	Signal stability	Endpointed measurement ¹
Auto endpoint	A	⌘	⌘
	A	Read ⇨	⌘
Manual endpoint	M	⌘ Read ⇨	⌘
	M	Read ⇨	⌘
Timed endpoint	T	⌘ ⌚ ⇨	⌘
	T	Read ⇨	⌘

¹The actual endpoint format (last column) and not the preselected is stored with the data.

6.9 Measurement limits

The upper and lower limits for measurement data can be defined. If a limit is either not reached or exceeded (in other words, less than or greater than a specific value), a warning is displayed on the screen and may be accompanied by an acoustic signal. The message "outside limits" also appears on the GLP printout.

6.10 Sensor ID/SN

Enter Sensor ID/SN

An alphanumeric sensor ID with up to 12 characters can be entered. The sensor ID will be assigned to each calibration and measurement value. This is valuable for tracing back data.

Up to 5 sensor IDs can be entered for each sensor type.

If a new sensor ID is entered, the theoretical calibration slope and offset for this type of electrode will be loaded. The sensor has to be newly calibrated.

If a sensor ID is entered, which is already in the memory of the meter and has been calibrated before, the specific calibration data for this sensor ID will be loaded.

When a new **ion sensor ID** is entered, the electrode type can be selected.

When connecting an **ISM® sensor** to the meter, the meter will:

- Automatically recognize the sensor when it's turned on (alternatively, when pressing **READ** or **CAL**)
- Load the stored sensor ID, sensor SN and sensor type as well as the latest calibration data of this sensor
- Use this calibration for the subsequent measurements

The sensor ID for analog ISM® sensors can be changed. However, this is not possible for digital ISM® sensors.

Select Sensor ID

Already entered sensor IDs can be selected from a list.

If a sensor ID is selected, which is already in the memory of the meter and has been calibrated before, the specific calibration data for this sensor ID will be loaded.

Note

- You can delete a sensor ID with its calibrations in the calibration data menu.

7 Data management

7.1 Menu structure of data menu

7.1.1 SG98

1.	Measurement Data		3.	ISM Data
	1. Review			1. pH
	2. Delete			1. Initial Calibration Data
2.	Calibration Data			2. Calibration History
	1. pH			3. Maximum Temperature
	1. Review			4. Reset ISM
	2. Delete			2. DO
	2. Ion			1. Calibration History
	1. Review			2. Maximum Temperature
	2. Delete			3. Sensor Cap Life
	3. DO			4. Reset ISM
	1. Review			
	2. Delete			

7.1.2 SG9

1.	Measurement Data		3.	ISM Data
	1. Review			1. Calibration History
	2. Delete			2. Maximum Temperature
2.	Calibration Data			3. Sensor Cap Life
	1. Review			4. Reset ISM
	2. Delete			

7.2 Measurement data

Review

All

All stored measurement data can be reviewed; the most recent data saved appears on the display.

Partial

The measurement data can be filtered according to 3 criteria.

- Memory number (MXXX)
- Sample ID
- Measurement mode

Memory number

- Enter the memory number of the data and press **Review**.
 - ➔ The measurement data is displayed.

Sample ID

- 1 Enter the sample ID and press **Review**.
 - ➔ The meter finds all stored measurements with this sample ID.
- 2 Scroll through the measurement data to review all measurements with the entered sample ID.

Measurement mode

- 1 Select a measurement mode from list and press **Review**. The meter finds all stored measurements of the selected measurement mode.
- 2 Scroll through the measurement data of the selected measurement mode.

Delete

All or partially stored measurement data can be deleted by filtering the measurement data. The filter works as described above in "Review".

Note

- Deletion is protected by a PIN. Upon delivery, the PIN is set to 000000. Change the PIN code to prevent unauthorized access.

7.3 Calibration data

Calibration data can be reviewed and deleted. Up to 5 calibrations per sensor ID are stored in the memory.

Review

- 1 Select between the sensor types: pH, ion or DO (only SG98).
- 2 Press **Review**.
 - ➔ A list of calibrated sensor IDs appears.
- 3 Select a sensor ID from the list and press **Review**.
- 4 Press **↑** and **↓** to navigate between the previous or next calibration data sets.
 - or —
 - Press and hold **CAL** for 3 seconds in the single channel measurement screen.
 - ➔ The current calibration data is displayed.

Delete

- 1 Select between the sensor types: pH, ion or DO (only SG98).
- 2 Press **Delete**.
 - ➔ A list of sensor IDs appears.
- 3 Select a sensor ID from the list and press **Delete**.
- 4 Press **Yes** when the message "This will delete all your selected data. Please confirm." appears
 - or —
 - Press **No** to cancel and exit.
 - ➔ After deletion, the sensor ID disappears from the list in the sensor ID menu.

Note

- An active sensor ID cannot be deleted.
- This menu is protected by a deletion PIN code. Upon delivery, the PIN code is set to 000000. Change the PIN code to prevent unauthorized access.

7.4 ISM data

SevenGo Duo pro™ SG98 and SevenGo pro™ SG9 incorporate Intelligent Sensor Management (ISM®) technology. This ingenious functionality provides extra security, safety and eliminates mistakes. The most important features are:

Extra security!

- After connecting the ISM® sensor, the sensor is automatically recognized and the sensor ID and serial number are transferred from the sensor chip to the meter.
- After calibration of the ISM® sensor, the calibration data is automatically stored from the meter to the sensor chip. The most recent data is always stored where it should be – on the sensor chip!

Extra safety!

After connecting the ISM® sensor, the five most recent calibrations are transferred to the meter. These can be reviewed to see the development of the sensor over time. This information provides an indication if the sensor should be cleaned or renewed.

Eliminate mistakes!

After connecting an ISM® sensor, the last set of calibration data is automatically used for measurements.

Additional features are described below.

Initial calibration data for pH (only SG98)

When an ISM® sensor is connected, the initial calibration data in the sensor can be reviewed or transferred. The following data is included:

- Response time
- Temperature tolerance
- Membrane resistance
- Slope (at pH 7) and offset
- Type (and name) of electrode (for example, InLab® Expert Pro ISM)
- Serial number (SN) and ordering (ME) number
- Production date

Calibration history

The last 5 calibrations data stored in ISM® sensor including current calibration can be reviewed.

Maximum temperature

The maximum temperature that the ISM® sensor has been exposed to during measurement is monitored automatically and can be reviewed for the evaluation of the electrode lifetime.

Sensor cap life

The date of first use of the optical DO sensor cap and the expiration date in UTC (universal time coordinated) can be reviewed. Once the first measurement is taken with the sensor cap, an internal clock starts counting down and the cap needs to be replaced after 12 months.

Reset ISM®

The calibration history in this menu can be deleted. This menu is protected by a deletion PIN. Upon delivery, the PIN for deletion is set to 000000. Change the PIN to prevent unauthorized access.

8 Maintenance

8.1 Meter maintenance

Never unscrew the two halves of the housing!

The meters do not require any maintenance other than an occasional wipe with a damp cloth. The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is sensitive to some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK).

Any spillage should be wiped off immediately.

8.2 pH electrode maintenance

Make sure the pH electrode is always kept filled with the appropriate filling solution.

For maximum accuracy, any filling solution that may have "crept" and encrusted the outside of the electrode should be removed with deionized water.

Always store the electrode according to the manufacturer's instructions and do not allow it to dry out.

If the electrode slope falls rapidly, or if the response becomes sluggish, the following procedures may help. Try one of the following, depending on your sample.

Problem	Action
Fat or oil build-up	Degrease the membrane with cotton wool soaked in either acetone or a soap solution.
pH sensor membrane has dried out	Soak the tip of the electrode overnight in 0.1M HCl
Protein build-up in the diaphragm of a pH sensor	Remove deposits by soaking the electrode in an HCl/pepsin solution.
Silver sulfide contamination of pH sensor	Remove deposits by soaking electrode in a thiourea solution.
Run a new calibration after treatment.	

Note

Cleaning and filling solutions should be handled with the same care as that given to toxic or corrosive substances.

8.3 InLab® OptiOx sensor maintenance

Storage

- 1 Do not remove the sensor cap.
- 2 The sensor can be stored dry. It is recommended to keep it in the white calibration tube to protect the sensor cap from mechanical damage.
- 3 Keep the sensor away from direct sunlight during storage.

Cleaning the sensor cap and body

- 1 Do not remove the sensor cap.
- 2 Rinse the sensor with distilled water.
- 3 Gently wipe with a soft-bristled brush or soft cloth if biofouling is present.
- 4 If extensive mineral build-up is present, soak the cap end in vinegar for 15 minutes.
- 5 Soak the sensor in deionized water for 15 minutes and blot it dry with a lint-free tissue.
- 6 After cleaning the sensor, 1-point calibration should be performed to check it.

Note

- Do not use organic solvents or soaps to clean the cap!
- Cleaning of the inner lens should only be performed when changing the cap.

8.4 Troubleshooting InLab® OptiOx

Issue	Recommended Action
Unable to calibrate	Verify the calibration setup and procedure. Make sure that no water droplets are on the surface of the cap. Verify that the cap lifetime has not expired yet.
Measurements are unstable	Measurements may take longer if the solution temperature is unstable.
Measurement is too low	Salt may be present in the sample. Set the salinity factor in the meter.
Wrong temperature displayed	Verify that the temperature sensor (metallic pin along the sensor shaft) is immersed in the solution.

- 1 Rinse the sensor thoroughly with distilled water, blot it dry with a lint-free tissue and examine the cap for scratches or discoloration.
- 2 Remove the cap from the sensor and make sure that there is no water inside the cap, the optical window is clean and clear, the O-rings are intact and have a thin coating of silicone grease and the spring contacts are clean and undamaged.
- 3 If readings continue to be erratic and unstable, a replacement of the cap or the whole sensor may be necessary.

8.5 Replacement of the OptiOx sensor cap


The sensor cap has a lifespan of 1 year once the first measurement is taken. The meter will display a message "sensor cap expired" when the cap needs to be replaced.

Note

- The OptiOx™ sensor contains an internal clock which counts down the 365 day life span of a new sensor cap. The countdown begins once the OptiOx™ cap has been attached, the sensor connected to the device and the first measurement or calibration carried out. This process cannot be undone once the first measurement has been carried out.
- 1 Pull the expired sensor cap off the sensor without twisting. See fig. A.
 - 2 Remove the existing O-rings from the sensor. See fig. B.
 - 3 Use a lint-free cloth to remove any moisture from the sensor body. Make sure there is no moisture in the O-ring grooves.
 - 4 Use your finger to apply a layer of lubricant around the O-ring grooves. Place the new O-rings on the sensor (included in the delivery of a new sensor cap). Apply another thin layer of lubricant to the O-rings and grooves. Do not transfer grease to the lens or sensor pins.
 - 5 Gently wipe the lens with a clean cloth and allow it to dry thoroughly. Do not wet the lens area with water or any other solution. Examine the lens for scratches or dirt.
 - 6 Install a new OptiOx sensor cap onto the optical sensor by aligning the arrow on the cap with the index mark on the sensor. See fig. C.
Without twisting, firmly press the cap onto the sensor until the cap is flushed with the sensor body. Make sure that the O-rings are not pinched or rolled between the cap and sensor. See fig. D.
 - 7 After replacing the OptiOx sensor cap, a vapor-saturated air calibration should be performed.

Note

- Do not remove the OptiOx cap after installation until a cap replacement is necessary.

8.6 Interfering substances of the InLab® OptiOx sensor

⚠ CAUTION

Damage of the sensor!

Do not use the InLab® OptiOx sensor in solutions that contain organic solvents, such as acetone, chloroform or methylene chloride.

Following substances may interfere with the dissolved oxygen measurement:

- Alcohols greater than 5%,
- Hydrogen peroxide (H₂O₂) greater than 3%,

- Sodium hypochlorite (NaClO_2) greater than 3%,
- Gaseous sulphur dioxide (SO_2) and gaseous chlorine (Cl_2).

Carbone dioxide (CO_2), ammonia (NH_3), pH, any ionic species like sulphide (S^{2-}), sulphate (SO_4^{2-}), chloride (Cl^-) and hydrogen sulphide (HS^-) do not interfere with the dissolved oxygen measurements.

DOMINIQUE DUTSCHER SAS

9 Error messages

Message	Description and Resolution
pH/mV/ion/temperature/dissolved oxygen exceeds max. limit	Measurement limits are activated in the menu settings and measured value is outside these limits.
pH/mV/ion/temperature/dissolved oxygen below min. limit	<ul style="list-style-type: none"> • Check the sample. • Check sample temperature. • Make sure that the pH electrode wetting cap has been removed and that the electrode is properly connected and placed in the sample solution.
Memory is full	<p>Max. 500 measurement data can be stored in the memory.</p> <ul style="list-style-type: none"> • Delete all or partial data in the memory, otherwise you will not be able to store new measurement data.
Please calibrate electrode	<p>Calibration reminder has been switched on in the menu settings and last calibration has expired.</p> <ul style="list-style-type: none"> • Calibrate the electrode.
Active sensor cannot be deleted	<p>Deleting the calibration data of the selected sensor ID is not possible, because it is currently the active sensor ID in the meter shown on the display.</p> <ul style="list-style-type: none"> • Enter new sensor ID in the menu settings. • Select another sensor ID from the list in the menu settings.
Wrong buffer	<p>Meter cannot recognize the buffer or standard/buffer has been used twice for calibration/two buffers differ less than 60 mV.</p> <ul style="list-style-type: none"> • Make sure that you have the correct buffer and that it is fresh. • Make sure that the buffer has not been used more than once during the calibration.
Slope out of range	<p>The calibration result is outside the following limits: Slope < 85% or > 105%, Offset < -35 mV or > + 35 mV.</p> <ul style="list-style-type: none"> • Make sure that you have the correct buffer and that it is fresh. • Check mV signal of electrode, clean or replace the electrode,
Offset out of range	
Buffer temp. out of range	<p>The ATC measured temperature is out of pH calibration buffer range: 5 ... 50 °C.</p> <ul style="list-style-type: none"> • Keep the buffer/standard temperature within the range. • Change the temperature setting.
Standard temp. out of range	
Temperature differs from setting	<p>ATC measured temperature differs by more than 0.5°C from the user-defined value/temperature range.</p> <ul style="list-style-type: none"> • Keep the buffer/standard temperature within the range. • Change the temperature setting.
ISM [®] sensor communication error	Data has not been transferred correctly between ISM [®] sensor and meter. Reconnect the ISM [®] sensor and try again.
Sensor cap expired	The lifetime of the optical DO sensor cap of 1 year has expired. Change the cap following the instructions of the optical DO sensor maintenance.
Sensor failure	<p>The cap of the optical DO sensor is missing or not installed properly</p> <ul style="list-style-type: none"> - Install a new DO cap onto the optical DO sensor following the instructions of the optical DO sensor maintenance - Remove the cap and reconnect it following the instructions of the optical DO sensor maintenance <p>There is no signal from the sensor. Switch off the instrument, reconnect the sensor and try again.</p>

Message	Description and Resolution
Over range	The measured oxygen value is out of the calibrated range. <ul style="list-style-type: none"> - Make sure that the sensor cap and the temperature sensor are covered with solution - Make sure that the sensor is free of salt or mineral contamination
Under range	The measured oxygen value is out of the calibrated range. <ul style="list-style-type: none"> - Make sure that the sensor cap and the temperature sensor are covered with solution - Make sure that the sensor is free of salt or mineral contamination
Self-test failure	Self-test has not been completed within 2 minutes or meter is defective. <ul style="list-style-type: none"> • Restart self-test and finish within 2 minutes. • Contact METTLER TOLEDO service if the problem persists.
Wrong settings	Entered value differs by less than 1 pH unit/5°C from other preset values. <ul style="list-style-type: none"> • Enter a higher/lower value in order to get a larger difference.
Out of range	Either entered value is out of range. <ul style="list-style-type: none"> • Enter a value, which is within the range shown on the display. or Measured value out of range. <ul style="list-style-type: none"> • Make sure the electrode wetting cap has been removed and that the electrode is properly connected and placed in the sample solution. • If no electrode is connected, put the shorting clip in the socket.
Wrong password	The entered PIN is not correct. <ul style="list-style-type: none"> • Re-enter the PIN. • Reset to factory settings, all data and settings will be lost.
Passwords do not match, try again	The confirmation PIN does not match with the entered PIN. <ul style="list-style-type: none"> • Reenter PIN.
Program memory error	Meter recognizes internal error during start-up. <ul style="list-style-type: none"> • Switch the meter off and back on. • Contact METTLER TOLEDO service if problem persists.
Data memory error	The data could not be stored in the memory. <ul style="list-style-type: none"> • Switch the meter off and back on. • Contact METTLER TOLEDO service if problem persists.
No matching data found in memory	The entered filter criterion does not exist. <ul style="list-style-type: none"> • Enter a new filter criterion.
Sensor ID already exists, previous SN will be overwritten	Two sensors with the same ID but different SN are not allowed in the meter. If a different SN has been entered for this sensor ID previously, the old SN will be overwritten. <ul style="list-style-type: none"> • Enter a different Sensor ID in order to keep the previous ID and SN.

9.1 Error limits

Message	Range not accepted	
Out of range	pH	<-2.000 or > 19.999
	mV	<-1999.9 or > 1999.9
	[O ₂]	< 0.1% or > 600 %
	[O ₂]	< 0.01 mg/L or > 80 mg/L
	[O ₂]	< 0.01 ppm or > 80 ppm
Buffer/standard temp. out of range	T (pH)	< 5 or > 50 °C
	T (DO)	< 0 °C or > 50 °C
Offset out of range	Eref1-Eb > 60 mV	
Slope out of range	Eref1-Eb > 60 mV	
Wrong buffer	ΔEref1 < 10 mV	
Invalid pH for user-defined buffer	ΔpH < 1 pH	
ATC measured temperature is different to the user-defined value	tATC-tbuffer > 1 °C	
Cal. 1 out of range	I	< 40 nA or > 110 nA
Cal. 2 out of range	T	< 0 nA or > 2 nA

10 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

DOMINIQUE DUTSCHER SAS

11 Sensors, solutions and accessories

Parts	Order No.
IP67 sensors with fixed cable	
InLab®Expert Go, rugged 3-in-1 pH sensor, IP67, PEEKshaft, ATC	51340288
Parts	Order No.
ISM®IP67 sensors with fixed cable	
InLab®Expert Go-ISM, rugged 3-in-1 pH sensor, IP67, PEEK shaft, ATC	51344102
InLab®Expert Go ISM-5m, rugged 3-in-1 pH sensor, PEEK shaft, ATC	51344103
InLab®Expert Go ISM-10m, rugged 3-in-1 pH sensor, IP67, PEEKshaft, ATC	51344104
InLab® OptiOx, optical DO sensor for use with SG9 and SG98 (1.8 m cable)	51344621
InLab® OptiOx-5m, optical DO sensor for use with SG9 and SG98 (5 m cable)	51344622
InLab® OptiOx-10m, optical DO sensor for use with SG9 and SG98 (10 m cable)	51344623
Parts	Order No.
ISM®IP67 sensors with multi-pin head	
InLab®Micro Pro ISM, 3-in-1 pH sensor, glass shaft, 5 mm shaft diameter, ATC, refillable	51344163
InLab®738-ISM, conductivity sensor, epoxy shaft, ATC, pressurized SteadyForce™ reference system	51344112
InLab®Pure Pro ISM, 3-in-1 pH sensor, glass shaft, immovable glass sleeve, ATC, refillable	51344172
InLab®Routine Pro ISM, 3-in-1 pH sensor, glass shaft, ATC, refillable	51344055
InLab®Science Pro ISM, 3-in-1 pH sensor, glass shaft, movable glass sleeve, ATC, refillable	51344072
InLab®Solids Pro ISM, 3-in-1 pH sensor, glass shaft, open junction, sharp membrane, ATC	51344155
Parts	Order No.
Solutions	
pH 4.01 buffer sachets, 30 x 20mL	51302069
pH 4.01 buffer solution, 6 x 250mL	51350018
pH 7.00 buffer sachets, 30 x 20mL	51302047
pH 7.00 buffer solution, 6 x 250mL	51350020
pH 9.21 buffer sachets, 30 x 20mL	51302070
pH 9.21 buffer solution, 6 x 250mL	51350022
pH 10.01 buffer sachets, 30 x 20mL	51302079
pH 10.00 buffer solution, 6 x 250mL	51350024
HCl/Pepsin solution (removes protein contamination)	51350100
Thiourea solution (remove silver sulfide contamination)	51350102
Zero oxygen standard tablets, 20 pieces	51300140
Parts	Order No.
Accessories	
Battery cover	51302328

Parts	Order No.
Bottles 50 mL	51300240
Bottom cap (blue)	51302324
Clip cover	51302327
Electrode weight	51303019
Replacement Cap OptiOx	51344630
Calibration Tube OptiOx	51344631
Protective Guard OptiOx (stainless steel)	51344632
BOD Adapter OptiOx	51344633
Rubber feet (2 pcs.)	51302335
SevenGo™ clip	51302325
SevenGo™ sealing kit	51302336
Wrist strap	51302331

12 Specifications

SevenGo Duo pro™ pH/ORP/Ion/DO meter SG98

Measurement range	pH	-2.000 ... 19.999
	mV	-1999.9 ... 1999.9 mV
	pH ATC	-5 ... 130°C
	pH MTC	-30 ... 130°C
	Ion	0.000 ... 999.9% 0.000...9999 ppm 1.00E-9...9.99E+9 mg/L 1.00E-9...9.99E+9 mmol/L
	DO	0.00 ... 50.00 mg/L 0.00 ... 50.00 ppm 0.0 ... 500.0 %
	Pressure	1100 mbar
	Temperature	0 ... 50 °C
	Resolution	0.1 / 0.01 / 0.001pH
1 / 0.1 mV		
pH Temperature		0.1°C
Ion		3 digits
DO		0.01 mg/L 0.01 ppm 0.1 %
Pressure		1 mbar
DO temperature		0.1 °C
Limits of error pH		± 0.002 pH
	± 0.2 mV	
	± 0.1°C	
Limits of error ion	± 0.5% (this limit only applies for meter)	
Limits of error DO	DO	±0.5 % of measured value
	Pressure	± 2 % of measured value
	Temperature	±0.1 °C
pH calibration	Up to 5 points	
Isopotential point	pH 7.00	
pH Calibration buffer	7 predefined groups	1 user-defined group of 5 buffers
DO calibration standard	Up to 2 points	Vapor-saturated air Zero oxygen solution
Power requirements	Ratings	6 V DC, 70 mA
	Batteries	4 x AA / LR6 1.5 V or NiMH 1.2 V rechargeable
Size / weight	220 x 90x 45 mm 368 g	
Display	Liquid crystal	

pH input	BNC (IP67), impedance > $3 \cdot 10^{12} \Omega$	
DO and DO temperature input	mini LTW(IP67), NTC 30k Ω	
pH T input	Cinch(IP67), NTC 30k Ω	
IP rating	IP67 with and without electrode	
Ambient conditions	Temperature	5 ... 40°C
	Relative humidity	5%...80% (non-condensing)
	Installation category	II
	Pollution degree	2
	Altitude	Up to 2000 m above sea level
Materials	Housing	ABS/PC reinforced
	Window	Polymethyl methacrylate (PMMA)
	Keypad	silicone rubber

SevenGo pro™ DO meter SG9

Measurement range	DO	0.00 ... 50.00 mg/L 0.00 ... 50.00 ppm 0.0 ... 500.0 %
	Pressure	1100 mbar
	Temperature	0 ... 50 °C
Resolution	DO	0.01 mg/L 0.01 ppm 0.1 %
	Pressure	1 mbar
	DO temperature	0.1 °C
Limits of error DO	DO	±0.5 % of measured value
	Pressure	± 2 % of measured value
	Temperature	±0.1 °C
DO calibration standard	Up to 2 points	Vapor-saturated air Zero oxygen solution
Power requirements	Ratings	6 V DC, 70 mA
	Batteries	4 x AA / LR6 1.5 V or NiMH 1.2 V rechargeable
Size / weight	220 x 90x 45 mm 368 g	
Display	Liquid crystal	
DO and DO temperature input	mini LTW(IP67), NTC 30k Ω	
IP rating	IP67 with and without electrode	

Ambient conditions	Temperature	5 ... 40°C
	Relative humidity	5% ... 80% (non-condensing)
	Installation category	II
	Pollution degree	2
	Altitude	Up to 2000 m above sea level
Materials	Housing	ABS/PC reinforced
	Window	Polymethyl methacrylate (PMMA)
	Keypad	silicone rubber

DOMINIQUE DUTSCHER SAS

13 Appendix

13.1 Buffer tables

METTLER TOLEDO USA (Ref. 25°C)

Temperature °C	1.68	4.01	7.00	10.01
5	1.67	4.00	7.09	10.25
10	1.67	4.00	7.06	10.18
15	1.67	4.00	7.04	10.12
20	1.68	4.00	7.02	10.06
25	1.68	4.01	7.00	10.01
30	1.68	4.01	6.99	9.97
35	1.69	4.02	6.98	9.93
40	1.69	4.03	6.97	9.89
45	1.70	4.04	6.97	9.86
50	1.71	4.06	6.97	9.83

METTLER TOLEDO Europe (Ref. 25°C)

Temperature °C	2.00	4.01	7.00	9.21	11.00
5	2.02	4.01	7.09	9.45	11.72
10	2.01	4.00	7.06	9.38	11.54
15	2.00	4.00	7.04	9.32	11.36
20	2.00	4.00	7.02	9.26	11.18
25	2.00	4.01	7.00	9.21	11.00
30	1.99	4.01	6.99	9.16	10.82
35	1.99	4.02	6.98	9.11	10.64
40	1.98	4.03	6.97	9.06	10.46
45	1.98	4.04	6.97	9.03	10.28
50	1.98	4.06	6.97	8.99	10.10

MERCK (Ref. 20°C)

Temperature °C	2.00	4.00	7.00	9.00	12.00
5	2.01	4.04	7.07	9.16	12.41
10	2.01	4.02	7.05	9.11	12.26
15	2.00	4.01	7.02	9.05	12.10
20	2.00	4.00	7.00	9.00	12.00
25	2.00	4.01	6.98	8.95	11.88
30	2.00	4.01	6.98	8.91	11.72
35	2.00	4.01	6.96	8.88	11.67
40	2.00	4.01	6.95	8.85	11.54
45	2.00	4.01	6.95	8.82	11.44
50	2.00	4.00	6.95	8.79	11.33

DIN(19266)/NIST (Ref. 25°C)

Temperature °C	1.68	4.008	6.865	9.183	12.454
5	1.668	4.004	6.950	9.392	13.207
10	1.670	4.001	6.922	9.331	13.003
15	1.672	4.001	6.900	9.277	12.810
20	1.676	4.003	6.880	9.228	12.627
25	1.680	4.008	6.865	9.183	12.454
30	1.685	4.015	6.853	9.144	12.289
35	1.691	4.026	6.845	9.110	12.133
40	1.697	4.036	6.837	9.076	11.984
45	1.704	4.049	6.834	9.046	11.841
50	1.712	4.064	6.833	9.018	11.705

DIN(19267) (Ref. 25°C)

Temperature °C	1.09	4.65	6.79	9.23	12.75
5	1.08	4.67	6.87	9.43	13.63
10	1.09	4.66	6.84	9.37	13.37
15	1.09	4.66	6.82	9.32	13.16
20	1.09	4.65	6.80	9.27	12.96
25	1.09	4.65	6.79	9.23	12.75
30	1.10	4.65	6.78	9.18	12.61
35	1.10	4.65	6.77	9.13	12.45
40	1.10	4.66	6.76	9.09	12.29
45	1.10	4.67	6.76	9.04	12.09
50	1.11	4.68	6.76	9.00	11.98

JJG119 (Ref. 25°C)

Temperature °C	1.680	4.003	6.864	9.182	12.460
5	1.669	3.999	6.949	9.391	13.210
10	1.671	3.996	6.921	9.330	13.011
15	1.673	3.996	6.898	9.276	12.820
20	1.676	3.998	6.879	9.226	12.637
25	1.680	4.003	6.864	9.182	12.460
30	1.684	4.010	6.852	9.142	12.292
35	1.688	4.019	6.844	9.105	12.130
40	1.694	4.029	6.838	9.072	11.975
45	1.700	4.042	6.834	9.042	11.828
50	1.706	4.055	6.833	9.015	11.697

JIS Z 8802 (Ref. 25°C)

Temperature °C	1.679	4.008	6.865	9.180
5	1.668	3.999	6.951	9.395
10	1.670	3.998	6.923	9.332

Temperature °C	1.679	4.008	6.865	9.180
15	1.672	3.999	6.900	9.276
20	1.675	4.002	6.881	9.225
25	1.679	4.008	6.865	9.180
30	1.683	4.015	6.853	9.139
35	1.688	4.024	6.844	9.102
40	1.694	4.035	6.838	9.068
45	1.700	4.047	6.834	9.038
50	1.707	4.060	6.833	9.011

13.2 Solubility of oxygen in water as a function of temperature and salinity

According to EN 25 814 and UNESCO tables (partly extrapolated)

Temperature (°C)	O ₂ solubility (mg/L)	Salinity correction factor F(T) (mg/L)
0	14.62	0.0875
1	14.22	0.0843
2	13.83	0.0818
3	13.46	0.0789
4	13.11	0.0760
5	12.77	0.0739
6	12.45	0.0714
7	12.14	0.0693
8	11.84	0.0671
9	11.56	0.0650
10	11.29	0.0632
11	11.03	0.0614
12	10.78	0.0593
13	10.54	0.0582
14	10.31	0.0561
15	10.08	0.0545
16	9.87	0.0532
17	9.66	0.0514
18	9.47	0.0500
19	9.28	0.0489
20	9.09	0.0475
21	8.91	0.0464
22	8.74	0.0453
23	8.58	0.0443
24	8.42	0.0432
25	8.26	0.0421
26	8.11	0.0407

27	7.97	0.0400
28	7.83	0.0389
29	7.69	0.0382
30	7.56	0.0371
31	7.43	0.0365
32	7.30	0.0353
33	7.18	0.0345
34	7.06	0.0339
35	6.95	0.0331
36	6.83	0.0323
37	6.72	0.0316
38	6.61	0.0309
39	6.51	0.0302
40	6.41	0.0296
41	6.32	0.0289
42	6.23	0.0283
43	6.14	0.0277
44	6.05	0.0272
45	5.96	0.0266
46	5.88	0.0261
47	5.79	0.0256
48	5.71	0.0251
49	5.63	0.0247
50	5.55	0.0242
51	5.47	0.0238
52	5.39	0.0234
53	5.31	0.0231
54	5.24	0.0228
55	5.16	0.0225
56	5.08	0.0222
57	5.00	0.0220
58	4.91	0.0218
59	4.83	0.0216
60	4.74	0.0215

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Dual and single channel portable meter
SG68 (SGx8-, SG2x and SGx - Series)
For additional types, see page type code

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
2006/95/EC Low voltage (LVD)
2004/108/EC Electromagnetic compatibility (EMC)

and that the standards have been applied.
und die Normen zur Anwendung gelangten.

Last two digits of the year in which the CE marking was affixed: **08**

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Chris Radloff
General Manager

Rolf Truttmann
Head SBU pH Lab

References of standards for this declaration of conformity, or parts thereof:
Harmonized standards of Europe and Switzerland:

Safety standards:

IEC/EN61010-1:2001

EMC standards:

EN61326-1:2006 (class B)
EN61326-1:2006 (Basic requirements)

Metrological standards:

--

IP standards:

--

Standards for Canada, USA and Australia:

CAN/CSA-C22.2 No. 61010-1-04
UL Std. No. 61010-1 (2nd Edition)

FCC, Part 15, class A (Declaration)
AS/NZS CISPR 11, AS/NZS 61000.4.3

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SG68	Tested type (pH / ion / dissolved oxygen)
SG78	Tested type (pH / ion / conductivity)
SG98	Tested type (pH / ion / optical dissolved oxygen)
SGx8	Series name SevenGo Duo pro
x	Function (software, outputs, display, ...)
SG2x	Series name SevenGo Duo
x	Function (software, outputs, display, ...) without backlight
SGx	Series name SevenGo pro (one channel)
x	Function (software, outputs, display, ...)
SGx	Series name SevenGo (one channel)
x	Function (software, outputs, display, ...) without backlight

Where x in the model designation may be up to 1 digit can be any number 2, 3, 6, 7, 8 or 9 denoting SELV/ELVEL secondary circuits or minor mechanical differences.

Remarks

Bemerkungen:

released: 9.12.2004 / ATL

Seite 2 von 2

printed: Okt. 2010

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

www.mt.com/phlab

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
Tel. +41 22 567 53 22
Fax +41 22 567 53 23
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 02/2022
51710917B en, fr, es

51710917