INTEGRA

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
CE
Swiss made


INTEGRA

This quick guide is intended to provide a brief overview of the main features and operation of PIPETBOY acu 2. For detailed information, please refer to the Operating Instructions of PIPETBOY acu 2 that can be found at www.integra-biosciences.com in several languages.

Battery Charge


Plug in

Charge the battery completely before first use.


The supply voltage must meet the requirements of the mains adapter: 100-240 VAC, 50/60 Hz.

Low Battery

Red Blink


Charging Battery

Red




~ 3.5 h

Charged Battery

Green

General Speed Setting

Plug in for Turbo Speed

50 %
Speed


100 %
Speed

125 %
Turbo Speed


INTEGRA

Manual Speed Control

A
S
P
I
R
A
T
E

Slow


Fast

Gravity Dispense


Slow

Fast


D
I
S
P
E
N
S
E

Pump not running

Pump running



Pump running

Filter Exchange

Press (1) and turn 1/8th counter clockwise (2)

Replace the filter

Colored side facing upwards

Intended Use

PIPETBOY acu 2 is a pipette controller designed for aspirating and dispensing aqueous solutions with plastic or glass pipettes in a volume range of 1 to 100 ml. It is intended for measurement, control and laboratory use.

Safety Notes

- 1) Do not use or charge PIPETBOY acu 2 in an atmosphere with danger of explosion. Do not pipette highly flammable liquids such as acetone or ether.
- 2) Comply with the material safety data sheet (MSDS) and with all safety guidelines when handling dangerous substances. Use protective clothing and safety goggles. Never point a pipette in anyone's direction.
- 3) Avoid pipetting of liquids whose vapours could attack the materials PA (polyamide), POM (polyoxymethylene), FPM (fluor-rubber), NBR (nitrile-rubber), CR (chloroprene) and silicone. Corrosive vapours could also damage metallic parts inside the device.
- 4) Use an original INTEGRA Biosciences mains adapter only. Protect it from moisture otherwise PIPETBOY acu 2 might be damaged.

Regardless of the listed safety notes, additionally applicable regulations and guidelines of trade associations, health authorities, trade supervisory offices, etc. must be observed.

Warranty

PIPETBOY pipette controllers are designed and manufactured following high standards in quality and durability. They are warranted for a period of 12 months.

To benefit from a 3 year extended warranty, register your PIPETBOY on www.pipetboy.info or www.integra-biosciences.com and click on Warranty registration. Your warranty extension will be activated after online submission of the registration form.

Equipment Disposal

PIPETBOY acu 2 device must not be disposed of with unsorted municipal waste. Dispose of PIPETBOY acu 2 in accordance with the laws and regulations in your area governing disposal of devices.

Manufacturer

Declaration of conformity | Konformitätserklärung |
Déclaration de conformité | Declaración de conformidad |
Dichiarazione di conformità

INTEGRA Biosciences AG – 7205 Zizers, Switzerland

declares on its own responsibility that the product | erklärt in alleiniger
Verantwortung, dass das Produkt | déclare sous sa responsabilité exclusive,
que le produit | declara bajo su propia responsabilidad que el producto |
dichiara sotto la propria responsabilità che il prodotto

PIPETBOY acu 2 Models: 155000, 155015, 155 016, 155 017, 155018, 155 019

in accordance with EC directives | gemäss der EU-Richtlinien |
est conforme au terme de la directives CE | de acuerdo con las directivas CE |
in conformità alle direttive CE

2006/95/EC	Low voltage equipment
2004/108/EC	Electromagnetic compatibility
2011/65/EC	Restriction of Hazardous Substances
2002/96/EC	Waste Electrical and Electronic Equipment
2009/129/EC	Eco Design Directive
2006/66/EC	Battery Directive

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EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements.

Standards for Canada and USA

CAN/CSA-C22.2 No. 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.
UL Std. No. 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use - General requirements.

FCC, Part 15, Class A Emission

Zizers, August 9, 2013

Elmar Morscher
CEO

Thomas Neher
Quality Manager

PIPETBOY acu 2 Launch Information – What’s new?

1 PIPETBOY acu 2

The new PIPETBOY acu 2 is a direct replacement of PIPETBOY acu.

The six part numbers and price of the instrument are maintained. As well, PIPETBOY acu 2 is delivered with the 3 year extended warranty program and a free Swiss pocket tool.

2 New battery technology

PIPETBOY acu 2 is equipped with the latest Lithium-Polymer technology, allowing the new generation of pipette controller to be lighter, faster and to have a longer operating time.

2.1 Battery

A new Lithium Polymer battery (also called Li-ion Polymer battery) for PIPETBOY acu 2 is replacing the NiMH battery of the earlier PIPETBOY acu generation.

Battery Compatibility Chart:

Instrument \ Battery	PIPETBOY acu 2	PIPETBOY acu
Lithium Polymer Part No. 155066	Yes	No¹
NiMH Part No. 153005	No²	Yes

1: Operation: PIPETBOY acu will operate.

Charging: the charging time is extended. Do not let the instrument charge for indefinite time.

2: Operation: PIPETBOY acu 2 will operate. The battery indicator is blinking alternatively red and green (battery error message) to inform the user that the wrong battery is used.

Charging: the battery will not be charged.

INTEGRA

The two batteries are based on different technologies. Their prices are different. Please contact your Area Sales Manager for further information.

2.2 Main board

To work with the new battery technology a new main board has been specifically developed for PIPETBOY acu 2. For service purposes, a new spare part has been introduced:

Instrument	PIPETBOY acu 2	PIPETBOY acu
Serial Number	From 1350001	Up to 1305558
Main board	Part No. 155309	Part No. 155300

The price for the main board has been maintained.

2.3 Mains adapter

PIPETBOY acu 2 is compatible with the latest generation of mains adapter for PIPETBOY acu.

It is not compatible with the previous AC mains adapters for PIPETBOY acu. PIPETBOY acu 2 will not be damaged when using AC mains adapter however it will not be charged.

Specification of the mains adapter for PIPETBOY acu 2:

- Input: 100-240 VAC, 50/60 Hz
- Output: 17 VDC $\pm 10\%$ / 180 mA $\pm 10\%$


2.4 Separate rapid charger

The separate rapid charger for two 9 V batteries - Part Number 153010 – is designed to only charge NiCd and NiMH batteries, as indicated on the charger itself. It must not be used to charge the new lithium polymer battery. The rapid charger 153010 is discontinued and will not be replaced.

3 Generation differentiation

To easily identify the two generations of PIPETBOY acu, several differentiation elements have been introduced for PIPETBOY acu 2:

- A new logo on the instrument (INTEGRA PIPETBOY 2)
- A gap in serial number (PIPETBOY acu 2: from SN 1350001)
- Two new labels on the packaging.
- A new label in the battery compartment.

4 New Quick Start Guide

PIPETBOY acu 2 is delivered with a Quick Start Guide that provides the user with clear information on the main features of the instrument. The document is available online on our website, under PIPETBOY acu 2 – Downloads. The Operating Instructions are available online in 5 different languages (EN, DE, FR, IT, ES).

DOMINIQUE DUTSCHER SAS