Gebrauchsanleitung TPP - Zellkulturflaschen / *Instructions for Use TPP Tissue Culture Flasks*

Anleitung

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Instructions

TPP tissue culture flasks are intended for cell culture in manual use. The growth area is opto-mechanically treated for an optimal cell adhesion and proliferation. The chanted neck decreases the risk of media contact at the inside of the screw cap. It also enables a better access for pipettes and scrapers to reach the corners of the flask.

- *Open product packaging and inspect for intact bag and product.*
- *Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the flask.*
- *Close flask:*
 1. *VENT screw cap for open culture systems, by opening the cap the gas exchange takes place*
 - *Introduce screw cap until it clicks in the VENT position.*
 - ✓ *Position for defined aeration of medium: rectangle on 12 o'clock.*
 - *Further turn screw cap over the click about 90° to close flask gas tight.*
 - ✓ *Position for gas tight: rectangle on 3 o'clock*
 2. *Filter screw cap with integrated hydrophobic membrane of pore size 0.22 µm for sterile gas exchange after the cap is screwed down tightly.*
 - *Screw cap down tightly.*
- *Incubate: Stapling ring and rim at the angled neck result in a stable and secure stand in the incubator.*
- *Harvest cells:*
 - *Harvest adherent cells by trypsinization or scrape with cell scraper.*
 - *Harvest suspension cells by exchange of medium.*

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Advice and Recommendation

For homogenous cultivation of the cells it is most important to level the incubator and the continuous supply of medium. Sources of vibrations must be eliminated to avoid circular growth pattern. The necessary gas exchange takes place through the screw cap.

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Materialien

Flasche PS
 Kappe PE
 Membrane PTFE,
 Porengrösse 0.22 µm

Technical Data

Materials

Flask PS
 Screw cap PE
 Membrane PTFE,
 pore size 0.22 µm

Abmessungen	Measurement	90025	90026	90075	90076	90150	90151	90300	90301
Kappe	Cap	VENT		VENT		VENT		VENT	
Kappe	Cap		Filter		Filter		Filter		Filter
Membrane µm	Membrane µm		0.22		0.22		0.22		0.22
Höhe mm	Heights mm	29		40		50		50	
Breite mm	Width mm	51		87		122		140	
Länge mm	Length mm	92		155		210		275	
Wachstumsfläche cm ²	Growth area cm ²	25		75		150		300	
Empf Volumen ml	Rec. volume ml	3 – 8		8 – 22		15 – 45		30 – 85	
Max Volumen ml	Max. volume ml	15		65		165		410	