
NOTICE RAMPE DE FILTRATION

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
DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

INSTRUCTIONS FOR USE FILTRATION MANIFOLD

Summary

1.

■ Introduction p.19

2.

■ Description of the equipment p.19

3.

■ Installation of the filtration manifold p.20

4.

■ Operation of the filtration manifold p.21

Filtration stations..... p.21

Usage of magnetized Funnels and membrane Supports..... p.22

Usage of the disposable filtration Units..... p.23

Filtration..... p.24

Revelation..... p.25

Rinsing and purging the filtration manifold..... p.26

5.

■ Malfunctions p.27

• IMPORTANT

This equipment is intended only for environmental analysis of hot domestic waters, drinking waters and

/ or cooling tower waters, and shall be used exclusively with Legio EZ reagents kits.

For technical assistance, please contact the technical support.

Contact information on the back of this leaflet.

Introduction

This filtration manifold allows to carry out the required filtration operations included within the **Legio EZ** analysis protocols.


You are required to carefully read the dedicated user manual before carrying out the **Legio EZ** analysis. If in doubt, contact the service support.

The equipment is safe in normal conditions of use.

Description of the equipment

The filtration manifold is equipped with 6 holes numbered from 1 to 6, allowing to conduct up to 6 analyzes simultaneously.

- Each hole is equipped with a **cap**, a **membrane Support** and a **magnetized Funnel**.
- On the rear part, the manifold is equipped with a **barbed fitting A** and a $\frac{1}{4}$ **turn valve B**.
- Each manifold is delivered with a length of **silicone tube D** and an **additional fitting apparatus E** allowing to connect 2 filtration manifolds in series, to simultaneously carry out a **maximum of 12 analyzes**.

Installation of the filtration manifold



The manifold(s) must be connected to a suction system respecting good laboratory practices. This suction system must include at a minimum an adequate **waste bottle**, an **antibacterial filter** and a **vacuum pump**.

The vacuum pump must have the following characteristics:

- Flow rate between 4 to 10 L / min (or 0.24 to 0.6m³ / h).
- Depression between 0.4 and 0.8 bar (or airspace from 200 to 400 mbar).

- It is strongly recommended to have a pressure gauge and an air regulator to respect the indicated ranges and thus guarantee the reliability of the analyzes results.
- When using a single manifold, connect the fitting A to the waste bottle with the silicone tube D supplied.
- When using two manifolds in series, connect the fitting A of the first manifold to the waste bottle using the silicone tube D supplied.
- Then connect the fitting A of the second manifold to the ¼ turn valve B of the first manifold using silicone tube D and additional fitting apparatus E supplied.

When using two manifolds in series, the ¼ turn valve of the first manifold must always remain in the open position (handle parallel to the valve body). Only the ¼ turn valve of the second manifold must be used during the analyzes.

Operation of the filtration manifold

Filtration stations

The **Legio EZ** analysis protocols require the use of reusable «**magnetized Funnels** and **membrane Supports**» supplied with the equipment (for the filtration procedure in the first part of the analysis protocol), and afterwards the use of «**disposable filtration Units**» supplied with the reagents kits (for the filtration steps in second part of the protocol).

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
Cap


Hole

Usage of magnetized funnels and membrane supports

- Unscrew and remove one or more caps according to the number of analyzes desired.
- Always use first the closest holes to the waste bottle as a priority.
- Unscrew the membrane Supports from the caps you just removed and screw them with their funnels attached on each hole, tightening without excess.
- Remove the magnetized Funnels by pulling them upwards.

- Place the membranes well centered and well oriented as indicated below. The tongue of the membrane should form an «L».
- Then place each funnel back on the membrane Support.

- Proceed in reverse to remove the membranes.

WARNING


After use, the membranes and funnels must be thoroughly **RINSED** with clean water and **DRIED**.

Usage of the disposable filtration Units



Disposable filtration Unit

- Unscrew and remove one or more membrane Supports with their Funnel, according to the number of analyzes desired, screw them onto their respective caps and set them aside.
- Always use first the closest holes to the waste bottle as a priority.
- Position one disposable filtration Unit per hole, center it and rotate by a quarter of a turn, while pressing down until it reaches the down-stop.
- Then, ensure that each disposable filtration Unit is equipped with a cap on its side.


Filtration


Before starting the filtration process, check that the $\frac{1}{4}$ turn valve is in closed position (handle perpendicular to the valve body). Ensure that the caps are in place on unused holes, screwing them without excess.

- After having poured the liquid to be filtered into each magnetized Funnel or disposable filtration Unit, start the filtration by starting the vacuum pump.
- The liquids are then vacuum filtered, and the filtrates are discharged into the waste bottle.


- When all the Funnels or disposable filtrations Units are empty, stop the pump. The suction circuit must then be returned to ambient pressure **by opening the valve $\frac{1}{4}$ turn.**
- Proceed the same way for the other filtrations.

Revelation

- The analysis protocol requires the installation of a collector in the filtration hole. This collector is delivered with each Legio EZ reagents kit.
- Take the collector and remove its cap. Then remove the disposable filtration Unit from its hole, place the collector inside the hole and replace the disposable filtration unit in its original position.
- During the next filtration, the revealing liquid will be gathered directly in the collector.
- Once this operation is complete, after stopping the filtration, remove the disposable filtration unit.
- Save the collector and screw it. At this point, the disposable filtration Unit can be discarded.

Rinsing and purging the filtration manifold

- The filtration manifold must be rinsed and purged after stopping for more than two hours.
- Turn the valve $\frac{1}{4}$ turn to the open position, then start the filtration. Then, pour about 25 ml of drinking water into each hole.
- Once the holes are empty, keep the vacuum pump running, and screw the caps one by one by screwing them without excess. Start with the one closest to the waste bottle.
- Once all caps are screwed, stop the pump.

Malfunctions

■ Filtration is insufficient or non-existent

POSSIBLE CAUSES	REASONS / SOLUTIONS
Air leak	Check that the ¼ turn valve is in the closed position (handle perpendicularly to the valve body).
Air leak	Check that the caps are correctly screwed on the unused holes.
Air leak	Check that the membrane Supports and their magnetized Funnels or the disposable filtration Units are correctly positioned on their respective holes.
Air leak	Check that each disposable filtration Unit has its cap correctly positioned on its side.
Air leak	Check that the tubes connecting the manifolds to each other and / or to the waste bottle are well connected and not damaged.
Clogged membrane(s)	<p>The event occurs at the beginning of the analysis</p> <p>The water sample is too heavily loaded with impurities, pre-filter your sample on a 10 µm filter.</p> <p>The event occurs after incubation</p> <p>The sample is too charged with microorganisms, contact technical support.</p>

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

Ref : 00-176-B

Rampe de filtration | Filtration Manifold

© Diamidex 2019

DOMINIQUE DUTSCHER SAS