

PRODUCT INFORMATION
Thermo Scientific
Thermo-Start Taq
DNA Polymerase

#AB-0908/B 10 x 250 U
Concentration: 5 U/μL
Lot _____ Expiry Date _____
Store at -20°C



www.thermoscientific.com/onebio

Ordering Information

Component	#AB-0908/B 10 x 250 U
Thermo-Start Taq DNA Polymerase	10 x 50 μL
10X Thermo-Start PCR Buffer	10 x 1.5 mL
25 mM MgCl ₂	10 x 1.5 mL

Rev.8

Estimation of primer melting temperature

For primers containing less than 25 nucleotides, the approx. melting temperature (T_m) can be calculated using the following equation:

$$T_m = 4 (G + C) + 2 (A + T),$$

where G, C, A, T represent the number of respective nucleotides in the primer.

If the primer contains more than 25 nucleotides, specialized computer programs, e.g. REViewer™ (www.thermoscientific.com/reviewer), are recommended to account for interactions of adjacent bases, effect of salt concentration, etc.

COMPONENTS OF THE REACTION MIXTURE

Template DNA

Optimal amount of template DNA in the 25 μL reaction volume is 0.5 – 125 ng. Higher amounts of template increases the risk of generation of non-specific PCR products. Lower amounts of template reduces the accuracy of the amplification.

All routine DNA purification methods are suitable for template preparation e.g., Thermo Scientific GeneJET Genomic DNA Purification Kit (#K0721) or GeneJET™ Plasmid Miniprep Kit (#K0502). Trace amounts of certain agents used for DNA purification, such as phenol, EDTA and proteinase K, can inhibit DNA polymerases. Ethanol precipitation and repeated washes of the DNA pellet with 70% ethanol normally removes trace contaminants from DNA samples.

MgCl₂ concentration

Due to the binding of Mg²⁺ to dNTPs, primers and DNA templates, Mg²⁺ concentration needs to be optimized for maximal PCR yield. The recommended concentration range is 1–4 mM. If the Mg²⁺ concentration is too low, the yield of PCR product could be reduced. On the contrary, non-specific PCR products may appear and the PCR fidelity may be reduced if the Mg²⁺ concentration is too high. If the DNA samples contain EDTA or other metal chelators, the Mg²⁺ ion concentration in the PCR mixture should be increased accordingly (1 molecule of EDTA binds one Mg²⁺).

dNTPs

The recommended final concentration of each dNTP is 0.2 mM. In certain PCR applications, higher dNTP concentrations may be necessary. Due to the binding of Mg²⁺ to dNTPs, the MgCl₂ concentration needs to be adjusted accordingly. It is essential to have equal concentrations of all four nucleotides (dATP, dCTP, dGTP and dTTP) present in the reaction mixture.

Description

Thermo Scientific Thermo-Start Taq DNA Polymerase is a chemically modified version of Thermo Scientific ThermoPrime Taq DNA Polymerase. It is completely inactive at room temperature, preventing the formation and subsequent amplification of non-specific products. The enzyme requires an activation step at 95 °C for 15 minutes. Thermo-Start has 5' to 3' polymerization and exonuclease activity but lacks 3' to 5' exonuclease activity (proofreading).

Applications

- Hot start PCR.
- RT-PCR.
- Highly specific amplification of complex genomic and cDNA templates up to 3 kb.
- Generation of PCR product for TA cloning.
- Amplification of low copy DNA targets.
- Real-time PCR.
- Multiplex PCR.

Source

E. coli cells with a cloned Taq DNA polymerase gene from *Thermus aquaticus* YT1.

Definition of Activity Unit

One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 min at 74°C.

Storage Buffer

The enzyme is supplied in: 20 mM Tris-HCl (pH 9.2 at 25 °C), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Tween® 20, 0.5% Nonidet® P40, 50% (v/v) glycerol.

Inhibition and Inactivation

- Inhibitors: ionic detergents (deoxycholate, sarkosyl and SDS) at concentrations higher than 0.06, 0.02 and 0.01%, respectively.
- Inactivated by phenol/chloroform extraction.

To achieve 0.2 mM concentration of each dNTP in the PCR mixture, use the following volumes of dNTP mixes:

Volume of PCR mixture	dNTP Mix, 2 mM each (#R0241)	dNTP Mix, 10 mM each (#R0191)	dNTP Mix, 25 mM each (#R1121)
50 μL	5 μL	1 μL	0.4 μL
25 μL	2.5 μL	0.5 μL	0.2 μL
20 μL	2 μL	0.4 μL	0.16 μL

Primers

The recommended concentration range of the PCR primers is 0.1–1 μM. Excessive primer concentrations increase the probability of mispriming and generation of non-specific PCR products.

For degenerate primers higher primer concentrations in the range of 0.3–1 μM are often favorable.

CYCLING PARAMETERS

Initial DNA denaturation and enzyme activation

Thermo-Start Taq DNA Polymerase requires an activation step at 95 °C for 15 min. For extra stringency, the enzyme can be activated gradually during the PCR in a series of steps. The initial activation step is replaced by longer (2 min) denaturation steps for the first 7–8 cycles of the reaction.

Denaturation

A DNA denaturation time of 20 seconds per cycle at 95 °C is normally sufficient. For GC-rich DNA templates, this step can be prolonged to 3–4 min. DNA denaturation can also be enhanced by the addition of either 10–15% glycerol or 10% DMSO, 5% formamide or 1–1.5 M betaine. The melting temperature of the primer-template complex decreases significantly in the presence of these reagents. Therefore, the annealing temperature has to be adjusted accordingly. In addition, 10% DMSO and 5% formamide inhibit DNA polymerases by 50%. Thus, the amount of the enzyme should be increased if these additives are used.

Primer annealing

The annealing temperature should be 5 °C lower than the melting temperature (T_m) of the primers. Annealing for 30 seconds is normally sufficient. If non-specific PCR products appear, the annealing temperature should be optimized stepwise in 1–2 °C increments. When additives which change the melting temperature of the primer-template complex are used (glycerol, DMSO, formamide and betaine), the annealing temperature must also be adjusted.

PROTOCOL

To prepare several parallel reactions and to minimize the possibility of pipetting errors, prepare a PCR master mix by mixing water, buffer, dNTPs, MgCl₂, primers and Thermo-Start Taq DNA Polymerase. Prepare sufficient master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and then add template DNA. Keep all reaction components on ice. Reaction set up can be performed at room temperature.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Add the following components for each 25 μL reaction:

	Volume	Final concentration
10X Thermo-Start PCR Buffer	2.5 μL	1X
dNTP Mix, 2 mM (#R0241)	2.5 μL	0.2 mM of each dNTP
Forward primer, 10 μM	1.25 μL	0.5 μM
Reverse primer, 10 μM	1.25 μL	0.5 μM
25 mM MgCl ₂	1.5 μL*	1.5 mM*
Template DNA	0.5–10 μL	0.5–125 ng
Thermo Start Taq DNA Polymerase (5 U/μL)	0.125 μL	0.625 U
Water, nuclease-free (#R0581)	to 25 μL	
Total volume	25 μL	

*Scale up or down the volumes and concentration as appropriate. Volumes of 25 mM MgCl₂, required for specific final MgCl₂ concentration:

Final concentration of MgCl ₂ , mM	1	1.5	2	2.5	3	4
Volume of 25 mM MgCl ₂ to be added for 25 μL reaction, μL	1	1.5	2	2.5	3	4

3. Gently vortex the samples and spin down.

4. If using a thermal cycler that does not use a heated lid, overlay the reaction mixture with half volume of mineral oil.

5. Perform PCR using recommended thermal cycling conditions:

Step	Temperature, °C	Time	Number of cycles
Initial Denaturation	95	15 min	1
Denaturation	95	20 s	30–40
Annealing	50–65	30 s	
Extension**	72	1 min	1
Final Extension	72	5 min	

**Increase length of time in proportion to size of amplicon. Thermo-Start Taq DNA Polymerase extends at approximately 1000 bp/min.

Extension

The optimal extension temperature for Thermo-Start Taq DNA Polymerase is 70–75 °C. The recommended extension step is 1 min/kb at 72 °C.

Number of cycles

The number of cycles may vary depending on the amount of template DNA in the PCR mixture and the expected PCR product yield.

If less than 10 copies of the template are present in the reaction, about 40 cycles are required. For higher template amounts, 30–35 cycles are sufficient.

Final extension

After the last cycle, it is recommended to incubate the PCR mixture at 72 °C for additional 5–15 min to fill-in any possible incomplete reaction products. If the PCR product will be cloned into TA vectors (for instance, using Thermo Scientific InstAclone PCR Cloning Kit (#K1213)), the final extension step may be prolonged to 30 min to ensure the highest efficiency of 3'-dA tailing of PCR product. If the PCR product will be used for cloning using Thermo Scientific CloneJET PCR Cloning Kit (#K1231), the final extension step can be omitted.

CERTIFICATE OF ANALYSIS

Deoxyribonuclease Assay

No degradation of single-stranded and double-stranded [³²P]labeled oligonucleotides was observed after incubation with 10 U of Thermo-Start Taq DNA Polymerase for 4 hours at 37 °C in 1X Thermo-Start PCR Buffer.

Ribonuclease Assay

No contaminating RNase activity was detected after incubation of 160 ng of 2 kb RNA transcript with 25 U of Thermo-Start Taq DNA Polymerase for 4 hours at 37°C.

Functional Assay (PCR)

Thermo-Start Taq DNA Polymerase was tested functionally for end-point PCR amplification of a 1988 bp DNA fragment from 50 ng of human genomic DNA in a 25 μL PCR mixture.

Functional Assay (qPCR)

Thermo-Start Taq DNA Polymerase was tested functionally for quantitative PCR amplification of DNA fragments from 100 pg – 100 ng of human genomic DNA in a 25 μL PCR mixture.

Quality authorized by: Jurgita Zilinskiene

GUIDELINES FOR PREVENTING CONTAMINATION OF PCR REACTION

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Use PCR-certified reagents, including high quality water (e.g., Water, nuclease-free, (#R0581)).
- Always perform "no template control" (NTC) reactions to check for contamination.

GUIDELINES FOR PRIMER DESIGN

Use the Thermo Scientific REViewer primer design software at www.thermoscientific.com/reviewer or follow general recommendations for PCR primer design as outlined below:

- PCR primers are generally 15–30 nucleotides long.
- Optimal GC content of the primer is 40–60%. Ideally, C and G nucleotides should be distributed uniformly along the primer.
- Avoid placing more than three G or C nucleotides at the 3'-end to lower the risk of non-specific priming.
- If possible, the primer should terminate with a G or C at the 3'-end.
- Avoid self-complementary primer regions, complementarities between the primers and direct primer repeats to prevent hairpin formation and primer dimerization.
- Check for possible sites of undesired complementary between primers and template DNA.
- When designing degenerate primers, place at least 3 conserved nucleotides at the 3'-end.
- Differences in melting temperatures (T_m) between the two primers should not exceed 5 °C.

(continued on reverse page)

NOTICE TO PURCHASER:

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,127,155, 5,677,152 (claims 1 to 23 only) and 5,773,258 (claims 1 and 6 only). The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.