**Operating Instructions
Bedienungsanleitung
Mode d'emploi
Manual de Instrucciones
Istruzioni d'uso**

METTLER TOLEDO

**METTLER TOLEDO
RS-P26 Printer
Drucker RS-P26
Imprimante RS-P26
Impresora RS-P26
Stampante RS-P26**

Page 3

English

Seite 15

Deutsch

Page 27

Français

Página 39

Español

Pagina 51

Italiano

DOMINIQUE DUTSCHER SAS

Contents

Overview of your RS-P26	4	5 Printer test	9
1 Introduction	5	6 Maintenance.....	9
1.1 Capabilities of the RS-P26 Printer.....	5	6.1 Changing the paper roll.....	10
1.2 Overview of functions.....	5	6.2 Replacing the ribbon.....	10
2 Startup procedure	5	6.3 Cleaning	10
2.1 Safety measures.....	5	6.4 Disposal	10
2.2 Power supply.....	6	7 Troubleshooting	11
2.3 Switching the RS-P26 on and off.....	6	8 Consumables and spare parts.....	12
2.4 Inserting paper.....	6	9 RS232C Interface.....	12
2.5 Connecting the RS-P26 to the balance...	6	10 Technical data RS-P26	13
3 Keystroke functions	7	11 Declaration of conformity	14
4 Menu.....	8		

Overview of your RS-P26

Controls:


- 1 On/Off Switch
- 2 Control panel
- 3 Pilot lamp – status indicator
- 4 RS232C Interface connector
- 5 AC power cord connection

Connections:


Functions:

- Zeroing / taring the balance
- Abort – exit

- Printout of date and time
- Scroll down menu option

- Paper feed
- Scroll up menu option

- Printing of stable weight values
- Opening the menu – select next menu option
- Save settings

1 Introduction

1.1 Capabilities of the RS-P26 Printer

The RS-P26 is a simple to operate dot matrix printer for METTLER TOLEDO balances with an RS232C interface. The printout on normal paper meets the requirements of modern quality assurance systems (GLP, GMP, ISO 9001, etc.). These norms define the principles and user requirements related to measuring and test instruments. The RS-P26 complies with the requirements of measurement traceability by precisely documenting the origin of data, as well as calibration procedures. It meets the requirement of archiving measurement results by providing a thermally stable and light resistant printout. When connected to METTLER TOLEDO balances, the RS-P26 prints all data sent by the balance without restriction, which are not only weight values and calibration procedures, but also the results of all applications, built into the balance.

1.2 Overview of functions

Printing of weight values

The RS-P26 prints out only the stable weight value

Date and time

Selected as an option with date and time in the record header (stored with battery back-up)

Zeroing/Taring of balance

Zero/Tare the connected balance at any time.

2 Startup procedure

2.1 Safety measures

The RS-P26 has been tested for the connections and intended purposes documented in these Operating Instructions. However, it is your responsibility to perform the necessary tests related to the methods and purposes of its intended use. Therefore, you should observe the following safety measures.

- The RS-P26 Printer must be used indoors only, and never in hazardous areas.
- The electronics of the RS-P26 Printer are protected against the ingress of splashing water and dust. However, since the printer housing is not water-tight, the printer should not be used in the presence of liquids.
- Comply with cleaning instructions (see section 6.3), protect power plug against wet conditions.
- Never open the printer housing. If you have any problems with your printer, please contact your responsible METTLER TOLEDO service.

2.2 Power supply

The power supply of the RS-P26 is integrated in the housing (100-240 V, 50-60Hz). This allows you to connect the RS-P26 to your local line voltage.

Attention


- Connect the RS-P26 only to AC Power receptacles with a connection to ground!
- Use an extension cable with an equipment grounding conductor!

2.3 Switching the RS-P26 on and off

Use the main switch <I/O> to power the RS-P26 on and off. As soon as the RS-P26 is powered on the unit will start with some line feeds. The pilot lamp will flash until a instrument is connected.

2.4 Inserting paper





- Remove paper cover (pull on the back) (1)
- Insert the paper rotary axis through the core of the paper roll (2).
- Feed paper through the slot in the printing unit in the direction of the arrow (3).
- Press and hold until enough paper has been fed (4).
- Lead paper through the paper cover slot
- Replace the paper cover (1).

2.5 Connecting the RS-P26 to the balance

- Attach the interconnect cable to the printer interface of the balance and to the printer.
- Make sure that there is a match between the communication parameters of the balance and the printer.
- The RS-P26 is ready for use as soon the pilot lamp stops flashing.

3 Keystroke functions

<p>Operation Level</p>  <p>Press briefly</p>  <p>zero / tare the balance</p>  <p>print Date and Time</p>  <p>one line feed</p>  <p>print</p>	 <p>Press and hold down</p>  <p>continuous feed (until release)</p>  <p>open menu</p>
<p>Menu Level</p>  <p>exit down up next</p>	 <p>store</p>

Example for setting the day:

RS-P26

Balance Display

* MENU

* rESEt

* DAY 18

* DAY 21

* StorEd

Open the menu:

- Hold down the «Menu» key until "reset" appears.
(Print out: --- See Balance Display ---)

- Select the menu option "Day" with the «Menu» key.

- Select the "Day" – scroll up with the «↑» key.
– scroll down with the «↓» key.

Save the settings:

- Hold down the «Menu» key until "Stored" appears.

5 Printer test

The RS-P26 is equipped with a self-test that automatically prints out the character set of the printer. (standard IBM/DOS).

Starting the test

- Press the key «↑» and power on the RS-P26.
The printer prints out the character set

Terminating the test

- Power off the RS-P26

6 Maintenance


Under normal conditions the RS-P26 requires practically no maintenance. The METTLER TOLEDO service department is available if servicing is required. Please ask your responsible METTLER TOLEDO dealer for details.

6.1 Changing the paper roll

Insertion of a new roll of paper is described in section 2.4. See section 8 (Consumables) for the paper roll's order number.

Note: You may obtain the paper rotary axis from the old paper roll.

6.2 Replacing the ribbon

- Remove the paper cover (1).
- Pull paper out of printing unit (2).

- To remove the ribbon cartridge, press in direction of arrow (3).
- Insert the new ribbon, and if necessary, take up slack with tensioning wheel (4). See section 8 (Consumables) for the ribbon cartridge order number.
- Insertion of a new roll of paper is described in section 2.4.

6.3 Cleaning

Since the printer housing is made of top grade, resistant materials, all commercially available mild cleaning agents may be used (see safety precautions in section 2.1).

6.4 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.


If this device is passed on to other parties (for private or professional use), compliance with this regulation will be required as well.

Thank you for your contribution to environmental protection.

7 Troubleshooting

Problem / Message	Possible cause	Correction
Printout unreadable	<ul style="list-style-type: none"> • Ribbon worn or entangled • Service life of printing unit at an end 	<ul style="list-style-type: none"> • Change ribbon cartridge or tense the ribbon properly • Please contact METTLER TOLEDO service
Pilot lamp does not light	<ul style="list-style-type: none"> • No power • Power supply fault 	<ul style="list-style-type: none"> • Switch printer on • Please contact METTLER TOLEDO service
Pilot lamp flashes ----- No Connection ----	<ul style="list-style-type: none"> • RS232C connection is disconnected / not attached • The interface parameters do not match the balance 	<ul style="list-style-type: none"> • Check the RS232C interconnect cable • Check interface parameter settings on printer and balance
-----Timeout -----	<ul style="list-style-type: none"> • Data transfer has not completed in the last 45 s 	<ul style="list-style-type: none"> • Move balance to more stable location • Check the balance
----- Over Load -----	<ul style="list-style-type: none"> • Balance in overload range 	<ul style="list-style-type: none"> • Check the weight of sample
----- Under Load -----	<ul style="list-style-type: none"> • Balance in underload range. Weighing pan is not in position 	<ul style="list-style-type: none"> • Check the balance
-- Transmission Error --	<ul style="list-style-type: none"> • Balance does not understand printer instruction 	<ul style="list-style-type: none"> • Check interface parameter settings on printer and balance
---- Not Executable ----	<ul style="list-style-type: none"> • Balance can not execute printer instruction 	<ul style="list-style-type: none"> • Check the balance
--- Date/Time Not Set --- ----- Memory Error -----	<ul style="list-style-type: none"> • Printer error messages 	<ul style="list-style-type: none"> • Check the balance • Reset the printer – set date/time • Back-up battery depleted Please contact METTLER TOLEDO service


8 Consumables and spare parts

Pos.	Item	Part number	Delivery quantity
1	Paper cover	12120734	1 pcs
2	Paper roll 57.5 x \varnothing 50 to 60 mm	12120799	2 pcs
3	Paper rotary axis	12120735	1 pcs
4	Ribbon cartridge, black	12120798	1 pcs
-	RS9 M/F Cable, 1.5 m	12120708	1 pcs

9 RS232C Interface

The RS-P26 is equipped with an RS232C interface to connect METTLER TOLEDO instruments.

- 9-pin male connector
- Matching to other device (transmission parameters) see section 4.

10 Technical data RS-P26

Power line voltage, frequency	100 – 240 VAC, 50/60 Hz
Power consumption	150 mA (100 V), 60 mA (240 V)
Altitude	up to 4000 m
Ambient temperature	0 °C – 40 °C
Relative humidity	max. 80 % at 31 °C, decreasing linearly to 50 % at 40 °C, noncondensing
Overvoltage category	Class II
Pollution degree	2
Printing unit	Dot matrix printer 5x7, 24-character line length
Printing speed	2.3 lines per second
Ribbon cartridge	Exchangeable, black
Paper roll	Standard paper 57.5 mm x ø 50 to 60 mm, integrated in housing, commercial size
Print quality	Light resistant and thermally stable printing (GLP, GMP, ISO 9001)
Interface	RS232C
Printer dimensions	(WxDxH) 120 mm x 203 mm x 73 mm
Packaging dimensions	(WxDxH) 255 mm x 205 mm x 150 mm
Net weight	740 g (paper roll incl.)
Compatibility limits	<ul style="list-style-type: none"> • Function Date / Time: not applicable on balances with built-in clock (built-in clock has priority) • Function Zeroing / Taring: not applicable on all models of balances

11 Declaration of conformity

The undersigned declare on behalf of

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that the Printer **METTLER TOLEDO RS-P26** to which this declaration relates (serial number specified on the product) is in compliance with the below mentioned EEC Directives (including all amendments)

73/23/EEC Low Voltage Directive

89/336/EEC Electromagnetic compatibility

and that the following standards have been applied

IEC/EN61010-1:2001,

EN61326:1997+ A1:98+A2:01+A3:03 (class B)

EN61326:1997+ A1:98+A2:01+A3:03 (Minimal requirements)

for Canada, USA and Australia

CAN/CSA-C22.2 No.61010-1-04, UL Std. No.61010A-1, FCC, Part 15, class A,
AS/NZS CISPR 22, AS/NZS 61000.4.3

Greifensee, 01.09.2006

Mettler-Toledo AG
Laboratory & Weighing Technologies

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