## LE Sensors Conductivity

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Any other use carries the risk of injury and damage. METTLER TOLEDO assumes that users are trained in the safe way to do their work, deal with potential hazards of their work area and take account of the necessary protective gear.

Any other use carries the risk of injury and damage. METTLER TOLEDO assumes that users are trained in the safe way to do their work, deal with potential hazards of their work area and take account of the necessary protective gear.

**NOTICE**  
Some chemical substances can damage the measuring cell and sensor.  
Do not expose the measuring cell to organic solvents, acids or highly concentrated chemicals neither when cleaning nor when being measurement.

**Commissioning**  
1. Visually check the sensor for damage.  
2. Connect the sensor to the measuring instrument.

**Operation**  
1. Rinse the sensor with deionized water.  
2. Carefully dry the sensor by with a tissue.  
3. Calibrate the sensor.  
4. The sensor is ready for measurement.

**Keep practice**  
• To reduce air bubbles, immerse the sensor of an angle and then move it to an upright position.  
• If an air bubble is trapped inside the measuring out, gently tap the bubble to remove it.  
• Make sure that the measuring cell is immersed completely.  
• Wash left the sensor and the sample have the same temperature before taking the measurement.

**Clean the sensor**  
1. Rinse the sensor with deionized water.  
2. If there is a build-up of solids on or in the measuring cell, immerse gently with a cotton swab soaked in deionized solution.  
3. Rinse the sensor with deionized water.  
4. Carefully dry the sensor by with a tissue.

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**Safety Information**  
METTLER TOLEDO conductivity sensors are intended to measure conductivity in aqueous solutions of atmospheric pressure.

**Further information**  
• Measurement and calibration: see Operating Instructions of the measuring instrument. <http://www.mt.com/Italy>  
• Technical data: see Data Sheet. <http://www.mt.com/Italy>

**Disposal**  
In accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of as domestic waste. This device belongs to countries outside the EU, per their specific requirements.

**WEEE**  
In accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of as domestic waste. This device belongs to countries outside the EU, per their specific requirements.

**WEEE**  
In accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of as domestic waste. This device belongs to countries outside the EU, per their specific requirements.

**WEEE**  
In accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of as domestic waste. This device belongs to countries outside the EU, per their specific requirements.

**WEEE**  
In accordance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this device may not be disposed of as domestic waste. This device belongs to countries outside the EU, per their specific requirements.

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© Mettler Toledo 2021 021018
Mettler Toledo GmbH