
# METTLER TOLEDO InLab® 73X Series InLab® 74X Series

METTLER TOLEDO

In accordance with the requirements of European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), this device must not be disposed of with household waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related. Thank you for your contribution to environmental protection.

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### **InLab® 73X series:**

4-electrode (4 pole) carbon sensors, ideal for general purpose use, general specifications: 0.01 to 1000mS/cm, 0.0 to 100°C, supplied with quality certificate

### **InLab® 74X series:**

2-electrode (2 pole) steel sensors, for low conductivity (pure water), general specifications: 0.001µS/cm to 500µS/cm, 0.0 to 100°C, supplied with calibration certificate stating the actual cell constant

### **Intended use and principles**

These products combine a state-of-the-art conductivity cell with a thermistor which is integrated in the sensor body enabling automatic temperature compensating (ATC) to provide an easy to use, reliable and accurate sensor. The electrode design minimizes the potential for external interferences from sample containers or other sensors. These sensors are equipped with cables in various lengths and two types of plugs to match METTLER TOLEDO conductivity meters, with and without Intelligent Sensor Management (ISM).

### **Intelligent Sensor Management (ISM®)**

The METTLER TOLEDO dual channel portable meters SevenGo Duo (pro) SG23, 68 and 78 have the ability to automatically recognize ISM sensors upon connection, transfer calibration data to and from the sensor etc. For more information regarding ISM® please read the respective chapter of the operating instructions supplied with the meter.

### **Handling of the sensor: tips and tricks**

For optimum performance, use the following operating procedures:

1. Remove all packaging material from sensor before use. Retain the quality certificate and the operating instructions on file for reference.
2. To prevent contamination from solution to solution, rinse the sensor with distilled or de-ionized water between measurements. If this quality of water is not available, rinse the sensor with the solution to be measured, and then discard the solution.
3. When measuring, make sure the solution is above the cell chamber so that the whole measuring cell is immersed.
4. Ensure the electrode measurement area is free of bubbles when measuring. To reduce air bubbles, immerse the sensor at an angle and then rise to the vertical position. Remaining small bubbles can be removed by gently tapping against the sensor shaft.
5. Allow sufficient time for the sensor to stabilize when measuring samples of different temperatures.
6. The InLab® 73X series is not recommended for low conductivity solutions (<10µS/cm); for these applications use the InLab® 74X series. A flow cell (ME-51302257) may also be needed for accurate work on low conductivity solutions to exclude the contamination by CO<sub>2</sub> or to save precious sample.
7. Clean the sensor with distilled or de-ionized water after use and keep it dry or wet when not in use. Keeping it wet speeds up the response time.
8. If there is a build up of solids inside the cell chamber, carefully remove it with a cotton wool bud soaked in detergent solution, and then rinse the sensor with distilled or de-ionized water.
9. Ensure that the cable connectors are kept free of moisture and dirt.
10. For calibration and measurement procedures, refer to the operating instructions of the conductivity meter. A booklet, "Guide to Conductivity and Dissolved Oxygen", is also available. See the operating instructions of the meter for details.
11. In unstirred solutions, measurement drift is possible. In general, measurement using stirred samples is preferred.

### **Reagents and supplies**

To achieve optimum performance of this conductivity sensor, the correct calibration solutions must be used. See the operating instructions supplied with the conductivity meter for details.

### **Precautions and limitations**

1. To prevent static damage to the meter, avoid touching the cell chamber area of the sensor. Always disconnect the sensor from the meter before cleaning.
2. The InLab® 73X series has limited chemical resistance. Do not use harsh chemicals or highly concentrated solvents when cleaning or when taking measurements.
3. Do not use the sensor outside the recommended temperature range. Irreparable damage might occur to the sensor.
4. The specifications and the performance of this sensor are guaranteed only when used with the appropriate conductivity meter.

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### Intelligent Sensor Management (ISM®)

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© 01/2016 Mettler-Toledo GmbH, 51710642C
**Mettler-Toledo GmbH**