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#### Storage

Membrane filters should be stored in cool and dry places, protected from sunlight. They should be stored separately.

#### Caution!

Membrane filters made of cellulose nitrate are inflammable. Flame point at 180° C.

#### Handling

Membrane filters are to be handled with forceps when placing them on the filter support of the filter holder device. Membrane filters of larger sizes (approx. 100 mm dia.) should be taken from the package covered by one interleaving paper using two forceps to hold the membrane at two spots. Membranes should never be folded or touched with bare fingers. Perspiration affects wetting properties of the membrane filters.

#### Sterilization

All membrane filters from Whatman can be autoclaved at 121° C and 1 bar overpressure. Alternatively, sterilization with ethylene oxide is recommended.

For steam sterilization filters made of cellulose acetate or mixed esters are to be placed wet into the filtration device.

All other membrane filters are placed dry into a **dry filtration device**.

All membrane filters can readily be tested for integrity by the bubble point test method. Minimum pressure for 0.2 µm removal rating pore size is about 3.5 bar and for 0.45 µm about 2.5 bar (cf. web catalog for exact values).

Polyamide membrane filters can be autoclaved more frequently; (steam sterilization (30 minutes) at 121° C either "in-line" or in autoclave).

For filtration of parenteralia pertinent regulations of the pharmacopoeia should strictly be observed e. g. prewashing.

Whatman also offers sterile membranes ready for use and individually packed.

#### Remark:

More detailed information on properties of membrane filters on the web site [www.gelifesciences.com/whatman](http://www.gelifesciences.com/whatman).

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GE Healthcare UK Limited
GE Healthcare Bio-Sciences AB
GE Healthcare Europe GmbH
GE Healthcare Bio-Sciences Corp.
Hyclone Laboratories, Inc.