

Featured Products

Amicon® Ultra Centrifugal Filters

Fast and Easy Protein Concentration

Amicon® Ultra Centrifugal filters provide fast sample processing and promote high sample recoveries, even in dilute samples, through ultrafiltration. The unique features of the Amicon® Ultra centrifugal filters give you the fastest, most efficient concentration for sensitive downstream applications.

Amicon® Ultra Centrifugal Filter Advantages:

Maximize Concentration with Highest Protein Recovery – True Engineered Dead Stop

- Avoids spinning to dryness
- Provides a predictable concentration factor
- No need to calibrate for several samples to run in parallel

Reverse Spin Recovery

- Reverse spin devices enable you to maximize protein recovery, especially with small dilute samples, without introducing pipetting errors
- Low binding membrane and polypropylene housing for > 90 % sample recovery



Ultracel® Low-binding Membranes

Fast and Efficient Concentration without Compromise


- Vertical membrane design aligns with filtrate rather than perpendicular for less clogging, less waste, and faster filtration
- Ultra-fast sample processing achieving concentration in as little as 10 minutes
- 25- to 80-fold concentration in a single step

Broad Chemical Compatibility

- Heat-sealed membrane eliminates adhesives and downstream extractables
- Large spectrum of compatibility
- Compatible with pH 1 to 9


Reliable Samples

- Spin precious samples with confidence in one robust, sleek unit that prevents leakage



Amicon® Ultra 4 mL Filters – Fast Spin Times with Excellent Recovery

Amicon® Ultra 4 mL filters were tested for percent recovery and spin time.



Consistently high recovery of diverse proteins with Amicon® Ultra filters

Concentration and percent recovery using Amicon® Ultra Filters: 4 different devices (Amicon® Ultra-0.5 mL, Amicon® Ultra-2 mL, Amicon® Ultra-4 mL, Amicon® Ultra-15 mL devices) were tested (3 kDa membrane with Cytochrome C, 10 kDa membrane with Cytochrome C, 30 kDa membrane with BSA and 100 kDa membrane with IgG) to determine percent recovery and concentration factor.

To select an Amicon® Ultra Centrifugal Filter, identify the starting volume, molecular weight of protein or nucleic acid being concentrated, final volume and concentration factor.





Then consult the product selection chart below to choose the Amicon® Ultra filter with the right nominal molecular weight cutoff (NMWCO).

	Amicon® Ultra-0.5 filter	Amicon® Ultra-2 filter	Amicon® Ultra-4 filter	Amicon® Ultra-15 filter	
Starting Volume	< 0.5 mL	< 2 mL	< 4 mL	< 15 mL	
Molecular Weight (MW)	Proteins				
	NMWCO (Da)				
	6 < MW < 20 k	3,000	3,000	3,000	3,000
	20 < MW < 60 k	10,000	10,000	10,000	10,000
	60 < MW < 100 k	30,000	30,000	30,000	30,000
	100 < MW < 200 k	50,000	50,000	50,000	50,000
200 k < MW	100,000	100,000	100,000	100,000	
Length	Single-Stranded and Double-Stranded Nucleic Acids				
	NMWCO (Da)				
137-1159 bp	30,000	30,000	30,000	30,000	
Particle Diameter (DIA)	Nanoparticles				
	NMWCO (Da)				
	1.5 < dia < 3 nm	3,000	3,000	3,000	3,000
	3 < dia < 5 nm	10,000	10,000	10,000	10,000
	5 < dia < 7 nm	30,000	30,000	30,000	30,000
	7 < dia < 10 nm	50,000	50,000	50,000	50,000
10 nm < dia	100,000	100,000	100,000	100,000	

NMWCO: Nominal Molecular Weight Cut Off

10,000 NMWCO Amicon® Ultra-4 and -15 filters are both CE marked and registered for *in vitro* diagnostic use.

Once you've chosen the right Amicon® Ultra filter for your needs, choose your rotor, G force and spinning time for concentrating your molecule. Designed as standard 1.5 mL, 15 mL conical or 50 mL conical tubes, Amicon® Ultra filters fit all standard rotor types.

	Amicon® Ultra-0.5 filter	Amicon® Ultra-2 filter	Amicon® Ultra-4 filter	Amicon® Ultra-15 filter	
					
Choose a rotor and G force	Starting Volume	< 0.5 mL	< 2 mL	< 4 mL	< 15 mL
	Final Volume	15–20 µL	15–70 µL	50 µL	200 µL
	Design of the Device	Standard 1.5 mL	Standard 15 mL	Standard 15 mL	Standard 50 mL
	Fixed-Angle (35 °) Rotor	14,000 g 1,000 g reverse spin	7,500 g 1,000 g reverse spin	5,000 g for 100,000 7,500 g for all other MWCO	5,000 g
	Swinging Bucket Rotor	N/A	4,000 g 1,000 g reverse spin	4,000 g	4,000 g
Concentration Factor	Final Volume	15–20 µL with reverse spin	15–70 µL with reverse spin	50 µL	200 µL
	Concentration Factor	X25–X30	X14–X67	X80	X75
Adjust spinning time	For Proteins and Nanoparticles				
	NMWCO (Da)				
	3,000	30 min.	60 min.	40 min.	40 min.
	10,000	15 min.	40 min.	15 min.	20 min.
	30,000	10 min.	20 min.	10 min.	20 min.
	50,000	10 min.	15 min.	10 min.	15 min.
100,000	10 min.	30 min.	10 min.	15 min.	
Single-Stranded and Double-Stranded Nucleic Acids					
30,000	10 min.	15 min., fixed angle 40 min., swinging rotor	10 min., 5,000 g, fixed angle	10 min., 5,000 g, fixed angle	

Amicon® Ultra Centrifugal Filters

	Product	Amicon® Ultra-0.5 filter	Amicon® Ultra-2 filter	Amicon® Ultra-4 filter	Amicon® Ultra-15 filter
	Maximum initial sample volume (mL)	0.5	2	4	15
	Final concentrate (retentate) volume (µL)	15-20	15-70	30-70	150-300
NMWCO (Da)	Qty/Pk				
3,000	8	UFC500308		UFC800308	UFC900308
	24	UFC500324	UFC200324	UFC800324	UFC900324
	96	UFC500396		UFC800396	UFC900396
	500	UFC5003BK			
10,000	8	UFC501008		UFC801008	UFC901008
	24	UFC501024	UFC201024	UFC801024	UFC901024
	96	UFC501096		UFC801096	UFC901096
	500	UFC5010BK			
30,000	8	UFC503008		UFC803008	UFC903008
	24	UFC503024	UFC203024	UFC803024	UFC903024
	96	UFC503096		UFC803096	UFC903096
	500	UFC5030BK			
50,000	8	UFC505008		UFC805008	UFC905008
	24	UFC505024	UFC205024	UFC805024	UFC905024
	96	UFC505096		UFC805096	UFC905096
	500	UFC5050BK			
100,000	8	UFC510008		UFC810008	UFC910008
	24	UFC510024	UFC210024	UFC810024	UFC910024
	96	UFC510096		UFC810096	UFC910096
	500	UFC5100BK			

Amicon® Ultra-4 and -15 Centrifugal Filters, registered for IVD use

Description	NMWCO	Qty/Pk	Catalog No.
Amicon® Ultra-4 Centrifugal Filter	10 KDa	8	UFC801008D
		24	UFC801024D
		96	UFC801096D
Amicon® Ultra-15 Centrifugal Filter	10 KDa	8	UFC901008D
		24	UFC901024D
		96	UFC901096D

