# MILLEX®-FG<sub>50</sub>

## Filter Unit

The Millex-FG<sub>50</sub> filter unit is a non-sterile non-pyrogenic filter device, manufactured of polypropylene with bi-directional filter support screens. An unlaminated, 50 mm diameter, 0.2 µm PTFE hydrophobic membrane is sealed into the unit between screen supports.

The filter unit is designed primarily for inline filtration of air or gases, for venting air from pressurized closed containers, the filtration of non-aqueous liquids such as HPLC solvents and as a moisture barrier between a vacuum source and the evacuated container. The Millex-FG<sub>50</sub> filter unit will withstand repeated autoclaving for at least 10 times.

### Directions for Use

#### With Syringes

Fill syringe with liquid, then insert tip firmly into one of the unit's connectors, and twist slightly clockwise to assure a tight connection. With filter unit positioned in upward direction, slowly expel air fro filter. When liquid just exits from the unit, filter into receiver, taking care not to use a pressure exceeding 2 bar.

#### With Hose Connections

Choose a silicone rubber or latex tubing of 4-9.5 mm I.D. Fit snugly over connector. To avoid blowing off tubing from unit, it is advisable to use hose clamps.

#### Autoclaving

If the device requires sterilization, it can be autoclaved. Remove the two connector plug caps, wrap and seal in Kraft paper, or other suitable autoclave paper. Autoclave at 121 °C for 20 minutes.

### Precautions

Do not apply pressures exceeding 2 bar.

If unit is used for filtering solvents or other non-aqueous solutions, check compatibility with polypropylene prior to filtration.

If unit is to be used in applications where pyrogens are of concern, do not remove connector plug caps until ready for autoclaving and use.

### Specifications

Temperature Limit: 130 °C (266 °F)  
Pressure Limit at 25 °C: 2 bar (30 psig) inlet/differentiation  
Filter Pore Size: 0.2 µm, FGUP filter  
Filtration Area: 19.6 cm<sup>2</sup> (3 in<sup>2</sup>)


### Typical Flow Rates

Air Flow Rate: ≥ 15 l/min at 0.3 bar (5 psig)  
Methanol Flow Rate: ≥ 150 ml/min at 0.3 bar (5 psig)

### Connections

Code	Connections
Code 0	Stepped hose 1/2"-1/4" (latex). Female Luer slip at both ends
Code 5	Stepped hose 1/2"-1/4" (latex) and 1/8" male NPT
Code 6	Stepped hose 1/2"-1/4" (silicone) and 1/8" male NPT
Code 7	1/8" male NPT at both ends
Code 8	Stepped hose 1/2"-1/4" (silicone). Female Luer slip at both ends

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### Millex-FG<sub>50</sub> Abmessungen / Code / Kat.-Nr.

## Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call **1-800-MILLIPORE** (1-800-645-5476). Outside the U.S., see your Millipore catalogue for the phone number of the office nearest you or go to our web site at [www.millipore.com/offices](http://www.millipore.com/offices) for up-to-date worldwide contact information. You can also visit the tech service page on our web site at [www.millipore.com/techservice](http://www.millipore.com/techservice).

Millipore provides innovative tools, services and application expertise that improve the scientist's productivity in the pharmaceutical, academic, life science, applied science, clinical and government research sectors. Our product range spans water purification, protein research, analytical sample preparation, cell-based assays and high throughput screening. Please visit [www.millipore.com/bio-science](http://www.millipore.com/bio-science) for more details.

## Standard Warranty

**Millipore Corporation** ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

**Millex®-FG<sub>50</sub>**

OPERATING INSTRUCTIONS  
MODE D'EMPLOI  
GEBRAUCHSANWEISUNG

**Millipore S.A.  
Zone Industrielle  
67120 Molsheim-France**

**MILLIPORE**