



User Guide

Swinny Stainless Steel 13 mm Filter Holder


Cat. No. XX3001200

Introduction

The Swinny Stainless Steel 13 mm Filter Holder is a syringe-operated filter holder designed to sterilize or clarify small volumes of liquids by pressure filtration. It can be used to sterilize tissue culture media or clarify HPLC samples, solvents, etc. The filter holder can be autoclaved with the filter in place for clarification/sterilization applications that do not require removing the filter intact for further use.

Filter Loading

1. Unscrew the inlet cover from the outlet base.
2. Remove the O-ring, support screen, and flat gasket.
3. Inspect the gasket, O-ring, and sealing surfaces for damage; imperfections can cause leaks during filtration. Replace damaged parts.
4. Place the flat gasket in the base and place the support screen on top of it.
5. Using smooth-tip forceps, place a 13 mm diameter filter on the support screen, then place the O-ring on top of the filter. If a depth-type prefilter (10 mm diameter) is required to extend the life of the final filter, center the prefilter on top of the final filter before adding the O-ring.
6. Screw the cover down into the base and finger tighten firmly so that the filter is securely sealed.



Operation

CAUTION: Do not exceed the maximum inlet and differential pressure of 6.9 bar (100 psi).

When using a syringe smaller than 10 mL, the syringe should have a Luer-Lok™ connection (not Luer-slip) to avoid leakage at the higher pressures generated by small syringes.

1. Draw the unfiltered sample into the syringe.
2. Attach the Luer-Lok™ inlet of the filter holder to the syringe. If required, attach a hypodermic needle to the filter holder outlet for dispensing the filtered liquid.
3. Press the syringe plunger to filter the sample.

Autoclaving (optional)

1. Make sure that all inner surfaces of the filter holder are completely dry, then assemble the holder with a filter in place. Tighten partially.
2. Wrap the filter holder with lint-free autoclave paper, and autoclave for 15 minutes at 121 °C and 1 bar (15 psi).
3. After the unit has cooled, fully tighten it, being careful not to touch the outlet tube.
4. Use aseptic technique when making outlet connections.

Maintenance

1. Disassemble the holder and inspect all components for signs of damage. Grit embedded in PTFE surfaces can sometimes be dislodged without damaging the surface, but it is preferable to replace the O-ring or gasket.
2. Wash all components thoroughly with a sponge, hot water, and a nonabrasive cleanser. Scrub threaded parts, recesses, and orifices with a stiff bristle brush and use a pipe cleaner to remove contaminants from inlet and outlet orifices.

CAUTION: To avoid component damage, never use steel wool or abrasive materials on any part of the filter holder.

Maintenance, continued

3. Rinse all parts thoroughly with hot running water, then rinse several times with cold, laboratory-grade water.
4. Allow the components to air dry while disassembled.

CAUTION: Do not wipe dry with paper or cloth, as this leaves fibers and lint on the surfaces and also generates electrostatic forces that can attract more dirt. Filtered compressed air can be used to facilitate drying.

Specifications

Maximum pressure	6.9 bar (100 psi) inlet and differential
Filter diameter	13 mm diameter
Pre-filter diameter	10 mm diameter (thick depth pre-filter)
Effective filter area	Approximately 0.8 cm ² (0.13 in ²)
Weight	23 g (0.8 oz)
Connections	Female Luer-Lok™ inlet, male Luer-slip outlet
Dimensions	Diameter: 16.0 mm (0.63 in.) Length: 36.6 mm (1.44 in.)
Materials of construction	Inlet cover: 316 stainless steel Outlet base: 316 stainless steel Support screen: 316 stainless steel Flat gasket: polytetrafluoroethylene (PTFE) O-ring: PTFE

Statement Regarding Compliance with the Pressure Equipment Directive, 97/23/EC

EMD Millipore Corporation certifies that this product complies with the European Pressure Equipment Directive, 97/23/EC of 29 May 1997. This product is classified under Article 3 § 3 of the Pressure Equipment Directive. It has been designed and manufactured in accordance with sound engineering practices to ensure safe use. The product is accompanied by user instructions and bears markings to permit identification of EMD Millipore Corporation as the manufacturer or authorized representative of this product within the European Community. In compliance with Article 3 § 3 of the Pressure Equipment Directive, this product does not bear the CE mark.

Ordering Information

This section lists the catalogue numbers for Swinny 13 mm Filter Holder, replacement parts, and accessories. See Technical Assistance for contact information. You can purchase these products on-line at www.millipore.com/products.

Description	Catalogue Number
Swinny Stainless Steel 13 mm Filter Holder (includes inlet cover [1], outlet base [1], O-ring [1], support screen [1], flat gasket [1])	XX3001200
Swinny 13 mm Filter Holder Replacement Parts Kit (includes O-ring [4], flat gasket [4], filter holder support screen [4])	XX30012RK
Support screen (1/pk)	XX3001210
Flat gasket, PTFE (10/pk)	XX3001202
Filter forceps, stainless, blunt end (3/pk)	XX6200006P
Plastic syringe, 20 mL capacity (12/pk)	XX1102012
Plastic syringe, 50 mL capacity (5/pk)	XX1105005

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by EMD Millipore Corporation ("Millipore") or an affiliate. Neither EMD Millipore Corporation nor any of its affiliates assumes responsibility for any errors that may appear in this document.

Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-221-1975. Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms ("Conditions of Sale").