Référence 053301

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With this type of thermometer you get minimum, maximum and present temperatures at a glance.

And this is how to use it:

Remove the magnetic stripe. By pressing the button the two blue markers will move down to the transparent column (Reset of maximum and minimum values). This enables you - at any given moment - to take the reading of the previous changes in temperature from the bottom parts of the markers; the reading of the maximum value can be taken on the right, the minimum value on the left side of the scale. You will find the current temperature at the top of each column.

Some practical advice:

The capillary is filled with a new non-polluting liquid. It is possible that the liquid column becomes separated or form bubbles, e.g. during transport. This is no manufacturing error! You can eliminate this by turning around the thermometer and pressing the button until the markers disappear. Then turn again and shake once strongly (similar to a fever thermometer). Repeat if necessary. Press the button to move down the markers to the column. Should one of the blue markers slip in the clear liquid, put it back in the red liquid by using a magnet (e.g. the magnetic stripe).

Avoid direct sunshine !