ANIOS CLEAN EXCEL D

DOMINIQUE DUTSCHER SAS

ANIOS'CLEAN EXCEL D

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Regulation and general information

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ANIOS CLEAN EXCEL D is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System

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ANIOS CLEAN EXCEL D is in compliance with the requirements of the European regulation related to biocides for disinfection of surfaces (Group 1 - PT 2.)

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ANIOS CLEAN EXCEL D formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments

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ANIOS CLEAN EXCEL D is labelled in accordance with the European regulation related to the classification and labelling of chemical products

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ANIOS CLEAN EXCEL D meets our voluntary eco-design commitment included in our ANIOSAFE charter

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Pure product ready to use :

- Storage between +5°C and +35°C
- Stability :
- Not opened bottle : 24 months from the manufacturing date indicated on the label.
- Opened bottle : maintenance of the 24 months in condition if closing the bottle after use.

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Diluted product

Dilution: 0.5% - 5mL per liter of water

Conservation of the diluted solution in a closed container (cap or spray gun) : 7 days.

Conservation of the diluted solution in an opened container : 24 hours (without use).

Any diluted preparation will be eliminated after use.

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ANIOS'CLEAN EXCEL D antimicrobial efficacy is tested against monobacterial biofilm with Pseudomonas aeruginosa, artificially developed inside a Tygon® tube, according to the ISO/TC 15883 – 5 : 2006 standard.

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After the passing of the 0.5% ANIOS'CLEAN EXCEL D diluted solution during 5 minutes, followed by rinsing with distilled water during 2 minutes, number of viable bacteria fixed by cm² of carrier and concentration of proteins are determined inside the Tygon® tube.

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Concentration in proteins determined in biofilm

	Concentration en protéines $\mu\text{g}/\text{cm}^2$ <i>Concentration in proteins $\mu\text{g}/\text{cm}^2$</i>	% de réduction <i>Reduction %</i>
Tube témoin <i>Control tube</i>	20,50	
ANIOS'CLEAN EXCEL D (0.5%)	2,80	<u>85%</u>

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
Present surfactants distinguish themselves by hydrophilic and hydrophobic properties, which ensures performances against proteic and glucidic soilings on one hand, and greasy soiling (lipidic) on another hand, by solubilisation and/or emulsification.

Dilution 0,5% - T° ambiante – Temps de contact de 15 min
Dilution 0,5% - room T° - Contact times of 15 min
Salissure : sang de Bœuf Citraté
Soil: Citrated bovine blood

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The pre-disinfectant detergent ANIOS'CLEAN EXCEL D, diluted at 0.5 % in tap water, presents a High Degreasing Cleaning Power against an organic soiling fixed on carriers made of brass (hydrophilic character).

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
It is demonstrated, by electrochemical study of pitting corrosion according to the NF S 94-402-1 standard (May 2004), that ANIOS'CLEAN EXCEL D does not present pitting corrosive character with stainless steel Z30 Cr13 in the conditions of use recommended by Laboratoires ANIOS.

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Not concerned by biodegradability notion

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*Readily biodegradable substances and/or inherently biodegradable*³

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Substances which not fulfil the biodegradability criteria, or substances without any biodegradability data.

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ANIOS'CLEAN EXCEL D contains at least 99.8 % of inorganic material and of biodegradable organic material.

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HAZARDS IDENTIFICATION AT DILUTION OF USE

Physicochemical hazards

Not classified at dilution of use in terms of physicochemical hazards.

Health

Not classified at dilution of use in terms of health hazards*.

Environment

Not classified at dilution of use in terms of environment hazards*.

* Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.

PERSONAL PROTECTION AT THE APPLICATION RATE UNDER THE CONDITIONS OF THE INTENDED USE

Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.

Respiratory protection

In case of insufficient ventilation, with the risk of exceeding the TWA/STEL**, wear suitable breathing apparatus (type A filtre).

Hand protection

When handling, wear suitable gloves (nitrile, latex or vinyle).
Gloves must be replaced immediately if they show signs of wear and tear.

Eye and face protection

Avoid all contact with eyes.
Water source nearby.

** please refer to SDS §8 for TWA (Time-Weighted Average)/STEL (Short-term exposure limit) of concerned substances present in the product

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The Safety Data Sheet (SDS) of this product is available on our Web site www.anios.com. In addition to visualization and/or impression, it is possible for you to download it.

The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.

In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)

* an access is also possible at the site <http://www.quickfds.com>

DOMINIQUE DUTSCHER SAS