# BAC MICROCHIRURGIE

## Protocole de traitement du bac

Bac de trempage autoclavable

### NETTOYAGE

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## Cleaning protocol of tank

### CLEANING

Bac Microchirurgie is designed so as to be cleaned in an **automatic washer machine** at **max temperature 90°C** (Machines to comply with ISO15883).

Bac Microchirurgie can be manually cleaned using suitable detergent or disinfectant (see protocol established in your department).

- 1 Separate the basket from the tank.
- 2 Position the components separately in the washer machine. Make sure the lid of the tank is open.

### STERILISATION

Bac Microchirurgie can be **sterilised in a steam autoclave** to a **temperature 134° during 18 minutes** (max temperature : 137°C) (Autoclaves to comply with ISO17665).

The packaging materials and methods as well as their positioning in the autoclave must be taken into consideration to avoid deformation.

- 1 Place the basket into the tank. Keep the lid of the box open.
- 2 Package up the whole thing with a sterilisation bag. Make sure the lid stays open.
- 3 Place the tank in the autoclave. Avoid all contact with other materials and products.

### STORAGE

Untreated tanks must be stored in their original packaging.

Secondary packaging is recommended, for example, a cardboard box.

The tanks must be stored **in a dry place** (minimum 4°C) **without direct contact** with the ground.

### REFERENCE

Bac Microchirurgie.....Ref. 404 010

#### RECOMMANDATIONS :

- On first use, perform a manual sterilisation with adapted disinfectant product.
  - A visual inspection is necessary to evaluate wear to the tank due to its use.
- If it is damaged, do not use the tank and report it.