



COMPLIANCE ENABLING ENVIRONMENTAL DISINFECTION

Incidin™ OxyFoam S

SPORICIDAL READY-TO-USE CLEANER AND DISINFECTANT WITH HI-SPEED H₂O₂™

The advanced formulation of Incidin™ OxyFoam S provides the benefits of a patented* hydrogen peroxide solution with broad spectrum efficacy.

They are effective against a wide range of pathogens and leave no active residues due to their rapid evaporation rate.

Incidin™ OxyFoam S provides:

- Hi-speed H₂O₂™
- Active ingredient decomposing into water and oxygen
- Sporicidal & virucidal offering
- Excellent material compatibility

*European Patent EP1987121



ECOLAB®

www.ecolab.com

LABORATOIRES
ANIOS
An Ecolab Company

www.anios.com

Incidin™ OxyFoam S

Sporicidal ready-to-use cleaner and disinfectant with Hi-speed H₂O₂™

EFFECTIVE

- Broad spectrum efficacy
- Short contact times

SAFE

- Active ingredient decomposes into water and oxygen
- Excellent material compatibility

CONVENIENT

- Cleans & disinfects in one step
- Ready to use formats

INSTRUCTIONS FOR USE



Wet surface



Wipe the surface ensuring it is completely wet



Let dry



Apply to clean wipe



Wipe the surface ensuring it is completely wet



Let dry

MICROBIOLOGICAL PROPERTIES

EFFICACY UNDILUTED	CONTACT TIME
Disinfection of surfaces acc. to VAH method 2015 (incl. EN 16615), dirty conditions	
Bactericidal, yeasticidal	5 min
Fungicidal	5 min
Mycobactericidal	15 min
Phase 2/Step 2 - 4 field-test	
EN 16615 (bactericidal)	1 min**
EN 16615 (yeasticidal)	1 min**
EN 16615 method (fungicidal)	5 min**
EN 16615 method (tuberculocidal, mycobactericidal)	5 min**
EN 16615 method acc. CEN TC216 prWI Oct 2020 / VAH method 19 (C. difficile)	15 min*
Phase 2/Step 1 - suspension tests	
EN 13727 (bactericidal)	1 min
EN 13624 (yeasticidal)	1 min
EN 13624 (fungicidal)	5 min
EN 14348 (tuberculocidal, mycobactericidal)	5 min
EN 14476 (efficacy against enveloped viruses)	30 sec
EN 14476 (limited spectrum virucidal)	1 min* / 15 min**
EN 14476 Adeno	30 sec
EN 14476 Noro	1 min* / 15 min**
EN 14476 virucidal	2 min* / 30 min**
EN 17126 (sporicidal)	60 min*
EN 17126 (C. difficile)	30min*/60 min**
Tested acc. DVV 2015	
Effective against enveloped viruses acc. to RKI recommendation 01/2004	15 sec
Limited spectrum virucidal & Norovirus (MNV)	30 min
Adeno	30 sec
Polyoma SV40	30 sec
Rota	15 sec

* clean conditions ** dirty conditions

PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Incidin™ OxyFoam S	6 x 750ml	2943E07

APPLICATION AREA



Fast acting, sporicidal ready-to-use surface cleaning and disinfection liquid

- acc. Directive 93 / 42 / EEC: for medical devices
- Biocide: for general surfaces

INSTRUCTIONS FOR USE

Incidin™ OxyFoam S: Wet surface, ensuring the surface is completely wet. Wipe the surface with a clean wipe and leave to dry.

Alternatively, apply the product to a clean wipe, wipe the surface and leave to dry. Do not rinse.

The use of a cleaning and disinfection product does not replace regular cleaning. Do not use on surfaces sensitive to oxidative agents such as marble, copper or brass. For reprocessing of medical devices always observe manufacturers' instructions for use, incl. material compatibility. Terminal disinfection of semi-critical medical devices (i.e. probes & TEE probes): If immersion is not possible, disinfect the device by wiping as described above. Rinse thoroughly with water of drinking quality in this case after contact time elapsed.

INGREDIENTS

100g contains the following active: 1.5g Hydrogen peroxide.

STORAGE

Product should be stored in the original container between 5°C and 25°C in a cool, dry location out of direct sunlight.

Shelf life after first opening: End of expiry date.

MEDICAL DEVICE AND BIOCIDES USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.

LABORATOIRES ANIOS

1, Rue de l'Espoir
59260 Lezennes
France



ECOLAB DEUTSCHLAND GmbH

Ecolab-Allee 1
40789 Monheim am Rhein
Germany



www.ecolab.com



www.anios.com