CATEGORY III

AQL 0.25, G1

User's Instructions

Catalogue number	Size	Scan to product page	Scan to chemical page
66 9251	Extra Small (XS/6)		
66 9252	Small (S/7)		
66 9253	Medium (M/8)		
66 9254	Large (L/9)		
66 9255	Extra Large (XL/10)		
66 9256	Extra Extra Large (XXL/11)		

Determination of material resistance to permeation by chemicals (EN 16523-1:2015) and degradation (EN 374-4:2013)

 ISO 374-1 Type A AJKLNPT	Code letter and chemical	CAS number	Permeation performance and breakthrough time	Degradation (DR) Mean in %
	A: Methanol 99%	67-56-1	Level 2, > 30 min.	59.6
	J: n-Heptane 99%	142-82-5	Level 2, > 30 min.	60.9
	K: Sodium hydroxide 40%	1310-73-2	Level 6, > 480 min.	-11.0
	L: Sulphuric acid 96%	7664-3-9	Level 2, > 30 min.	78.2
	N: Acetic acid 99%	64-19-7	Level 3, > 60 min.	61.0
	P: Hydrogen peroxide 30%	7722-84-1	Level 6, > 480 min.	-9.1
	T: Formaldehyde 37%	50-00-0	Level 6, > 480 min.	-5.0

Marking of gloves protecting against viruses, bacteria and fungi, as defined in ISO 374-5:2016 (Protective gloves against dangerous chemicals and micro-organisms – Part 5: Terminology and performance requirements for micro-organisms risks). Protection against viruses: tested according to ISO 16604:2004 (Procedure B) Phi-X, Herpes, HIV.

 ISO 374-5:2016 Level 3 VIRUS	Performance level	Acceptable quality level unit	Inspection level
	Level 3	< 0.65	G1
	Level 2	< 1.5	G1
	Level 1	< 4.0	S4

(en) Storage recommendations: Store in a cool and dry place. Gloves should be protected from exposure to direct sunlight, intense artificial light, X-ray machines and other sources of ozone. Disposal: Dispose of in accordance with local authority regulations. This is an MDR Class 1 examination glove according to Regulation (EU) 2017/745 and a PPE Category III protective glove according to Regulation (EU) 2016/425. Gloves provide protection against dangerous chemicals, micro-organisms and mechanical risks. Gloves shall not be worn when there is a risk of entanglement by moving parts of machines. Caution: These gloves should not be used for thermal protection. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Declaration of conformity accessible from: <http://shieldscientific.com/resource-center/glove-declaration-of-conformity>.

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Website: www.shieldscientific.com • E-mail: info@shieldscientific.com • ™: trademark of SHIELD Scientific.
All rights reserved. © Copyright 2020.
EC REP SHIELD Scientific B.V.
SGS Fimko Oy (Notified Body No: 0598) Takomotie 8, FI-00380 Helsinki, Finland
SHIELD Scientific B.V. • Dr. Willem Dreeslaan 1 • 6721 ND Bennekom • The Netherlands • Tel: +31(0)317700202