

Omni Sonic Ruptor 4000

Ultrasonic Homogenizer

The **Omni Sonic Ruptor 4000 Ultrasonic Homogenizer** offers precision engineering with all the features necessary to create a total system for ultrasonic disruption. It can disintegrate most cells, bacteria, spore or tissue. It can prepare an emulsion down to 1/100 of a micron, homogenize "immiscible" liquids, accelerate enzymatic and chemical reactions, stimulate bacterial activity, disperse solids in liquids and degas liquids. The Omni Sonic Ruptor 4000 includes an integrated **Sound Abating Chamber** to reduce cavitation sound emitted during processing.

How it works


The ultrasonic power supply transforms line voltage to high frequency 20 kHz of electrical energy. This electrical energy is transmitted to the probe where it is changed to mechanical energy. The vibrations from the probe are coupled to and intensified by the titanium tip. The probe vibrates in a longitudinal direction and transmits this motion to the titanium tip immersed in the solution. Cavitation results, in which microscopic vapor bubbles are formed momentarily and implode, causing powerful infinitesimal shock waves to radiate throughout the solution in proximity to the radiating face of the tip.

Automatic Tuning

The Omni Sonic Ruptor employs a proprietary feedback system, insuring that the Ultrasonic Homogenizer is always working at its maximum efficiency regardless of the application. There is no need for constant adjustment by the operator. The percentage of ultrasonic power emitted is indicated by the output meter, enabling accurate, reproducible results.

Performance

The high powered instrument provides continuous duty operation on difficult or large scale applications. It delivers up to 300 watts to the probe and tip.



> Specifications

Variable Power: 0 - 400 watts

Timer: 0 - 15 minutes

Automatic Tuning: Allows easy interchange of probe tips

Pulse Mode Power: Allows sample to be processed at full power while limiting temperature rise. Essential for heat sensitive samples.

Constant Power Mode: Selectable for non-heat critical sample applications.

Dimensions:

Control Unit: 9 5/8" (l) x 11 1/8" (w) x 6" (h)
24.45 cm (l) x 28.26 cm (w) x 15.24 cm (h)

Transducer: 3" (w) x 5 5/8" (h)
7.62 cm (w) x 14.29 cm (h)

Temperature/Humidity Range:
0-40°C / 28% - 80%

Power Consumption: 400 watts

Omni Sonic Ruptor

Ultrasonic Homogenizer Accessories

Accessory Tips



PRODUCT	INTENSITY	TIP SIZE	VOLUME RANGE	ORDER NO.
Stepped Micro Processing Tip	Very High	Diameter: 5/32" (3.8 mm) Length: 10.1" (25.6 cm)	250 µL - 10 mL	OR-T-156
Intermediate Processing Tip	High	Diameter: 3/8" (9.5 mm) Length: 8.6" (21.8 cm)	10 mL - 250 mL	OR-T-375
1/2" Processing Tip	Medium-High	Diameter: 1/2" (12.7 mm) Length: 5.38" (13.65 cm)	10 mL - 300 mL	OR-T-500
Standard Processing Tip	Medium	Diameter: 3/4" (19 mm) Length: 4.1" (10.5 cm)	25 mL - 500 mL	OR-T-750
Full Size Processing Tip	Low	Diameter: 1" (25.4 mm) Length: 4.85" (12.3 cm)	50 mL - 1 L	OR-T-1000



ORDER NO. OR-A-01

Sound Abatement Chamber

The **Omni Sonic Ruptor 4000** has an **Integrated Sound Abatement Chamber** which reduces cavitation sound emitted during processing. Harmonics are produced by the vessel walls and fluid surface and can be disturbing with extended operation. The clear Plexiglas door permits viewing of the sample while protecting the operator against accidental splashing. An access port for tubing is also provided for use with the accessory cooling jacket or continuous flow chamber.



Continuous Flow Chamber

ORDER NO. OR-F-01

This chamber permits continuous processing of liquids and is designed for emulsifying and homogenizing applications. The continuous flow chamber may be sealed in a closed system when mounted so that infectious materials can be used. It is equipped with a water flow through jacket to enable the temperature of the processed solution to be maintained at a desired level.

Cup Tip, 250 mL

ORDER NO. OR-C-250

Used to process contents of sealed tubes, vials or ampoules. Samples can also be processed by placing the sample in a beaker, and the beaker is suspended in the water filled cup tip.



> For more information —
omni-inc.com

Omni International, Inc.
935C Cobb Place Blvd.
Kennesaw, GA 30144
800.776.4431 • 770.421.0058

