Thermo
S C I E N T I F I C

**Cimarec+ stirrers, hotplates,
and stirring hotplates**

OPERATING MANUAL
AND PARTS LIST

CIC0000809 V17 02/22/17

This manual covers the models shown below

Model	Voltage	Description
HP88854100	100-120V	4x4 Ceramic Hot Plate
HP88854105	220-240V	4x4 Ceramic Hot Plate
S88854100	100-120V	4x4 Ceramic Stirrer
S88854105	220-240V	4x4 Ceramic Stirrer
SP88854100	100-120V	4x4 Ceramic Stirring Hotplate
SP88854105	220-240V	4x4 Ceramic Stirring Hotplate
HP88857100	100-120V	7x7 Ceramic Hot Plate
HP88857105	220-240V	7x7 Ceramic Hot Plate
S88857100	100-120V	7x7 Ceramic Stirrer
S88857105	220-240V	7x7 Ceramic Stirrer
SP88857100	100-120V	7x7 Ceramic Stirring Hotplate
SP88857105	220-240V-EU	7x7 Ceramic Stirring Hotplate
SP88857106	220-240V-AUS/CHN	7x7 Ceramic Stirring Hotplate
SP88857107	220-240V-UK	7x7 Ceramic Stirring Hotplate
HP88857104	100-120V	7x7 Aluminum Hot Plate
HP88857108	220-240V	7x7 Aluminum Hot Plate
S88857104	100-120V	7x7 Aluminum Stirrer
S88857108	220-240V	7x7 Aluminum Stirrer
SP88857104	100-120V	7x7 Aluminum Stirring Hotplate
SP88857108	220-240V	7x7 Aluminum Stirring Hotplate
HP88850100	100-120V	10x10 Ceramic Hot Plate
HP88850105	220-240V	10x10 Ceramic Hot Plate
S88850100	100-120V	10x10 Ceramic Stirrer
S88850105	220-240V	10x10 Ceramic Stirrer
SP88850100	100-120V	10x10 Ceramic Stirring Hotplate
SP88850105	220-240V-EU	10x10 Ceramic Stirring Hotplate
SP88850106	220-240V-AUS/CHN	10x10 Ceramic Stirring Hotplate
SP88850107	220-240V-UK	10x10 Ceramic Stirring Hotplate

Do You Need Information or Assistance on
Thermo Scientific Products

If you do, please contact us 8:00 a.m. to 6:00 p.m. (Eastern Time) at:

1-740-373-4763

Direct

1-800-438-4851

Toll Free, U.S. and Canada

1-877-213-8051

FAX

<http://www.thermoscientific.com>

Internet Worldwide Web Home Page

service.led.marietta@thermofisher.com

Tech Support Email Address

www.unitylabservices.com

Certified Service Web Page

Our Sales Support staff can provide information on pricing and give you quotations. We can take your order and provide delivery information on major equipment items or make arrangements to have your local sales representative contact you. Our products are listed on the Internet and we can be contacted through our Internet home page.

Our Service Support staff can supply technical information about proper setup, operation or troubleshooting of your equipment. We can fill your needs for spare or replacement parts or provide you with on-site service. We can also provide you with a quotation on our Extended Warranty for your Thermo Scientific products.

Whatever Thermo Scientific products you need or use, we will be happy to discuss your applications. If you are experiencing technical problems, working together, we will help you locate the problem and, chances are, correct it yourself...over the telephone without a service call.

When more extensive service is necessary, we will assist you with direct factory trained technicians or a qualified service organization for on-the-spot repair. If your service need is covered by the warranty, we will arrange for the unit to be repaired at our expense and to your satisfaction.

Regardless of your needs, our professional telephone technicians are available to assist you Monday through Friday from 8:00 a.m. to 6:00 p.m. Eastern Time. Please contact us by telephone or fax. If you wish to write, our mailing address is:

Thermo Fisher Scientific
401 Millcreek Road, Box 649
Marietta, OH 45750

International customers, please contact your local Thermo Scientific distributor.

Table of Contents

Safety Information.....	5
Alert Signals	5
Warnings	5
Información de seguridad.....	7
Señales de alerta	7
Advertencias.....	8
Sicherheitshinweise	9
Gefahrensymbole	9
Warnungen.....	9
Informations de sécurité.....	11
Signaux d'alerte.....	11
Avertissements.....	11
Informação de Segurança.....	13
Sinais de Aviso	13
Avisos.....	13
Информация по безопасности	15
Сигналы тревоги.....	15
Предупреждения.....	15
Specifications	17
Ceramic Hot Plate Specifications	17
Ceramic Stirrer Specifications	17
Ceramic Stirring Hotplate Specifications.....	18
Aluminum Top Unit Specifications	19
Heating Specifications*	20
Stirring Speed Specifications	20
Environmental Conditions.....	21
Declaration of Conformity	21
Introduction.....	22
General Usage	22
Principles of Operation	23
General Cleaning Instructions	23
Unpacking and Installation.....	24
Unpacking	24
Installation	24
Operation.....	25
Power Switch	25
Setting the Stirring Speed	25
Setting the Temperature	26
Controlling Solution Temperature Using External Probe.....	27
Heating Small Volumes.....	27
Heating Metal Vessels and Sand Baths.....	27
Troubleshooting	28
Error Codes.....	28
Accessories	29
THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY	30

Safety Information

Alert Signals

Warning

Warnings alert you to a possibility of personal injury.

Caution

Cautions alert you to a possibility of damage to the equipment.

Note

Notes alert you to pertinent facts and conditions.

Hot Surface

Hot surfaces alert you to a possibility of personal injury if you come in contact with a surface during use or for a period of time after use.

Your Thermo Scientific Cimarec[®] Hot Plate, Stirrer or Stirring Hot Plate has been designed with function, reliability, and safety in mind. It is your responsibility to install it in conformance with local electrical codes. For safe operation, please pay attention to the alert signals throughout the manual.

Warning: These products should be used only under the operating conditions specified in the Operating Manual.

- If the equipment is used in a manner not specified by the manufacturer, protection provided by the equipment may be impaired
- Always use safe laboratory practices and do not leave the hotplate in operation while unattended as product functionality or laboratory practice failures could occur that might lead to uncontrolled or excessive heating of the top surface.
- Safety procedures (including, but not limited to, unplugging when not in use) and emergency response plans should be put in place to address the worst case possibility.
- If an over-temperature failure occurs, the top surface temperature could rise to the maximum temperature (300-540°C depending on your model's specification) and remain at that temperature indefinitely. Under these conditions, the material being heated on the surface of the hotplate could reach levels in excess of the maximum temperature.

This manual contains important operating and safety information. The user must carefully read and understand the contents of this manual prior to the use of this equipment.

Warnings

To avoid electrical shock, always:

1. Use a properly grounded electrical outlet of correct voltage and current handling capacity.

SAFETY INFORMATION

2. Disconnect from the power supply prior to maintenance and servicing.

To avoid personal injury:

1. Do not use in the presence of flammable vapors outside of a fume hood or near combustible materials where heating could result in a fire or explosion. This device contains components which may ignite such materials. This device is not rated for use in designated hazardous locations where flammable vapors may reach the lower flammability limit.
2. Use caution when heating volatile materials; top surface and element can reach the "Flash Point Temperature" of many chemicals. These hot plates are not explosion proof. Fire or explosion may result. Unit contains components which may ignite such materials.
3. Keep top surface clean. Use a non-abrasive cleaner. Alkali spills, hydrofluoric acid spills or phosphoric acid spills may damage top and lead to thermal failure. Unplug unit and remove spills promptly. Do not immerse unit for cleaning.
4. Replace the top immediately if damaged by etching, scratching or chipping. A damaged top can break in use.
5. Do not use metal foil on hot plate which may block air flow. Overheating will result.
6. Check and tighten the removable cord periodically making sure it is secure. If loosened, the cord could become hot and/or spark and be a potential fire hazard. If cord appears damaged, replace immediately. If cord is repeatedly loosened it is recommended to purchase the power cord.
7. Do not remove or modify grounded power plug. Use only properly grounded outlets to avoid shock hazard.
8. Use appropriate hand and eye protection when handling hazardous chemicals.
9. Gross weight of items placed on top of hot plates should not exceed 35 lbs. (15.9 kg.) on the 10" x 10" models, 25 lbs. (11.3 kg) on the 7" x 7" models and 15 lbs. (6.8 kg) on the 4" x 4" models.
10. "Caution: Hot Top. Avoid Contact." The top plate of the unit can remain hot for some time after use. A "CAUTION - HOT TOP" light will remain on until top plate temperature cools to below 50°C.
11. Do not leave an active probe out of the fluid. This may cause uncontrolled heating of the fluid on the hot plate and unintentional boiling or an explosion could occur.
12. Do not heat flammable liquids outside of a fume hood. When heating flammable liquids, heater element temperature can be significantly higher than the temperature indicated on the display. If flammable vapors concentrations reach the internal heating element a fire or explosion may result.
13. Note that the exterior housing will be hot during and for a period of time after use.
14. Refer servicing to qualified personnel.

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Specifications

Ceramic Top Hot Plate Specifications

Overall Dimensions	HP88854100	HP88854105	HP88857100	HP88857105	HP88850100	HP88850105
Width - in.(cm)	5.1 (13)	5.1 (13)	8.2 (20.8)	8.2 (20.8)	11.3 (28.8)	11.3 (28.8)
Height - in.(cm)	3.9 (9.8)	3.9 (9.8)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)
Depth - in.(cm)	10.2 (25.9)	10.2 (25.9)	14.2(36)	14.2(36)	17.2(43.8)	17.2(43.8)
Weight - lbs.(kg)	3.3 (1.5)	3.3 (1.5)	7.5 (3.4)	7.5 (3.4)	11.5 (5.2)	11.5 (5.2)
Top Plate Dimensions						
Width - in.(cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Height - in.(cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in.(cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Volts	100-120	220-240	100-120	220-240	100-120	220-240
Amps	4.4	1.6	12	4.4	13	6.0
Watts	540	384	1440	1056	1560	1440
Freq.	50/60	50/60	50/60	50/60	50/60	50/60
Phase	1	1	1	1	1	1
Max. Temp °C (°F)	540°C(1004°F)	540°C(1004°F)	540°C(1004°F)	540°C(1004°F)	400°C(752°F)	400°C(752°F)

Ceramic Top Stirrer Specifications

Overall Dimensions	S88854100	S88854105	S88857100	S88857105	S88850100	S88850105
Width - in. (cm)	5.1 (13)	5.1 (13)	8.2 (20.8)	8.2 (20.8)	5.1 (13)	5.1 (13)
Height - in. (cm)	3.9 (9.8)	3.9 (9.8)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)
Depth - in. (cm)	10.2 (25.9)	10.2 (25.9)	14.2(36)	14.2(36)	17.2(43.8)	17.2(43.8)
Weight - lbs. (kg)	5.0 (2.2)	5.0 (2.2)	9.0 (4.1)	9.0 (4.1)	13.0 (5.9)	13.0 (5.9)
Top Plate Dimensions						
Width - in. (cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Volts	100-120	220-240	100-120	220-240	100-120	220-240
Amps	0.2	0.1	0.2	0.1	0.2	0.1
Watts	24	24	24	24	24	24
Freq.	50/60	50/60	50/60	50/60	50/60	50/60
Phase	1	1	1	1	1	1

SPECIFICATIONS

Ceramic Top Stirring Hotplate Specifications

Overall Dimensions	SP88854100	SP88854105	SP88857100	SP88857105	SP88850100	SP88850105
Width - in. (cm)	5.1 (13)	5.1 (13)	8.2 (20.8)	8.2 (20.8)	11.3 (28.8)	11.3 (28.8)
Height - in. (cm)	3.9 (9.8)	3.9 (9.8)	4.1(10.5)	4.1(10.5)	4.1(10.5)	4.1(10.5)
Depth - in. (cm)	10.2 (25.9)	10.2 (25.9)	14.2 (36)	14.2 (36)	17.2 (43.8)	17.2 (43.8)
Weight - lbs. (kg)	5.0 (2.2)	5.0 (2.2)	9.0 (4.1)	9.0 (4.1)	13.0 (5.9)	13.0 (5.9)
Top Plate Dimensions						
Width - in. (cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	4.25 (10.8)	4.25 (10.8)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Volts	100-120	220-240	100-120	220-240	100-120	220-240
Amps	4.4	1.6	12	4.4	13	6.0
Watts	540	384	1440	1056	1560	1440
Freq.	50/60	50/60	50/60	50/60	50/60	50/60
Phase	1	1	1	1	1	1
Max. Temp °C (°F)	540°C(1004°F)	540°C(1004°F)	540°C(1004°F)	540°C(1004°F)	400°C(752°F)	400°C(752°F)

Overall Dimensions	SP88857106	SP88857107	SP88850106	SP88850107
Width - in. (cm)	8.2 (20.8)	8.2 (20.8)	11.3 (28.8)	11.3 (28.8)
Height - in. (cm)	4.1(10.5)	4.1(10.5)	4.1(10.5)	4.1(10.5)
Depth - in. (cm)	14.2 (36)	14.2 (36)	17.2 (43.8)	17.2 (43.8)
Weight - lbs. (kg)	9.0 (4.1)	9.0 (4.1)	13.0 (5.9)	13.0 (5.9)
Top Plate Dimensions				
Width - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Volts	100-120	220-240	100-120	220-240
Amps	12	4.4	13	6.0
Watts	1440	1056	1560	1440
Freq.	50/60	50/60	50/60	50/60
Phase	1	1	1	1
Max. Temp °C (°F)	540°C(1004°F)	540°C(1004°F)	400°C(752°F)	400°C(752°F)

Aluminum Top Unit Specifications

Overall Dimensions	HP88857104	HP88857108	S88857104	S88857108	SP88857104	SP88857108
Width - in. (cm)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)
Height - in. (cm)	4.1(10.5)	4.1(10.5)	4.1(10.5)	4.1(10.5)	4.1(10.5)	4.1(10.5)
Depth - in. (cm)	14.2 (36)	14.2 (36)	14.2 (36)	14.2 (36)	14.2 (36)	14.2 (36)
Weight - lbs. (kg)	7.5 (3.4)	7.5 (3.4)	9.0 (4.1)	9.0 (4.1)	9.0 (4.1)	9.0 (4.1)
Top Plate Dimensions						
Width - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Volts	100-120	220-240	100-120	220-240	100-120	220-240
Amps	6.0	2.6	0.2	0.1	6.0	2.6
Watts	720	624	24	24	720	624
Freq.	50/60	50/60	50/60	50/60	50/60	50/60
Phase	1	1	1	1	1	1
Max. Temp °C (°F)	300°C(572°F)	300° C(572° F)	-	-	300°C(572°F)	300° C(572° F)

GENERAL SPECIFICATIONS

Heating Specifications*

Top Plate Surface - Solid Ceramic

Models	4" x 4"	7" x 7"	10" x 10"
Temperature range	30°C - 540°C (86°F - 1004°F)**		30°C - 400°C (86°F - 752°F)**
Heat-up time to maximum temperature (unloaded top plate)	8 minutes		
Accuracy of the temperature display vs. the actual average temperature of a 2" diameter of setting area at the center of the top plate (setpoint 100°C unloaded)	± 5.0°C		
Temperature stability at the center of the top plate surface (@ 100°C unloaded)	± 2.0°C		
Typical accuracy of remote probe at user selected calibration temperature after calibration procedure	±1.8°C		
Temperature stability using remote probe (500 ml of water in a 1000 ml flask at 70°C)	±0.5°C		

Top Plate Surface - Solid Aluminum

Models	7" x 7"
Temperature range	30°C - 300°C (86°F - 572°F)**
Heat-up time to maximum temperature (unloaded top plate)	10 minutes
Accuracy of the temperature display vs. the actual average temperature of a 2" diameter of setting area at the center of the top plate (setpoint 100°C unloaded).	± 10.0°C
Temperature stability at the center of the top plate surface (@ 100°C unloaded)	± 2.0°C
Typical accuracy of remote probe at user selected calibration temperature after calibration procedure	± 1.8°C
Temperature stability using remote probe (500ml of water in a 1000 ml flask at 70°C)	± 0.5°C

Stirring Speed Specifications

Speed Range 50 to 1500 RPM(Maximum speed is dependent on the viscosity of the solution).

Stability of the stirring speed setpoint (600 ml of water in a 1000 ml glass flask) ± 5.0% at 1000 RPM

Top Plate Size	Max Recommended Flask Size	Max Weight on Top Plate
4" x 4"	1 liter	15 lbs
7" x 7"	4 liters	25 lbs
10" x 10"	6 liters	35 lbs

* The specification defined in this section is derived from testing data statistics of Thermo Fisher Scientific lab. We reserve the right to make changes in design and specification without prior notice.

** As the top plate becomes dirty, the maximum temperature will decrease. To return the unit to its maximum temperature performance, use a mild abrasive to remove stained areas.

Environmental Conditions

Operating:	5°C to 35°C; 20% to 80% relative humidity, non-condensing. Installation category II (overvoltage) in accordance with IEC 664. Pollution degree 2 in accordance with IEC 664. Altitude Limit: 2,000 meters.
Storage:	-25°C to 65°C 10% to 85% relative humidity

Declaration of Conformity

We hereby declare under our sole responsibility that this product conforms with the technical requirements of the following standards:

EMC:	EN 61000-3-2	Limits for harmonic current emissions
	EN 61000-3-3	Limits for voltage fluctuations and flicker
	EN 61326-1	Electrical equipment for measurement, control, and laboratory use; Part I: General Requirements
Safety:	EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use; Part I: General Requirements
	EN 61010-2-010	Part II: Particular requirements for laboratory equipment for the heating of materials
	EN 61010-2-051	Part II: Particular requirements for laboratory equipment for mixing and stirring

per the provisions of the Electromagnetic Compatibility Directive 2014/30/EU, and per the provisions of the Low Voltage Directive 2014/35/EU.

Copies of the Declaration of Conformity are available upon request.

Introduction

Cimarec⁺ Stirring Hotplate

Please read all the information in this manual before operating the unit.

Your Cimarec⁺ hot plate, stirrer or stirring hot plate is a heating and/or stirring plate designed for laboratory procedures requiring precise control of temperature and/or stirring speed. Each Cimarec⁺ model includes a digital temp display and/or a digital speed display for monitoring actual temperature and/or stirring speed.

The hot plate is capable of producing accurately controlled top plate temperatures from 30°C through 300°C, 400°C or 540°C (depending on the model). The temperature is controlled at the plate surface by an internal sensor.

The stirrer will accurately maintain stirring speeds from 50 rpm up to 1500 rpm. The top plate on the Cimarec⁺ unit is solid ceramic or solid aluminum, and is suitable for use with glass or metal vessels.

Your Cimarec⁺ model may be used for general purpose heating applications and/or general laboratory mixing of solutions, including sample preparation, heating reagents, melting paraffin, warming resinous chemicals, content analysis, solvent evaporations, digestions, media preparation and sterilization, titrations, sand baths, and microscale chemistry application.

General Usage

Do not use this product for anything other than its intended usage.

Principles of Operation

Each Cimarec⁺ unit utilizes the latest microprocessor technology to deliver a reliable, controlled, ceramic top stirring hot plate.

Your Cimarec⁺ stirring hot plate has an electronic closed-loop feedback control which will accurately maintain temperature setpoints from 30°C through 300°C, 400°C or 540°C, depending on the model. Top plate surface temperatures are tested using the latest infrared temperature measurement. The measurements are made with an infrared camera mounted approximately 26" above the top surface of the stirring hot plate. If the temperature measurement of the ceramic top is made with measurement devices other than infrared, the error of the measuring technique may be greater than the error of the unit.

The electronic stirring speed control will maintain the speed setpoint when the unit is loaded. The motor in the Cimarec⁺ stirring hot plates produces maximum stirring torque under normal laboratory load conditions, and is combined with a powerful magnet to provide exceptional magnetic coupling with a stir bar.

General Cleaning Instructions

Wipe exterior surfaces with lightly dampened cloth containing mild soap solution.

Unpacking and Installation

Warning

Use a properly grounded electrical outlet of correct voltage and current handling capacity.

Unpacking

Remove your Cimarec⁺ unit from the carton. Inspect to ensure that the unit has not been damaged during shipment. If the unit appears to have sustained shipping damage contact the distributor from whom you purchased this product or Customer Service at 800-553-0039. Check for stir bar and thumbscrew prior to discarding packaging.

The following items are included in the shipment:

Cimarec⁺ Unit
Cord
Stir Bar (Stirring models only)
Knob
Operating Manual

If any of these items are missing from the carton, contact Customer Service.

Installation

Set the unit on a flat stable surface at least 12" away from combustible materials, and plug the cord set into a properly grounded electrical outlet of correct voltage and current handling capacity.

Operation

Warning

Do not heat volatile materials outside of a fume hood. The top surface and heating element can reach the autoignition temperature of some materials. A fire or explosion may result if this warning is ignored. This stirring hot plate is not explosion proof.


Use appropriate hand and eye protection when handling hazardous chemicals.

“Caution: Hot Top. Avoid Contact.” The top plate of the unit can remain hot for some time after use. A “CAUTION - HOT TOP” light will remain on until top plate temperature cools to below 50°C.

Caution

To avoid damage to the top plate or heating element, always keep a vessel filled with liquid on the top plate of a stirring hot plate when the unit is heating or cooling.

Note

For 4" models - SP88854100 and SP88854105, the speed is indicated by the marker on the membrane.

Power Switch

A power switch, located on the left side, is used to power the equipment on and off.

When the power switch is turned on, the unit briefly displays the following according to this order: 1) type of unit “SP”, “HP” or “S” *, 2) software version and 3) current temperature and RPM.

It is recommended to unplug the unit when not in use for 24 hours or longer.

Setting the Stirring Speed

Your Cimarec+ stirring hot plate has an electronic feedback speed control which will maintain a speed set-point from 50 rpm through 1500 rpm at $\pm 5.0\%$. (Maximum speed is dependent on viscosity of the solution.)

To set the speed, turn the STIR knob to your desired setting. Turn the knob clockwise to increase the speed or counter-clockwise to decrease the speed with 10 rpm increments. To discontinue stirring, turn knob counter-clockwise to OFF.

* "SP" stands for *Stirring Hotplate models*.

"HP" stands for *Hot Plate models*.

"S" stands for *Stirrer models*.

Note

The solution temperature is approximately 25% cooler than the hot plate surface temperature.

Note

If you allow the top plate to reach the maximum temperature of 540°C while preheating and then turn the control down to a setpoint less than 200°C, the temperature of the top will drop rapidly to 200°C. Because of the natural cooling characteristics of ceramic, the temperature of the top will drop much more gradually after the top plate temperature reaches 200°C.

Note

Boiling times are dependent on solution volume and the surface area of the flask that is exposed to the hot plate. For example, when heating the same amount of solution in a 2L flask vs. a 1L flask, the solution will heat about 20% faster.

Setting the Temperature

Your Cimarec⁺ stirring hot plate has an electronic closed-loop feedback control which will accurately maintain temperature setpoints in 1° increments from 30°C through 300°C, 400°C or 540°C, depending on the model. An unloaded hot plate will heat to maximum temperature in just 8 minutes. The temperature is controlled at the top plate by the internal sensor. A “CAUTION - HOT TOP” light on the front panel will illuminate whenever the top surface temperature exceeds 50°C.

Your Cimarec⁺ stirring hot plate will display the temperature in °C. When choosing a setpoint, the display will indicate the setpoint for 2 seconds, after which the display will flash until the desired setpoint is reached.

When plugging in the remote probe, turn the temperature knob to “OFF”, then adjust new temperature setpoint Value. The probe light will illuminate when the remote probe is plugged in. “CAUTION - HOT TOP” light will flash until the top surface has cooled to below 50°C.

Controlling Solution Temperature Using External Probe

To control the solution temperature plug an external probe into the probe receptacle located on the back side of the unit. When the probe is connected, the adjustable temperature range is 30-250 °C. Place the probe into the solution. The display will indicate the actual temperature of the solution as measured by the probe. When the probe is connected, the adjustable temperature range is 30-250 °C.

The external probe offers more exact temperature control than regulating the top plate by the internal sensor. If you need to maintain a precise setpoint it is recommended to use a probe to control the solution temperature instead of controlling by the top plate temperature.

When using a probe with the Cimarec⁺ unit it is recommended that a clamp on a support rod be used to hold the probe in the solution.

To ensure accurate probe readings, as much of the probe sheath as possible should be immersed in the solution. Make sure the probe is immersed in the liquid and is not located in air or outside of the solution. If the probe is plugged into the hot plate, but is not in solution while the heat control is operating, the temperature display will continue to indicate an ambient temperature, and a probe out of solution error will occur. Because the set point cannot be reached the heating element will continue to supply heat to the top plate, and the maximum top plate temperature (300°C, 400°C or 540°C) may be reached. If the remote probe does not

sense a temperature change in three minutes the unit will display E03 and shut down.

Heating Small Volumes

Preheating small volumes is not necessary as it may cause the temperature to overshoot the desired setpoint.

Heating Metal Vessels and Sand Baths

Cimarec⁺ is capable of precisely regulating the top plate temperature, metal vessels and sand baths may be heated safely without the danger of the ceramic top breaking. Use the lowest temperature setting possible for applications to limit thermal stress to the ceramic top. Using a metal vessel or sand bath may reduce the life of the stirring hot plate.

Note

It is not allowed to plug in the external probe while the temperature knob is ON. Please turn the temperature knob to OFF state while plugging in the external knob. Otherwise, the display will blink and the heating element will be turned off until reset.

Troubleshooting

Error Codes

Errors E01 through E04 are heating errors. Error Handler will lock out heating functions if heating error is detected. Stirring functionality is unaffected. Error E05 is stirring error. Error Handler will lock out stirring functions if stirring error is detected. Heating functionality is unaffected. Error E08 is failure of temperature of speed knob, Error Handler will lock out both heating and stirring functions if it is detected. If the condition that caused the error is no longer present, for E01, E03 (out of range) and E05, the unit will clear the errors automatically, for E03 (out of solution), rotating temperature knob to OFF position or turning off the switch on the left side of the unit will clear the error, and for E02, E04 and E08, only unplugging the unit will clear the errors.

Displayed Message	Intended to Detect	Cause	Solution
E01	Internal thermocouple out of range.	Internal thermocouple not connected. Thermocouple open. Thermocouple connected backwards (reversed polarity).	Call technical service
E02	Excessive top heat-up time.	Internal thermocouple short circuit. Failure in internal thermocouple. Failure in heating element. Failure in optocoupler/triac circuit.	Call technical service
E03	External probe left out of solution or out of range.	External probe left out of solution. Failure in external probe.	Place external probe into solution. Replace external probe.
E04	Over temperature condition detected by electronic circuit.	Internal thermocouple temperature is above the Over Temperature Protection triggering value .	Call technical service
E05	Motor system failure.	Locked rotor condition. Failure in motor. Failure in motor circuit. Failure in motor sensor.	Free locked rotor.
E08	Knob failure	Failure in temperature knob or speed knob.	Call technical service

Accessories

Order No	Description
CIC0001444	Splash guard protection shield 4", for 4" x 4" Stirrers
CIC0001445	Splash guard protection shield 7", for 7" x 7" Stirrers
CIC0001446	Splash guard protection shield 10", for 10" x 10" Stirrers
CIC0001447	PT100 External probe, Stainless steel, for all sizes of Hotplates and Hotplate Stirrers
CIC0001448	PT1000 External probe, Stainless steel, for all sizes of Hotplates and Hotplate Stirrers
CIC0001449	Hastelloy PT1000 External probe, Hastelloy, for all sizes of Hotplates and Hotplate Stirrers
CIC0001830	Supporting rod, Aluminum Alloy, for Hotplates and Hotplate Stirrers
CIC0001831	Clamp holder, Stainless steel, for Hotplates and Hotplate Stirrers
CIC0001468	Power cord assembly China plug 10A 250V
CIC0001469	Power cord assembly EURO plug 10A 250V
CIC0001470	Power cord assembly British BS 1363 plug 10A 250V
CIC0001471	Power cord assembly US&JP plug 16A 125V

THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first three (3) years, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo Fisher Scientific postage paid and replacement parts are shipped FOB destination.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY. Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Department at 1-800-438-4851 or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.

WEEE Compliance

WEEE Compliance. This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96EC. It is marked with the following symbol. Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, and this product should be disposed of or recycled through them. Further information on our compliance with these Directives, the recyclers in your country, and information on Thermo Scientific products which may assist the detection of substances subject to the RoHS Directive are available at www.thermo.com/

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DOMINIQUE DUTSCHER SAS