Thermo
S C I E N T I F I C

SUPER-NUOVA⁺
Hot Plates, Stirrers and
Stirring Hotplates

OPERATION MANUAL
AND PARTS LIST

CIC0000810 -V14 ·10/29/15

This manual covers the models shown below

Model	Voltage	Description
HP88857190	100-120V	7x7 Ceramic Hot Plate
HP88857195	220-240V	7x7 Ceramic Hot Plate
S88857190	100-120V	7x7 Ceramic Stirrer
S88857195	220-240V	7x7 Ceramic Stirrer
SP88857190	100-120V	7x7 Ceramic Stirring Hotplate
SP88857195	220-240V	7x7 Ceramic Stirring Hotplate
HP88857194	100-120V	7x7 Aluminum Hot Plate
HP88857198	220-240V	7x7 Aluminum Hot Plate
SP88857194	100-120V	7x7 Aluminum Stirring Hotplate
SP88857198	220-240V	7x7 Aluminum Stirring Hotplate
HP88850190	100-120V	10x10 Ceramic Hot Plate
HP88850195	220-240V	10x10 Ceramic Hot Plate
S88850190	100-120V	10x10 Ceramic Stirrer
S88850195	220-240V	10x10 Ceramic Stirrer
SP88850190	100-120V	10x10 Ceramic Stirring Hotplate
SP88850195	220-240V	10x10 Ceramic Stirring Hotplate

Do You Need Information or Assistance on
Thermo Scientific Products?

If you do, please contact us 8:00 a.m. to 6:00 p.m. (Eastern Time) at:

1-740-373-4763

Direct

1-800-438-4851

Toll Free, U.S. and Canada

1-877-213-8051

FAX

<http://www.thermoscientific.com>

Internet Worldwide Web Home Page

service.led.marietta@thermofisher.com

Tech Support Email Address

www.unitylabservices.com

Certified Service Web Page

Our Sales Support staff can provide information on pricing and give you quotations. We can take your order and provide delivery information on major equipment items or make arrangements to have your local sales representative contact you. Our products are listed on the Internet and we can be contacted through our Internet home page.

Our Service Support staff can supply technical information about proper setup, operation or troubleshooting of your equipment. We can fill your needs for spare or replacement parts or provide you with on-site service. We can also provide you with a quotation on our Extended Warranty for your Thermo Scientific products.

Whatever Thermo Scientific products you need or use, we will be happy to discuss your applications. If you are experiencing technical problems, working together, we will help you locate the problem and, chances are, correct it yourself...over the telephone without a service call.

When more extensive service is necessary, we will assist you with direct factory trained technicians or a qualified service organization for on-the-spot repair. If your service need is covered by the warranty, we will arrange for the unit to be repaired at our expense and to your satisfaction.

Regardless of your needs, our professional telephone technicians are available to assist you Monday through Friday from 8:00 a.m. to 6:00 p.m. Eastern Time. Please contact us by telephone or fax. If you wish to write, our mailing address is:

Thermo Fisher Scientific
401 Millcreek Road, Box 649
Marietta, OH 45750

International customers, please contact your local Thermo Scientific distributor.

Table of Contents

Safety Information	6
Alert Signals	6
Warnings	6
Información de seguridad.....	8
Señales de alerta.....	8
Advertencias	8
Sicherheitshinweise	10
Gefahrensymbole.....	10
Warnungen	10
Informations de sécurité.....	12
Signaux d'alerte	12
Avertissements	12
Informação de Segurança.....	14
Sinais de Aviso.....	14
Avisos	14
Информация по безопасности	16
Сигналы тревоги	16
Предупреждения	16
General Specification.....	18
Ceramic Hot Plate Specifications.....	18
Ceramic Stirrer Specifications.....	18
Ceramic Stirring Hotplate Specifications	19
Aluminum Hot Plate Specifications	19
Heating Specifications	20
Stirring Speed Specifications	20
Environmental Conditions	21
Declaration of Conformity.....	21
Introduction.....	22
General Cleaning Instructions.....	22
Unpacking and Installation.....	23
Unpacking.....	23
Installation	23
SUPER-NUOVA ⁺ Hot Plate Control Panel	24
SUPER-NUOVA ⁺ Stirrer Control Panel	24
SUPER-NUOVA ⁺ Stirring Hot Plate Control Panel.....	24
Control Panel.....	25
Principles of Operation.....	26

Operation.....	27
Hard Switch	27
Setting the Stirring Speed.....	27
Setting the Temperature.....	28
Controlling Solution Temperature Using External Probe	29
Using Preset Programs	29
Heating Metal Vessels and Sand Baths	30
Setting the Time	30
Hardware Over-Temperature-Protection.....	30
System Configuration	31
Temp Unit.....	32
Over-temperature Set Point –OTP.....	32
Preset Continue	33
Temperature Probe Calibration	34
Troubleshooting	35
Error Codes	35
Accessories	36
THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY	37

DOMINIQUE DUTSCHER SAS

Safety Information

Alert Signals

Warning

Warnings alert you to a possibility of personal injury.

Caution

Cautions alert you to a possibility of damage to the equipment.

Note

Notes alert you to pertinent facts and conditions.

Hot Surface

Hot surfaces alert you to a possibility of personal injury if you come in contact with a surface during use or for a period of time after use.

Note

The SUPER-NUOVA⁺ hot plates are not explosion proof. If explosion proof models are required, contact Customer Service at 1-800-553-0039 for more information.

Warning

Refer servicing to qualified personnel.

Your Thermo Scientific SUPER-NUOVA⁺ hot plate, stirrer or stirring hot plate has been designed with function, reliability, and safety in mind. It is your responsibility to install it in conformance with local electrical codes. For safe operation, please pay attention to the alert signals throughout the manual.

Warning: These products should be used only under the operating conditions specified in the Operating Manual. Always use safe laboratory practices and do not leave the unit in operation while unattended as product functionality or laboratory practice failures could occur that might lead to uncontrolled or excessive heating of the top surface. Safety procedures (including, but not limited to, unplugging when not in use) and response plans should be put in place to address the worst case possibility. If an over-temperature failure occurs, the top surface temperature could rise to the maximum temperature (300-450°C depending on your model's specification) and remain at that temperature indefinitely. Under these conditions, the material being heated on the surface of the hotplate could reach levels in excess of the maximum temperature.

This manual contains important operating and safety information. The user must carefully read and understand the contents of this manual prior to the use of this equipment.

Warnings

To avoid electrical shock, always:

1. Use a properly grounded electrical outlet of correct voltage and current handling capacity.
2. Disconnect from the power supply prior to maintenance and servicing.

To avoid personal injury:

1. Do not use in the presence of flammable or combustible materials — fire or explosion may result. This device contains components which may ignite such materials. Not rated for use in hazardous atmospheres.
2. Use caution when heating volatile materials; top surface and element can reach the “Flash Point Temperature” of many chemicals. These hot plates are not explosion proof. Fire or explosion may result. Unit contains components which may ignite such materials.
3. Keep top surface clean. Use a non-abrasive cleaner. Alkali spills, hydrofluoric acid spills or phosphoric acid spills may damage top and lead to thermal failure. Unplug unit and remove spills promptly. Do not immerse unit for cleaning.
4. Replace the top immediately if damaged by etching, scratching or chipping. A damaged top can break in use.
5. Do not use metal foil on hot plate which may block air flow. Overheating will result.
6. Check and tighten the removable cord periodically making sure it is secure. If loosened, the cord could become hot and/or spark and be a potential fire hazard. If cord appears damaged, replace immediately. If cord is repeatedly loosened it is recommended to purchase the power cord.
7. Do not remove or modify grounded power plug. Use only properly grounded outlets to avoid shock hazard.
8. Use appropriate hand and eye protection when handling hazardous chemicals.
9. Gross weight of items placed on top of hot plates should not exceed 35 lbs. (15.9 kg.) on the 10” x 10” models and 25 lbs. (11.3 kg) on the 7” x 7” models.
10. “Caution: Hot Top. Avoid Contact.” The top plate of the unit can remain hot for some time after use. A “CAUTION-HOT TOP” light will remain on until top plate temperature cools to below 50°C.
11. Do not leave an active probe out of the fluid. This may cause uncontrolled heating of the fluid on the hot plate and unintentional boiling or an explosion could occur.
12. Localized heater element temperature can be significantly higher than the temperature indicated on the display. If flammable concentrations reach internal element, a fire could result.
13. Note that the exterior housing will be hot during and for a period of time after use.
14. Refer servicing to qualified personnel.

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Ceramic Stirring Hotplate Specifications

Overall Dimensions	SP88857190	SP88857195	SP88850190	SP88850195
Width - in. (cm)	8.2 (20.8)	8.2 (20.8)	5.1 (13)	5.1 (13)
Height - in. (cm)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)
Depth - in. (cm)	14.2 (36)	14.2 (36)	17.2 (43.8)	17.2 (43.8)
Weight - lbs. (kg)	9.0 (4.1)	9.0 (4.1)	13.0 (5.9)	13.0 (5.9)
Top Plate Dimensions				
Width - in. (cm)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	7.25 (18.4)	7.25 (18.4)	10.5 (26.7)	10.5 (26.7)
Volts	100-120	220-240	100-120	220-240
Amps	12.3	4.4	13.2	6.0
Watts	1480	1056	1580	1440
Freq.	50/60	50/60	50/60	50/60
Phase	1	1	1	1
Max. Temp °C (°F)	450°C(842°F)	450°C(842°F)	400°C(752°F)	400°C(752°F)

Aluminum Hot Plate Specifications

Overall Dimensions	HP88857194	HP88857198	SP88857194	SP88857198
Width - in. (cm)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)	8.2 (20.8)
Height - in. (cm)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)	4.1 (10.5)
Depth - in. (cm)	14.2 (36)	14.2 (36)	14.2 (36)	14.2 (36)
Weight - lbs. (kg)	7.5 (3.4)	7.5 (3.4)	9.0 (4.1)	9.0 (4.1)
Top Plate Dimensions				
Width - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Height - in. (cm)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)	1.0 (2.5)
Depth - in. (cm)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)	7.25 (18.4)
Volts	100-120	220-240	100-120	220-240
Amps	6.0	2.6	6.0	2.6
Watts	720	624	720	624
Freq.	50/60	50/60	50/60	50/60
Phase	1	1	1	1
Max. Temp °C (°F)	300°C(572°F)	300°C(572°F)	300°C(572°F)	300°C(572°F)

GENERAL SPECIFICATIONS

Heating Specifications*

Top Plate Surface - Solid Ceramic

Temperature range: 30°C - 450°C (86°F - 842°F)* 7" x 7"

30°C - 400°C (86°F - 752°F)* 10" x 10"

Heat-up time to maximum temperature (unloaded top plate).

Accuracy of the temperature display vs. the actual average temperature of a 2" diameter of setting area at the center of the top plate (setpoint 100°C unloaded).

- Temperature stability at the center of the top plate surface (@ 100°C unloaded).

Accuracy of remote probe at user selected calibration temperature after calibration procedure.

- Temperature stability using remote probe (500 ml of water in a 1000 ml flask at 70°C)

7" x 7" 10" x 10"

5 minutes

7 minutes

± 5.0°C

± 5.0°C

± 1.0°C

± 1.0°C

±0.5 typical

±0.5 typical

±0.5°C

±0.5°C

Top Plate Surface - Solid Aluminum

Temperature range: 30°C - 300°C (86°F - 572°F)* 7" x 7"

Heat-up time to maximum temperature (unloaded top plate).

Accuracy of the temperature display vs. the actual average temperature of a 2" diameter of setting area at the center of the top plate (setpoint 100°C unloaded).

- Temperature stability at the center of the top plate surface (@ 100°C unloaded).

Accuracy of remote probe at user selected calibration temperature after calibration procedure.

- Temperature stability using remote probe (500 ml of water in a 1000 ml flask at 70°C)

7" x 7"

5 minutes

± 5.0°C

± 1.0°C

±0.5 typical

±0.5°C

** As the top plate becomes dirty, the maximum temperature will decrease. To return the unit to its maximum temperature performance, use a mild abrasive to remove stained areas.

Stirring Speed Specifications

Speed Range 50 to 1500 RPM (Maximum speed is dependent on the viscosity of the solution).

Stability of the stirring speed setpoint (600 ml of water in a 1000 ml glass flask) ± 2.0% at 1000 RPM

Top Plate Size

7 x 7

10 x 10

Max Recommended Flask Size

4 liters

6 liters

Max Weight on Top Plate

25 lbs

35 lbs

* The specification defined in this section is derived from testing data statistics of Thermo Fisher Scientific lab. We reserve the right to make changes in design and specification without prior notice.

Environmental Conditions

Operating:	5°C to 35°C; 20% to 80% relative humidity, non-condensing. Installation category II (overvoltage) in accordance with IEC 664. Pollution degree 2 in accordance with IEC 664. Altitude Limit: 2,000 meters.
Storage:	-25°C to 65°C 10% to 85% relative humidity

Declaration of Conformity

We hereby declare under our sole responsibility that this product conforms with the technical requirements of the following standards:

EMC:	EN 61000-3-2	Limits for harmonic current emissions
	EN 61000-3-3	Limits for voltage fluctuations and flicker
	EN 61326-1	Electrical equipment for measurement, control, and laboratory use; Part I: General Requirements
Safety:	EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use; Part I: General Requirements
	EN 61010-2-010	Part II: Particular requirements for laboratory equipment for the heating of materials
	EN 61010-2-051	Part II: Particular requirements for laboratory equipment for mixing and stirring

per the provisions of the Electromagnetic Compatibility Directive 2014/30/EU, and per the provisions of the Low Voltage Directive 2014/35/EU.

Copies of the Declaration of Conformity are available upon request.

Introduction

SUPER-NUOVA⁺ Stirring Hotplate

Please read all the information in this manual before operating the unit.

Your SUPER-NUOVA⁺ hot plate, stirrer or stirring hot plate is a heating and/or stirring plate designed for laboratory procedures requiring precise control of temperature and/or stirring speed. Each SUPER-NUOVA⁺ model includes a liquid crystal display for monitoring actual and setting temperature and/or stirring speed. The hot plate is capable of producing accurately controlled top plate temperatures from 30°C through 300°C, 400°C or 450°C depending on models. The temperature is controlled at the plate surface by an internal Type K thermocouple sensor, or the solution temperature may be controlled by utilizing the included PT100 probe. A PT1000 probe may be ordered separately. The 7 x 7 or 10 x 10 top plate on the SUPER-NUOVA⁺ units are solid ceramic or solid aluminum, and are suitable for use with glass or metal vessels and sand baths.

Other features of the SUPER-NUOVA⁺ include a time feature that can be used to shut off heating. Specific temperature setpoints may be set using the PROGRAM buttons on the SUPER-NUOVA⁺. An over temperature protection (OTP) can be displayed in the system configuration menus. See details in "System Configuration".

Your SUPER-NUOVA⁺ model may be used for general purpose heating applications including sample preparation, heating reagents, melting paraffin, warming resinous chemicals, content analysis, solvent evaporations, digestions, media preparation and sterilization, titrations, sand baths, and micro-scale chemistry applications.

General Cleaning Instructions

Keep top surface clean. Use a non-abrasive cleaner. Alkali spills, hydrofluoric acid spills or phosphoric acid spills may damage top and lead to thermal failure. Unplug unit and remove spills promptly. Do not immerse unit for cleaning. Wipe exterior housing with lightly dampened cloth containing mild soap solution.

Unpacking and Installation

Warning

Use a properly grounded electrical outlet of correct voltage and current handling capacity.

Do not remove or modify grounded power plug. Use only properly grounded outlets to avoid shock hazard. Not rated for use in hazardous atmospheres.

Do not use in the presence of flammable or combustible materials; fire or explosion may result. This device contains components which may ignite such materials. If explosion proof models are required, contact Customer Service at 1-800-553-0039 for more information.

Do not use in highly corrosive atmospheres; corrosive fumes and spills may damage top and internal components, creating shock hazard.

Note

Some misalignment of the motor bearings in this product may have occurred during shipping. Prior to using this product, run the stirrer at maximum speed for ten minutes to realign the bearings.

Unpacking

Remove your SUPER-NUOVA⁺ unit from the carton. Inspect to ensure that the unit has not been damaged during shipment. If the unit appears to have sustained shipping damage contact the distributor from whom you purchased this product or Customer Service at 800-553-0039. Check for remote probe and thumbscrew prior to discarding packaging.

The following items are included in the shipment:

SUPER-NUOVA⁺ unit

Cord

Remote Probe – PT100

Stir Bar (Stirring models only)

Knob

Operation's Manual


If any of these items are missing from the carton, contact Customer Service.

Installation


Set the unit on a flat stable surface at least 12" away from combustible materials, and plug the cord set into a properly grounded electrical outlet of correct voltage and current handling capacity.

UNPACKING AND INSTALLATION


SUPER-NUOVA⁺ Hot Plate Control Panel

SUPER-NUOVA⁺ Stirrer Control Panel

SUPER-NUOVA⁺ Stirring Hot Plate Control Panel

Control Panel

Key	Description
1	CAUTION! HOT TOP Indicator
2	Display
3	Knob
4	Temperature Set Button
5	Speed Set Button
6	Time Set Button
7	Program Buttons

DOMINIQUE DUTSCHER SAS

Principles of Operation

Each SUPER-NUOVA⁺ unit utilizes the latest microprocessor technology to deliver the most reliable, accurately controlled, ceramic top or aluminum top hot plate on the market.

Your SUPER-NUOVA⁺ unit with heating functionality has an electronic, closed-loop feedback control which will accurately maintain temperature setpoints from 30°C through 300°C, 400°C or 450°C depending on models.

The SUPER-NUOVA⁺ units use a Type K thermocouple for measuring the top plate temperature, OTP and the external probe. Both top plate surface and OTP temperature are measured with independent analog circuits. This provides a significant increase in safety due to the redundancy of the circuits. It is possible to heat certain types of metal vessels and sand baths on the SUPER-NUOVA⁺ without damaging the ceramic top or aluminum top.

The electronic stirring speed control will maintain the speed setpoint when the unit is loaded. The motor is combined with a powerful magnet to provide exceptional magnetic coupling with a stir bar.

The time feature on the SUPER-NUOVA⁺ can be programmed to turn off the stirring, heating or both. See “Setting the Time” section of this manual.

Operation

Warning

Use caution when heating volatile materials; top surface and element can reach the “Flash Point Temperature” of many chemicals. These hot plates are not explosion proof. Fire or explosion may result. Unit contains components which may ignite such materials.

Use appropriate hand and eye protection when handling hazardous chemicals.

“Caution: Hot Top. Avoid Contact.” The top plate of the unit can remain hot for some time after use. A “CAUTION - HOT TOP” light will remain on until top plate temperature cools to below 50°C.

Note

Boiling times are dependent on solution volume and the surface area of the flask that is exposed to the hot plate. For example, when heating the same amount of solution in a 2L flask vs. a 1L flask, the solution will heat about 20% faster.

Note

The top plate temperature will be higher than the solution temperature.

Hard Switch

The SUPER-NUOVA⁺ unit has a hard switch on the left side of the unit. The unit may be powered on by pressing the hard switch. When the unit is turned on, the unit will initialize and then model name will be displayed on the first row of display screen while software version in second row. To turn off power to the unit, press the hard switch off.


Setting the Stirring Speed

The SUPER-NUOVA⁺ stirring unit has an electronic feedback speed control which will maintain a precise speed setpoint from 50 rpm through 1500 rpm.

(Maximum speed is dependent on viscosity of the solution.) The SUPER-NUOVA⁺ stirring hot plates are equipped with a strong magnet and high torque motor which will draw a vortex in fluids up to 700 cp.

To set the stirring speed, press “SPEED” key located under display. Then the speed setpoint will flash. Use the knob on the front panel to set a speed. Turn the knob clockwise to increase the speed or counterclockwise to decrease the speed. When desired setpoint has been reached, press the “SPEED” key or wait 10 seconds, the speed setpoint will no longer flash indicating that it has been set.

Once the unit begins stirring, the display will indicate both actual speed and setpoint. To stop the stirring action without changing your setpoint press the knob once. To resume stirring at the set speed, press the knob again.

To turn stirring off, press “SPEED” key under the display and turn the knob counterclockwise until the display reads “OFF.”

Note

Boiling times are dependent on solution volume and the surface area of the flask that is exposed to the hot plate. For example, when heating the same amount of solution in a 2L flask vs. a 1L flask, the solution will heat about 20% faster.


Note


The top plate temperature will be higher than your solution temperature.

Setting the Temperature

The SUPER-NUOVA⁺ heating unit has an electronic closed-loop feedback control which will accurately maintain temperature setpoints from 30°C through 300°C, 400°C or 450°C (depending on the model). An unloaded hot plate will heat to maximum temperature in just 5 minutes. A “CAUTION HOT TOP” light on the front panel will illuminate whenever the top surface temperature exceeds 50°C. The temperature may be controlled either at the top plate by the internal Type K thermocouple sensor, or in the solution using an accessory ungrounded PT100 probe. If a probe is used, the display will indicate “R” after actual temperature. The SUPER-NUOVA⁺ heating units will accept any ungrounded PT100(PT1000) probe, however the accuracy of probes other than those tested and specified may not reach the probe accuracy rating designated.

To set the temperature, press “TEMP” key located under display. Then the temp setpoint will flash. Use the knob on the front panel to select a temp. Turn the knob clockwise to increase the temp or counterclockwise to decrease the temp. When desired setpoint has been reached, press the “TEMP” key or wait 10 seconds, the temp setpoint will no longer flash indicating that it has been set.

Controlling Solution Temperature Using External Probe

To control the solution temperature plug a PT100 probe into the probe receptacle located on the back side of the unit. When the probe is connected, the adjustable temperature range is 30-250 °C. Place the probe into the solution. The display will indicate the actual temperature of the solution as measured by the probe.

The external probe offers more exact temperature control than regulating the top plate by the internal sensor. If you need to maintain a precise setpoint it is recommend to use a probe to control the solution temperature instead of controlling by the top plate temperature.

When using a probe with the SUPER-NUOVA⁺ unit it is recommended that a clamp on a support rod be used to hold the probe in the solution.

To ensure accurate probe readings, as much of the probe sheath as possible should be immersed in the solution. Make sure the probe is immersed in the liquid and is not located in air or outside of the solution. If the probe is plugged into the hot plate, but is not in solution while the heat control is operating, the temperature display will continue to indicate an ambient temperature, and a probe out of solution error will occur. Because the set point cannot be reached the heating element will continue to supply heat to the top plate, and the maximum top plate temperature(300°C, 400°C or 450°C) may be reached. If the remote probe does not sense a temperature change in three minutes the unit will display ExProbe Err and shut down.

Using Preset Programs

The SUPER-NUOVA⁺ unit features the option of presetting 3 setpoints into its memory. To program a preset, press any preset button, and the indicator LED on the number button will be illuminated. Then adjust the time, speed and heat features to your desired setpoints. Press the number button again to exit the preset program. The setpoints are stored in memory automatically. The presets will be held in memory, even if the unit is unplugged, until you choose to override it with another preset.

To activate one of the presets, press the corresponding number then press the knob, the time, speed and heat features will automatically adjust itself.

It may be helpful to set one of the presets at "OFF." This way, the heating can be shut off quickly.

The SUPER-NUOVA⁺ unit has a Preset Continue function. It means different presets can run continuously when Preset Continue is being utilized. The default setting is to turn off the function. By pressing 1 and then pressing the knob the unit will run by continuous presets by order 1, 2, 3. By selecting 2 and then pressing the knob the unit will run by continuous presets by order 2, 3.

See "System Configuration" section of the manual to activate the Preset Continue function.


OPERATION

Heating Metal Vessels and Sand Baths


Metal vessels and sand baths cannot be heated on most solid ceramic tops because of the tendency for the metal and the sand to reflect heat back into the top, eventually exceeding the maximum temperature rating for a ceramic top, which causes it to break. Because of its advanced electronic control, the SUPER-NUOVA⁺ is capable of precisely regulating the top plate temperature. Metal vessels and sand baths may be heated safely without the danger of the ceramic top breaking.

Setting the Time

The time function on the SUPER-NUOVA⁺ can be set to shut off heating. To set the time, press "TIME" key located under display. Hold will begin to flash. Turn the knob until time shows then press "TIME" key again. Use the knob on the front panel to select a time. Turn the knob clockwise to increase the time or counterclockwise to decrease the time. When desired setpoint has been reached, press the "TIME" key or wait 10 seconds, the time setpoint will no longer flash indicating that it has been set. In the temperature rising mode, when the temperature is more than 2 degree below the temperature setpoint, the timer will start to count. In the temperature declining mode, when the temperature is less than 2 degree above the temperature setpoint, the timer will start to count.


When the time setpoint display shows "HOLD" it means there is no setting for time function.

Hardware Over-Temperature-Protection

A hardware Over-Temperature-Protection(OTP) analog circuit is designed separately to protect unexpected heating. A small straight blade screwdriver can be used to adjust the scale at right side of the unit. The protection value can be set from 60°C to 500°C(or 300°C/450°C depending on models).

Please notice that this function is only purposed to avoid abnormal heating if unit is out of control. Users should not use this function to limit the temperature in normal applications. Refer to the section 'Over-temperature Set Point – OTP' for accurate temperature protection.

Default set is at maximum scale.

System Configuration

The System Configuration has three features that will allow a user to customize their unit.

A table of the features available in the System Configuration is given below. The features are given in order of their appearance in the System Configuration when rotating the knob in a clockwise rotation. The Feature column is the name of the feature, the Display column is the designation of that feature as shown on the display of the unit, the Availability column gives which type of unit the feature is applicable (HP = hot plate SP = stirring hot plate S=stir plate), and the last column gives a brief description of the purpose of the feature.


Feature	Display	Availability	Brief Description
Temperature Unit Conversion Over-temperature Set Point	Temp Unit	HP , SP	Allows the user to converse temperature unit between °C and F.
Over-temperature Set Point	User OTP	HP , SP	Allows the user to see the setting of the Over Temperature Protection control.
Preset Continue Setting	Preset Continue	HP , S, SP	Allows the user to set the program presetting function.
Temperature Probe Calibration	Calibrate Temp	HP , SP	Allows calibration of the temperature probe and associated circuitry.

System Configuration


Temp Unit

This option is only available for units that have heating functionality.

Plug the unit into the appropriate power, enter the System Configuration by pressing and holding the knob. After about 5 seconds you may remove your finger from the knob. The display indicate “System Config” on the first row, the “Temp Unit” will show on the second row of the display.

Press the knob to enter the Temp Unit setting interface.(Long press will return to operation menu.) The “Temp Unit” will display on the first row, the Actual temperature and unit “c” or “F” flash on the second row. Turn the knob clockwise or counterclockwise to change the temperature unit. When the desired temperature unit has reached, wait 10 seconds, the display will return to operation menu indicating that it has been set.


Over-temperature Set Point –OTP

This option is only available for units that have heating functionality. The Over-temperature Set Point is the set point of the independent temperature control. The purpose of the control is to limit the top temperature of the element independently of the standard heating control. This is helpful if there is a failure of the main control that would cause the top to heat excessively. Typically the set point will be adjustable from 60 to 300°C, 400°C or 540°C(depending on the model).


The Over-temperature can be set or viewed with this menu option. To do so, follow the instructions given below:

Enter the System Configuration by pressing and holding the knob. After about 5 seconds a single you may remove your finger from the knob. The display indicate “system config” on the first row, the “Temp Unit” will show on the second row of the display.

Turn the knob clockwise until the “User OTP” show on the second row of the display. If you do no action, the display will return to the operation status.


Press the knob to enter the Over-temperature setting interface.(Long press will return to operation menu.) The “User OTP” will display on the first row, the Over-temperature setpoint flash on the second row. Turn the knob clockwise to increase the temp or counterclockwise to decrease the temp. When the desired Over-temperature setpoint has been reached, wait 10 seconds, the display will return to operation menu indicating that it has been set.

Preset Continue

Enter the System Configuration by pressing and holding the knob. After about 5 seconds a single you may remove your finger from the knob. The display indicate “system config” on the first row, the “Temp Unit” will show on the second row of the display. Turn the knob clockwise until the “Preset Continue” show on the second row of the display. If you do no action, the display will return to the operation status.


Press the knob to enter the Preset Continue setting interface.(Long press will return to operation menu.) The “User OTP” show on the first row, “OFF” flash on the second row. Turn the knob clockwise to change the status. When “On” has reached, wait 10 seconds, the display will return to operation menu indicating that it has been set.


Temperature Probe Calibration

This section applies only to the calibration of the temperature probe system. The hot plate surface and the motor speed cannot be calibrated. This method is only applicable to units that have heating functionality.

Enter the System Configuration by pressing and holding the knob. After about 5 seconds a single you may remove your finger from the knob. The display indicate “system config” on the first row, the “Temp Unit” will show on the second row of the display. Turn the knob clockwise until the “Calibrate Temp” show on the second row of the display. If you do no action, the display will return to the operation status. The option “Calibrate Temp” will only show when there is a remote probe PT100 or PT1000 adopted.

Press the knob to enter the Calibrate Temp setting interface. The “Calibrate Temp” show on the first row, the measured temperature flash on the second row. If the data is not matched the temperature you measure by external thermo detector, you can turn the knob to calibrate. When it’s matched, press the knob to complete the calibration.

When setting and/or viewing is complete, the knob can be pressed to return to the previous operation menu.

Troubleshooting

Error Codes

Error Handler will lock out heating functions if heating error is detected. Stirring functionality is unaffected. Error Handler will lock out stirring functions if stirring error is detected. Heating functionality is unaffected. If the condition that caused the error is no longer present, the unit will clear Errors automatically.

Displayed Message	Intended to Detect	Cause	Solution
InProbe Err	Internal thermocouple out of range.	Internal thermocouple not connected. Thermocouple open. Thermocouple connected backwards (reversed polarity).	Ensure proper connection and polarity of thermocouple. Replace thermocouple (attached to element). Ensure proper connection and polarity.
Check Heater	Excessive top heat-up time	Internal thermocouple short circuit. Failure in Internal thermocouple. Failure in Element. Failure in optocoupler/triac circuit.	Remove short. Replace thermocouple (attached to element). Replace Element. Replace Control Board.
ExProbe Err	External Probe left out of solution.	External Probe left out of solution. External probe connected backwards.	Place external probe into solution. Correct orientation of external probe.
Over Temp	OTP thermocouple out of range.	OTP Circuit failure. OTP potentiometer set too low.	Replace Control Board. Increase OTP setting.
Check Motor	Motor system failure.	Locked rotor condition. Failure of motor. Failure of motor circuit. Failure of motor sensor.	Free locked rotor. Replace motor. Replace Control Board. Replace motor sensor.

Accessories

Order No	Description
CIC0001444	Splash guard protection shield 4", for 4" x 4" Stirrers and Hotplate Stirrers
CIC0001445	Splash guard protection shield 7", for 7" x 7" Stirrers and Hotplate Stirrers
CIC0001446	Splash guard protection shield 10", for 10" x 10" Stirrers and Hotplate Stirrers
CIC0001447	PT100 External probe, 0-400°C, Stainless steel, for all sizes of Hotplates and Hotplate Stirrers
CIC0001448	PT1000 External probe, 0-400°C, Stainless steel, for all sizes of Hotplates and Hotplate Stirrers
CIC0001449	Hastelloy PT1000 External probe, 0-400°C, Hastelloy, for all sizes of Hotplates and Hotplate Stirrers
CIC0001830	Supporting rod, Aluminum Alloy, for Hotplates and Hotplate Stirrers
CIC0001831	Clamp holder, Stainless steel, for Hotplates and Hotplate Stirrers
CIC0001468	Thermo ASSY CORD CHINA 10A 250V PLUG
CIC0001469	Thermo ASSY CORD 18-3 HSJO 10A EURO 86 PE
CIC0001470	Thermo ASSY CORD 1MM-2 BRIT 240 10A PE
CIC0001471	Thermo ASSY CORD US&JP 120V PLUG

THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first three (3) years, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo Fisher Scientific postage paid and replacement parts are shipped FOB destination.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY. Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Department at 1-800-438-4851 (USA and Canada) or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.

DOMINIQUE DUTSCHER SAS