P.O. Box 3201  
314 Ella T. Grasso Ave.  
Torrington, CT 06790  
[www.waringproducts.com](http://www.waringproducts.com)

018800 RV09 0802

**WARING**  
**COMMERCIAL**

# 4 Litre Blender

## Operating Manual

Blender 4 Litre 120 volt - 240 volt

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Licuada de 4 litros, 120 voltios - 240 voltios

## Mode d'emploi

Mixer 4 Litres 120 volt - 240 volt

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Mixer 4 Litre 120 v - 240 v

**DOMINIQUE DUTSCHER SAS**

## Wartung

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## IMPORTANT SAFEGUARDS

When using electrical appliances, basic safety precautions should always be followed including the following:

1. **READ ALL INSTRUCTIONS.**
2. To protect against electrical hazards, do not immerse the Blender Base in water or other liquid.
3. Close supervision is necessary when any appliance is used by or near children.
4. Unplug from outlet when not in use, before putting on or taking off parts, and before cleaning.
5. Avoid contacting moving parts.
6. Do not operate any appliance with a damaged cord or plug, or after the appliance malfunctions or is dropped or damaged in any manner. Return appliance to the nearest authorized service facility for examination, repair or electrical or mechanical adjustment.
7. The use of attachments including canning jars not recommended by the appliance manufacturer may cause hazards.
8. Do not use outdoors.
9. Do not let cord hang over edge of table or counter, or touch hot surfaces.
10. Keep hands and utensils out of container while blending to prevent the possibility of severe personal injury and/or damage to the blender. A rubber scraper may be used but must be used only when the blender is not running.
11. Blades are sharp. Handle carefully.
12. Always operate blender with cover in place.
13. When blending hot liquids, completely remove the stopper in the cover before turning on the blender. See Hot Blending Instructions.

## SAVE THESE INSTRUCTIONS

## CIRCUIT BREAKER OPERATION

You will find a small, plastic button located on the underside of the motor base of some models.

This is the reset button for a circuit breaker which prevents the motor from overheating in case of an overload.

If an overload occurs, the circuit breaker will shut off the motor and cause the button to protrude slightly from its normal position.

To reset the circuit breaker and resume operation, do the following:

- 1) SET THE POWER SWITCH TO ITS OFF (O) POSITION.
- 2) UNPLUG THE POWER CORD FROM ELECTRICAL OUTLET.
- 3) DETERMINE THE REASON FOR THE OVERLOAD, AND ELIMINATE THE CAUSE (FOR EXAMPLE: REDUCE THE AMOUNT OF MATERIAL BEING PROCESSED).
- 4) PUSH THE RESET BUTTON BACK INTO ITS NORMAL POSITION.
- 5) PLUG THE POWER CORD BACK INTO ELECTRICAL OUTLET.
- 6) RESUME NORMAL OPERATION.

## Specifications

- Capacity: 4 Litre Working Capacity
- Switch Setting: (3 speeds) OFF (0) LOW (1) MED (2) HIGH (3) PULSE (PULS)
- Revolutions Per Minute: 15,800 18,000 20,800 20,800  
(free running-approximate)
- Rating: 120 Volts-60HZ-Single Phase AC (~)-15.0 Amps.  
220-240 Volts-50-60HZ-Single Phase AC (~)-1500 Watts.
- Base: Die-cast – gray coating
- Container: Stainless steel with stainless steel blades - clear polycarbonate clamp-on lid, with gasket and stopper.
- Bearing: High performance ball bearing.
- Cord: 7 foot (2.44 m) - three conductor ground cord with strain relief in base.
- Overall Height: 26 inches (66 cm).
- Weight: 34 pounds (15.4 kg) net - Shipping Weight 38 pounds (17.2 kg).
- UL, NSF, CSA Listed (120 volts ONLY).
- Timer: Some models have a timer, adjustable from 1 to 120 seconds for automatic blending.
- Operating Cycle: 3 min ON/1 min OFF

## GROUNDING INSTRUCTIONS

For your protection, this appliance is equipped with a 3-conductor cordset.

120 Volt units are supplied with a molded 3-prong grounding-type plug (NEMA 5-15P) and should be used in combination with a properly connected grounding type outlet as shown in Figure A.

If a grounding type outlet is not available, an adapter, shown in Figure B, may be utilized to enable a 2-slot wall outlet to be used with a 3-prong plug except in Canada and Europe.

Referring to Figure C, the adapter must be grounded by attaching its grounding lug under the screw of the outlet cover plate.


**Caution:** Before using an adapter, it must be determined that the outlet cover plate screw is properly grounded. If in doubt, consult a licensed electrician. Never use an adapter unless you are sure it is properly grounded.

**Note:** Use of an adapter is not permitted in Canada and Europe.

220-240 volt "GS Safety-Approved" units are supplied with a molded European grounding-type plug (CEE 7/7), for use in Europe.

220-240 volt "Standard" units are supplied without a plug, and an appropriate one must be installed by user in accordance with local electrical codes observing the following cordset conductor color coding:

Green/Yellow: Earth Ground • Blue: Neutral • Brown: Live  
Note: User installed cordset plugs are not permitted in Europe.

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
## Operation

1. Before starting, make sure that the blender cord is unplugged and that the switch is set to the OFF (0) position.
2. Place container assembly on motor base so that it is firmly seated and not tilted to one side.
3. Be sure the switch is still in the OFF (0) position, and plug the power cord into a convenient outlet.
4. Put ingredients into container and cover it. Always add liquid ingredients first, unless otherwise specified in recipe. Do not fill above MAX line.  
For continuous operation, set switch to the LOW (1), MED (2), or HIGH (3) speed position. For momentary operation, push the PULSE ("PULS") button; motor will run on HIGH (3) speed until the PULSE ("PULS") button is released.
5. Always hold container with one hand when starting, and while blending. Never operate without jar lid in place.

### WHEN FINISHED BLENDING, SET SWITCH TO OFF (0) POSITION.

**NOTE:** For blenders with timers, turn the timer control knob clockwise to the desired time setting before starting the motor. If time setting is less than 20 seconds, turn knob past "20", then back to desired setting.

When switch is set to the desired speed, motor will run. Motor will automatically stop at the end of the set time.

6. Allow motor to come to a complete stop before removing container. Never place the container on the motor base or remove it from the motor base while the motor is running.
7. Remove the container from the motor base and empty the contents.
8. Wash and dry container after each use.
9. **Not For Waterjet Cleaning.**
10. The use of ear protection is recommended during extended operating periods.

### To realize maximum service life, observe the following precautions:

- Do not lubricate container bearing, drive shaft, or shaft seals. These parts are lubricated and sealed at the factory and do not require any further lubrication.
- Reduce amounts of material processed to less than the full 4 litre capacity when working with extremely viscous or fibrous materials.
- When cleaning exterior of container, keep it in an upright normal operating position to prevent liquids from running into bearing through its bottom end.

## Hot Blending Instructions: Important!

To minimize the possibility of scalding when blending hot materials, the following precautions and procedures must be observed:

Hot materials are those which are 100°F (38°C) or hotter.

If you are using the standard container lid that originally came with your blender, limit the amount of material processed to one (1) litre and completely remove the stopper in the lid before turning on blender. This will allow venting of the rapidly heated and expanding air present inside container when blender is running.

Always start blending on lowest speed setting, then progress to desired speed.

Never stand directly over blender when it is being started or while it is running.

Constant monitoring of the condition of the container lid gasket and the container lid clamps is mandatory. Replace these parts if they become worn or loose, or are damaged in any way.

**Note:** A special accessory Hot Blending lid assembly is available which allows safe processing of larger loads. Order: One-Gallon Hot Blending Lid, Catalog No. CBL-10. If you are using the special accessory Hot Blending Lid, follow the instructions packed with the lid. **BEFORE PROCEEDING WITH HOT BLENDING, PLEASE REREAD ALL INSTRUCTIONS IN THIS OPERATING MANUAL.**

## Maintenance

### Motor - Not For Waterjet Cleaning

The Waring Commercial Blender ball bearing motor requires no greasing or oiling. Sealed-in lubrication provides smooth unit operation.

Container - The stainless steel container can be sterilized without injury to the bearing.

Depending on the operating conditions, service life of the main components of the container blade assembly, socket head shaft and shaft seals is 80 to 100 hours actual running time.

When any one of these main components has worn to the extent that replacement is required, it is likely that the other components will need to be replaced also. It is therefore recommended that container blade assemblies be replaced as a unit.

Die-cast lock nut must be tightened securely (65-70 footpounds torque) when reassembling container.

Power Supply Cord - The power supply cord should not be replaced by the end user. The power supply cord should only be replaced by the manufacturer's authorized importer (or, in the United States and Canada, by an authorized service center).

**Note:** All parts require constant monitoring for conditions requiring replacement. For assistance, contact your nearest Service Center.

## Limited One Year Warranty

Waring warrants every new Waring Commercial Blender to be free from defects in material and workmanship for a period of one year from the date of purchase when used with foodstuffs, nonabrasive liquids (other than detergents) and nonabrasive semiliquids, providing it has not been subject to loads in excess of the indicated rating. In addition, the motor in this appliance has a full three year warranty.

Under this warranty, the Waring Commercial Division will repair or replace any part which, upon our examination, is defective in materials or workmanship, provided the product is delivered prepaid to the Waring Service Center at P.O. Box 3201, 314 Ella T. Grasso Ave., Torrington, CT 06790, or any Factory Approved Service Center.

This warranty does not: a) apply to any product which has become worn, defective, damaged or broken due to use or due to repairs or servicing by other than the Waring Service Center or a Factory Approved Service Center or due to abuse, misuse, overloading or tampering; or b) cover incidental or consequential damages of any kind.

The one year warranty is applicable only to appliances used in the United States or Canada; this supersedes all other express product warranty or warranty statements. For Waring products sold outside of the U.S. and Canada, the warranty is the responsibility of the local importer or distributor. This warranty may vary according to local regulations.

**Warning:** This warranty is void if appliance is used on Direct Current (D.C.)

## Service

Should service on this product in USA or Canada be required at any time, call or write to:  
Waring Service Center, P.O. Box 3201, 314 Ella T. Grasso Ave., Torrington, CT 06790  
Tel. 1-800-492-7464 • FAX (860) 496-9017  
or  
Nearest Authorized Service Center

Replacement parts for Waring "Out of Warranty" repairs are obtainable at the above address.

**Note:** If your application involves the processing of biologically active materials, any items returned for service must be sterilized before being shipped. Such items must be clearly identified as having been sterilized both by a tag on the item itself and by a statement in the accompanying correspondence. Waring will not accept delivery of any returns unless they are so identified.

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