Gebrauchsanleitung TPP - Kryoröhrchen / *Instruction for Use TPP Cryo Tube*

Anleitung

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Instruction

Cryo tubes are designed to store samples in vapor phase of liquid nitrogen (LN₂) or in freezers or similar low-temperature equipment only.

- *Take product packaging and look for intact bag and product.*
- *Handling: The cryo rack # 99016 enables extra stability and single handling during filling. Due to the star sight locking system the tubes will not turn in the holder.*
- *Filling: Aseptically fill sample no higher than 90 % of the maximum graduated volume. For freezing, an expansion of the content in the closed tube must be allowed for.*
- *Identification: Use colored insert and/or white writing area*
- *Closure: Both threads of the tube must be absolutely dry, any dampness could impair tightness. Check correct and tight closure.*
- *Cooling: Rate of 1 °C/min.*
- *Storage: According to your SOPs.*
- *Thawing: Immerse the tube in a water bath 37 – 40 °C and keep it in motion. Remove the tube as soon as the content has melted and continue with your SOP.*

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Advice and Recommendation

1. If the tube is used in the liquid phase LN₂ may penetrate over time into the tube for example through an **improperly** sealed cap.
2. When retrieving the tube from liquid LN₂ to ambient temperature, take precautions for the risks of the immediate conversion of liquid nitrogen to vapor inside the tube, which could lead to tube cracking /exploding and dissemination of bio hazardous material. However, during storage in liquid LN₂ the tube shall be enclosed in an extra covering such as welded tubing. Always use appropriate safety equipment and procedures.
3. When undertaking cryogenic preservations, the tubes must be evenly exposed to freezing temperatures. Uneven temperature exposure can cause formation of ice plugs (i.e. at tube top) that inhibit the expansion of freezing liquid (i.e. at tube bottom), resulting in dangerous high pressure and subsequent harm or damage of tubes.
4. Never exceed maximum working volumes as specified in our product catalogue.
5. Centrifugation: perform test run with water at designated speed to assure performance, check for proper fit in the rotor. TPP covers no guarantee for centrifugation.

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Instruction

Cryo box L, # 99014 / Cryo box S, # 99015

Designed for storage of cryo tubes and most other similarly designed vials in ultralow temperature.

- *Working temperature: Suitable temperature range of +121 to -196 °C, autoclavable.*
- *Identification: Permanent alpha-numerical scale placed on the base as well as on the outer side of box. Clear lid for quick identification of vials.*
- *Storage: Fits most existing stainless steel racks.*
- *Size of tubes held: Two lid sizes available, S fits up to 2.0 ml, L fits tubes up to 5.0 ml.*

Cryo Rack # 99016

Designed for single handed manipulation of TPP cryo tubes.

- *The locking system in the grooves of the rack in combination with the star shaped foot of the TPP tubes stop twisting of the tube and enables single handed operation.*

Abmessungen	Measurements	99014	99015	99016
Typ	Type	Box	Box	Rack
Anordnung Reihe	Array	9 x 9	9 x 9	4 x 10
Total Röhrchen Stk.	Total tubes held pcs	81	81	40
Länge mm	Length mm	133	133	200
Breite	Width mm	133	133	100
Höhe inkl. Deckel	Height incl. lid mm	95	45	25
Deckelhöhe	Lid height mm	65	26	---
Passend für Röhrchen	Size of tubes	>5 ml	>2 ml	all
Material	Materials		PP	
Einsatz-Temperatur	Working temperature	+121 bis/to -196 °C		