

Panasonic

CO₂ incubators

Optimising cell culture productivity

MCO-230AIC



Optimum cell growth

Panasonic MCO-230AIC CO₂ incubators offer outstanding quality in performance for successful cell growth.



Discovery powered by
precision

Optimising cell culture productivity

We understand that creating successful cell cultures requires a CO₂ incubator that offers the highest levels of precision, security and ease of use. Our latest MCO-230AIC CO₂ incubators with Panasonic's innovative technologies offer outstanding quality in performance, maximise cell culture productivity and provide optimum results and reproducibility.



Providing a precisely controlled environment for sensitive cell cultures

Delivering long-term performance, optimal cell viability and successful experiments, each Panasonic incubator provides precise control of CO₂ concentration and temperature, while remaining easy to operate and maintain. The MCO-230AIC CO₂ incubators support a reliable, stable cell culture environment across all shelf positions, meaning each and every cell is safely maintained under ideal conditions.

DHA Direct heat and air jacket system

- Provides high precision temperature control for advanced uniformity and rapid recovery after door opening.

NEW PID Control of CO₂ and temperature

- Precise control guarantees exceptional performance and optimum results.

NEW Dual IR CO₂ sensor

- Minimises the effect of temperature and humidity changes during and after door openings for outstanding CO₂ control and fast recovery.

Optimum protection for your cell cultures

We know how valuable your cultures can be. Our advanced contamination control systems are designed to prevent the loss of your irreplaceable cell cultures with continuous background contamination control, all supported by a new security system.

InCu saFe® copper-enriched stainless steel interior

- Protects cell cultures by eliminating surface contamination sources and mitigating the effect of airborne contaminants.

Optional Safe Cell UV® with NEW increased UV lamp life

- Sterilises airborne and water pan contamination with 2.5 times longer lamp life.¹

NEW Optional electric door lock


- Provide enhanced security for your cultures and their conditions with an optional automatic electric door lock with password. Multiple registered users can have access with their individual ID and password, ensuring that valuable cell lines are only accessed by authorized persons.

1) Compared to the previous Panasonic MCO-20AIC CO₂ incubator.



NEW

Increased capacity
230 litre



Increasing
work efficiency



We have designed our incubators with ease of use and efficiency in mind. By delivering a user friendly cell culture incubator with rapid systems and processes, Panasonic can help make your work as simple as possible.

NEW More space for more cultures

- In a laboratory environment it is important to make the most of all the space available. With new integrated shelf supports, the Panasonic MCO-230AIC CO₂ incubators provide space for up to 25% more culture vessels².

NEW Integrated shelf supports

- Save valuable time and reduce the risk of contamination with an easy to clean incubator interior featuring fully rounded corners and integrated shelf supports.

NEW USB port

- Optimise cell culture protocols and adhere to standard operating procedures by conveniently transferring data to a USB memory stick to pass on to a PC. Logged parameters include chamber temperature, CO₂ level, door open status and alarms.




NEW Full colour LCD touch screen

- The user friendly full colour LCD touch screen provides clear information of incubator status with easy access to controller functions for fast, convenient set-up.


NEW Rapid sterilisation cycle

- Panasonic's H₂O₂ vapour sterilisation cycle reduces downtime to less than 3 hours for complete, validatable decontamination for increased productivity.



²) The MCO-230AIC CO₂ incubators can accommodate 25% more BD Falcon 225 cm² flasks per shelf compared to the previous Panasonic MCO-20AIC CO₂ incubator.

Specifications and options



- The MCO-230AIC is designed for stacking, allowing one unit to be positioned on top of another, doubling interior volume without additional floor space.
- An optional roller base is available for stacked installations for easier mobility.

See table below for details.

Double stacking table

Spacer for double-stacking		Lower unit
		MCO-230AIC
Upper unit	MCO-230AIC	MCO-170PS-PW
	MCO-19AIC(IM)/MCO-170AIC	MCO-230SB-PW
	MCO-18AC	
	MCO-20AIC	MCO-230SB-PW

The MCO-230AIC can also be stacked on top of a MCO-20AIC by using the stacking plate MCO-230SB-PW.



The MCO-230AIC series are certified as a Class IIa Medical Device (93/42/EEC and 2007/47/EC) for medical purposes of culturing cells, tissues, organs and embryos.

Notes:

- 1) Exterior dimensions of main cabinet only, excluding handle and other external projections
- See dimensions drawings on website for full details
- 2) Ambient temp 23°C, setting 37°C, CO₂ 5%, no load
- 3) Nominal value

MCO-230AIC CO₂ incubators

Characteristics

- V = Visual alarm
- B = Buzzer alarm
- R = Remote alarm



MODEL	MCO-230AIC-PE (standard)	MCO-230AICUV-PE (standard + UV)	MCO-230AICUVH-PE (standard + UV + H ₂ O ₂)
Dimensions			
External dimensions (WxDxH) ¹⁾	mm	770 x 730 x 905	
Internal dimensions (WxDxH)	mm	643 x 523 x 700	
Volume	litres	230	
Net weight (approx)	kg	90	
Performance			
Temperature control range and fluctuation	°C	AT +5 ~ +50, ±0.1	
Temperature uniformity ²⁾	°C	±0.25	
CO ₂ control range and fluctuation	%	0 ~ 20, ±0.15	
Humidity level and fluctuation	%RH	95, ±5	
Control			
Temperature sensor		Thermistor	
CO ₂ sensor		Dual IR	
Display		LCD touch screen	
Construction			
Exterior material		Painted steel (rear cover not painted)	
Interior material		SS copper alloyed	
Insulation material		Extruded polystyrene	
DHA heating system		Y	
Outer door	qty		1
Outer door lock		option	option
Reversible door			Y
Inner door	qty		1
Shelves	qty		4
Max. load per shelf	kg		7
Max. total load	kg		20
Max. shelf capacity	qty		10
Access port	qty		1
- position			Rear
- diameter	∅ mm		30
Alarms			
Power failure			R
Out of temperature setting			V-B-R
High temperature			V-B-R
Out of CO ₂ setting			V-B-R
Door open			V-B
Electrical and noise level			
Power supply	V		230
Frequency	Hz		50
Noise level ³⁾	dB		25
Options			
SafeCell UV [®] system		MCO-170UVS-PE ⁴⁾	standard
H ₂ O ₂ decontamination board		MCO-170HB-PE ⁴⁾	MCO-170HB-PE ⁴⁾
Electric door lock with password		MCO-170EL-PW ⁴⁾	MCO-170EL-PW ⁴⁾
H ₂ O ₂ vapour generator			MCO-HP-PW ⁴⁾
H ₂ O ₂ reagent, pack of 6 bottles			MCO-H202-PE
CO ₂ gas pressure regulator			MCO-100L-PW
Automatic CO ₂ cylinder changeover system			MCO-21GC-PW
Semi-automatic one point gas calibration kit			MCO-SG-PW
InCu saFe [®] shelf			MCO-230ST-PW
InCu saFe [®] half tray system			MCO-35ST-PW
Double stacking bracket *			MCO-170PS-PW
Stacking plate *			MCO-230SB-PW
Roller base			MCO-230RB-PW
Optional communication systems⁵⁾			
Ethernet interface (LAN)			MTR-L03-PW
Digital interface (RS232C/RS485)			MTR-480-PW

- 4) MCO-230AIC series requires MCO-170HB-PE, MCO-170EL-PW, MCO-HP-PW and SafeCell UV option for H₂O₂ decontamination
- 5) MCO-230AIC series can only be fitted with one communications interface

* If stacking two incubators, make sure the double-stacking dedicated securing hardware and spacer are used. (please see double stacking table)

Appearance and specification are subject to change without notice.

for more information:

www.biomedical.panasonic.eu