# PHCbi

## Operating Instructions

Ultra-Low Temperature Freezer

### MDF-DU300H

Please read the operating instructions carefully before using this product, and keep the operating instructions for future use.

See page 36 for model number.

# CONTENTS

INTRODUCTION .....	3
PRECAUTIONS FOR SAFE OPERATION .....	4
LABELS ON UNIT .....	9
SYMBOLS ON UNIT .....	10
ENVIRONMENTAL CONDITIONS .....	10
INTENDED USE AND PRECAUTIONS .....	11
FREEZER COMPONENTS	
Main body .....	12
Control panel .....	14
INSTALLATION SITE .....	16
INSTALLATION .....	17
START-UP PROCEDURE .....	19
DURING/AFTER POWER FAILURE .....	20
SETTING THE CHAMBER TEMPERATURE .....	21
LOCKING THE CHAMBER TEMPERATURE SETTING .....	22
SETTING THE ALARM TEMPERATURE	
Setting the high-temperature alarm .....	23
Setting the low-temperature alarm .....	24
SETTING THE BUZZER SUSPENDED PERIOD .....	25
SETTING THE COMPRESSOR DELAY .....	26
REMOTE ALARM TERMINAL .....	27
ALARM FUNCTIONS .....	28
SAFETY FUNCTIONS .....	28
SELF-DIAGNOSTIC FUNCTIONS .....	29
ROUTINE MAINTENANCE	
Cleaning the exterior, interior, and accessories .....	30
Cleaning the filter .....	30
Defrosting the chamber .....	31
Checking the alarm operation .....	31
REPLACEMENT OF WORN-OUT PARTS	
Replacement of the battery for the power-failure alarm .....	32
Replacement of the fan motor .....	32
CALIBRATION .....	32
TROUBLESHOOTING .....	33
OPTIONAL COMPONENTS	
Temperature recorder .....	34
Back-up cooling kit .....	34
Inventory rack .....	34
SPECIFICATIONS .....	35
PERFORMANCE .....	36
EMC PERFORMANCE .....	36
SAFETY CHECK SHEET .....	37

# INTRODUCTION

---

- Read the operating instructions carefully before using the product and follow the instructions for safe operation.
- PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.
- Keep the operating instructions in a suitable place so that they can be referred to as necessary.
- The operating instructions are subject to change without notice for improvement of performance or function.
- Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.
- No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

## **IMPORTANT NOTICE**

PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for any loss or damage to the contents of the product.

# PRECAUTIONS FOR SAFE OPERATION

**It is imperative that the user complies with the operating instructions as they contains important safety advice.**

Items and procedures are described so that you can use this unit correctly and safely. Following these precautions will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:

## **WARNING**

Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

## **CAUTION**

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

Symbols have the following meanings:

-  This symbol means caution.
-  This symbol means an action is prohibited.
-  This symbol means an instruction must be followed.

Be sure to keep the operating instructions in a place that is accessible to users of this unit.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Do not use the unit outdoors.** Exposure to rain may cause leakage and/or electric shock.
-  **Only qualified engineers or service personnel should install the unit.** Installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit in a location capable of bearing the total combined weight (product + optional accessories + stored items). After installing the unit, be absolutely sure to take precautions to prevent the unit from falling over.** If the unit is installed in a location which is not strong enough or if the proper precautions are not taken, the unit may fall over and cause injuries.
-  **Do not install the unit where there are high levels of moisture or where it may be splashed with water.** Installing the unit where there are high levels of moisture or where it may be splashed with water may cause the insulation to deteriorate and give rise to leakage and/or electric shock.
-  **Do not install the unit in a location where flammable or volatile substances are present.** Installing the unit in a location where flammable or volatile substances are present may cause explosions and/or a fire.
-  **Do not install the unit in a location where corrosive gases such as acids are present.** Installing the unit in a location where corrosive substances are present may cause electrical components to corrode, leading to leakage and/or electric shock due to the deterioration of insulation resulting from corroded electrical components.
-  **Do not place this unit in a location where it is difficult to disconnect the power supply plug.** Failure to disconnect the power supply plug may cause fire in the event of a problem or malfunction.
-  **Be absolutely sure to earth (ground) the unit to prevent electric shock.** Failure to earth the product may give rise to electric shock. If necessary, ask a qualified contractor to do this work.
-  **Do not connect the earth wire to a gas pipe, water pipe or lightning rod when earthing the unit.** Earthing the unit improperly may give rise to electric shock.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **Never store volatile or flammable substances in this unit except in a sealed container.** Such substances may cause explosion or fire if they leak.
-  **Never insert metal objects such as pins and wires into any vent, gap, or outlet on the unit.** This may cause electric shock or injury by accidental contact with moving parts.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **When handling harmful samples (for example, those which consist of toxic, pathogenic or radioactive substances), install the unit inside a designated isolation facility.** If the unit is installed in a location which is not an isolation facility, there may be detrimental effects on both people and the natural environment.
-  **Before proceeding with maintenance or checking the unit, set the power switch to OFF, and disconnect the power supply plug.** Performing the work while power is still flowing to the product or while the power supply plug is still connected may give rise to electric shock and/or injury.
-  **Do not touch any electrical parts (such as power supply plug) or operate switches with a wet hand.** This may cause electric shock.
-  **Wear protective gloves and mask during maintenance.** Touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.
-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never put containers with liquid on top of the unit** as this may cause electric shock or short circuit if the liquid is spilled.
-  **Never damage the power supply cord or power supply plug (by breaking, adapting, placing near a source of heat, bending with force, twisting, pulling, adding weight, or binding).** A damaged power supply cord or power supply plug may cause electric shock, short circuit, or fire
-  **Never disassemble, repair, or modify the unit yourself.** A high-voltage area is located inside the unit. Any work carried out by unauthorized personnel may result in electric shock. Contact our sales representative or agent for maintenance or repair.
-  **Make sure the power supply plug is pushed fully in.** Faulty insertion of the power supply plug may cause electric shock or fire due to generation of heat. Never use a damaged power supply plug or loose power outlet
-  **Disconnect the power supply plug if there is anything wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  **Grip the power supply plug when disconnecting the power supply cord from the outlet.** Pulling the power supply cord may cause electric shock or short circuit.
-  **Remove dust from the power supply plug periodically.** Dust on the power supply plug may cause insulation failure due to moisture and thus cause a fire. Disconnect the power supply plug and wipe it with a dry cloth

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

 **Disconnect the power supply plug before moving the unit.** Take care not to damage the power supply cord. A damaged power supply cord may cause electric shock or fire.

 **Disconnect the power supply cord when the unit is not in use for long periods.** Keeping the unit connected may cause electric shock, leakage, or fire due to the deterioration of insulation.

 If the unit is to be stored unused in an unsupervised area for a long period, **ensure that children do not have access and that doors cannot be closed completely.**

 **Ask a qualified contractor to carry out disassembly and disposal of the unit.** Leaving the unit in a location that can be accessed by third parties may result in unexpected accidents (e.g. the unit may be used for unintended purposes).

 **Do not leave the plastic bags used for packing in a place where they can be reached by small children** as this may result in unexpected accidents such as suffocation.

 **Never replace the battery for the power-failure alarm yourself.** Only qualified engineers or service personnel should replace the battery.

 **When moving the unit, be sure to take precautions to prevent it from falling over.** Moving the unit with too much force may cause it to fall over, possibly resulting in injury. A qualified individual must be assigned to supervise the safe movement and relocation of the unit.

 **Install the unit in a well-ventilated (airy) location to prevent the accumulation of flammable refrigerant.** The flammable refrigerant may cause fire if it leaks.

 **Never damage the chamber wall or pipework in the chamber when removing frost.** The refrigerant is flammable and may cause a fire if it leaks.

 **Flammable and explosive product.** The unit contains flammable refrigerant. Be sure to follow the below instructions when servicing or recycling.

- Well ventilate the room to prevent refrigerant accumulation.
- Keep fire away when the refrigerant is contained in the product.
- Do not damage or break the pipework.

As with any equipment that uses CO<sub>2</sub> gas, there is a likelihood of oxygen depletion in the vicinity of the equipment. It is important that you assess the work site to ensure there is suitable and sufficient ventilation. If lack of ventilation is suspected, then other methods of ensuring a safe environment must be considered. These may include atmosphere monitoring systems and warning devices with alarms.

 **Do not touch the condenser directly when the filter is removed for cleaning.** Touching the condenser may cause injury due to its hot surface.

# PRECAUTIONS FOR SAFE OPERATION


## CAUTION

-  **Never install the unit in a location where corrosive materials such as sulphur compounds are likely to be generated (e.g. near a drainage facility).** Corrosion of the copper pipes may result in the deterioration and consequently the failure of the cooling unit.
-  **This unit must be plugged into a dedicated circuit protected by a branch circuit-breaker.**
-  **Use a dedicated power source as indicated on the rating label attached to the unit.** A multiple-tap source may cause fire due to excessive heating.
-  **Do not climb on top of the unit or put any objects on the unit.** Falling from the unit may cause injury; falling objects may cause damage to the unit
-  **Never store corrosive substances such as acids or alkalis in this unit except in a sealed container.** These may be harmful to your health and may cause corrosion of internal components or electrical parts.
-  **Check the settings when restarting operation after a power failure or after turning the power off.** The settings may have changed as a result of stopping the unit. Stored items inside the unit may be adversely affected when operation is resumed if the settings have changed.
-  **To ensure the safety of the service engineer, submit a safety check sheet with the required items filled out.** This is provided as the photocopyable "Safety Check Sheet" at the end of these operating instructions

# LABELS ON UNIT

<Labels applied to the unit>

To avoid accidents, users are advised to read carefully the hazard labels found at key locations on the interior and exterior of the unit.

Possible Danger	Warning/Caution Type Location of Danger	Warning/Caution Label	Description of Danger
Personal injury	<b>Flammable and explosive product</b> Exterior, Interior		This product contains flammable refrigerant. Be sure to follow the instructions on the warning label when servicing or recycling.
Personal injury	<b>Outer door latch</b> Exterior		Take care not to trap your hand in the outer door latch.
Personal injury Sample damage	<b>Frostbite Rise in chamber temperature</b> Interior		To prevent frostbite, wear protective gloves when handling frozen items in the chamber. Too much frost may cause chamber temperature rise resulting from incomplete door close.
Sample damage	<b>Chamber temperature</b> Interior		Always use the outer door latch when closing the outer door. Improper closure may cause the temperature to rise in the chamber.
Sample damage	<b>Chamber temperature</b> Interior		Clean the filter once a month. A dusty filter may cause poor cooling performance.

# SYMBOLS ON UNIT

The following symbols are attached to the unit. The table describes the meaning of the symbols.

	This symbol is attached to covers that access high-voltage electrical components to prevent electric shock. Only a qualified engineer or service personnel should be allowed to open these covers.
	This symbol indicates that caution is required. Refer to product documentation for details.
	This symbol indicates Incorrect usage could lead to a fire hazard.
	This symbol indicates an earth.
	This symbol means "ON" for a power switch.
○	This symbol means "OFF" for a power switch.

# ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):

- Indoor use;
- Altitude up to 2000 m;
- Ambient temperature 5 °C to 40 °C;
- Maximum relative humidity 80 % for temperature up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C;
- Mains supply voltage fluctuations up to  $\pm 10$  % of the nominal voltage;
- Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;
- Temporary OVERVOLTAGES occurring on the mains supply;
- Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases);

# INTENDED USE AND PRECAUTIONS

---

This equipment is designed for low temperature storage of human cells, organs, plasma and DNAs.

Temperature and duration of storage:

cells: 1month - 1year at -80 °C

organs: 11months at -80 °C

DNA: Long term (8 years) at -80 °C - -70 °C


plasma: 2-3 months at -80 °C

- The effective storage period depends on the sample condition and storage temperature. It is necessary to determine the storage temperature and period suitable for the purpose.
- For the live cells, a lower storage temperature should be required for long-term storage. It is recommended to store the live cells at -130 °C or lower.

DOMINIQUE DUTSCHER SAS

# FREEZER COMPONENTS

## Main body

# FREEZER COMPONENTS

**1. Inner door:** This prevents cold air from escaping when the outer door is opened. Always be sure to close the inner door securely before closing the outer door. The inner door can be removed for cleaning or defrosting [page 31].

**2. Outer door latch:** Always lock the outer door latch when the outer door is closed. Provision has been made for use of an additional padlock (not included).

**3. Inner door latch:** Always lock the inner door latch when the inner door is closed.

**4. Access port <rear and bottom>:** These ports are used to pass the sensor or cable of measuring equipment, the sensor of a temperature recorder (optional), or the nozzle of a back-up cooling kit (optional) to the chamber.

✧ Replace the insulation and the rubber caps when the access port is not in use. Improper replacement may cause an increase in chamber temperature or condensation around the access port.

**5. Air intake vent:** This is an air intake vent for circulating the air in the chamber. Do not block this vent.

✧ Blocking this vent may cause unstable chamber temperature.

✧ Do not insert fingers or other objects into this vent.

**6. Levelling foot <bottom>:** These are screw bolts used to install and fix the unit. Adjust the height of the levelling feet by turning the screw bolts until the two front casters are away from the floor [page 17].

**7. Filter <behind the grille>:** This filter prevents dust from accumulating on the condenser. A dusty filter may cause failure of the refrigerating device. Clean the filter once a month [page 30].

**8. Space for temperature recorder:** A temperature recorder (optional) can be mounted here so that the chamber temperature can be recorded automatically [page 34].

✧ Contact our sales representative or agent for the installation.

**9. Control panel:** The chamber temperature and other alarms/functions can be set using the keys on the control panel. The operational status can be checked on the temperature display and indicators [pages 14 and 15].

**10. Fixture:** Use the fixtures and secure the unit to a wall with a strong rope or chain [page 17].

**11. Power switch:** This is the power switch of the unit. (ON=" | ", OFF="O")


**12. Remote alarm terminal:** A remote alarm device (separately available) can be connected to this terminal. The remote alarm relays the alarm to an operator in a remote location if the unit is unattended. Refer to page 27 for details.

✧ Contact our sales representative or agent to arrange the installation.

**13. Battery switch:** This is an ON-OFF switch for the battery for the power-failure alarm. Always turn this switch on when the unit is operating to ensure that the power-failure alarm is working. Turn this switch off when the unit is not used for a long period in order to protect the battery.

# FREEZER COMPONENTS

## Control panel

**1. Temperature display:** In normal operation, this displays the current chamber temperature. During temperature alarms, the display blinks [page 28]. When a self-diagnostic warning is triggered, a code is displayed here [page 29].

**2. Battery indicator (BATTERY):** The orange LED lamp lights up to indicate that the battery needs to be replaced (every 3 years). The orange LED lamp blinks to indicate that the fan motor needs to be replaced (every 6 years). Refer to page 29 for details.

To arrange replacement, contact our sales representative or agent.


**3. Alarm indicator (ALARM):** The red LED lamp blinks in the event of an alarm [page 28].

**4. Filter check indicator (FILTER):** The orange LED lamp lights up and the buzzer sounds when excessive dust has accumulated on the filter. When this indicator lights up and the buzzer sounds, clean the filter following the procedure described on page 28.

◇“temperature display mode”: the status when the current chamber temperature is displayed on the temperature display.

◇“setting mode”: the status when the input on the temperature display can be changed by pressing the up arrow key for more than 5 seconds.

# FREEZER COMPONENTS

**5. Up arrow key (▲):** During “temperature display mode”: pressing this key for more than 5 seconds activates the setting mode. During “setting mode”: pressing this key allows the (blinking) figure on the temperature display to be changed [pages 21 - 26]. The lock for the chamber temperature setting can also be set to ON or OFF by pressing this key [page 22].

**6. Scroll key (▶▶):** During “temperature display mode”: pressing this key for more than 5 seconds enables “locking of the chamber temperature setting” [page 22]. During “setting mode”: pressing this key moves the user between the different (blinking) digits to be entered on the temperature display [pages 21 - 26].

**7. Set key (SET):** Pressing this key activates the chamber temperature setting mode [page 21]. Press this key to enter the required value during the setting procedure [pages 21 - 26].

**8. Alarm test key (ALARM TEST):** This key is used to check the operation of buzzer, alarm indicator and remote alarm terminal [page 31].

**9. Buzzer stop key (BUZZER):** The buzzer is silenced by pressing this key when the alarm indicator blinks and the buzzer sounds. (The remote alarm cannot be cancelled.)

◇“temperature display mode”: the status when the current chamber temperature is displayed on the temperature display.

◇“setting mode”: the status when the input on the temperature display can be changed by pressing the up arrow key for more than 5 seconds.

# INSTALLATION SITE

This unit must be installed in a location which meets all the conditions described below.

◇ If the unit is installed in a location which does not meet the conditions, its specified performance may not be achieved or malfunctions and accidents may occur.

## ■ A location not exposed to direct sunlight

Avoid any location which is exposed to direct sunlight. Installing the unit in a location exposed to direct sunlight may reduce its cooling performance.

## ■ A well-ventilated (airy) location

In order to ensure ventilation, leave clearances of at least 10 cm around the unit (at the left, right, back and top). Blocking the ventilation may reduce the unit's cooling performance or cause malfunctions.

## ■ A location away from sources of heat

Avoid any location which is close to a major source of heat (such as a heater or boiler). Installing the unit near a major source of heat may reduce the unit's cooling performance.

## ■ A location with minimal changes in temperature

Avoid any location where the ambient temperature is subject to sudden changes. If the unit is installed in a location where the ambient temperature is subject to sudden changes, it will not be possible to achieve a stable cooling performance.

## ■ A flat surface where the floor is also capable of bearing the total combined weight (product + optional accessories + stored items)

Install the unit on a flat surface which is even and which is capable of bearing the total combined weight (product + optional accessories + stored items). If the unit is installed where the surface is uneven or where the unit will be inclined at an angle, the unit will be unstable, and accidents or injuries may occur and/or unnecessary vibration or noise may be generated.

## ■ A location with minimal humidity

Install the unit in a location where the relative humidity is less than 80 % R.H. Installing the unit in a very humid location may cause earth faults and/or electric shock.

## ■ A location free of flammable or corrosive gases

Avoid any location exposed to flammable or corrosive gases. Flammable or corrosive gases can cause explosions and/or a fire. Furthermore, corrosion of the electrical parts may cause the insulation to be reduced and result in earth faults and/or electric shock.

## ■ A location where nothing can fall onto the unit

Avoid locations where objects may fall onto the unit. Objects falling and hitting the unit may cause it to break down or fail.

# INSTALLATION

When installing the unit, follow the steps below to secure the unit properly, and also be absolutely sure to earth the unit.

✧In addition, install an earth leakage circuit breaker (on the unit's power supply side), which is mandatory under the applicable laws and regulations.

## 1. Preparations after unpacking

Remove all the tape used to secure the doors and interior parts, and leave the doors open for a short while for ventilation.

If any surfaces of the outer cabinet are dirty, wipe the surface using a cloth moistened with a diluted neutral dish-washing detergent.

✧Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.

✧After wiping the unit using the diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth, allowing the surfaces of the outer cabinet to dry out completely, and then proceed with the installation.

### Note:

Remove the cable tie that bands the power supply cord. Prolonged contact with the tie may cause corrosion of the cord coating.



## 2. Securing and levelling the unit using the levelling feet

Rotate the front levelling feet clockwise until the casters are raised 5 mm to 10 mm above the floor surface [Fig. 1]. In addition, rotate the levelling feet slightly clockwise or anticlockwise, and adjust them so that the unit is completely level.

✧When the casters are raised from the floor surface, the unit will be secured. If they are left touching the floor, the unit may accidentally move when its door is opened or closed.

## 3. Securing the unit by using the fixtures

Use the fixtures on the rear panel of the unit, and secure the unit to a wall with a strong rope or chain [Fig. 2].

# INSTALLATION

## 4. Preventing electric shock by earthing the unit

When installing the unit, be absolutely sure to earth (ground) it. Earthing is necessary to prevent electric shocks resulting from deterioration of electrical insulations

✧ This unit comes with a 3-pin plug having one earth pin. Earthing work is not required in the case of a 3-pin power outlet equipped with an earth contact.

✧ If the power outlet is not a 3-pin outlet equipped with an earth contact, ask a qualified contractor to do the earthing work.

## 5. Setting up the shelves

Three shelves are packaged at the bottom of the chamber. Set the shelves firmly in place on the shelf stoppers at the standard locations [Fig. 3].

✧ The three lips of each shelf should be facing upwards at the back and sides of the chamber.

## 6. Installing an earth leakage circuit breaker

Install an earth leakage circuit breaker (on the unit's power supply side), which is mandatory under the applicable laws and regulations.

✧ Contact our sales representative or agent to arrange the installation of an earth leakage circuit breaker.


Fig. 3

# START-UP PROCEDURE

Follow this procedure for the initial operation of the unit and for consequent operations (after temporary stoppage for cleaning, maintenance or moving) .

◇After a power failure, the unit will restart operation automatically with the same settings as before the power failure [page 20]

1. Check that the following switches are turned off: [power switch, battery switch, switch of the optional back-up cooling kit (if installed)] .

2. Connect the power supply cord to the dedicated power source with the appropriate rating with the chamber empty.

3. Turn on the power switch .

▶ When the power is switched on, the high-temperature alarm is activated. The alarm indicator and the temperature display blink. (This is not a malfunction.)

◇The high-temperature alarm will be cancelled when the chamber temperature reaches the approximate set temperature [(factory setting): -80 °C].

◇When operation of the unit is started in high ambient temperature, the filter check indicator may lights up. The indicator will turn off automatically when the ambient temperature decreases.

4. Turn on the battery switch.

#### <Important>

Always turn the battery switch on when the unit is operating to ensure that the power-failure alarm is working.

5. Set the desired chamber temperature [page 21].

◇The factory setting of chamber temperature is -80 °C.

6. Using the temperature display, check that the chamber temperature has cooled to the set temperature.

◇Check that the chamber temperature falls to the set temperature when the start-up after cleaning, maintenance or moving.

7. Turn on the switch of the optional back-up cooling kit (if installed).

8. Press the alarm test key (ALARM TEST) to check the alarm operation.

▶ The alarm indicator blinks and the buzzer sounds. The remote alarm is also operated.

◇The alarm indicator (ALARM) blinks and "E09" appears on the temperature display if the battery switch is off.

9. Press the test switch of the optional back-up cooling kit (if installed) to check it is working.

10. Gradually place the material inside the chamber.

◇Putting a large amount of material into the chamber at one time causes the temperature to rise.

11. Set the alarm temperature [pages 23 - 24] and the buzzer suspended period [page 25], lock the chamber temperature setting [page 22], and the compressor delay [page 26] as required.

# DURING/AFTER POWER FAILURE

## ■ Display of chamber temperature during power failure

By pressing the buzzer stop key (BUZZER) during a “power-failure alarm”, the buzzer stops and the chamber temperature is displayed on the temperature display for 5 seconds.

◇ Check the chamber temperature as appropriate by pressing the buzzer stop key (BUZZER) during power failure.

## <Important>

The battery for the power-failure alarm has a limited lifespan. Replace the battery every 3 years. The alarm may not be activated if the power fails and the battery has not been replaced regularly. Contact our sales representative or agent to arrange replacement of the battery.

## ■ Operation check after recovery from power failure

After recovery from a power failure, the operation will resume automatically with the same settings as before the power failure. There is no need to re-set. However, always check the operational status after recovery.

◇ The settings are memorized during power failure by nonvolatile memory.

## <Important>

Upon recovery from a power failure, check that the unit restarts operation properly. Also check that there are no changes to the settings.

# SETTING THE CHAMBER TEMPERATURE

Set the chamber temperature as required to keep the stored material at an appropriate temperature for a long period of time.


- Setting range of chamber temperature: between -90 °C and -50 °C

- Initial setting (factory setting): -80 °C

◇The guaranteed minimum chamber temperature is -85 °C at an ambient temperature of 30 °C with no load, although the set value for the chamber temperature can range from -90 °C to -50 °C.

For example, to change the setting of chamber temperature to -75 °C

(The following is only an example. Choose the actual setting according to the desired chamber temperature.)

	Operation	Key operated	Display after operation
1	Plug-in the power supply plug and turn on the power switch. (only when start-up of unit)	----	The current chamber temperature blinks. 
2	Press the "Set" key once.	SET	The current setting (-80) is displayed and the first figure (8) blinks. 
3	Press the up arrow key (as many times as needed to reach the desired figure).	▲	The display changes from -80 to -70. 
4	Press the scroll key once.	▶▶	The second figure (0) starts to blink. 
5	Press the up arrow key (as many times as needed to reach the desired figure).	▲	The display changes from -70 to -75. 
6	Press the "Set" key once.	SET	The new setting is memorized and the current chamber temperature blinks. 

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed.

# LOCKING THE CHAMBER TEMPERATURE SETTING

The chamber temperature setting can be locked to avoid accidental changes. When the lock is ON, a change of chamber temperature setting is not accepted even if the key on the control panel is operated.

- Initial setting (factory setting): lock OFF

Display	Lock status	Setting of chamber temperature
L 0	Lock OFF	Enabled to change
L 1	Lock ON	Disabled from change

For example, to change the lock for chamber temperature setting from OFF (factory setting) to ON  
(The following is only an example. Operate the key according to the desired lock status.)

	Operation	Key operated	Display after operation
1		----	The current chamber temperature is displayed. 
2	Press the scroll key for 5 seconds.		The current chamber temperature changes to L 0 and the figure (0) blinks. 
3	Press the up arrow key once.		The display changes from L0 to L1. 
4	Press the "Set" key once.	SET	The lock is set to ON and the current chamber temperature is displayed. 

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed.

# SETTING THE ALARM TEMPERATURE

## Setting the high temperature alarm

An abnormality (chamber temperature rise) is indicated by the alarm indicator, and the chamber temperature display blinking, and the buzzer sounding (15 minutes after blinking) if the chamber temperature exceeds the set value of the high-temperature alarm. Always set the high-temperature alarm to protect stored items from damage resulting from an increase in temperature.

- Setting range for high-temperature alarm:  
between 5 °C and 40 °C higher than the chamber set temperature
- Initial setting (factory setting): 10 °C higher than the chamber set temperature

For example, to change the alarm temperature from chamber set temperature plus 10 °C (i.e. the high-temperature alarm works at 10 °C higher than the chamber set temperature) to chamber set temperature plus 5 °C (i.e. the high-temperature alarm works at 5 °C higher than the chamber set temperature)

(The following is only an example. Choose the actual setting according to the desired alarm temperature.)

	Operation	Key operated	Display after operation
1		----	The current chamber temperature is displayed.
2	Press the up arrow key for 5 seconds.	▲▲	The current chamber temperature changes to function code F00 and the right-hand figure (0) blinks.
3	Press the up arrow key once.	▲	The display changes from F00 to F01.
4	Press the "Set" key once.	SET	The current setting (010) is displayed and the right-hand figure (0) blinks.
5	Press the up arrow key (as many times as needed to reach the desired figure).	▲	The display changes from 010 to 015.
6	Press the scroll key once.	▶▶	The middle figure (1) starts to blink.
7	Press the up arrow key (as many times as needed to reach the desired figure).	▲	The display changes from 015 to 005.
8	Press the "Set" key once.	SET	The new setting of the high-temperature alarm is memorized and the current chamber temperature is displayed.

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed.

# SETTING THE ALARM TEMPERATURE

## Setting the low temperature alarm

An abnormality (chamber temperature drop) is indicated by the alarm indicator, and the chamber temperature display blinking, and the buzzer sounding (15 minutes after blinking) if the chamber temperature falls below the set value of low-temperature alarm. Always set the low-temperature alarm to protect stored items from damage resulting from a drop in temperature.

- Setting range for low-temperature alarm:  
between 5 °C and 40 °C lower than the chamber set temperature
- Initial setting (factory setting): 10 °C lower than the chamber set temperature

For example, to change the alarm temperature from chamber set temperature minus 10 °C (i.e. the low-temperature alarm works at 10 °C lower than the chamber set temperature) to chamber set temperature minus 5 °C (i.e. the high-temperature alarm works at 5 °C lower than the chamber set temperature)

(The following is only an example. Choose the actual setting according to the desired alarm temperature.)

	Operation	Key operated	Display after operation
1		----	The current chamber temperature is displayed.
2	Press the up arrow key for 5 seconds.	▲▲	The current chamber temperature changes to function code F00 and the right-hand figure (0) blinks.
3	Press the up arrow key twice.	▲▲	The display changes from F00 to F02.
4	Press the "Set" key once.	SET	The current setting (-10) is displayed and the right-hand figure (0) blinks.
5	Press the up arrow key (as many times as needed to reach the desired figure).	▲▲	The display changes from -10 to -15.
6	Press the scroll key once.	▶▶	The middle figure (1) starts to blink.
7	Press the up arrow key (as many times as needed to reach the desired figure).	▲▲	The display changes from -15 to -05.
8	Press the "Set" key once.	SET	The new setting of the low-temperature alarm is memorized and the current chamber temperature is displayed.

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed.

# SETTING THE BUZZER SUSPENDED PERIOD

An abnormality is indicated again by the buzzer sounding after a certain period has elapsed (the buzzer recovery time) if the alarm mode continues even after the buzzer has already been silenced by pressing the buzzer stop key (BUZZER). Always set the "buzzer suspended period" to avoid incorrect recognition of alarm mode.

■ Setting range of buzzer suspended period:

between 10 and 60 minutes in 10-minute intervals (display of setting: between 010 and 060) or no recovery (display of setting; 000)

■ Initial setting (factory setting): 30 minutes (display of setting: 030)

◇The buzzer does not sound again when the setting is 000. The setting should be between 010 and 060 to ensure the safety of stored items.

For example, to change the setting of the buzzer suspended period from 30 minutes (initial setting) to 20 minutes

(The following is only an example. Choose the actual setting according to the desired buzzer suspended period.)

	Operation	Key operated	Display after operation
1		---	The current chamber temperature is displayed.
2	Press the up arrow key for 5 seconds.		The current chamber temperature changes to function code F00 and the right-hand figure 0 blinks.
3	Press the up arrow key 5 times.		The display changes from F00 to F05.
4	Press the scroll key once.		The middle figure (0) starts to blink.
5	Press the up arrow key (as many times as needed to reach the desired figure).		The display changes from F05 to F25.
6	Press the "Set" key once.	SET	The current setting (030) is displayed and the middle figure (3) blinks.
7	Press the up arrow key (as many times as needed to reach the desired figure).		The display changes from 030 to 020.
8	Press the "Set" key once.	SET	The new setting is memorized and the current chamber temperature is displayed.

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed.

◇It is recommended that the buzzer suspended period is set when the freezer is not under alarm conditions.

# SETTING THE COMPRESSOR DELAY


The delay time before the compressor starts can be changed to reduce the load on the power line and to facilitate the start-up (reset) of the freezer during recovery after power failure.

- Setting range of delay time: between 3 and 15 minutes, in 1-minute intervals
- Initial setting (factory setting): 3 minutes

**Note:**

There is no need to change the delay time when the capacity of the power source is adequate.

For example, to change the setting of the compressor delay from 3 minutes (initial setting) to 4 minutes (The following is only an example. Choose the setting according to the desired start delay.)

	Operation	Key operated	Display after operation
1		---	The current chamber temperature is displayed. 
2	Press the up arrow key for 5 seconds.	▲	The current chamber temperature changes to function code F00 and the right-hand figure (0) blinks. 
3	Press the up arrow key 5 times.	▲	The display changes from F00 to F05. 
4	Press the "Set" key once.	SET	The current setting (003) is displayed and the right-hand figure (3) blinks. 
5	Press the up arrow key (as many times as needed to reach the desired figure).	▲	The display changes from 003 to 004. 
6	Press the "Set" key once.	SET	The new setting is memorized and the current chamber temperature is displayed. 

◇The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). In this case, the setting is not changed..

# REMOTE ALARM TERMINAL


The alarm is relayed to a remote location when a remote alarm device (commercially available) is connected to the remote alarm terminal. Installation of a remote alarm device is recommended when the unit is installed in an unattended location so that the operator of the unit is notified of the alarm.

✧ Contact our sales representative or agent to arrange the installation of a remote alarm device.

- Location of remote alarm terminal: lower right side
- Allowable contact capacity: DC 30 V•2 A

Contact output:

Connection	Normal status	Abnormal status
Between COM. and N.O.	open	closed
Between COM. and N.C.	closed	open

Use a twisted shielded wire for the connection.

Type: UL 2343, UL 2448, UL 2464, UL 2552, UL 2623

Length: 30 m max.

The remote alarm terminal will be in alarm mode when the power supply cord is unplugged or the power switch is turned off because this will be interpreted as a power failure.

✧ The operation of the remote alarm terminal is not conjunction with the alarm buzzer. Therefore, the operation of the remote alarm terminal cannot be canceled by pressing the alarm buzzer stop key (BUZZER).

# ALARM FUNCTIONS

This unit has the alarm functions listed below.

When an alarm goes off, check the cause on the basis of the type of alarm generated, and take the necessary actions to deal with it without delay.

◇Contact our sales representative or agent if the alarm is not resolved by eliminating the following causes.

◇Before calling for service, take some precautions to prevent the stored items (for example by transferring the items to another freezer).

Alarms	Cause	Indication	Buzzer	Remote alarm
High- temperature alarm	The chamber temperature is higher than the setting of the high-temperature alarm.	The alarm indicator blinks. The chamber temp. blinks.	Intermittent tone with a 15-minute delay.	Alarm mode with a 15-minute delay.
Low- temperature alarm	The chamber temperature is lower than the setting of the low-temperature alarm.	The alarm indicator blinks. The chamber temp. blinks.	Intermittent tone with a 15-minute delay.	Alarm mode with a 15-minute delay.
Power-failure alarm	A power failure has occurred. The power switch is OFF. The power supply cord is disconnected.	The alarm indicator blinks.	Intermittent tone.	----
Filter alarm	Excessive dust has accumulated on the filter.	The filter check indicator lights up.	Intermittent tone	----

# SAFETY FUNCTIONS

This unit has the safety functions listed below.

Safety	Cause	Indication, Buzzer	Safety operation
Auto-return	During setting mode, there is no key operation for 90 seconds.	----	Exit from setting mode and return to temperature display mode.
Key lock	The key lock is ON (L1).	----	Change of chamber temperature setting is disabled.
Battery switch check	The battery switch is OFF when pressing the alarm test key (ALARM TEST).	The alarm indicator blinks. E09 blinks. The buzzer sounds.	◇The battery switch must be ON to activate the alarm.

# SELF-DIAGNOSTIC FUNCTIONS

Alarms involving an error code indicated on the temperature display are generated by the unit's self-diagnostic function.

If any of the error codes are displayed, contact our sales representative or agent without delay, report the error code, and request repairs.

Self-diagnostic	Cause	Indication	Buzzer	Remote alarm
Sensor abnormality	The thermal sensor is disconnected.	The alarm indicator blinks. E01 and 50 are displayed alternately.	Intermittent tone.	Alarm mode
	The thermal sensor is short-circuited.	The alarm indicator blinks. E02 and -180 are displayed alternately.		
	The cascade sensor is disconnected.	The alarm indicator blinks.. E03 and chamber temperature are displayed alternately.		
	The cascade sensor is short-circuited.	The alarm indicator blinks. E04 and chamber temperature are displayed alternately.		
	The filter sensor is disconnected.	The alarm indicator blinks. E05 and chamber temperature are displayed alternately.		
	The filter sensor is short-circuited.	The alarm indicator blinks. E06 and chamber temperature are displayed alternately.		
	The ambient temperature sensor is disconnected.	The alarm indicator blinks. E07 and chamber temperature are displayed alternately.		
	The ambient temperature sensor is short-circuited.	The alarm indicator blinks. E08 and chamber temperature are displayed alternately.		
Condenser temperature abnormality	The fan motor for cooling the compressor is faulty. The ambient temperature exceeds the usable environmental condition, etc.	E10 and chamber temperature are displayed alternately.	Intermittent tone.	Alarm mode
Battery check	When about 3 years have passed with the power switch on (time to replace the battery)	The battery indicator lights up.	-----	-----
Fan motor check	When about 6 years have passed with the power switch on (time to replace the fan motor)	The battery indicator blinks.	-----	-----

✧Contact our sales representative or agent to replace the battery for the power-failure alarm when the battery indicator lights up.

✧Contact our sales representative or agent to replace the fan motor when the battery indicator blinks.

# ROUTINE MAINTENANCE

## Cleaning the exterior, interior, and accessories

Use a dry cloth to wipe down the outside and inside of the unit and all accessories. If the outside panels are dirty, clean them with a diluted neutral dish-washing detergent.

Wipe off condensation from the exterior of the cabinet with a dry, soft cloth.

◇Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.

◇After the wiping the cabinet or accessories with a diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth.

### <Important>

■ Do not use a brush, an acid, a thinner, laundry soap, a powder detergent, or boiling water for cleaning. These may cause damage to painted surfaces or cause perishing of plastic and rubber components. Moreover, do not wipe plastic and rubber components with a volatile material.

■ In order to maintain the unit's intended level of performance, always replace accessories that have been removed for cleaning.

## Cleaning the filter

Clean the filter when the filter check indicator lights up and the buzzer sounds. Clean the filter once a month even if the filter check indicator is not lit since a clogged filter may reduce life of the compressor as well as causing reduced cooling performance.

Clean the filter by the procedure below.

1. Mute the buzzer by pressing the buzzer stop key (BUZZER) when the buzzer sounds.


2. Open the grille by pulling it toward you, as shown in the figure.

3. Take out the filter.

4. Wash the filter with water.

5. Replace the filter and the grille. (The handle of the filter should be on the top.)

6. If the filter check indicator was previously on, check that this indicator is now off.

# ROUTINE MAINTENANCE


## Defrosting the chamber

Frost may build up in the upper portion of the chamber and the inner door. Excessive accumulation of frost is likely to create gaps between the cabinet and the magnetic door gasket, which can reduce the cooling performance. Remove the frost from the chamber and the inner door with the scraper enclosed with the unit. Use the following procedure for defrosting when excessive frost builds up in the chamber.

**Note:**

Do not use tool with a sharp edge (such as a knife or screw-driver) to remove the frost.

1. Turn off the switch for the back-up cooling kit (if installed).
2. Take out all contents from the freezer and transfer them to another freezer or a container which is refrigerated by liquid carbon dioxide or dry ice.
3. Turn off the power switch and the battery switch of the freezer.
4. Open the outer door and inner door. Remove the inner door by lifting it up as shown in the figure.
5. Leave the freezer in this state until the frost in the chamber melts.
6. Wipe up the water that accumulates at the bottom of the chamber with a dry cloth.
7. After cleaning the chamber, replace the inner door and start up the unit according to the procedure on page 19.
8. Check that the chamber temperature reaches the set temperature and then replace the contents.
9. Turn on the switch for the back-up cooling kit (if installed).

## Checking the alarm operation

Check the alarm operation regularly (once every 6 months) to ensure that the alarm is activated in case of a power failure. The alarm operation can be checked by the procedure below.

1. Press the alarm test key (ALARM TEST).
2. The alarm indicator blinks, the buzzer sounds intermittently, and the remote alarm is in alarm mode. This continues for 90 seconds and then finishes automatically. During the alarm operation, the temperature display turns off.  
◇During the 90 seconds, the buzzer cannot be stopped by pressing the buzzer stop key (BUZZER).
3. The alarm operation can be cancelled by pressing the alarm test key (ALARM TEST) again during the 90 seconds.

---

# REPLACEMENT OF WORN-OUT PARTS

---

## Replacement of the battery for the power failure alarm

Replace the battery for the power-failure alarm every 3 years (when the battery indicator lights up) to ensure the alarm will operate in the event of a power failure. Contact our sales representative or agent to arrange the replacement of the battery when the battery indicator lights up.

◇The alarm function (blinking of alarm indicator, sounding of buzzer) will not operate when the battery for the power-failure alarm is flat.

◇The alarm indicator blinks and the buzzer sounds as a result of power from the battery for the power-failure alarm. The regular replacement of the battery for the power-failure alarm is important to prevent an increase in chamber temperature in the event of unexpected circumstances.

### <Important>

The used battery is a recyclable resource. Do not dispose of the battery. Always follow the procedure for recycling.

## Replacement of the fan motor

Replace the fan motor for the cooling circuit every 6 years (when the battery indicator blinks) to ensure the appropriate operation of the freezer. Contact our sales representative or agent to arrange the replacement of the fan motor for the cooling circuit when the battery indicator blinks.

◇The appropriate operation of the unit is maintained by cooling the compressor with a fan. The motor to operate the fan (the fan motor for the cooling circuit) has a limited lifespan. A fan motor with diminished performance may lead to reduced cooling performance.

---

# CALIBRATION

---

During continuous operation, the following service tasks must be performed:

- Perform a temperature calibration at least once a year.

For temperature calibration, contact our sales representative or agent.

# TROUBLESHOOTING

If the unit malfunctions, check out the following before calling for service.

**<Attention>**

If the malfunction is not resolved after checking the following items or if the malfunction is not shown in the table below, contact our sales representative or agent.

Malfunction	Check/Remedy
Nothing operates even when the power supply plug is plugged in	<ul style="list-style-type: none"> <li><input type="checkbox"/> The unit is not connected to the power supply properly.</li> <li><input type="checkbox"/> The capacity and voltage of the power supply is not sufficient.</li> <li><input type="checkbox"/> There is a power failure.</li> <li><input type="checkbox"/> The circuit breaker on the supply circuit is activated.</li> <li><input type="checkbox"/> The fuse on the supply circuit is blown.</li> </ul>
The alarm is activated during operation	<ul style="list-style-type: none"> <li><input type="checkbox"/> The unit is not connected to the power supply properly.</li> <li><input type="checkbox"/> The capacity and voltage of the power supply is not sufficient.</li> <li><input type="checkbox"/> There is a power failure.</li> <li><input type="checkbox"/> The circuit breaker on the supply circuit is activated.</li> <li><input type="checkbox"/> The fuse on the supply circuit is blown.</li> <li><input type="checkbox"/> The chamber temperature setting has been changed.</li> <li><input type="checkbox"/> The door has been kept open for a long time.</li> <li><input type="checkbox"/> Containers with a high temperature (load) have been put in the chamber.</li> <li><input type="checkbox"/> The door is open.</li> </ul>
No key operation is available	<ul style="list-style-type: none"> <li><input type="checkbox"/> The key lock is set in ON (L 1).</li> <li>→ Set the key lock in OFF (L 0).</li> </ul>
During setting mode, the mode returns to temperature display mode	<ul style="list-style-type: none"> <li>◇ The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function).</li> </ul>
Excessive noise	<ul style="list-style-type: none"> <li><input type="checkbox"/> The floor is not stable.</li> <li><input type="checkbox"/> The installation site is not level.</li> <li><input type="checkbox"/> The freezer is tilted.</li> <li><input type="checkbox"/> The cabinet is touching the surrounding wall.</li> </ul>
The chamber does not get cold enough	<ul style="list-style-type: none"> <li><input type="checkbox"/> Warm material has been put in the chamber.</li> <li><input type="checkbox"/> The door is frequently opened.</li> <li><input type="checkbox"/> The chamber temperature setting is too high.</li> <li><input type="checkbox"/> The unit is in direct sunlight.</li> <li><input type="checkbox"/> The unit is not installed in suitable location as specified in these operating instructions.</li> <li><input type="checkbox"/> The ventilation around the unit is blocked.</li> <li><input type="checkbox"/> There is a nearby heat source.</li> <li><input type="checkbox"/> The ambient temperature is too high.</li> <li><input type="checkbox"/> There are too many items stored inside the chamber.</li> <li><input type="checkbox"/> The access port is not covered.</li> <li>→ The access port should be covered with insulation and rubber caps when not in use.</li> <li><input type="checkbox"/> The door seal is damaged.</li> <li>→ If it is damaged, contact our sales representative or agent for replacement.</li> <li><input type="checkbox"/> A foreign substance is located between door gaskets.</li> </ul>

# OPTIONAL COMPONENTS

## Temperature recorder

The chamber temperature can be recorded and checked by installing an optional temperature recorder (MTR-G85C or MTR-85H).


✦Contact our sales representative or agent to arrange purchase of the temperature recorder.

Main specifications of temperature recorder


	MTR-G85C	MTR-85H
Recording range	Between -100 °C and +40 °C	Between -100 °C and +50 °C
Feed speed of recording paper	1-day/turn, 7-day/turn, 32-day/turn	60-day/batch
Recording paper	Circular type	Strip type
Power source	Supplied from the freezer	Dry cell

✦For the installation of the temperature recorder MTR-G85C, the recorder mounting bracket MTR-C8 is necessary.

✦For the installation of the temperature recorder MTR-85H, the recorder mounting bracket MTR-C8 and recorder sensor cover MDF-S3085 are necessary.

MTR-G85C

MTR-85H

## Back-up cooling kit

The back-up cooling kit (CVK-UB2) is available as an optional component. Installation of the back-up cooling kit is recommended to protect stored items from damage resulting from any unexpected rise in temperature.

✦Contact our sales representative or agent to arrange purchase of the back-up cooling kit.

## Inventory rack

Optional inventory racks (IR-A216U, IR-220U, IR-312UN) are useful to store important items in the chamber effectively. When these racks are used, it is necessary to change the location of the shelves.

✦Contact our sales representative or agent to arrange purchase of an inventory rack.

# SPECIFICATIONS

Product name	Ultra-Low Temperature Freezer MDF-DU300H
External dimensions	W750 mm x D870 mm x H1830 mm
Internal dimensions	W490 mm x D600 mm x H1140 mm
Effective capacity	333 L
Exterior	Painted steel
Interior	Painted steel
Outer door	Painted steel
Inner door	ABS resin panel with stainless frame, 2 doors
Shelf	Stainless steel, 3 shelves (adjustable) W464 mm x D535 mm, Load: 50 kg/shelf
Access port	17 mm diameter, 3 locations (back x 1, bottom x 2)
Insulation	Rigid polyurethane foamed-in place
Compressor	High stage side; Hermetic type, Output: 450 W Low stage side; Hermetic type, Output: 450 W
Evaporator	Tube on sheet type
Condenser	High stage side; Fin and tube type, Low stage side; Shell and tube type
Refrigerant	High stage side; R-290, Low stage side; R-170
Temperature controller	Microcomputer control system
Temperature display	Digital display
Thermal sensor	Platinum resistance (Pt 1000Ω)
Alarm	High-temperature alarm, Low-temperature alarm, Power-failure alarm
Remote alarm contact	Allowable contact capacity: DC 30 V, 2 A ※1
Battery	Nickel-metal-hydrate battery, DC 6 V, 1100 mAh, Auto-recharge
Accessories	2 keys, 1 scraper
Weight	241 kg
Optional components	Temperature recorder (MTR-G85C) ※2, (MTR-85H) Inventory rack (IR-A216U, IR-220U, IR-312UN) Interface board (MTR-480), LAN interface board (MTR-L03) ※1, ※3 Backup cooling kit (CVK-UB2): LCO <sub>2</sub>

※1: It is recommended to use standard signal and interface cables with a maximum length of 30 meters.

※2: Power source of the temperature recorder must be 220 V.

※3: For the data acquisition system MTR-5000 user only. Contact our sales representative or agent for purchase.

**Note:**

◇Design or specifications are subject to change without notice.

◇Refer to the updated catalogue when ordering an optional component.

# PERFORMANCE

Product name	Ultra-Low Temperature Freezer MDF-DU300H
Model number	MDF-DU300H-PE
Cooling performance	-86 °C at the center of the chamber (ambient temperature: 30 °C, no load) ※1
Temperature control range	-86 °C to -50 °C (ambient temperature: 30 °C, no load)
Power source	AC 230 V/240 V, 50 Hz
Rated power consumption	405 W/435 W
Noise level	52 dB [A] (background noise: 20 dB)
Maximum pressure	2.6 MPa
Usable environment condition	Temperature; 5 °C to 30 °C Humidity; equal or less than 80 %R.H.

※1: Maximum cooling performance.

The chamber temperature can be reached at -86 °C at ambient temperature 30 °C with no load.

Note:

- ◇ The above data is measured based on our internal basis.
- ◇ The unit with CE mark complies with EC directives.

# EMC PERFORMANCE

Emission: EN 61326-1

Immunity: EN 61326-1

This product is intended for use in a basic electromagnetic environment.

## CAUTION

Please fill in this form before servicing.  
Hand over this form to the service engineer to keep for his and your safety.

### Safety check sheet

1. Units contents :

- Risk of infection:  Yes  No  
Risk of toxicity:  Yes  No  
Risk from radioactive sources:  Yes  No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit

Unit interior

- No contamination  Yes  No  
Decontaminated  Yes  No  
Contaminated  Yes  No

Others:

3. Instructions for safe repair/maintenance/disposal of the unit

- a) The unit is safe to work on  Yes  No  
b) There is some danger (see below)  Yes  No

Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

Product name: Ultra-low Temperature Freezer	Model: MDF-	Serial number:	Date of installation:
---	----------------	----------------	-----------------------

Please decontaminate the unit yourself before calling the service engineer.

# DISPOSAL OF UNIT

## Recycle of battery

Ni-MH

- Label indication is obliged to comply with Japanese battery regulation.

- Label indication is obliged to comply with Taiwanese battery regulation.

## Decontamination of unit

Before disposal of unit with biological hazards, decontaminate the unit as much as possible.

## Disposal of Old Equipment and Batteries

### Only for European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

#### Note for the battery symbol (bottom symbol):

This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

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**PHC Europe B.V.**