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DOMINIQUE DUTSCHER SAS

Table of contents	Page
Short-Form Instructions	5
Extent of delivery	7
Warranty	7
schuett count	7
Installation instructions	9
Putting into operation	9
Cleaning	9
Mounting the accessories	9
Controls	11
Operation (counting)	13
Average count calculation	17
USB-interface	17
Optional accessories	19
Troubleshooting	21
Technical Data	23
Optional accessories	25
Spare-parts	25

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

## Short-Form Instructions

**Fig. 1: Back of the device**


- 1) On/Off LED top light lamp
- 2) USB-connector
- 3) 5 V DC-connector (power supply)
- 4) Jack for magnifying glass (optional accessory)
- 5) Jack for top light lamp (optional accessory)
- 6) Contrast area with counting sensors
- 7) Foil keypad

**Fig. 2: Top view on the device**

**Fig. 3: Foil keypad with 4-digit LED-display**

- a) ON/OFF-switch
- b) Key for resetting to display "0"
- c) "Up" and "Down" arrows for manual adding or deleting of colonies and for adjusting the settings of d) e) and f)
- d) Key for setting the light intensity
- e) Key for setting the sensor sensitivity
- f) Key for setting the volume of the counting buzzer
- g) Key for entering a value for average count calculation. At the same time the value of the display is transferred via the USB-interface, see fig. 1 (2).
- h) Key for displaying the average values
- i) Key for deleting the memory for average calculation
- j) 4-digit LED-display (0-9999)

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## **Extent of delivery**

- schuett count colony counter
- Universal power supply 12V DC
- Black contrast disk
- Adapter made of Plexiglas for counting colonies on 90 mm Ø Petri dishes on clear agar (illumination from below)
- Adapter (black side/white side) for counting colonies on 90 mm Ø Petri dishes with non-transparent agar or filter disks (illumination from the side)
- Marker pen and holder
- USB-cable for data transfer
- Instruction manual
- Warranty card with reply card

## **Warranty**

The manufacturer grants a one year guarantee on the schuett count.

To establish your right for warranty coverage, send in the reply card accompanying the warranty card to the manufacturer within 14 days.

Device faults and malfunctions resulting from improper use, e.g. when opening the device, are not covered by the warranty.

Only original accessories and spare-parts supplied by the manufacturer may be used.

## **schuett count**

The schuett count colony counter has been designed for reliable and efficient counting of bacterial colonies and bacteriophage plaques growing on nutrient agar or nutrient disks in Petri dishes.

The circular illumination from the side with an adjustable light intensity allows for working without straylight or eyestrain. For easy counting of white and transparent colonies which hardly differ from the colour of the agar, place the black contrast disk below the Petri dish.

In order to reach optimum illumination for all applications, the schuett count operates on two levels: clear agar is illuminated from below, while non-transparent agar and filter disks are illuminated from top. The Petri dishes are placed with the corresponding adapter into the schuett count. Optionally, there is a gooseneck LED-top light lamp available for connection to the schuett count.

The electronic counting is activated by tapping the marker pen onto the Petri dish, an audible signal is released. The sensitivity of the pressure sensor may be reproduced and adjusted in a range of 1-30.

The volume of the buzzer which is confirming each count may be reproduced and adjusted in a range of 0-25. The buzzer may also be switched off completely.

The schuett count colony counter may be individually completed with our optional accessories.

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## Installation instructions

Take the device out of its shipping package and place it in its intended location.

## Putting into operation

Connect the universal power supply 5V DC of the schuett count to the line/mains supply by using a shockproof socket. Connect the power supply to the 5V DC-connector at the back of the schuett count, fig. 1 (3).

If applicable, connect the USB-connector at the device with the optional USB-connecting cable to your PC / Laptop in order to transfer the counting results, see fig. 1 (2).

## Cleaning

### Cleaning the colony counter

Clean the schuett count only with commercial cleaning solutions and soft cloths.

### Cleaning of magnifying glasses made of Acrylglas

Disinfection with 70% ethanol or isopropanol is possible with short-term exposure (wiping off).

**Important:** Please do not immerse the magnifiers in the cleaning solution - it must not remain on the surface of the magnifier for longer than a few seconds. Dry with a lint-free cloth, if necessary. Sharp chemicals will damage the surfaces.

## Mounting the accessories

### **Magnifying glass (gooseneck):**

Remove the cover cap of the jack at the left-hand side of the device.  
Insert the gooseneck magnifying glass into the jack; see fig. 1 (4).

### **LED-top light lamp (gooseneck):**

Remove the cover cap of the jack at the right-hand side of the device.  
Insert the gooseneck LED-top light lamp into the jack; see fig. 1 (5). The LED-top light lamp is now energized. Bend and turn the LED-top light lamp so that the Petri dish to be counted is optimally lighted. Please hold on to the gooseneck's basis to reduce pressure on the holder. Make sure that no reflections of the LED-top light lamp are visible on the Petri dish and that the gooseneck does not constrain easy reading of the LED-display.

By pushing the switch on the backside of the device, you may switch the lamp on and off; see fig. 1 (1).

### **Attention:**

**Please first start the count colony counter, then the LED-top light lamp.**

**Magnifiers and incident light lamps with goosenecks should ONLY be held by the gooseneck for adjustment to find the correct position for the work, and not by the magnifier, otherwise it may break off.**

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## Operation (counting)

Choose the optimum settings for counting your colonies with the schuett count:

Configuration	Contrast disk	90 mm adapter (Plexiglas)	90 mm adapter (aluminium)	Top light lamp (optional accessory)
Light colonies on clear agar	+	+	-	-
Dark colonies on clear agar	-	+	-	-
Colonies on non-transparent agar	-	-	+	+
Colonies on filter disks	-	-	+	+

- If applicable, place the black contrast disk on the contrast area, see fig. 2 (6).
- Put the three legs of the corresponding counting adapter into the holes in the contrast area, see fig. 2 (6), so that the two short pins are in one row. The feet are asymmetrical, so that they only fit the holes in this position.
- If applicable, mount the gooseneck magnifying glass (optional), see fig. 2 (4) and the gooseneck LED-top light lamp (optional), see fig. 2 (5).
- Please hold on to the gooseneck's basis to reduce pressure on the holder.
- Insert the Petri dish to be counted into the counting adapter. The two short pins centre the Petri dish. Reset the display to "0" by pressing "Reset", see fig 3 (b).
- For starting the schuett count colony counter, push the ON/OFF-switch on the foil keypad, see fig. 3 (a).
- Adjust the light intensity of the LED illumination by using the key "Dimmer", see. fig. 3 (d). The display, see fig. 3 (j), shows "HE:" and the actual level of light intensity. For setting the light intensity push the arrows to increase „↑“ or decrease “↓“, see fig. 3 (c). The light intensity is reproducible and adjustable in a range from 1-30 shown in the display, see fig. 3 (j). "1" is the lowest light intensity and "30" is the strongest light intensity. Press the button "Dimmer" again to confirm the settings.
- Adjust the volume of the buzzer by pressing the key "Buzzer", see fig 3 (f). The display, see fig. 3 (j) shows "to:" and the actual volume. For setting the volume, push the arrows to increase "↑" or decrease "↓", see fig. 3 (c). The volume is reproducible and adjustable in a range from 0-30 shown in the display, see fig. 3 (j). At "0" the tone is switched off, "30" is the display of the highest volume. Press the button "Buzzer" again to confirm the settings.

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- Adjust the sensitivity of the counting plate by pressing the key "Sensitivity", see fig. 3 (e). The display, see fig. 3 (j) shows "SE:" and the actual sensitivity. For setting the sensitivity press the arrows to increase "↑" or decrease "↓", see fig. 3 (c). The sensitivity is reproducible and adjustable in a range from 1-30 shown in the display, see fig. 3 (j). "1" is the lowest sensitivity and "30" is the highest sensitivity. Press the key "Sensitivity" again to confirm the settings.
- For resetting the count display to "0", press the key "Reset", see fig. 3 (b).
- Start counting the colonies by lightly tapping the marker pen on the Petri dish. The sensor records each pressure pulse and the result is shown on the 4-digit LED-display, see, fig. 3 (j). At the same time, the marker pen marks the counted colonies in order to avoid double counts.
- Inadvertently counted colonies add or delete manually by press the arrows for increase "↑" or decrease "↓".

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### Average count calculation

The schuett count may calculate the average count of several counting results:

- Count the first Petri dish as usual.
- Press the key “Enter count“, see. fig. 3 (g) to save the first value.
- Count the second Petri dish and press again the key “Enter count“ and so on. The schuett count may save up to 99 results with a maximum of 32,500 colonies.
- For showing the average value, press the key “Average count“, see fig. 3 (h).
- For deleting the saved values, press the key “Clear count“, see fig. 3 (i).
- The schuett count may calculate the average value only with 3-digit numbers, Petri dishes with more than 1,000 colonies may not be regarded for average count calculation.

### USB-interface

schuett count provides an USB-interface, which may transfer the counting result via software directly to a Microsoft Excel<sup>®</sup>-chart or to a field of the analysis report.

A corresponding USB cable for data transfer is already included in the scope of delivery.

The software can be acquired optionally with license code and matrix description, see accessories.

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## Optional accessories

### **Adapter:**

For counting colonies in Petri dishes with 50 to 60 mm diameter, there is an optional adapter available.

As well for Petri dishes with 140-150 mm diameter.

### **Magnifying glasses\*:**

In general, we recommend the gooseneck magnifying glass (3-fold) for counting.

Alternatively, we offer a gooseneck magnifying glass with 6-8-fold magnification.

### **Counting needle:**

For direct marking of the colonies, there is a counting needle available. The pressure of the needle is registered by the counting sensors.

### **Wolffhügel disk:**

For segmenting the Petri dish surface, place the optional Wolffhügel counting disk onto the contrast area, see fig. 2 (6). For high colony numbers, simply count one or more of the square segments (1 cm<sup>2</sup> each) and calculate the total number of colonies per area of the Petri dish, e.g. Petri dishes with a diameter of 90 mm Ø have an area of 63,62 cm<sup>2</sup>. The Wolffhügel disk has one black and one white background and is suitable for light as well as for dark colonies/agar types.

### **Spiral-Plater disk:**

For counting Petri dishes according to spiral plating (counting of ring segments in 1/8 circular sections), place the optional spiral plater disk on the contrast area, see fig. 2 (6). For maximum contrast, it offers a black and a white background.

\*Please pay attention to our cleaning recommendations for the working magnifiers.

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### Malfunction

### Measures to eliminate the fault

A simple fault may often be the reason for a malfunction, and as such can easily be corrected by the user himself.

schuett count may not be activated by pressing the ON/OFF-switch, see fig. 3 (a).  
The control lamp of the universal power supply 5 VDC does not light up.

- Check whether there is any working voltage from the line/mains supply, e.g. by plugging in another device.
- The power supply is overheated, disconnect from line/mains supply and let it cool down for 10 minutes and try again. Check the cable of the power supply for possible damage.

No counts registered when tapping the marker pen on the Petri dish.

- Make sure that the adapter is fully inserted with all legs into the holes at the contrast area and that neither the Petri dish nor the bottom of the adapter is in contact with the device.
- Increase the sensitivity by pressing key "Sensitivity" see fig. 3 (e) and then pressing the increase arrow "↑" see fig. 3 (c), until figure 30 is shown in the display. Confirm the setting by pressing the key "Sensitivity" see fig. 3 (e) again.

Counts registered due to shocks of the device.

- Decrease the sensitivity by pressing the key "Sensitivity" see fig. 3 (e) and then pressing the decrease arrow "↓" see fig. 3 (c), until the figure 1 is shown in the display. Confirm the setting by pressing the key "Sensitivity" see fig. 3 (e) again.

It is not possible to insert the gooseneck LED-top light lamp into the jack; see fig. 1 (5).

- Carefully use some oil (e.g. WD40 Ballistol) to restore ease of movement.

If following these instructions the equipment does not return to its proper functioning mode, it should be returned to the manufacturer for repair.

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DOMINIQUE DUTSCHER SIS

## Technical Data

Dimensions (w x h x d)	approx. 260 x 130 x 250 mm
Height with magnifying glass/ top light lamp	approx. 300 mm
Weight	approx. 1.5 kg
Safety class	III
IP	30
Universal power supply	100 V - 240 V, 50 / 60 Hz
Power input	7,5 W
Fuse	Overheating protection of the power supply, which resets automatically when cooled down.

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Optional accessories	Cat.-No.
<b>Magnifying glass, 3-fold</b> for 3-fold magnification of the complete Petri dish incl. gooseneck support, individually adjustable, Ø 100 mm <b>Caution with direct sunlight! The magnifying glass may work like a burning glass.</b>	3.081 602
<b>Magnifying glass, 6-8-fold</b> for 6-8-fold magnification of selected parts of the Petri dish incl. gooseneck support, individually adjustable, Ø 60 mm <b>Caution with direct sunlight! The magnifying glass may work like a burning glass.</b>	3.081 612
<b>LED-top light lamp</b> for optimum additional illumination of the colonies incl. gooseneck support	3.081 702
<b>Adapter for Petri dishes</b> with Ø 50-60 mm for center positioning	3.081 802
<b>Adapter for Petri dishes</b> with Ø 140-150 mm for center positioning	3.081 812
<b>Spiral plater disk, black/white</b> for counting with the spiral plater method	3.081 902
<b>Wolffhügel disk, black/white</b> for sectorial counting of the Petri dish	3.081 912
<b>Counting needle</b> for counting by touching the colony, works as replacement for puncture electrodes	3.081 922
<b>Felt-tip marker pen (waterproof)</b>	3.081 932
<b>USB-connection cable for PC</b>	3.081 992
<b>Software for transferring counting in Excel®, per USB, Deutsch/English, license code</b>	3.881 900
Spare-parts	Cat.-No.
<b>Contrast disk</b> black, 93 mm diameter	3.908 140

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

**schuett-biotec GmbH**  
Rudolf-Wissell-Straße 13  
D-37079 Göttingen  
Germany

Fon +49 (0) 551/5 04 10-0  
Fax +49 (0) 551/5 04 10-99

*info@schuett-biotec.de*  
**www.schuett-biotec.de**