**Operating Instructions | Betriebsanleitung | Mode d'emploi |
Instrucciones de funcionamiento | Istruzioni per l'uso**

Sartorius YTM05MA


Temperature Adjustment Set for MA45, MA50, MA100 and MA150

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English – page 3

In cases involving questions of interpretation, the German-language version shall prevail.

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
Overview


Pos.	Description
1	Carrying case
2	Adjustment disk
3	Disk holder
4	Marking
5	Serial number (ref. ID no.)
6	Operating instruction manual
7	Factory calibration certificate

Equipment Supplied

Immediately after you unpack the temperature adjustment set, check to make sure that there is no visible damage to the equipment.

- Adjustment disk
- Disk holder
- Carrying case
- Operating instruction manual
- Factory calibration certificate

Intended Use

The disk for temperature adjustment is used for calibrating and adjusting the temperature of Sartorius moisture analyzers MA45, MA50, MA100 or MA150. Equipment that has been calibrated and adjusted with the temperature adjustment disk can be used as reference devices or as inspection, measurement and testing equipment. The printed calibration record can be used as proof of compliance with the standards for inspection and measuring equipment. Every temperature adjustment disk has a serial number (ref. ID no.) and a factory calibration certificate. The quality of the temperature adjustment disk should be similarly tested. We recommend that you have the temperature adjustment disk tested once a year. For this purpose, you must send the temperature adjustment set back to Sartorius.

The moisture analyzer should be adjusted:

- at room temperature (approx. 20°C)
- under the same ambient conditions
- at regular intervals defined by the user
- using the same complete temperature adjustment set.

The moisture analyzer should be set up on a clean workbench in surroundings that are supplied with clean air. A change in any of the following ambient conditions will affect the results obtained:

- Sunlight
- Drafts caused by open windows or doors
- Air conditioning or fume hoods
- Gases or dust particles.

Leave a space of at least 20 cm around the moisture analyzer and any other objects in the direct vicinity of it, and at least 1 m above the moisture analyzer. After you have started the analyzer, it will carry out calibration and adjustment fully automatically.

Important Note:

Read these operating instructions carefully before calibrating and adjusting the temperature.

Contents

3	Overview
3	Equipment Supplied
4	Intended Use
5	Contents
6	Warning and Safety Information
7	Adjusting the Temperature
7	2-Point Adjustment (MA50 and MA100)
11	2-Point Adjustment (MA45 and MA150)
13	1-Point Adjustment
14	Care and Maintenance
15	Error Codes
16	Factory Calibration Certificate

DOMINIQUE DUTSCHER SAS


Warning and Safety Information

- ⚠ Do not insert or remove the adjustment disk, pan retainer or draft shield immediately after a drying routine. Warning: Danger of severe burns! Do not remove the disposable pan and pan retainer until they have cooled down enough to be touched.
- ⚠ Always use the holder to remove the adjustment disk from the dryer after completion of a measurement or to transport the disk.
- ⚠ Allow the adjustment disk to cool on a clean and fire-proof pad after completion of a measurement.
- ⚠ After cooling, place the adjustment disk back into its packaging and store it in a dust-free and protected place.


Adjusting the Temperature

2-Point Adjustment (MA50 and MA100)

Select the 2-point adjustment for samples that are usually measured at different temperatures. This adjustment option should also be used when components such as the temperature sensor or the sample chamber have to be serviced or replaced.


- Open the sample chamber by pressing the  key.


- Remove the disposable pan and pan retainer and place them on a clean pad.

- > ERR 54 or "L" appears on the display (star-shaped retainer missing).

- Use the disk holder to remove the adjustment disk from its carrying case.
- Insert the adjustment disk under slight pressure. The marking should be facing towards the user.

- Close the sample chamber by pressing the \updownarrow key.

SETUP	isoTEST
Heater adjustment	

- Press the TEST key; the following is displayed: “Heater adjustment”.

SETUP	isoTEST
Heater adjustment	



- Start “Heater adjustment” by pressing the \rightarrow soft key.

SETUP	isoTEST	HEATER
2-point temperature adjustment		
1-point temperature adjustment		


- Select “2-point temperature adjustment”; start by pressing the \rightarrow soft key.

80/160°C	20.0/20.0min	123
		Ref ID

- If required, enter the serial number (ref. ID number) (4). Confirm by pressing the soft key [Ref ID]. The reference ID number is a serial number; for example: 11901617 (on the bottom side of the temperature adjustment disk).

SETUP		isoTEST	HEATER
Nom.,	Actual	Temperatures	
88	82	°C	
168	159	°C	



- Press the [Start] soft key to activate the temperature adjustment routine.
- > The heating symbol appears in the display; the heating time and the previous drying time are displayed continuously.
- After you have finished the measurement for calibrating and adjusting the temperature, you will hear 3 beeps. Display: Analysis ended
- Press [Info] to display the measured values.
- Go back by pressing the ⏪ soft key
- > Return to the main menu without performing an adjustment by pressing the [End] soft key. A record is printed out automatically.
- Activate the calibration/adjustment and temperature adjustment routine on the unit by pressing the [Adjust.] soft key. This concludes the adjustment routine. If the final temperature shows 163°C, for example, the temperature is 1°C outside the tolerance limit of 160°C ± 2°C. This adjustment automatically corrects the deviation. At the same time, a data record is printed out (for example, on the Sartorius YDP20 data printer if you are using the MA45 or MA150).

27.05.2009 10:23
Mod. MA100C000230V
Ser. no. 90805350
Ver. no. 01-38-02
ID

2-PT.ADJ.:

Ref ID 90801234
L.temp. 80°C: 78°C
U.temp. 140°C: 163°C
Adjustment
completed

27.05.2009 10:23
Name:

Important Note:

If the deviation is too high, the analyzer is not adjusted. In this case, please inform your local Sartorius office, dealer or service center.

- Open the sample chamber by pressing the  key
> Allow the moisture analyzer to cool down.


- Snap the disk holder into the adjustment disk and remove the disk. Caution: Equipment is extremely hot after drying!

To avoid severe burns, do not touch until sufficiently cooled.


- ⚠ Allow the adjustment disk to cool on a clean, fire-proof pad.

2-Point Temperature Adjustment (MA45 and MA150)


Select the 2-point temperature adjustment for samples that are usually measured at different temperatures. This adjustment option should be also used when components such as the temperature sensor or the sample chamber have to be serviced or replaced.


- Open the sample chamber.


- Remove the disposable pan and/or pan retainer and place them on a clean pad.

- > ERR 54 or "L" appears on the display (star-shaped retainer missing).

- Use the disk holder to remove the adjustment disk from its carrying case.
- Insert the adjustment disk under slight pressure. The marking should be facing towards the user.

- Close the sample chamber.


- Press the **Enter** key; select “isoTEST” SETUP by pressing the **↗** **↙** keys. Confirm by pressing the **Enter** key. The following appears on the display “Heater adjustment”. Press the **Enter** key again.

- Select “2-point temperature adjustment”; confirm by pressing the **Enter** key.

- Press the **Enter** key again to activate the temperature adjustment routine.

- > The heating symbol appears on the display; the heating time and the previous drying time are displayed continuously.

SETUP	isoTEST	HEATER	S
Non..	Actual	Temperatures	
80	82	°C	
160	159	°C	

2-PT.ADJ. 80/160°C Analysis end. S			
+	153	°C	
END	INFO	ADJUST	

2-PT.ADJ. 80/160°C Adjustment OK S			
+	148	°C	
			END

- After you have finished the measurement for calibrating and adjusting the temperature, an acoustic signal (3 beep tones) is emitted.
Display: Analysis end(ed)
- Press [Info] to display the measured values.
- > Any values outside the tolerance limit are displayed in “Info” (calibration). Press the (CF) key to exit the display.
- Activate the calibration/adjustment and temperature adjustment routine on the unit by pressing the (Enter) key. This concludes the temperature adjustment routine. If the final temperature shows 163°C, for example, the temperature is correct.
- Return to the main menu by pressing the (CF) key; this concludes the adjustment routine. A record is printed out automatically if a printer is connected to the unit.

1-Point Adjustment

You should select the 1-point adjustment for samples that are usually measured at one temperature. In the 1-point adjustment mode, several analyzers can be adjusted to the same drying temperature. (↗) (↖) In the weighing mode, select one of the three programs (P1, P2, P3). In the heating program (standard drying), select and save the desired temperature. The heating time is a factory setting that cannot be changed.

Care and Maintenance

Storage and Shipping Conditions

- During transport, the adjustment set is sufficiently protected by the packaging. Keep all parts of the packaging for storing the adjustment set or in case you have to return it to Sartorius.
- Store at: -20°C ... $+75^{\circ}\text{C}$

Cleaning

- Do not use aggressive cleaning agents (solvents etc.). Do not scratch the adjustment disk. The results can be distorted if the paint finish is damaged.
- Keep the adjustment disk clean. Use a clean cloth that has been wet with water or conventional glass cleaner to wipe off any contamination. Dust and dirt can be cleaned off with a small brush or a dry, soft and lint-free cloth. Do not touch the adjustment disk with oily or greasy fingers.

Error Codes

“Error 53”	The temperature adjustment disk has not been inserted properly.
“L”	Star-shaped retainer missing. The temperature adjustment routine has not been selected yet. Press the TEST key.
“Canceled”	The sample chamber was opened during the measurement.
“Outliers”	The temperature adjustment disk is contaminated. The moisture analyzer is not ready for operation. The moisture analyzer has been running too long (overheated). The moisture analyzer is defective; please inform your local Sartorius office, dealer or service center.
“too many charac.”	You have entered more than 14 characters for the reference ID no.
“no entry”	You have not entered a temperature for the 1-point adjustment.
“too high”	The temperature entered is too high.
“too low”	The temperature entered is too low.

Sartorius Factory Calibration Certificate

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*Test equipment used: Reference temperature adjustment disk, test equipment no.: TEM 0118
 Thermometer, test equipment no.: TEM0121
 Oven for calibration, test equipment no.: TSY 0196*

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*Traceable to: - Physikalisch-Technische Bundesanstalt (PTB) [German Federal Institute of Physics and Metrology]
 - DKD calibration laboratory (DKD= German calibration Service) for temperature measuring equipment*

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The user shall be responsible for meeting reasonable repeat calibration deadlines. One-year regular calibration intervals are recommended; calculated from the buying date for the first recalibration.

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The Temperature Adjustment Set was in perfect condition on the date of initial calibration.

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*Test equipment used: Reference temperature adjustment disk, test equipment no.: TEM 0118
Thermometer, test equipment no.: TEM0121
Oven for calibration, test equipment no.: TSY 0196*

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*Traceable to: - Physikalisch-Technische Bundesanstalt (PTB) [German Federal Institute of Physics and Metrology]
- DKD calibration laboratory (DKD= German calibration Service) for temperature measuring equipment*

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
Sartorius AG
Weender Landstrasse 94–108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-mechatronics.com

Copyright by Sartorius AG,
Goettingen, Germany.
All rights reserved. No part
of this publication may
be reprinted or translated in
any form or by any means
without the prior written
permission of Sartorius AG.
The status of the information,
specifications and illustrations
in this manual is indicated
by the date given below.
Sartorius AG reserves the
right to make changes to
the technology, features,
specifications and design of
the equipment without notice.

Status:
Februar 2010, Sartorius AG,
Goettingen, Germany

Printed in Germany on paper that has
been bleached without any use of chlorine
W4A000 · KT
Publication No.: WYT6005-p10021