

**INSTRUCTION MANUAL**

**BlueVertical™ PRIME™**

**Vertical Electrophoresis Chamber**



SERVA Electrophoresis GmbH □ Carl-Benz-Str. 7 □ D-69115 Heidelberg  
Phone +49-6221-138400, Fax +49-6221-1384010  
e-mail: [info@serva.de](mailto:info@serva.de) □ <http://www.serva.de>

## WARNING

These units are capable of delivering potentially lethal voltage when connected to a power supply and are to be operated only by qualified technically trained personal.

Please read the entire operator's manual thoroughly before operating this unit.

The BlueLine vertical electrophoresis systems are designed to give long service and reproducible results in your laboratory. A few moments spent reading these instructions will ensure that your expectations are reflected in the successful use of the apparatus.

First check that the apparatus has been received complete and undamaged following shipment. Any faults or losses must be notified to responsible **SERVA Electrophoresis GmbH** distributor immediately. **SERVA Electrophoresis GmbH** Heidelberg cannot accept responsibility for goods returned without prior notification.

Refer to the packing list and check that all components and accessories are present.

**Warranty is 12 months from the date of delivery.**

**Please retain all packaging materials until the warranty period has expired.**

## Contents

|   |           |
|---|-----------|
| <b>1. Packing List</b>  | <b>4</b>  |
| <b>2. Specifications</b>  | <b>4</b>  |
| <b>3. Operational conditions</b>  | <b>5</b>  |
| <b>4. Using the vertical gel electrophoresis unit</b>                       | <b>5</b>  |
| 4.1. Safety precautions   | 5         |
| 4.2. General care and maintenance   | 6         |
| 4.3. Gel and buffer volumes/running conditions                              | 6         |
| 4.4. Performing electrophoresis   | 7         |
| 4.4.1. Insertion of pre-cast gel cassettes into the inner core running unit | 7         |
| 4.4.2. Sample loading   | 8         |
| 4.4.3. Start of electrophoresis   | 8         |
| 4.4.4. End of run   | 8         |
| <b>5. Recommended reagents for vertical electrophoresis</b>                 | <b>10</b> |

## 1. PACKING LIST

**BlueVertical™ PRiME™**

**Cat. no.: BV 104**

| No. of Items | Description | Cat. no. |
|--------------|-------------|----------|
| 1            | Main Unit   | BV 104   |
| 1            | Dummy Plate | BV 104-7 |

## 2. Specifications

- Rugged acrylic construction
- All acrylic joints chemically bonded
- Doubly insulated cables, rated safe up to 1000 volt
- Gold plated electrical connectors, corrosion-free and rated safe up to 1000 volt
- Safety lid with integrated power connectors
- Pure platinum electrodes (0.2 mm diameter)
- Size: 16 cm x 15.6 cm x 9.5 cm (w x h x d)
- Weight: 1.2 kg



Table 1: Specifications

|                             |                                      |
|-----------------------------|--------------------------------------|
| Volume inner buffer chamber | 200 ml                               |
| Volume outer buffer chamber | 450 ml                               |
| Voltage (max)               | 500 V                                |
| Current (max)               | 250 mA                               |
| Operating temperature       | 4 °C – 65 °C                         |
| Electrodes                  | Platinum wire (0.2 mm, 99.99 % pure) |
| Dimensions                  | 16 x 15.6 x 9.5 cm (WxHxD)           |
| Weight                      | 1.2 kg                               |

### 3. Operational conditions

This apparatus is intended for indoor use only. Maximum relative humidity up to 80 % (for temperatures up to 31 °C) decreasing linearly to 50 % relative humidity (for temperatures up to 40 °C), at maximum altitude of 2000 m (MSL).

### 4. Using the vertical gel electrophoresis unit

#### 4.1. Safety precautions

- **Read** the instructions carefully before using the apparatus.
- Always disconnect the electrophoresis unit from the power supply before removing the safety cover.
- **Do not** attach the safety lid to the internal gel running unit while it is out of the bottom buffer chamber.
- **Do not** exceed the maximum operating voltage or current (see Table 1).
- **Do not** operate the electrophoresis units in metal trays.
- Polymerised gels may contain some unpolymerised acrylamide monomer, which is a volatile, cumulative neurotoxin and suspected carcinogen. Wear effective protective clothing and gloves. Follow recommended handling and disposal procedures.
- **Do not** fill the unit with running buffer above the maximum fill lines.
- **Do not** move the unit when it is running.
- **Caution:** During electrophoresis very low quantities of various gases are produced at the electrodes. The type of gas produced depends on the composition of the buffer employed. To disperse these gases make sure that the apparatus is run in a well ventilated area.

## 4.2. General care and maintenance

- Before use clean and dry the apparatus with **distilled water only**.  
**Important: Acrylic plastic is not resistant to aromatic or halogenated hydrocarbons, ketones, esters, alcohols (over 30 %) and acids (over 25 %).**
- Before use, and then on a monthly basis, check the unit for any leaks at the bonded joints. Place the unit on a sheet of dry tissue and then fill with **distilled water only** to the maximum fill line. Any leakage will cause dampening of the tissue paper. If any leakage is seen **do not attempt to repair or use the apparatus**, but notify **SERVA Electrophoresis GmbH Heidelberg** resp. the **SERVA Electrophoresis GmbH** distributor immediately.
- The platinum electrodes are partially shrouded for protection. However, when cleaning the main tank **do not** use cleaning brushes in the electrode area. Usually a thorough rinse with distilled water is all that is required.
- Ensure that the connectors are clean and dry before usage or storage.

## 4.3. Buffer volume/running conditions

The buffer volume can be seen in Table 1. The running conditions (Table 2) vary according to the number of gels and their composition. The current required will increase in proportion to the number of gels providing that the voltage is not limiting, e. g. two gels require twice the current of one but the same voltage.

By increasing the gel concentration the electrical resistance is increased and the rate of migration decreases. Higher voltages can be applied but be careful not to overheat the gel.

Table 2: Running conditions using BlueVertical™ in different electrophoresis applications


| Type of electrophoresis     | Type of SERVAGel™ Precast Gel | Running conditions   |
|-----------------------------|-------------------------------|--|
| <b>SDS PAGE</b>             | SERVAGel™ HSE                 | <ul style="list-style-type: none"><li>• Volt<sub>const.</sub>: 400 Volt</li><li>• Time: 20 minutes</li></ul>   |
| <b>SDS PAGE</b>             | SERVAGel™ TG PRIME™           | <ul style="list-style-type: none"><li>• Volt<sub>const.</sub>: 300 Volt</li><li>• Time: 35 minutes</li></ul>   |
| <b>SDS PAGE</b>             | SERVAGel™ TG Neutral          | <ul style="list-style-type: none"><li>• Volt<sub>const.</sub>: 150 Volt</li><li>• Time: 70 minutes</li></ul>   |
| <b>Native PAGE</b>          | SERVAGel™ N                   | <ul style="list-style-type: none"><li>• Volt<sub>const.</sub>: 130 Volt</li><li>• Time: 60 – 120 minutes</li></ul>   |
| <b>Isoelectric Focusing</b> | SERVAGel™ IEF                 | <ul style="list-style-type: none"><li>• Volt<sub>const.</sub>: 100 Volt for 60 minutes</li><li>• Volt<sub>const.</sub>: 200 Volt for 60 minutes</li><li>• Volt<sub>const.</sub>: 500 Volt for 30 minutes</li><li>• Time: 150 minutes</li></ul> |

For detailed information about the running conditions of SERVAGel™ precast vertical mini gels, please see the manuals supplied with the gels.

## 4.4. Performing electrophoresis

### 4.4.1. Insertion of precast gel cassettes into the inner core running unit

- Clean the silicone gasket, located in the inner core unit. If the gasket becomes unseated from its groove simply press it back into place.
- Insert the precast gel cassette into the inner core running unit.



- When running only one gel, a dummy plate is required on the opposite side to retain the required buffer level.
- Place the gel running unit into the buffer tank using the locating pegs as guides.
- Add the appropriate volumes of running buffer to the inner and outer buffer chambers (see Table 1).
- **Important: Do not** overfill (please see marking “max fill”).

#### 4.4.2. Sample loading


- Load the samples using a gel loading pipette tip. During sample loading the pipette tip should be 1 - 2 mm above the bottom of the well to minimise dilution of the sample and to keep the sample as a tight layer.
- Avoid cross-contamination by overloading the sample wells.
- Fill unused wells with the equivalent volume of sample buffer to maintain uniform electrical resistance across the gel.

#### 4.4.3. Start of electrophoresis

- Close the safety lid firmly making sure that the electrical connectors form a good contact.
- Connect the electrophoresis apparatus to the power pack and connect the power pack to the mains supply, e. g. SERVA BluePower 500x4 (BP-500x4). Turn all settings to zero before turning on the mains supply. Adjust the controls to the desired settings. Follow manufacturers instructions.

#### 4.4.4. End of run

- Turn the power supply settings to zero, turn off mains supply and disconnect the power leads.
- Remove the safety lid, take out the internal gel running unit and discard the buffer.





- Discard the buffer from the outer buffer tank. Rinse the chambers twice with distilled water then dry the electrode connectors with tissue. Ensure that the connectors are clean and dry before usage or storage (see chapter 4.2).

**For further information about the SERVA BlueLine, please contact the technical service of SERVA Electrophoresis GmbH in Heidelberg, phone Tel.: +49 (0)6221 13840-44.**

DOMINIQUE DUTSCHER SIS

## 5. Recommended reagents for vertical electrophoresis

SERVA reagents for electrophoresis undergo stringent quality and application control to ensure best performance and results. We recommend the usage of SERVA electrophoresis reagents especially along with BlueLine electrophoresis instruments as the quality of consumables is fine-tuned to the equipment (application tests).

| <b>Precast gels and Starter Kits</b>   |                     |                 |
|--|---------------------|-----------------|
| <b>SERVAGel™ Gels for SDS PAGE</b>   |                     |                 |
| <b>Product</b>   | <b>Sample wells</b> | <b>Cat. No.</b> |
| SERVAGel™ TG PRiME™ Starter Kit  |                     | 43206           |
| SERVAGel™ TG PRiME™ 8 % Tris-Glycine Gel 8 %   | 12                  | 43260           |
|  | 10                  | 43261           |
| SERVAGel™ TG PRiME™ 10 % Tris-Glycine Gel 10 %   | 12                  | 43263           |
|  | 10                  | 43264           |
| SERVAGel™ TG PRiME™ 12 % Tris-Glycine Gel 12 %   | 12                  | 43266           |
|  | 10                  | 43267           |
| SERVAGel™ TG PRiME™ 14 % Tris-Glycine Gel 14 %   | 2D                  | 43268           |
|  | 12                  | 43269           |
| SERVAGel™ TG PRiME™ 4-12 % Tris-Glycine Gel 4-12 %                                       | 10                  | 43270           |
|  | 2D                  | 43271           |
| SERVAGel™ TG PRiME™ 4-20 % Tris-Glycine Gel 4-20 %                                       | 12                  | 43273           |
|  | 10                  | 43274           |
| SERVAGel™ TG PRiME™ 4-20 % Tris-Glycine Gel 4-20 %                                       | 12                  | 43276           |
|  | 10                  | 43277           |
| SERVAGel™ TG PRiME™ 8-16 % Tris-Glycine Gel 8-16 %                                       | 12                  | 43279           |
|  | 10                  | 43280           |
| SERVAGel™ TG 8 % Tris-Glycine Gel 8 %<br>(Available upon request, please inquire.)       | 2D                  | 43281           |
|  | 12                  | 43208           |
| SERVAGel™ TG 10 % Tris-Glycine Gel 10 %<br>(Available upon request, please inquire.)     | 10                  | 43209           |
|  | 12                  | 43210           |
| SERVAGel™ TG 12 % Tris-Glycine Gel 12 %<br>(Available upon request, please inquire.)     | 10                  | 43211           |
|  | 12                  | 43212           |
| SERVAGel™ TG 14 % Tris-Glycine Gel 14 %<br>(Available upon request, please inquire.)     | 10                  | 43213           |
|  | 2D                  | 43226           |
| SERVAGel™ TG 16 % Tris-Glycine Gel 16 %<br>(Available upon request, please inquire.)     | 12                  | 43214           |
|  | 10                  | 43215           |
| SERVAGel™ TG 4-12 % Tris-Glycine Gel 4-12 %<br>(Available upon request, please inquire.) | 2D                  | 43227           |
|  | 12                  | 43216           |
| SERVAGel™ TG 4-20 % Tris-Glycine Gel 4-20 %<br>(Available upon request, please inquire.) | 10                  | 43217           |
|  | 12                  | 43232           |
| SERVAGel™ TG 8-16 % Tris-Glycine Gel 8-16 %<br>(Available upon request, please inquire.) | 10                  | 43238           |
|  | 12                  | 43230           |
| SERVAGel™ Neutral pH 7.4   | 10                  | 43236           |
|  | 12                  | 43231           |
| SERVAGel™ Neutral pH 7.4 Gradient  | 10                  | 43237           |
|  | 2D                  | 43228           |
| SERVAGel™ Neutral HSE Starter Kit  | 12                  | 43220           |
|  | 10                  | 43222           |
| SERVAGel™ Neutral HSE  | 12                  | 43221           |
|  | 10                  | 43223           |
| SERVAGel™ Neutral HSE  | 12                  | 43207           |
|  | 10                  | 43245           |
| SERVAGel™ Neutral HSE  | 10                  | 43246           |
|  | 2D                  | 43247           |

| <b>SERVAGe™ Gels for native PAGE (Blue Native BN/Clear Native CN)</b> |                     |                 |
|---|---------------------|-----------------|
| <b>Product</b>  | <b>Sample wells</b> | <b>Cat. No.</b> |
| SERVAGe™ N Native Starter Kit   | 12/10               | 43204           |
| SERVAGe™ N 3-12, Vertical Native Gel 3-12 %                           | 12                  | 43250           |
|   | 10                  | 43251           |
| SERVAGe™ N 4-16, Vertical Native Gel 4-16 %                           | 12                  | 43252           |
|   | 10                  | 43253           |
| <b>SERVAGe™ Gels for vertical IEF</b>                                 |                     |                 |
| SERVAGe™ IEF Starter Kit  | 12                  | 43205           |
| SERVAGe™ IEF 3-10, Vertical IEF Gel pH 3-10                           | 12                  | 43240           |
|   | 10                  | 43242           |
| <b>Buffers for SDS PAGE</b>   |                     |                 |
| <b>Product</b>  | <b>Cat. No.</b>     |                 |
| SERVA Tris-Glycine/SDS Sample Buffer (2x)                             | 42527               |                 |
| Laemmli Buffer 10x, for SDS PAGE                                      | 42556               |                 |
| SERVA Tris-Glycine/SDS Electrophoresis Buffer (10x)                   | 42529               |                 |
| Laemmli Sample Buffer (2x), for SDS PAGE                              | 42526               |                 |
| SERVA Tris-Tricine/SDS Sample Buffer (2x)                             | 42551               |                 |
| SERVA Tris-Tricine/SDS Electrophoresis Buffer (10x)                   | 42552               |                 |
| SERVA Tris-Tricine/SDS Electrophoresis Buffer (20x)                   | 42560               |                 |
| <b>Buffers for native PAGE</b>  |                     |                 |
| SERVA Tris-Glycine Native Sample Buffer (2x)                          | 42528               |                 |
| SERVA Tris-Glycine Native Electrophoresis Buffer (10x)                | 42530               |                 |
| Native Anode Buffer for BN/CN (10x)                                   | 42535               |                 |
| Native Cathode Buffer for BN/CN (10x)                                 | 42536               |                 |
| Sample Buffer for BN (2x)   | 42533               |                 |
| Sample Buffer for CN (2x)   | 42534               |                 |
| SERVA Blue G solution for BN , 1 %                                    | 42538               |                 |
| <b>Buffers for vertical IEF</b>                                       |                     |                 |
| IEF sample buffer (2x), sterile filtered                              | 42537               |                 |
| SERVAGe™ IEF Running Buffer Kit                                       | 42539               |                 |
| <b>Reagents</b>   |                     |                 |
| Glycine   | 23390               |                 |
| Tris(hydroxymethyl)aminomethane (TRIS)                                | 37190               |                 |
| Dodecylsulfate-Na-salt (SDS)  | 20763 / 20770       |                 |
| 2-Mercaptoethanol   | 28625               |                 |
| Dithiothreitol (DTT)  | 20710               |                 |
| Dithioerythritol (DTE)  | 20697               |                 |
| SERVA Blue G  | 35050               |                 |
| SERVA Blue R  | 35051               |                 |
| Bromophenol Blue-Na-salt  | 15375               |                 |
| SDS Solution, 20 %  | 20767               |                 |

A comprehensive range of SERVA products for electrophoresis is available listed in the Main Catalogue available from **SERVA Electrophoresis**. Please inquire.