

# Microflow $\alpha$

## MICROBIOLOGICAL AIR SAMPLER



**Microflow  $\alpha$  60**

**Microflow  $\alpha$  90-C**

**Microflow  $\alpha$  60/90-C**

<b>INDEX</b>	pag. 2
<b>GENERAL INFORMATION</b>	pag. 4
Working principle	pag. 4
<b>DESCRIPTION OF THE UNIT</b>	pag. 5
Case content	pag. 5
Main features	pag. 5
Technical data	pag. 5
Control and command panel	pag. 6
<b>OPERATIVE INSTRUCTIONS</b>	pag. 7
Operative conditions	pag. 7
Preparing the unit	pag. 7
Placing the plate	pag. 7
Sampling program setting	pag. 8
Function: “Sampling”	pag. 8
Manual sampling	pag. 9
Programmed sampling	pag. 10
Sequential sampling	pag. 13
Sampling end	pag. 16
<b>OTHER FUNCTIONS</b>	pag. 16
Battery level	pag. 17
Setup date and clock	pag. 18
Stored samplings	pag. 19
Download data on PC	pag. 20
Delete memory	pag. 21
Flow rate calibration	pag. 22
Language Setup	pag. 22
Choice sampling head	pag. 23
Reset system	pag. 24
<b>POWER SUPPLY</b>	pag. 24
Low battery	pag. 24
Battery charging	pag. 25
Battery replacement	pag. 26
<b>UNIT DISINFECTION</b>	pag. 26
Sampling head fire disinfecting	pag. 27
<b>MAINTENANCE</b>	pag. 28
<b>WARRANTY</b>	pag. 28
<b>TROUBLE SHOOTING</b>	pag. 29
<b>Microflow <math>\alpha</math> 60 CONVERSION TABLE</b>	pag. 30
<b>Microflow <math>\alpha</math> 90 CONVERSION TABLE</b>	pag. 31
<b>CE CERTIFICATE</b>	pag. 32
<b>ISO CERTIFICATE</b>	pag. 33

## **USER DESTINATION**

---

Airborne bio-pollutants sampler

## **MANUFACTURER and TECHNICAL ASSISTANCE**

---

### **AQUARIA srl**

Via della Levata 14  
20084 Lacchiarella (MI) - Italia  
Tel. (+39) 02 90091399  
Fax (+39) 02 9054861  
E-mail: info@aquariasrl.com  
Web: www.aquariasrl.com

DOMINIQUE DUTSCHER SAS

## GENERAL INFORMATION


This instrument has been developed for airborne biological contaminants sampling in order to access the level of bio-contamination in critical places (e.g.: clean rooms, operating rooms, pharmaceutical production plants, food and cosmetic industries, water purification plants, landfills, etc.) and in normal working life situations and not (air conditioning and ventilation controls).

MICROFLOW  $\alpha$  allows a wide range of airborne bio-contaminant sampling by impacting on agar plates (60 mm  $\pm$  2 mm contact plate or 90 mm  $\pm$  2 mm Petri dish).

The instrument complies with the recommendations indicated in the **Method UNICHIM n° 1962-2/2016** regarding the determination of microbiological contamination of air in the workplaces by active samplers with orthogonal impact and with **ISPESL guidelines** for sampling regarding sequential samplings in operating rooms.

### Working principle

The sampling air is drawn through a perforated head, at constant velocity and for a period of time depending on the condition of the area to be controlled. The inflow air impacts on the surface of the culture media, chosen on the basis of the microbiological investigation to be carried out. At the end of the sampling period the plate is removed and placed into an incubator. At the end of the incubation period it is possible to count the number of colonies developed per cm<sup>2</sup> (CFU / cm<sup>2</sup>) and extrapolate the data to the volume in cubic meters expressed



## UNIT DESCRIPTION

### Case content

MICROFLOW  $\alpha$  kit includes:

- Carrying case
- Microbiological sampler "Microflow  $\alpha$ "
- Battery charger
- IR remot control
- Pencil torch
- Additional 60 mm aluminium sampling head (for **Microflow 60 kit, cod. G.1004**)
- Additional 90 mm aluminium sampling head (for **Microflow 90-C kit, cod. G.1015**)
- 60 or 90 mm aluminium head assembly (for **Microflow 60/90-C kit, cod. G.1027**)
- Base for sampler vertical positioning
- English user manual
- Calibration Report
- Warranty card



### Main features

- Used with: 60  $\pm$  2 mm contact plates (*Microflow  $\alpha$  60*)  
 90  $\pm$  2 mm Petri dish (*Microflow  $\alpha$  90*)  
 Both with contact and Petri dish (*Microflow  $\alpha$  60/90-C*)
- NI-Mh rechargeable battery, without memory effect
- Silent microprocessor controlled fan to ensure maximum reliability on the sampled air volume
- Compensation airflow system (in compensated models)
- Low battery control and alarm
- User friendly. Prompts in 5 languages: English, French, German, Spanish, Italian. All functions are "step-by-step" visualised
- Delay start
- Total volume of sampled air programmable from 1 up to 2000 litres, with 1 litre steps
- Manual sampling for air volumes greater than 2000 litres
- Sequential sampling also as required by **ISPESL guidelines**
- Storing and displayed up to 99 samplings, with record of progressive number, date, hour, sampling volume or time, flow rate, type of sampling (manual, sequential, programmable).
- USB port to download data to a PC (optional)
- Calibration control function (with optional calibration kit)
- Remote switch with infra-red rays control
- Vertical or horizontal sampling (90° and 180°)
- Prepared for mounting on photographic tripod

### Technical data

*Dimensions:* 310 x 130 x 170 mm (W x D x H)

*Weight:* 1,9 Kg

*Chassis:* disinfected polyurethane, joint-free,

*Sampling Head:*

60 mm: autoclavable anodised aluminium, with 219 holes ( $\varnothing$  1 mm)

90 mm: autoclavable anodised aluminium, with 380 holes ( $\varnothing$  1 mm)

*Plate support:* autoclavable anodised aluminium for plate with nominal diameter of 60 mm and 90 mm (adaptable to  $\pm$ 2 mm).

Flow rates: 30 – 60 – 90 - 100 -120 l/min  
 Sampled Air Volume: 1 - 2000 litres (and more with sampling manual)

## Control pannel



Fig. 1 – Control pannel

- **Display**  
High contrast graphic display matrix
- **Battery charger connector**
- infra-red sensor **Remote control**
- **ON/OFF button**  
This is the power switch.  
Put in the ON position to power the unit.  
Place on OFF at the end of the use to save the battery charge.
- **Prog ; Set ; Up ; Down ; Start / Stop ; Clear ; Enter**  
Programming and operating multifunctional keys
- **Data Output**  
USB port to transfer on the PC the stored data (the port is always present but you need to buy the optional software linked to download and manage data in PC)

## OPERATIVE INSTRUCTIONS

### Operation conditions

Microflow sampler needs to be used in the following operating conditions:

- Temperature: 0 – 50°C
- Umidity: 20 – 80% HR

### Preparing the unit


Before to start the sampling, MICROFLOW  $\alpha$  must be disinfected, programmed, loaded with an agar plate and put on a plane surface or on a suitable photographic tripod using the appropriate threaded insert.

**WARNING:** *don't screw too hard the tripod to avoid damage on the threaded insert*

MICROFLOW  $\alpha$ , thanks to spring steel flexible support-plate, allows to operate with plates with nominal diameter of 60 mm (+ 2 mm) or 90 mm (+ 2 mm) depending on the mechanical group chosen. To change the mechanical group simply unscrew the two screws placed on the support plates.

Disinfect properly aluminum plates support prior to sampling.

**WARNING:** *If you use the pencil torch unscrew the support plates and remove it from the body of the Microflow  $\alpha$  in order to avoid damage to the chassie and the fan motor.*




### Placing the plate

Wear sterile gloves. Remove the sampling head rotating clockwise, handling with care to avoid any contamination of its surfaces. Deposit the head on a clean surface and covered with a sterile gauze.

Disinfect the aluminium plate support with a suitable disinfectant, avoiding liquid spillage inside the instrument.

Place a sterile filled plate on one of the three flexible supports. Press gently in outer direction and lay the plate on the two other supports (see figure). Remove the cover plate and put it in a place protected from contamination






Close the sampling sterile unit (or flamed with the accessory provided in the kit) screwing clockwise.

## Sampling program setting

Switch on the unit by pressing [On/Off] button.

On the display will appear sequentially the following **start information**:




**NOTE:** In the second screen instead of Aquaria's address it is likely to occur the customer's address

## Function: "Sampling"

Press [Start] key to start a sampling using the *last stored parameters*.

Press [Prog] key to enter a new sampling program.

Press [Up] and [Down] key to display the following screenshots:




When the needed sampling is shown, press [Enter].



## Manual sampling


The operator can carry out sampling with timings or volumes different from the pre-set ones by manual instrument control

To program a manual sampling, view on display:



and press **[Enter]** to accept.


It will appear :



Where:  
 xxx = the last selected flow rate.

If the flow rate is correct, press **[Enter]**. If not, select the needed flow rate using the **[Up]** and **[Down]** keys. When the flow rate needed is displayed press **[Enter]**.

It will appear:




Where:  
**Manual Sampampling** = manual sampling selected  
 xx x = flow rate selected

If it is necessary to change some data, press the **[Clear]** key and the first display of Manual sampling will appear.

If the programmed data sampling are correct, press **[Enter]**: it will appear the sequence of the three initial displays.

Press **[Start]** key to start the sampling. The two red LED will switch on and, on the display, it will appear:



Where:

- T** = sampling time in **hour:minutes:seconds** (during sampling you will see an increase in the elapsed time)
- L** = sampling **litres** (during sampling you will see an increase in the volume of air sampled)
- Fr** = **selected flow** rate in l/min

Press **[Stop]** key to end the sampling.

The instrument will start to sound a beep, the red LEDs will flash and the following message appears:

```

Sampling end
Fr yyy t h:mm:ss
ENTER ESC
CLEAR Delete
  
```

The beep will continue to sound intermittently until the **[Enter]** or **[Clear]** key is pressed

**N.B.** Pressing **[Enter]** the sampling will be stored.  
 Pressing **[Clear]** the sampling will NOT be stored

Pressed the keys, the LEDs will turn off, will appear "*please wait*" and then the initial screenshots will appear.

## Programmed sampling

To run a programmed sampling, once the display shows:

```

Sampling type
Programmed
↑↓ change
ENTER select
  
```

Press **[Enter]** to accept.

It will appear:

```

Volume
L xxxx
↑↓ change
ENTER select
  
```

where  
**Xxxx** = the last selected volume.

If the volume is correct, press **[Enter]**. If not, select the needed volume using the **[Up]** and **[Down]** keys, to increase or decrease the value, and the **[Set]** key to move the cursor.

When the needed volume is displayed, press **[Enter]**.

**WARNING:** *"Programmed" mode allows to sample an air volume between 1 and 2000 Lt. To sample greater air volumes it is necessary to work in "Manual" mode.*

If you set a volume greater than 2000 liters, for a few seconds the following screen will appear

```

Volume
L 2000
MAX 2000
ENTER select
  
```

It will be memorized a sampling volume of 2000 liters.

It will appear:

```

Flow rate
  yy I/min
  ↑↓ chsnge
  ENTER select
  
```

Where: **yy** = the last flow rate selected.

If the flow rate is correct, press **[Enter]**. If not, select the needed flow rate using the **[Up]** and **[Down]** keys. When the needed flow rate is displayed, press **[Enter]**.

It will appear:

```

Delay Start
  00 : 00
  ↑↓ chsnge
  ENTER select
  
```

Set the needed delay using the **[Up]** and **[Down]** keys, to increase or decrease the value, and the **[Set]** key to move the cursor.

When the delay needed is displayed, press **[Enter]**.

**WARNING:** Set *EXACTLY* how many hour and minutes from the beginning of the sampling.

If the delay start is not needed, leave 00:00

After the **[Enter]** key has been pressed, all the parameters will be displayed (Volume, Flow rate, the eventual delay start and sampling time calculated by volume and flow rate set)

```

Fr yy I/min
V vvvv
Ds 00:00 t mm:ss
  ENTER select
  
```

Where:

**yy** = set flow rate (L/min.)  
**vvvv** = sampling volume (Litres)  
**Ds 00:00** = delay start in hours and minutes  
**t mm:ss** = minutes and seconds of sampling

If the set data are correct, press **[Enter]** to save the program. The sequence of the first displays will appear.


If it is necessary to change one or more data, press the **[Clear]** key and it will be shown the initial display with the programmed sampling settings. Press the **[Enter]** key up to the parameter that it is necessary to change. Modify the parameter and press **[Enter]** to save the setting.

Then return to the screen with list of data set by repeatedly pressing the **[Enter]** and accept all the programming by pressing **[Enter]**.

The three initial displays will appear and the Microflow α is ready for the sampling. Press **[Start]** key to start the sampling.

a) **Delay start not active:**


Press **[Start]** to sample. The two red LED will brighten and, on the display, will be shown:



Where:

- t = sampling time in hours: minutes: seconds (during sampling you will see the countdown of the set time)
- L = sampling **liters** (it will see an increase in the volume of air sampled)
- P = **flow rate se** in liters/min

At the end of the sampling will begin to emit a beep sound, the red LED will begin to flash and the following message appears:



Where:

- Fr is the selected flow rate
- L indicates the air liters sampled


The beep will continue to sound intermittently until the **[Enter]** or **[Clear]** key is pressed

**N.B.** Pressing **[Enter]** the sampling will be stored.  
 Pressing **[Clear]** the sampling will NOT be stored

Pressed the keys, the LEDs will turn off, will appear "please wait" and then the initial screenshots will appear.

b) **Delay start activated:**

After programming e pressing **[Start]**, the two red LED will lighten and, on the display, will be shown:



Where the countdown shows the hours / minutes / seconds missing at the beginning of the sampling.

At the end of the delay start the LED will brighten and, on the display, will be shown:

- Program sampling
- Countdown to the end of sampling
- Increase liters sampled
- Flow rate set

At the end of the sampling period the instrument will operate as already indicated for the programmed sampling without Delay Start

## Sequential sampling

To carry out a sequential sampling, press **[Enter]** after seeing on display:

```

Sampling type
Sequential
↑↓ change
ENTER select
  
```

It will appear:

```

Volume
L xxxx
↑↓ change
ENTER seleziona
  
```

Dove  
xxxx = ultimo volume selezionato.

If the volume is correct, press **[Enter]**. If not, select the needed volume using the **[Up]** and **[Down]** keys, to increase or decrease the value, and the **[Set]** key to move the cursor.

When the needed volume is displayed, press **[Enter]**.

**WARNING:** *“Programmed Seq ” mode allows to sample ONLY an air volume between 1 and 2000 Lt.*

If you set a volume greater than 2000 liters, for a few seconds the following screen will appear

```

Volume
L 2000
MAX 2000
ENTER select
  
```

It will be memorized a sampling volume of 2000 liters.  
 It will appear:

```

Flow rate
yyy l/min
↑↓ chsng
ENTER select
  
```

Where:  
yyy = the last flow rate selected.

If the flow rate is correct, press **[Enter]**. If not, select the needed flow rate using the **[Up]** and **[Down]** keys. When the needed flow rate is displayed, press **[Enter]**.

It will appear:

```

NUM SQ 02
BreaksSQ 001 min
↑↓ change
ENTER select
  
```

where:

**NUM SQ** = number of sequential samplings  
**BreaksSQ** = break in minutes between two samplings

Set the needed data using the **[Up]** and **[Down]** keys to increase or decrease the value, the **[Set]** key to move the cursor and **[Enter]** to move from the first to the second line

**WARNING:** *You can set from 1 to 99 sampling cycles.  
 The time interval between samplings is between 1 and 120 minutes (2 hours). If a larger range is set, the 120 minute interval value will automatically be stored*

When the needed sequence is displayed press **[Enter]**.

It will appear:

```

  Delay Start
    00 : 00
  ↑↓  change
  ENTER select
  
```

Set the needed delay using the **[Up]** and **[Down]** keys, to increase or decrease the value, and the **[Set]** key to move the cursor.

When the delay needed is displayed, press **[Enter]**.

**WARNING:** *Set EXACTLY how many hour and minutes from the beginning of the sampling.*

If the delay start is not needed, leave 00:00

After the **[Enter]** key has been pressed, all the parameters will be displayed (Volume, Flow rate, the eventual delay start and sampling time calculated by volume and flow rate set)

```

  Fr yyy l/min
  V vvvv
  Ds 00:00  t mm:ss
  Nseq nn   Bseq
  
```

Where:

- yyy** = set flow rate (L/min.)
- vvvv** = sampling volume (Litres)
- Ds 00:00** = delay start in hours and minutes
- t mm:ss** = minutes and seconds of sampling
- Nseq nn** = number of sequences set
- Bseq sss** = break in seconds between a sampling and the following

If the set data are correct, press **[Enter]** to save the program. The sequence of the first displays will appear.

If it is necessary to change one or more data, press the **[Clear]** key and it will be shown the initial display with the programmed sampling settings. Press the **[Enter]** key up to the parameter that it is necessary to change. Modify the parameter and press **[Enter]** to save the setting.

Press **[Enter]** up to the summing up screen and save the set parameters pressing the **[Enter]** key

The sequence of the first displays will appear and the unit is ready for the sampling.

Press **[Start]** and the Microflow  $\alpha$  will start the sampling as shown for the programmed sampling, with the following screenshots

**During sampling:**

<b>Sequen. sampling</b>	
t hh:mm:ss	
L xxxx	Sq 01/NN
Battery	www%

Where:

- t = countdown hours/minutes/seconds to end sampling
- L = air sampling volume (Litres increasing)
- Sq = sampling number on the cycles number set (es : 01/02 says that we are doing the first sampling of two)

In this moment the LED are switch on.

**During the interval between a sampling and the next in the same sequence:**

<b>Sequen. sampling</b>	
t hh:mm:ss	
<b>Next sampling</b>	
Sq 02/NN	

Where a countdown will indicate the hours/minutes/seconds before the beginning of the next sample

In this moment the LED are flashing.

**During the second sampling:**

<b>Sequen. sampling</b>	
t hh:mm:ss	
L xxxx	Sq 02/NN
Battery	www%

Where the screen is identical to that of the first sampling but only the indication of the No. corresponding to the sequence varies

In this moment the LED are switch on.

The screens will continue to interleave with the same sequence until the end of the last sampling set. At this point the red LED will begin to flash and the following message appears:

<b>Sampling end</b>	
Fr xxx	Sq 02/02
<b>ENTER esc</b>	
<b>CLEAR delete</b>	

The beep will continue to sound intermittently until the **[Enter]** or **[Clear]** key is pressed

<b>N.B.</b> Pressing <b>[Enter]</b> the sampling will be stored. Pressing <b>[Clear]</b> the sampling will NOT be stored
---

Pressed the keys, the LEDs will turn off, will appear "please wait" and then the initial screenshots will appear.



## Sampling end

Switch off the Microflow  $\alpha$  pressing the **[Off]** key.  
 Remove the sampling head unscrewing anticlockwise.  
 Close the plate with the lid, and after the identification, bring it to incubation.  
 Prepare the unit to the next sampling.


**WARNING:** *In order to save energy, after 10 minutes of sampler inactivity, the LCD backlighting system will automatically turn off (only "Stanby" will appear on the display) to immediately turn on pressing the ENTER key.*

## OTHER FUNCTIONS

Press **[on/off]** to switch on the instrument.

Press **[Set]** to enter in the setting menu of the MICROFLOW  $\alpha$ .

Press the **[Up]** key to scroll the following menus on the Display::



In the **Microflow α 60/90-C (cod. G.1069)**, and ONLY in this instrument, after the display

```

Language setup

↑↓  change
ENTER  select
  
```

It will appear:

```

Setup head

↑↓  change
ENTER  select
  
```

In models in which the optional download data to the PC is enabled, and ONLY in these models, after the display

```

Stored sampling

↑↓  change
ENTER  select
  
```

It will appear:

```

Downloaded data

↑↓  change
ENTER  select
  
```

When the required submenu is displayed, press the **[Enter]** key to accept.

## FUNCTION: "Battery level"

Displayed

```

Battery level

↑↓  change
ENTER  select
  
```

Press **[Enter]** and it will appear:

```

Level: HIGH
  xxx%
Charge current
  000%
  
```

Where:

xxx% battery level

000% means that the instrument isn't under recharge

Pressing the **[Enter]** key the fun will simultaneously switch on so that the battery percentage is the actual one

The battery level is shown by a number (%) and by the words HIGH, WARNING and LOW .

WARNING word will appear when the level is around 40% and LOW word will appear when the level arrives to 5%.

Shortly after the LOW level appears, the the following alarm signal will be displayed:

<b>Level</b>
<b>Xx5%</b>
<b>Charge current</b>
<b>000%</b>

At this point all the keys are disabled, the instrument can no longer work and the only possible thing is to turn off the instrument with the **[On / off]**

Press the **[Clear]** key to exit from the "Battery level" menu.

Press again the **[Clear]** key to visualize the first three displays

**WARNING:** *When WARNING level is reached, the battery charge decreases very rapidly (eg at 40%, you can sample for no more than 20 minutes). Therefore we suggest to recharge the sampler when you reach this level.*

## FUNCTION: "Setup date and clock"

Displayed:

<b>Setup</b>
<b>Date and clock</b>
↑↓ <b>change</b>
<b>ENTER select</b>

Press **[Enter]** and it will appear:

<b>Data</b>	<b>22/07/17</b>
<b>Hour</b>	<b>12:42:03</b>
↑↓ <b>change</b>	
<b>ENTER select</b>	

Set data and time using the **[Up]** and **[Down]** keys, to increase or decrease the value, and the **[Set]** key to move the cursor.

Press the **[Clear]** key to exit from the "Setup date and clock" menu.

Press one time the **[Clear]** key to stay in the Set menu. Press two times the **[Clear]** key to visualize the first three displays.

## FUNCTION: "Stored samplings"

Displayed:

```

Stored sampling
↑↓  change
ENTER  select
  
```

Press the [Enter] key: the sampling progressive sequence will be displayed.

The screen will be different depending on the sampling performed:

### a) Manual sampling

```

N. xx ↑
dd/mm/yy hh:mm
Fr yyy l/min H 90
t www sec ↓
  
```

Where:

**N. xx** = sampling sequential number  
**dd/mm/yy** = sampling date (day/month/year)  
**hh:mm** = sampling start time (hour/minute)  
**Fr = yyy** = flow rate (l/min)  
**H** = sampling head (it can be indicated either 60 or 90)  
**t = www** = sampling duration (seconds)

### b) Programmed sampling

```

N. xx ↑
dd/mm/yy hh:mm
Fr yyy l/min H 90
V zzzz l ↓
  
```

Where:

**N. xx** = sampling sequential number  
**dd/mm/yy** = sampling date (day/month/year)  
**hh:mm** = start sampling time (hour/minute)  
**Fr = yyy** = flow rate (l/min)  
**H** = sampling head (it can be indicated either 60 or 90)  
**V = zzzz** = sampled volume (litres)

### c) Sequential sampling

```

N.xx Seq 01/02 ↑
dd/mm/yy hh:mm
Fr yyy l/min H 90
V www l IS mmm ↓
  
```

```

Nxx+1 Seq 02/02 ↑
dd/mm/yy hh:mm
Fr yyy l/min H 90
V www l IS mmm ↓
  
```

In this case, the sequence sampling will be displayed on so many subsequent screens as are the sampling that has been performed sequentially.

Where:

- N°xx** = number of the first sequential sampling
- N°xx + 1** = identifies the number of the second sequential sampling  
(eg. if the first sequential sampling is identified with n°1, the second sequential sampling will be identified with n°2)
- dd/mm/yy** = sampling date (day/month/year)
- hh:mm** = sampling start time (hour/minute)
- Seq** = sequential sampling
- 01/02** = the first sampling of a sequence of 2 samplings
- 02/02** = the second sampling of a sequence of 2 samplings
- Fr yyy** = flow rate (l/min)
- H** = sampling head (it can be indicated either 60 or 90)
- V wwww** = sampled volume (litres)
- IS** = delay start between sequences of sampling

To switch from one sampling to the next, press the **[Up]** key while pressing the **[Down]** key to switch to the previous sampling


Press the **[Clear]** key to exit from the "Stored samplings" menu.

Press one time the **[Clear]** key to stay in the Set menu. Press two times the **[Clear]** key to visualise the first two displays.


### FUNCTION: "Download data on PC"

(only for Microflow  $\alpha$  with PC interface is enabled)

Displayed:



Press **[Enter]** and it will appear:



Where **Tx** is the number of sampling that will be downloaded to PC:

**NNN** is the total number of stored sampling

**nnn** indicates the sampling number that is downloading to a PC (a progressive increase of nnn until it is equal to NNN will be observed)

Press the **[Clear]** key to exit from the "Stored samplings" menu.

Press one time the **[Clear]** key to stay in the Set menu. Press two times the **[Clear]** key to visualise the first two displays.

**NOTE:** *during this step, the Microflow a has to be connected to a PC with the proper cable.  
We suggest to download only after the connection with the PC.*

**WARNING:** *Every time the download is done, the Microflow a transfer all the data to the PC.  
We suggest to erase the data stored on the Microflow a after the download.*

## FUNCTION: "Delete memory"

Displayed:

```
Delete memory
↑↓  change
ENTER select
```

Press **[Enter]** and it will appear:

```
Delete
samplings
in memory
from 001 to NNN
```

Where:

**NNN** is the total number of samplings stored in the memory

In this case use the **[Set]** key to move the cursor and the **[Up]** and **[Down]** keys to select the required sampling number. After selecting the first sampling number to be deleted, use the **[Enter]** key to switch to the last sampling number to be deleted.

Impostato l'intervallo prescelto premere **[Enter]**. Comparirà:

Set the selected range, press **[Enter]**. and it will appear:

```
Are you sure ?
No
```


If it is not necessary to delete the stored data press the **[Enter]** key to exit from this menu.

If it is necessary to delete the data press the **[Down]** key and it will appear:

```
Are you sure ?
Yes
```

Pressing the **[Enter]** key all the data in the selected range will be erased (if you do not want to erase the data, pressing **[Up]** you will return to the previous screen and you can exit the "Deleted Data Erase" menu).

It will appear the following sequence:




and the instrument automatically returns to the setup menu.

Press the **[Clear]** key to return to the initial screens

## FUNCTION: " Flow rate calibration"

Displayed:



Press **[Enter]** and it will appear:




**WARNING:** *In order to avoid an undesired deletion of the calibration data, the password will be supplied ONLY to the customers who buy the MICROCONTROL calibration kit.*

Press the **[enter]** key to exit from the "Flow rate calibration" menu.

Press the **[Clear]** key to to return to the initial screens

## FUNCTION: " Language Setup"

Displayed:





Press **[Enter]** and it will appear:

```
Language Setup
English
↑↓ change
ENTER select
```

Press the **[Up]** and **[Down]** keys to scroll the available languages (Italian, English, Francais, Deutsch, Espanol).

Press **[Enter]** to confirm. All displayed messages will be now in the selected language

Press the **[Clear]** key to return to the initial screens.

### FUNCTION: " Choice sampling head" (only for Microflow α 60/90-C - cod. G.1027)

Displayed:

```
Setup head
↑↓ change
ENTER select
```

Press **[Enter]** and it will appear:

```
Head
90
↑↓ change
ENTER select
```

Press the **[Up]** and **[Down]** keys to scroll the available heads: 60 and 90.

Selected the desired head, press **[Enter]** and the following display will appear:


```
Sure?
NO
```

Press the **[Up]** and **[Down]** keys to choose yes or no and confirm with **[Enter]** key.

The "Setup head" menu will be displayed again. Press the **[Clear]** key to come back to the main menu.

## FUNCTION: "Reset system"

Displayed:



Press [Enter] and it will appear:



**WARNING:** *This menu can be only used by the manufacturer.*

Press the [enter] key to exit from the "Flow rate calibration" menu.

Press the [Clear] key to return to the initial screens

## POWER SUPPLY

The unit is powered by a 8 elements Ni/Mh battery pack (each element: 1,2V) protect against short-circuits and termic overcharges.

**NOTE:** *These batteries have no memory effect and can be recharged at any time without damaging the elements*

The chip installed on the instrument avoid to run a sampling if the battery level is not enough to maintain a constant flow rate as long as necessary

## Low Battery

a) **Microflow  $\alpha$  switch on but not operative**

If the MICROFLOW $\alpha$  has been left switched on, although not operating, the battery level will slowly decrease. If the battery level decrease over the limit "WARNING", it will appear:



b) **Microflow  $\alpha$  Operative**

If MICROFLOW  $\alpha$  is sampling and during the sampling the battery charge decreases under the WARNING level (and it is not possible to maintain the constant flow rate), the display will show:

<b>Smpling End</b>	
Fr xxx	T hh:mm:ss
ENTER	select
CLEAR	esc

Where: **Fr** is the selected flow rate  
**T** is Past sampling time (in hour, minute, second)

Push the button [ **Enter** ] the sampling will be stored  
 Push the button [ **Clear** ] the sampling will not be stored

**WARNING:** *The sampling will not be stored on the basis of the setted values but on the basis of the real sampling time*

After pressing the button [ **Enter / Clear** ] it will appear :

<p><b>Battery          Level          Very low          WARNING</b></p>
---

and it's not possible proceed with other samplings until the instrument will not be charged.

## Battery charging

**WARNING:** *to avoid any damages to batteries or power adapter we suggest to follow the following sequence:*

- Connect recharger jack to Microflow  $\alpha$
- Connect the charger to the power supply

*After the recharge time, disconnect the charger to the power supply then disconnect the jack from the instrument.*

The complete recharging time is about 8 hours but it depends on the residual battery charge.

**WARNING:** *we recommend to recharge the batteries when the instrument is off. Nevertheless, the instrument can be recharged even if in operation. In this case the charging time could be longer and and the total battery charge is not guaranteed.*

a) **Instrument Off**

The charging of the instrument may be followed by observing the LED on the instrument

**LED flash rapid:** the instrument is in recharge

**LED flash slow :** the instrument is fully charged (battery level 100%) and in maintenance phase

**LED off :** the instrument is fully charged (battery level 100%)

#### b) **Instrument On**

If the instrument is switched on, but not active, the LED will be switched off and the display will show the word "Battery charging"

**WARNING:** *If necessary, the Microflow  $\alpha$  can also be used connected to the power supply through the battery charger.*

*However, we recommend not to use the sampler connected to the network for long periods to avoid damaging the batteries.*

During the recharging, you can enter the set menu and select "level battery charge": on the display it will appear both the charge percentage and the current recharge percentage.

If the charger is connected to the instrument, the battery charge level showed on the display will be 100% (this means that the battery is charging till 100%) while the current recharge will be 100% at the beginning but it will decrease during the battery charge.

If the charger is not connected to the power supply the current recharge will be equal to 0% (any current recharge is issued) and the battery real percentage value will be displayed.

**WARNING:** *To avoid the total batteries discharge, under the threshold value, that would cause a very long recharge battery time (also over 10 hours), WE RECOMMEND to recharge the batteries at least once a month even if the instrument has not been used and in any case, to recharge ALWAYS the instrument after using it.*

## **Battery replacement**

To change the batteries, simply unscrew the metal panel at the bottom of the instrument. When lifting the lid, be careful not to disconnect the connection cable to the IR sensor (if it is disconnected, re-add it simply by making sure that the color wires are coincident) and pay attention to that the batteries do not protrude from the battery compartment by tearing the cable (the batteries are simply housed inside the compartment without any fixing).

Disconnect the batteries from the power cord and replace them with the new battery pack.

## **UNIT DISINFECTION**

MICROFLOW  $\alpha$  must be cleaned and disinfected after each sampling cycle to avoid cross contamination which could modify the reliability of the sampling.

Switch off the unit.

Clean the surface of MICROFLOW  $\alpha$  using a sterile cloth wet with a suitable disinfectant. After drying, the unit is ready for a new sampling cycle.

**WARNING:** **DON'T SPRAY disinfectant solution directly on the instrument or the keyboard: these solutions could enter inside the instrument and they cause damage as example and oxidation of the circuits.**

***AQUARIA will not recognize the warranty caused of unsuitable disinfectant solution.***


The sampling head and the support plate can be sterilized in autoclave, following the standard procedure of sterilisation ( i.e. 120 °C, 18 minutes).

Between two samplings is necessary to disinfect the sampling head with pencil torch supplied with the starting kit. In the absence of this, it is possible to use a sterile cloth with disinfectant over this. It is very important to avoid the liquid penetration inside the instrument.

Let dry at room temperature. After drying the sampling head is ready to be used.

## Sampling head fire disinfecting

The MICROFLOW kit is equipped with a pencil to disinfect the head. Disinfect always the head between two samplings to avoid a cross-contamination.



To switch on the pencil torch rotate clockwise the control-valve. When the gas goes out from the head, bring the fire to light the pencil.

**WARNING: Pass quickly the fire on the interior part of the sampling head.  
Do not stop the fire on a specific point of the head, it can be damaged.**

To switch off the pencil torch rotate counter-clockwise the control-valve.

To recharge the pencil torch it is possible to use the Butan gas used for cigarette-lighters.

Before to recharge the pencil torch be sure it is switched off. The complete recharge needs about 10 seconds. It is necessary to wait some minutes so that the gas is stable.

For correct recharging it is recommended to keep horizontal the pencil torch and the and the recharge cartridge

## MAINTENANCE

---

**All operations shall be performed by qualified personnel only!**

There are no parts replaceable by the user himself. In case of problems please call your dealer for service.

In the event of problems or malfunctions contact your dealer or AQUARIA technical service.

**WARNING: It is recommended to calibrate the instrument at least once a year**

For further information please contact



AQUARIA S.R.L.  
Via delle Levata 14 - 20084 Lacchiarella (MI) - ITALY  
Tel: +39 - 02 90091399  
Fax +39 - 02 9054861  
E-mail: [assistenza@aquariasrl.com](mailto:assistenza@aquariasrl.com)

## WARRANTY

---

AQUARIA guarantees the instrumentation for **12 months**

During the warranty period, any faults resulting from defective components by repair or, in the opinion of its technicians, will be eliminated without replacement by the defective component or the whole appliance.

Damage caused by improper use or improper use of the appliance, parts subject to normal wear, defects that do not affect the normal operation of the appliance and consumables (batteries, etc.) are excluded from the warranty.

The warranty is void if repair work is carried out by unauthorized personnel from AQUARIA and when non-original spare parts are used.

The manufacturer of the appliance can not be held liable, in any way, against the customer, the user and / or third parties for damage caused by persons or things (including the same appliance) in connection with the installation, installation and / or operation of the appliance due to:

- Improper execution, according to the rules of good construction technique
- Works or modifications of any kind made on the device, or its accessories, without the prior approval of the manufacturer and / or without its technical direction, even if such work is carried out by employees dependent on the manufacturer without prior authorization.
- Use of the appliance and / or its accessories for purposes other than that for which it was designed, constructed and sold
- Operation of the appliance and / or its accessories with the use of chemicals, and / or temperatures, and / or concentrations other than those specified in the manufacturer's instructions

## TROUBLE SHOOTING

<i>Nothing lights up</i>	Verify that: <ol style="list-style-type: none"> <li>1. the appliance is switched on</li> <li>2. the battery is charged</li> <li>3. If the problem persists, contact technical support.</li> </ol>
<i>The engine does not work</i>	contact technical support
<i>The remote does not activate the "START"</i>	Verify that: <ol style="list-style-type: none"> <li>1. The remote control battery is charged</li> <li>2. If the problem persists, contact technical support.</li> </ol>
<i>The LCD display does not work properly</i>	contact technical support.
<i>Questions or other abnormal symbols appear on the LCD</i>	contact technical support
<i>Red LEDs do not turn on during sampling</i>	contact technical support
<i>The "LOW BATTERY" or "BATTERY 30%" warning appears with the charged battery</i>	contact technical support
<i>The engine is noisy</i>	To verify: <ol style="list-style-type: none"> <li>1. the presence of foreign bodies in the engine fan.</li> <li>2. If the problem persists, contact Technical Support.</li> </ol>
<i>The engine is out of control</i>	<ol style="list-style-type: none"> <li>1. If anomaly occurs at a single flow rate, correct the relative flow calibration by using the calibration kit.</li> <li>2. If anomaly occurs at all flow rates, check the calibration of all the flow rates.</li> <li>3. If the problem persists, contact technical support.</li> </ol>
<i>Data can not be set via the membrane keyboard</i>	contact technical support



## CONVERSION TABLES Microflow α 60

To correlate the CFUs (Colony Forming Units) present on the agar plate to the most probable number (MPN) of micro-organisms per cubic meter of air sampled, the following formula is used:

$$P_r = F [1/F + 1/(F-1) + 1/(F-2) + \dots + 1/F-(n+1)]$$

Where:

$P_r$  = most probable number of micro-organisms in the volume of air sampled

F = number of holes on Microflow α sampling head

n = number of CFUs on the agar plates after incubation

Table published on UNICHIM method  
n ° 1962-2, edition 2006

n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$
1	1	30	32	59	69	88	112	117	167	146	240	175	349	204	580
2	2	31	33	60	70	89	114	118	169	147	243	176	351	205	595
3	3	32	34	61	71	90	116	119	171	148	246	177	359	206	611
4	4	33	36	62	73	91	117	120	173	149	249	178	365	207	627
5	5	34	37	63	74	92	119	121	175	150	252	179	370	208	646
6	6	35	38	64	76	93	121	122	178	151	255	180	375	209	666
7	7	36	39	65	77	94	122	123	180	152	258	181	381	210	688
8	8	37	40	66	78	95	124	124	182	153	261	182	387	211	712
9	9	38	42	67	80	96	126	125	185	154	265	183	393	212	739
10	10	39	43	68	81	97	128	126	187	155	268	184	399	213	771
11	11	40	44	69	83	98	130	127	189	156	271	185	405	214	807
12	12	41	45	70	84	99	131	128	192	157	275	186	412	215	851
13	13	42	46	71	86	100	133	129	194	158	278	187	416	216	906
14	14	43	48	72	87	101	135	130	196	159	282	188	425	217	979
15	15	44	49	73	88	102	137	131	199	160	286	189	432	218	1088
16	17	45	50	74	90	103	139	132	201	161	289	190	439	219	1307
17	18	46	51	75	92	104	141	133	204	162	293	191	447		
18	19	47	53	76	93	105	142	134	206	163	297	192	455		
19	20	48	54	77	95	106	144	135	209	164	301	193	463		
20	21	49	55	78	96	107	146	136	212	165	305	194	471		
21	22	50	57	79	98	108	148	137	214	166	309	195	480		
22	23	51	58	80	99	109	150	138	217	167	313	196	489		
23	24	52	59	81	101	110	152	139	220	168	317	197	499		
24	25	53	60	82	102	111	154	140	222	169	322	198	508		
25	26	54	62	83	104	112	156	141	225	170	325	199	519		
26	28	55	63	84	106	113	158	142	228	171	331	200	530		
27	29	56	64	85	107	114	160	143	231	172	335	201	542		
28	30	57	66	86	109	115	162	144	234	173	340	202	554		
29	31	58	67	87	110	116	165	145	237	174	344	203	567		

## CONVERSION TABLES Microflow α 90

To correlate the CFUs (Colony Forming Units) present on the agar plate to the most probable number (MPN) of micro-organisms per cubic meter of air sampled, the following formula is used:

$$P_r = F [1/F + 1/(F-1) + 1/(F-2) + \dots + 1/F-(n+1)]$$

Where:

$P_r$  = most probable number of micro-organisms in the volume of air sampled

F = number of holes on Microflow α sampling head

n = number of CFUs on the agar plates after incubation

Table published on UNICHIM  
 method n° 1962-2, edition 2006

n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$	n	$P_r$
1	1	46	49	91	104	136	168	181	245	226	342	271	473	316	674	361	1129
2	2	47	50	92	105	137	170	182	247	227	345	272	477	317	680	362	1149
3	3	48	51	93	106	138	171	183	249	228	347	273	480	318	686	363	1170
4	4	49	52	94	108	139	173	184	251	229	350	274	484	319	693	364	1192
5	5	50	54	95	109	140	174	185	253	230	352	275	487	320	699	365	1216
6	6	51	55	96	110	141	176	186	255	231	355	276	491	321	705	366	1242
7	7	52	56	97	112	142	178	187	257	232	358	277	495	322	712	367	1269
8	8	53	57	98	113	143	179	188	259	233	360	278	498	323	718	368	1298
9	9	54	58	99	115	144	181	189	261	234	363	279	502	324	725	369	1330
10	10	55	59	100	116	145	182	190	263	235	365	280	506	325	732	370	1364
11	11	56	60	101	117	146	184	191	265	236	368	281	510	326	738	371	1402
12	12	57	62	102	119	147	186	192	267	237	371	282	514	327	745	372	1444
13	13	58	63	103	120	148	187	193	269	238	373	283	517	328	753	373	1492
14	14	59	64	104	121	149	189	194	271	239	376	284	521	329	760	374	1546
15	15	60	65	105	123	150	190	195	273	240	379	285	525	330	767	375	1609
16	16	61	66	106	124	151	192	196	275	241	381	286	529	331	775	376	1685
17	17	62	68	107	125	152	194	197	277	242	384	287	533	332	783	377	1780
18	18	63	69	108	127	153	195	198	279	243	387	288	537	333	791	378	1907
19	19	64	70	109	128	154	197	199	281	244	390	289	542	334	799	379	2097
20	21	65	71	110	130	155	199	200	283	245	392	290	546	335	807	380	2477
21	22	66	72	111	131	156	200	201	285	246	395	291	550	336	815		
22	23	67	74	112	132	157	202	202	288	247	398	292	554	337	824		
23	24	68	75	113	134	158	204	203	290	248	401	293	559	338	833		
24	25	69	76	114	135	159	206	204	292	249	404	294	563	339	842		
25	26	70	77	115	137	160	207	205	294	250	407	295	567	340	851		
26	27	71	78	116	138	161	209	206	296	251	410	296	572	341	861		
27	28	72	80	117	140	162	211	207	298	252	413	297	576	342	871		
28	29	73	81	118	141	163	213	208	301	253	415	298	581	343	881		
29	30	74	82	119	143	164	214	209	303	254	418	299	586	344	891		
30	31	75	83	120	144	165	216	210	305	255	421	300	590	345	901		
31	32	76	85	121	145	166	218	211	307	256	425	301	595	346	912		
32	33	77	86	122	147	167	220	212	310	257	428	302	600	347	923		
33	34	78	87	123	148	168	221	213	312	258	431	303	605	348	935		
34	36	79	88	124	150	169	223	214	314	259	434	304	610	349	947		
35	37	80	90	125	151	170	225	215	316	260	437	305	615	350	959		
36	38	81	91	126	153	171	227	216	319	261	440	306	620	351	972		
37	39	82	92	127	154	172	229	217	321	262	443	307	625	352	985		
38	40	83	94	128	156	173	230	218	323	263	447	308	630	353	998		
39	41	84	95	129	157	174	232	219	326	264	450	309	635	354	1012		
40	42	85	96	130	159	175	234	220	328	265	453	310	641	355	1027		
41	43	86	97	131	160	176	236	221	330	266	456	311	646	356	1042		
42	44	87	99	132	162	177	238	222	333	267	460	312	652	357	1058		
43	46	88	100	133	163	178	240	223	335	268	463	313	657	358	1075		
44	47	89	101	134	165	179	242	224	338	269	466	314	663	359	1092		
45	48	90	103	135	167	180	243	225	340	270	470	315	669	360	1110		



## dichiarazione di conformità

*E.C. declaration of conformity*



**AQUARIA srl**  
Via della Levata n°14 - 20084 Lacchiarella (MI) - Italy

represented by its President  
**Emilio Dadati**

### DECLARE that

the MICROFLOW microbiological air sampler, "**MICROFLOW  $\alpha$** "  
serial number N°

conforms to the following E.C. standard:

<i>Immunity:</i>	IEC 1000 - 4- 3 IEC 1000 - 4 - 2	(R.F.E.M. Field) (Electrostatic Discharge)
<i>Emission:</i>	EN 55011 EN 55011	(Conducted Emission) (Radiated Emission - Pre-Compl.)
<i>Others:</i>	CEI EN 50082 - 2 CEI EN 50081 - 1	(Generic Immunity Std) (Generic Emission Std)

Lacchiarella,

AQUARIA srl.

*Emilio Dadati*  
President



THE INTERNATIONAL CERTIFICATION NETWORK

# CERTIFICATE

CISQ/IMQ has issued an IQNet recognized certificate that the organization:

**AQUARIA SRL**  
 VIA DELLA LEVATA 14 - 20084 LACCHIARELLA (MI)

has implemented and maintains a  
**Quality Management System**  
 for the following scope:

**Design, assembly, reselling and repairing of instruments and accessories  
 for sampling and analysis of environmental pollution**

Further clarifications regarding the applicability of ISO 9001:2015 requirements may be obtained by consulting the organization  
 which fulfills the requirements of the following standard:

**ISO 9001:2015**  
 Issued on: **2018 - 09 - 07**  
 Expires on: **2020 - 07 - 14**

This attestation is directly linked to the IQNet Partner's original certificate  
 and shall not be used as a stand-alone document

Registration Number: IT - 12627




Alex Stoichitou  
President of IQNET




Ing. Claudio Provetti  
President of CISQ

IQNet Partners\*:  
 AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy  
 CQC China CQM China CQS Czech Republic Cto Cert Croatia DQS Holding GmbH Germany FCAV Brazil  
 FONDONORMA Venezuela ICONTEC Colombia Inspecia Seritfloat Oy Finland INTECO Costa Rica  
 IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland  
 NYCE-SIGE Mexico PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia  
 SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia  
 IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

\* The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under [www.iqnet-certification.com](http://www.iqnet-certification.com)



www.iqnet.it

CISQ is a member of



THE INTERNATIONAL CERTIFICATION NETWORK  
[www.iqnet-certification.com](http://www.iqnet-certification.com)

IQNet, the association of the world's first class  
 certification bodies, is the largest provider of management  
 System Certification in the world.  
 IQNet is composed of more than 20 bodies and counts  
 over 100 subsidiaries all over the globe.

**CERTIFICATO N. 9190.AQU2**  
**CERTIFICATE N. 9190.AQU2**

SI CERTIFICA CHE IL SISTEMA QUALITA' DI  
 WE HEREBY CERTIFY THAT THE QUALITY SYSTEM OPERATED BY

**AQUARIA SRL**  
 VIA DELLA LEVATA 14 - 20084 LACCHIARELLA (MI)  
 UNITA' OPERATIVE / OPERATIVE UNITS

VIA DELLA LEVATA 14 - 20084 LACCHIARELLA (MI)

E' CONFORME ALLA NORMA / IS IN COMPLIANCE WITH THE STANDARD

**ISO 9001:2015**

PER LE SEGUENTI ATTIVITA' / FOR THE FOLLOWING ACTIVITIES

Progettazione, assemblaggio, commercializzazione e riparazione di strumenti  
 ed accessori per il campionamento e l'analisi degli inquinanti ambientali  
*Design, assembly, reselling and repairing of instruments and accessories  
 for sampling and analysis of environmental pollution*

Ulteriori informazioni riguardanti l'applicabilità dei requisiti ISO 9001:2015 possono essere ottenute consultando l'organizzazione  
 Further clarifications regarding the applicability of ISO 9001:2015 requirements may be obtained by consulting the organization

IL PRESENTE CERTIFICATO E' SOGGETTO AL RISPETTO DEL  
 REGOLAMENTO PER LA CERTIFICAZIONE DEI SISTEMI DI GESTIONE  
 THE USE AND THE VALIDITY OF THE CERTIFICATE SHALL SATISFY THE  
 REQUIREMENTS OF THE RULES FOR CERTIFICATION OF MANAGEMENT SYSTEMS

DATE:	PRIMA CERTIFICAZIONE FIRST CERTIFICATION	EMISSIONE CORRENTE CURRENT ISSUE	SCADENZA EXPIRY
	1999-12-22	2018-09-07	2020-07-14



IMQ S.p.A. - VIA QUINTILIANO, 43 - 20138 MILANO ITALY  
 Management Systems Division - Fluido Orange



SGQ N° 005 A

Numero dell'Attestato di  
 Accreditazione: 005/19/005  
 Data di emissione: 02/04/2019  
 Data di scadenza: 02/04/2021  
 Modalità di accreditazione:  
 ISO 9001:2015 (SISTEMI DI GESTIONE)  
 ISO 14001:2015 (SISTEMI DI GESTIONE)  
 ISO 45001:2018 (SISTEMI DI GESTIONE)

IAF: 19, 29



Organismo di Certificazione Federato CISQ  
[www.imq.it](http://www.imq.it)



CISQ è la Federazione Italiana di Organismi di  
 Certificazione dei Sistemi di gestione aziendale.  
 CISQ is the Italian Federation of management  
 system Certification Bodies.