BagMixer[®] 400

Lab blenders

User's manual

GB

FR

CN

interscience

BagMixer® 400 Range

BagMixer® 400 P

The essential

- Fixed speed (8 strokes/sec.)
- Variable blending time (30 - 210 sec. or open-running)
- Also available in CoolMix® version (refrigerated/heated door)

Ref 021 230

BagMixer® 400 W

Best-seller

- Window door
- Fixed speed (8 strokes/sec.)
- Variable blending time (30 - 210 sec. or open-running)
- Security drip tray

Ref 022 230

BagMixer® 400 VW

Comfort and precision

- Window door
- Multifunction digital screen
- Variable speed (4 - 10 strokes/sec.)
- Variable blending time (1sec-1h or open-running/count-down)
- Adjustable blending pressure with on-screen indication
- Q-Tight® & side-by-side stopping paddle systems
- Silent brushless motor
- Security drip tray

Ref 023 230

BagMixer® 400 CC®

Click & Clean®

- Window door
- Multifunction digital screen
- Variable speed (4 - 10 strokes/sec.)
- Variable blending time (1sec-1h or open-running/count-down)
- Adjustable blending pressure with on-screen indication
- Click & Clean® System (removable paddles)
- Q-Tight® & side-by-side stopping paddle systems
- Silent brushless motor
- Security drip tray

Ref 024 230

BagMixer[®] 400

User's manual

interscience INTERNATIONAL

30, ch. Bois Arpents - 78860 St Nom - FRANCE
T: +33(0) 1 34 62 62 61 - F: +33(0) 1 34 62 43 03
info@interscience.fr - www.interscience.fr

interscience USA & CANADA

199 Weymouth ST. - ROCKLAND - MA 02370 - USA
T: + 1 781 792 2133 - F: +1 781 792 2134
info@intersciencelab.com - www.intersciencelab.com

interscience CHINA

上海市漕溪路222号航天大厦北楼1007室
电话: 021-52541629 - 手机: 13040616422
sales.china@interscience.fr-www.interscience.fr/china

interscience AFTER SALES

T: +33 (0)1 30 56 64 01 - tech@interscience.fr

Thank you for choosing **BagMixer®**.

Your new lab blender:

- ensures an excellent preparation of solid samples before their analysis.
- avoids the use of mixers or high speed blenders which damage the cells and bacteria contained in the sample, and which require fastidious cleaning and sterilization.

Using the **BagMixer® 400** together with interscience filter-bags⁽¹⁾ will guarantee optimal analysis results in ideal conditions of reproductibility and hygiene.

To get the best use from your **BagMixer® 400** and high-quality sample preparations in compliance with ISO 7218 / NF V 08-010 / ISO 6887-1983 specifications, please read this manual carefully.

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CONTENT

I.	SETTING-UP	6
II.	OPERATING	6
III.	CAUTION	8
IV.	STERILE BAGS	8
V.	BAGMIXER® 400 ACCESSORIES	8
VI.	MAINTENANCE AND CLEANING	9
VII.	CHANGING THE SPACE BETWEEN THE DOOR AND PADDLES	10
VIII.	REMOVING THE HOUSING	10
IX.	TROUBLE-SHOOTING	11
X.	TECHNICAL DATA	11
XI.	GUARANTEE	12
	GUARANTEE CARD	29
	CERTIFICATE	31


I. SETTING-UP

1. Take the **BagMixer® 400** out of its packaging and place it on a flat, stable surface.
2. For **BagMixer® 400 W**, **BagMixer® 400 VW** and **BagMixer® 400 CC®** models: slide the drip tray under the blending chamber; this will collect any drops or leaks during blending.

II. OPERATING

BagMixer® P and BagMixer® W models

1. Plug the power cord into the 220-240 V~ mains (100-127 V~ option available). Check that the main switch at the back of the unit is in position "I".
2. Adjust the "Timer" button on top of the unit to the required setting.

3. Open the front door by pulling the lever up. Place a bag containing the sample in the unit and close the door while holding the bag. Push the lever down until you feel a click.
BagMixer® 400 starts automatically.
4. **BagMixer® 400 stops automatically at the end of cycle.**
You may monitor the blending process through the window. (model **BagMixer® 400 W**)
5. Open the door and lift out the bag.
BagMixer® 400 is ready for the next cycle.

BagMixer® 400 VW, BagMixer® 400 CC® models

1. Plug the power cord into the 100-240 V~ mains. Check that the main switch at the back of the unit is in position "I".
2. Set up the time using the "time" keys (+ or -) on top of keyboard. You can also set up open running with the "∞" key.
3. Set up the blending speed of the paddles (1, 2, 3 or 4) using the "speed" keys (+ or -).

4. Set the blending pressure by using the crank at the back of the unit. (chapter VI) Check the settings with the "gap" screen on the keyboard⁽¹⁾.
5. Open the front door by pulling the lever up. Place a bag containing the sample in the unit and close the door while holding the bag. Push the lever down until you feel a click.

BagMixer® 400 starts automatically.

6. **BagMixer® 400 stops automatically at the end of cycle.** You may monitor the blending process through the window.

7. Open the door and lift out the bag.

BagMixer® 400 is ready for the next cycle.

⁽¹⁾ "gap" screen only available on BagMixer® 400 CC® model

III. CAUTION

The unit must be installed upright, on a stable, horizontal surface; never run the unit in any other position! Always use an earthed plug.

The housing should only be removed by qualified and duly authorized staff (see Chapter IX). Never open the housing of the machine without first unplugging the unit from the electric mains.

In normal use, do not run the **BagMixer® 400** for more than 15 minutes continuously (depending on sample size and nature, paddle position and speed selection).

IF YOU ARE NOT SURE, PLEASE CONSULT US FOR ADVICE BEFORE USING!

IV. STERILE BAGS

You may use all good quality standard blender bags with **BagMixer® 400**. However, we strongly recommend **interscience** bags:

BagFilter®: bags with non-woven side filter. Highly robust. Excellent for general use.

BagPage®: bags with a full-width microperforated filtering membrane.

BagLight®: plain blender bags without filter, made in PolySilk®.

V. BAGMIXER® 400 ACCESSORIES⁽²⁾

As well as the drip tray, **interscience** proposes a full range of accessories to improve user working conditions and to increase safety and analysis quality:

BagOpen®: bag holder for use when filling the bag

BagRack®: storage racks for bags

BagClip®: closing clips for blender bags

BagPipet®: fixed-volume pipette (0.1 - 0.9 - 1 ml)

BagTips®: sterile straws for use with the **BagPipet®**

BagSeal®: heat-sealing unit for bags


Ideal complements: **Gravimat®** and **BabyGravimat®** gravimetric dilutors

⁽²⁾ See *BagTools®* and *BagSystem®* on www.interscience.fr

VI. MAINTENANCE AND CLEANING

Maintenance and cleaning are very easy, thanks to its fully opening door system.


1. The housing may be cleaned with standard household products for cleaning stainless steel.
2. To clean the blending chamber, proceed as follows:

- **Unplug the power cord from the mains.** (Fig. 1)

- Lift the lever to open the door, then unhook the lower edge of the door from the axis bar.

- Swing the lower edge of the door upwards and lay the door flat on top of the unit.


- **Click the paddles and remove them.**⁽³⁾ (Fig.2)

- **Clean the blending chamber.** (Fig.3)


- **Click the paddles to reposition them.**⁽³⁾ (Fig.4)

- Reposition the door, **ensuring that it is correctly centered on the lower axis bar.**

- Plug the unit into the power socket.

- Put the switch at the back to position "I".


- Close the door gently.

- Start and test the unit.

⁽³⁾ Only on BagMixer® 400 CC® model


VII. CHANGING THE SPACE BETWEEN THE DOOR AND THE PADDLES (BAGMIXER® 400 VW and BAGMIXER® 400 CC®)

You may reduce or increase the space between the doors and the paddles for special samples from -5 mm to +10 mm. With a negative value, the paddles will knock on the door and produce an audible blending sound. It is normal.

1 . Release the crank at the back of the unit by pulling the handle outwards. (Fig.5)

2 . **Turn the crank to push the paddles backwards or forwards.** (Fig.6)

-  Turn the crank 360° clockwise to push the paddles **forward**.

-  Turn the crank 360° anticlockwise to push the paddles **backwards**.

VIII. REMOVING THE HOUSING

This operation should only be carried out exceptionally, by authorized and qualified staff.

Disconnect the unit from the mains then:

1. Lift the lever to open the door, then unhook the lower edge of the door from the axis bar.
2. Swing the lower edge of the door upwards and lay the door flat on top of the unit.
3. Lift off the green cap from the "timer" knob with a Phillips screwdriver, tap gently but firmly on the brass screw and lift off the black button.
4. Unscrew the screws with a Torx screwdriver.
5. Gently remove the housing, disconnect the earth wire and the flat cable. (Only on **BagMixer® 400 VW and BagMixer® 400 CC® models for the flat cable.**)

To close the housing, follow the instructions in reverse order.

IX. TROUBLESHOOTING

PROBLEMS	CHECK
Unit does not start	<ul style="list-style-type: none"> - Power supply is correctly connected - Main switch (at the back of the unit) is in position "I" (ON) - Door is closed correctly - Fuse is OK - Voltage of power supply is correct
Door closes incorrectly	Check that it is correctly positioned on the axis bar at the bottom of the door
Blending time is wrong	Adjust the timer
Blending speed is wrong	Adjust the speed function
Bags leak during the cycle, or at the end of blending	Use good quality bags: BagFilter® , BagPage® or BagLight® bags from interscience
Bags often leak during the cycle or at the end of the blending, when blending hard (wood, fibers...) or heavy material	If necessary, increase the dead space between door and paddles from and/or reduce blending speed

X. TECHNICAL DATA

BagMixer® 400	P	W	VW	CC®
Useful volume	50 - 400 ml			
Power	100-127/220-240V - 50/60Hz		100/240V - 50/60Hz	
Size (L x W x H)	39 x 26 x 29 cm		42 x 26 x 26.5 cm	
Weight	16.5 kg		16 kg	
Blending time	30-210 sec. & open running	30-210 sec. & open running	1sec.-1h. and open running/countdown	1sec.-1h. and open running/countdown
Blending speed	8 strokes/sec.		4-10 stroke/sec.	
Drip tray	option		✓	
Reference	021 230	022 230	023 230	024 230

**BagMixer® 400 is delivered with:
FREE bag samples, power cord and instruction manual.**

XI. GUARANTEE

BagMixer® 400 is a laboratory unit. Industrial use, over-loaded and daily non-stop use can damage your **BagMixer® 400** and void your guarantee.

See the guarantee card at the end of this handbook.

3 ways to validate your guarantee:

- **online:** www.interscience.fr
- **by post:** **return this card**
- **by fax to:** **+33(0)1 34 62 43 03**

interscience will repair the unit **free-of-charge** according to the guarantee conditions on the guarantee card. The unit must be returned carriage-paid to our laboratory, and **only after** prior agreement with the interscience after-sales department.

Technical support online: +33 (0)1 30 56 64 01
and email: tech@interscience.fr

interscience will extend the guarantee period to **3 years** for users who return the guarantee card to interscience in the 8 days following the purchase.

Our guarantee becomes void in case of incorrect use or opening of the unit by unqualified or unauthorized persons.

GENERAL GUARANTEE CONDITIONS

BagMixer® 400 is guaranteed for a period of 12 months starting on the shipment date indicated on the delivery note. The guarantee is extended to 3 years if the end-user returns the card within 8 days following the purchase. This guarantee covering parts and workmanship will be ensured by calling your dealer, or faxing interscience on +33 (0) 1 34 62 43 03. As malfunctions are extremely rare, please read the chapter 10 "Trouble-shooting" and check your **BagMixer® 400** before any call.

The return of a unit will only be accepted with a repair reference number (issued by the interscience after-sales services). The unit must be returned to us freight pre-paid, in its original packaging or in a similarly protected case after our written agreement.

It is hereby stipulated that this unit must be used for laboratory applications only, and must be maintained according to the instructions in the present user's handbook.

In the case of damage caused by use of the unit for which it was not designed or by unauthorized repairs, the liability of interscience to honor the guarantee will be void.

UNIT SUBJECT TO ARTICLE R 233-84 OF THE WORK CODE

CERTIFICATE OF COMPLIANCE WITH THE RULES OF HYGIENE AND WORK SECURITY
(Application of Article R.233.68 of the Work Code)

interscience certifies that the "BagMixer® 400" complies with the articles R.233-85 to R.233-106 of the Work Code and with the Decree N° 80-544 of 15 July 1980, 1/2.

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3 ways to validate your guarantee:

- **online: www.interscience.fr**
- **by post: return this card**
- **by fax to: +33(0)1 34 62 43 03**

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GUARANTEE

3 WAYS TO VALIDATE YOUR GUARANTEE:

- Online: www.interscience.fr
- By fax: +33 (0)1 34 62 43 03
- By postal mail

Thank you for choosing **BagMixer® 400**.

PLEASE VALIDATE YOUR GUARANTEE WITHIN 8 DAYS AND GET THE FREE GUARANTEE EXTENSION OF 3 YEARS. (see general conditions Chapter XI)

To validate your guarantee by fax or by postal mail, please send us the card below. Keep this part as proof.

Company (user):
First name:
Family name:
Address:
Country:
Phone:
Fax:
Email:
Applications:

BagMixer® 400 Serial number:
Date of purchase:
Distributor's name:
Any comments? Please inform us, we appreciate your advice.

interscience

BagMixer® 400

GUARANTEE

3 WAYS TO VALIDATE YOUR GUARANTEE:

- Online: www.interscience.fr
- By fax: +33 (0)1 34 62 43 03
- By postal mail

Thank you for choosing **BagMixer® 400**.
PLEASE VALIDATE YOUR GUARANTEE WITHIN 8 DAYS AND GET THE FREE GUARANTEE EXTENSION OF 3 YEARS. (see general conditions Chapter XI)

To validate your guarantee by fax or by postal mail, please send us the card below. Keep this part as proof.

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CERTIFICATE OF APPROVAL

This is to certify that the Quality Management System of:

**INTERLAB
Mourjou
France**

has been approved by Lloyd's Register Quality Assurance
to the following Quality Management System Standards:

ISO 9001:2008

The Quality Management System is applicable to:

**Manufacture of laboratory equipment, sterile bags
and filter bags for laboratories and of mixers using these bags.
Processing of laboratory products
and equipment for INTERSCIENCE.**

Approval
Certificate No: FQA 0950488

Original Approval: 6 June 1997

Current Certificate: 1 July 2009

Certificate Expiry: 30 June 2012

A. Juscas

Issued by: Lloyd's Register Quality Assurance France SAS

001

This document is subject to the provision on the reverse

1, boulevard Vivier Merle, 69443 Lyon cedex 03

This approval is carried out in accordance with the LRQA assessment and certification procedures and monitored by LRQA.
The use of the UKAS Accreditation Mark indicates Accreditation in respect of those activities covered by the Accreditation Certificate Number 001

March 2009 Issue 13

CERTIFICATE OF COMPLIANCE

We, interscience, guarantee that the following products:

- MiniMix® 100 VP
- MiniMix® 100 VW
- MiniMix® 100 CC
- CoolMix® 100 VP
- JumboMix® 3500 VP
- JumboMix® 3500 VW
- JumboMix® 3500 VW WarmMix®
- JumboMix® 3500 CoolMix®
- Scan 100®
- Scan 300®
- Scan 500®
- Scan 1200®
- BagMixer® 400 P
- BagMixer® 400 W
- BagMixer® 400 VW
- BagMixer® 400 CC
- CoolMix® 400
- IMix®
- EasySpiral®
- EasySpiral® Pro
- Scan 100®
- Scan 300®
- Scan 500®
- Scan 1200®
- Baby Gravimat®
- Gravimat®

➤ Are compliance with the following directives:

- **2002/95/EC** (the RoHS Directive) restricting the use of hazardous substances in electrical and electronic equipment (based on information provided by our suppliers)
- **2002/96/EC** (the WEEE Directive) imposing selective recovery of Electrical and Electronic Components (return cartilage costs to interscience to be borne by the sender)
- **94/62/CE** relating to packaging and packaging waste
- **SJ/T11363-2006** restricting the use of hazardous substances in electrical and electronic equipment (RoHS China Concentration Standard)
- **SJ/T11364-2006** imposing the marking of the products concerned with the standard above and their packaging (RoHS China Labeling Standard)
- **GB18455-2001** describing the imposed marking of packaging of the RoHS products.

St Nom-ia-Breteche,
March 24th, 2011

Emmanuel JALENQUES
Director

(« CE ») CERTIFICATE OF COMPLIANCE

We, interscience, guarantee that the following products:

- MiniMix® 100 VP
- MiniMix® 100 VW
- MiniMix® 100 CC
- CoolMix® 100 VP
- JumboMix® 3500 VP
- JumboMix® 3500 VW
- JumboMix® 3500 VW WarmMix®
- JumboMix® 3500 CoolMix®
- Scan 100®
- Scan 300®
- Scan 500®
- Scan 1200®
- BagMixer® 400 P
- BagMixer® 400 W
- BagMixer® 400 VW
- BagMixer® 400 CC
- CoolMix® 400
- IMix®
- EasySpiral®
- EasySpiral® Pro
- Scan 100®
- Scan 300®
- Scan 500®
- Scan 1200®
- Baby Gravimat®
- Gravimat®

➤ Are made in CE compliance with the following directives:

- European Directive 2006/42/EC relating to machines.
- European Directive 2004/108/EC about electromagnetic compatibility.
- European Directive 2006/95/EC relating to low voltage electrical goods

St Nom-ia-Breteche,
March 22th, 2011

Emmanuel JALENQUES
Director

NOTES

A series of horizontal dotted lines for writing notes.

Get the best results with...

Gravimat®

Gravimetric dilutor

Automatically dilutes your sample with the appropriate mass of solvent.

BagFilter® 400

Blender bags with side filter

Sterile blender bag with a side non-woven filter.

BagPage® 400

Full page sterile filter bags


Robust & food compatible, multilayer.

BagLight® 400

Plain sterile bags

Made of PolySilk®, extremely robust. For samples that do not require any filtration.

RollBag®

Sampling bags

For collecting, sampling & preservation of any sample. Stands up alone.

BagOpen® 400

Bag opener

Easy opening, extremely stable.

BagRack® 400

Storage rack

All stainless steel.

BagClip®

Closing clips

Water- and air-tight. Close any kind of bag.