## Minisart® High Flow | NML | Air | HY | Acticosart

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# Minisart® High Flow | NML | Air | HY | Acticosart

Sartorius Syringe Filters for Liquid Filtration and Venting

Sartorius Minisart® filters with MBS housing and a choice of hydrophilic and hydrophobic membrane filters are suitable for clarification (0.8–1.2 µm), particle removal (0.45–0.65 µm) and sterilization of liquids (0.1–0.2 µm), air and gases. The housing material MBS could be glued with tubing by using suitable glue or residue-free solvents.

Minisart® High Flow with hydrophilic PES membrane offers highest flow rates and is optimized for filtration of additives and cell culture buffers. Article 16553-K with 0.1 µm is suitable for Mycoplasma removal.

Minisart® NML with surfactant-free cellulose acetate is virtually free of extractables and is optimized for clarification of aqueous solutions and drugs.


Minisart® Air and HY with hydrophobic PTFE membrane are designed for the sterile filtration and cleaning of air and gases and are optimized for the venting of bags and containers, e.g. fermenters.

Minisart® Acticosart is designed for the ultracleaning of air and gases and contains active carbon for binding of aerosols and volatile substances.

## Applications

- Minisart® High Flow for clarification and sterilization of liquids with highest flow rates and total throughput
- Minisart® NML for surfactant-free clarification and particle removal of aqueous liquids
- Minisart® Air | HY for sterilization and cleaning of air and gases and venting of containers
- Minisart® Acticosart for ultracleaning of air and gases

## Instructions for Use

If possible, draw a slight amount of air (at least 1 ml) into a syringe (a). Then fill the syringe with the liquid to be filtered: pull the plunger upwards to draw this liquid from a suitable container into the syringe (b). If required, remove any remaining liquid from the tip of the syringe and attach a Minisart® to the luer lock or luer slip connector on the filled syringe. Sterile blister units are opened by peeling off the protective backing. Hold the opened blister part on the outside and connect the inside sterile syringe filter to the syringe (c). Apply consistent pressure to press in the plunger of the syringe in order to filter the liquid through the Minisart® filter into a suitable vial for collection (d). Afterwards, press the plunger all the way in so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter. As a result, this will reduce the hold-up volume.

## Caution!

1. If sterile Minisart® blister units are used please check that the blister packing is not damaged!
2. Take care when using syringes with a volume of less than 10 ml, as they can easily generate a pressure greater than the maximum recommended pressure resistance of 6 bar (87 psi) for Minisart®. Therefore, slowly press in the plunger on a syringe with a volume of less than 10 ml. As soon as you detect considerable resistance (= max. filter pressure capacity), do not use force while continuing to press the plunger exceeding 6 bar. Otherwise, you can damage the filter, and filtration will no longer guarantee an ultrapure filtrate!
3. Minisart® is designed for bidirectional use. However, once you have selected one direction of filtration, be sure to maintain this direction. Never use the same syringe filter for both directions!

## Note

1. If you need to spike a closure first before injecting the filtrate, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile needle (for a male luer slip) to the outlet of the Minisart® and remove the protective cap before injection.
2. If you would like to achieve maximum recovery of your sample to be filtered, detach the syringe from the syringe filter and draw air once more into the syringe. After reconnecting the syringe and filter, use the air cushion created in this way to press out the remaining liquid. To do so, you will have to apply pressure beyond the bubble point of the particular Minisart® membrane incorporated (see Table 1). This will reduce the hold-up volume to the lowest possible amount.
3. Our filters are leading in terms of purity, so prerinsing is usually not necessary. For very sensitive methods or samples, you may need to additionally prerinse the particular syringe filter you use with 1–3 ml ultrapure water or buffer. This way, you can be sure that you have removed even the slightest traces of extractables.
4. Use a Minisart® for only one sample to reliably prevent carrying over residues of one sample to the next.
5. You will find an overview of chemical compatibilities in the 'Minisart® Chemical Compatibility Guide' on our website.

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Table 1

Specifications for Minisart® High Flow   NML with 28 mm filtration area Ø, Air with 15 mm Ø, HY   Acticosart with 26 mm Ø										
Housing material	MBS (Methacrylate butadiene styrene)									
Membranes	High Flow: – PES = Polyethersulfone – NML: (SF)CA = (Surfactant-free) Cellulose Acetate – Air HY: PTFE = Polytetrafluoroethylene – Acticosart: PTFE and a dome reservoir filled with active carbon									
Application limits High Flow	Max. recommended operating pressure 6.0 bar   87 psi									
Application limits NML HY Air	Max. recommended operating pressure 4.5 bar   65 psi									
Application limits Acticosart	Max. recommended operating pressure 1 bar   14.5 psi									
Housing burst pressure	>7 bar   102 psi (not determined for Acticosart)									
Max. temperature	60°C not autoclavable!									
Sterilization	Non-sterile Minisart® High Flow   NML can be sterilized by ethylene oxide (EO) or Gamma sterilization Non-sterile Minisart® Air   HY   Acticosart can be sterilized by ethylene oxide (EO) sterilization*									
Minisart <sup>3</sup> type	PES 0.1 µm	PES 0.2 µm	PES 0.45 µm	CA 0.65 µm	CA 0.8 µm	CA 1.2 µm	CA 5.0 µm	PTFE 0.2 µm	Actico-sart	PTFE 1.0 µm
Bubble point (≥)	with water 5.0 bar   73 psi	with water 3.2 bar   46 psi	with water 2.0 bar   29 psi	with water 1.3 bar   19 psi	with water 0.8 bar   12 psi	with water 0.7 bar   10 psi	with water 0.4 bar   6 psi	with ethanol 1.1 bar   16 psi	with ethanol 0.9 bar   13 psi	with ethanol 0.5 bar   7 psi
<b>Flow rate:</b> High Flow   NML, 28 mm Ø = 6.2 cm <sup>2</sup> filter area   Air, 15 mm Ø = 1.7 cm <sup>2</sup> filter area   HY   Acticosart, 26 mm Ø = 5.3 cm <sup>2</sup> filter area										
<b>Hold-up volume<sup>1</sup>:</b> High Flow   NML: 100-150 µl   Air   HY   Acticosart: n.a.										
with water at 1 bar:	40 ml/ min	140 ml/ min	220 ml/ min	250 ml/ min	400 ml/ min	500 ml/ min	600 ml/ min	– <sup>3</sup>	– <sup>3</sup>	– <sup>3</sup>
with air at 0.1 bar:	– <sup>2</sup>	– <sup>2</sup>	– <sup>2</sup>	– <sup>2</sup>	– <sup>2</sup>	– <sup>2</sup>	– <sup>2</sup>	2.0 l/min	2.3 l/min	4 l/min
Water penetration point	–	–	–	–	–	–	–	>3.0 bar   44 psi <sup>3</sup>	n.a.	>1.5 bar   22 psi
Sterile filtration capability <sup>4</sup> acc. to BCT	yes	yes	no	no	no	no	no	yes	n.a.	no
Main applications	hydrophilic filters with PES membrane for fast sterilization or clarification of liquids			hydrophilic filters for particle removal or clarification of aqueous solutions			hydrophobic filters for cleaning or ultracleaning (Actico-sart) of air and gases and venting of containers			

<sup>1</sup> Hold-up volume after air purge. For minimum hold-up please see "Note | Hinweis | Remarques | Nota | Note 1"

<sup>2</sup> Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

<sup>3</sup> Hydrophobic membranes cannot be wetted with aqueous solutions unless you overcome their water penetration point

<sup>4</sup> According to bacterial challenge test (BCT) with 10<sup>7</sup> *Brevundimonas diminuta*.

All non-sterile Minisart<sup>3</sup> types can be sterilized according to the sterilization processes listed above.  
n.a. = not applicable

\* Minisart<sup>6</sup> Air can be sterilized by Gamma sterilization according to the following parameters:  
Range 25–40 kGy (validated with 50 kGy).

**Table 2**

**Sartorius Syringe Filters for Liquid Filtration and Venting**

Minisart® High Flow | NML | Air | HY | Acticosart

Ø mm	Pore Size	Outlet	Qty.   Pk	Order Number	Sterility Status
<b>Minisart® High Flow (PES)</b>					
28 mm	0.1 µm	Male Luer Lock	50	16553-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Lock	50	16532-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Lock	50	16532-----GUK	single-blister-packed, Gamma sterile
28 mm	0.2 µm	Male Luer Lock	500	16532-----Q	non-sterile
28 mm	0.2 µm	Male Luer Slip	50	16541-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Slip	500	16541-----Q	non-sterile
28 mm	0.45 µm	Male Luer Slip	50	16533-----K	single-blister-packed, EO sterile
28 mm	0.45 µm	Male Luer Slip	50	16533-----GUK	single-blister-packed, Gamma sterile
28 mm	0.45 µm	Male Luer Slip	500	16533-----Q	non-sterile
28 mm	0.45 µm	Male Luer Lock	50	16537-----K	single-blister-packed, EO sterile
28 mm	0.45 µm	Male Luer Lock	500	16537-----Q	non-sterile
<b>Minisart® NML (SFCA)</b>					
28 mm	0.65 µm	Male Luer Slip	50	16569-----K	single-blister-packed, EO sterile
28 mm	0.8 µm	Male Luer Lock	50	16592-----K	single-blister-packed, EO sterile
28 mm	0.8 µm	Male Luer Lock	50	16592-----GUK	single-blister-packed, Gamma sterile
28 mm	0.8 µm	Male Luer Lock	500	16592-----Q	non-sterile
28 mm	1.2 µm	Male Luer Lock	50	17593-----K	single-blister-packed, EO sterile
28 mm	1.2 µm	Male Luer Lock	500	17593-----Q	non-sterile
28 mm	5.0 µm	Male Luer Lock	500	17594-----Q	non-sterile
<b>Minisart® Air (hydrophobic PTFE)</b>					
15 mm	0.2 µm	Male Luer Slip	500	1751A-----Q	non-sterile
15 mm	0.2 µm	Needle	50	16596-----HMK	single-blister-packed, Gamma sterile
<b>Minisart® HY (hydrophobic PTFE)</b>					
26 mm	1.0 µm	Male Luer Lock	500	1659A-----HYQ	non-sterile
26 mm	1.0 µm	Male Luer Slip	500	1659B-----HYQ	non-sterile
<b>Minisart® Acticosart with Dome Reservoir + hydrophobic PTFE</b>					
26 mm	0.45 µm	Male Luer Slip	500	17840-----Q	non-sterile, with active carbon
26 mm	0.45 µm	Male Luer Slip	500	40080-----Q	Non-sterile, empty reservoir

- Would you like to filter solvents, acids or bases?
- Do you need to filter small volumes?
- Would you like to use PP housings and other membranes?  
Please refer to Minisart® RC, Minisart® NY, or Minisart® SRP for highest chemical compatibility;  
also available in 4 mm or 15 mm filtration area
- Do you need Minisart® with pre-filters for filtration of highly particle laden samples?  
Please refer to Minisart® NML Plus and Minisart® NY Plus which incorporate a glass fiber pre-filter
- Are you looking for CE-certified Minisart®?  
Please request CE-certified article numbers for Minisart® NML, Minisart® HY and Minisart® SRP

## Quality Assurance Certificate

Minisart is manufactured by Sartorius in accordance with the Applicable Good Manufacturing Practice Standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Sterile, individually packed products were sterilized with ethylene oxide gas using a validated process following DIN/EN ISO 11135 regulations, or were sterilized with gamma irradiation using a validated process following DIN/EN ISO 11137 regulations.

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no. alone or lot no. with one of the following: pore size, ID no., material type, order no. are printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biosafe according to the USP Class VI Biological tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.72 for non-fiber-releasing filters.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

Non-Sterile and Sterile Products  
 Burst Pressure Test  
 Bubble Point Test  
 Endotoxine Test  
 Pressure Hold Test  
 Flow Rate Performance

Additional Tests for 0.2 µm Pore Size  
 Sterile Filtration Capability (Bacterial Challenge Test)

Additional Tests for Sterile Products  
 Sterility

09.03.2015

Date

*Hartmut Hennig*

Dr. Hartmut Hennig,  
 Site Manager QA Filtration Products

Manufactured by  
 Sartorius Stedim Biotech GmbH  
 August-Spindler-Strasse 11  
 37079 Goettingen, Germany  
 www.sartorius.com

## Symbols

- |  |   |  |
|--|---|--|
| <p> Nicht zur Wiederverwendung<br/>                 Single Use<br/>                 Ne pas réutiliser<br/>                 No reutilizable<br/>                 Monouso</p>  | <p><b>LOT</b> Chargenbezeichnung<br/>                 Lot Number   Code du lot<br/>                 Numéro de lote<br/>                 Numero di lotto</p>   | <p> pyrogenfrei<br/>                 non-pyrogenic<br/>                 non pyrogène<br/>                 Apirógeno<br/>                 Apirogeno</p>   |
| <p><b>REF</b> Bestellnummer<br/>                 Order Number<br/>                 Référence du catalogue<br/>                 Número de referencia<br/>                 Codice d'ordine</p>   | <p><b>STERILE ETO</b> Sterilisation mit Ethylenoxid<br/>                 sterilized, ETO<br/>                 Méthode de stérilisation<br/>                 utilisant de l'oxyde d'éthylène<br/>                 Esterilización por ETO<br/>                 Sterilizzazione con ETO</p>  | <p> Porengröße<br/>                 Pore size<br/>                 Taille de pore<br/>                 Tamaño de poro<br/>                 Porosità</p>  |
| <p> Achtung!   Warning!<br/>                 Attention!   Atención!<br/>                 Attenzione!</p>   | <p><b>STERILE R</b> Sterilisation durch Bestrahlung<br/>                 sterilized, irradiated<br/>                 Méthode de stérilisation<br/>                 utilisant l'irradiation<br/>                 Esterilización por radiación<br/>                 Sterilizzazione per irradiazione</p>                                  | <p> Temperaturbegrenzung<br/>                 maximum and<br/>                 minimum temperature<br/>                 Limite de température<br/>                 Temperatura máxima<br/>                 y mínima<br/>                 Temperatura massima<br/>                 e minima</p> |
| <p> Hersteller   Manufacturer<br/>                 Fabricant   Fabricante<br/>                 Produttore</p>  | <p><b>STERILE I</b> Sterilisation mit Dampf<br/>                 sterilized, hot steam<br/>                 Méthode de stérilisation<br/>                 utilisant la vapeur ou<br/>                 la chaleur sèche<br/>                 Esterilización por vapor caliente<br/>                 Sterilizzazione con vapore caldo</p> |  |
| <p><b>IVD</b> In-Vitro-Diagnostika<br/>                 In-Vitro-Diagnostics<br/>                 Dispositif médical<br/>                 de diagnostic in vitro<br/>                 Para uso en diagnóstico<br/>                 in vitro<br/>                 Per uso diagnostico<br/>                 in vitro</p> | <p> Gebrauchsanweisung beachten<br/>                 Consult instructions for use<br/>                 Consulter le mode d'emploi<br/>                 Consulte las instrucciones de uso<br/>                 Consultare le istruzioni per l'uso</p>  |  |
| <p> verwendbar bis<br/>                 Use before   Utiliser jusqu'au<br/>                 Fecha caducidad<br/>                 Data di scadenza</p>  |   |  |