Minisart® High Flow | NML | Air | HY | Acticosart

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Sartorius Minisart® filters with MBS housing and a choice of hydrophilic and hydrophobic membrane filters are suitable for clarification (0.8–1.2 µm), particle removal (0.45–0.65 µm) and sterilization of liquids (0.1–0.2 µm), air and gases. The housing material MBS could be glued with tubing by using suitable glue or residue-free solvents.

Minisart® High Flow with hydrophilic PES membrane offers highest flow rates and is optimized for filtration of additives and cell culture buffers. Article 16553-K with 0.1 µm is suitable for Mycoplasma removal.

Minisart® NML with surfactant-free cellulose acetate is virtually free of extractables and is optimized for clarification of aqueous solutions and drugs.


Minisart® Air and HY with hydrophobic PTFE membrane are designed for the sterile filtration and cleaning of air and gases and are optimized for the venting of bags and containers, e.g. fermenters.

Minisart® Acticosart is designed for the ultracleaning of air and gases and contains active carbon for binding of aerosols and volatile substances.

Applications

- Minisart® High Flow for clarification and sterilization of liquids with highest flow rates and total throughput
- Minisart® NML for surfactant-free clarification and particle removal of aqueous liquids
- Minisart® Air | HY for sterilization and cleaning of air and gases and venting of containers
- Minisart® Acticosart for ultracleaning of air and gases

Instructions for Use

If possible, draw a slight amount of air (at least 1 ml) into a syringe (a). Then fill the syringe with the liquid to be filtered: pull the plunger upwards to draw this liquid from a suitable container into the syringe (b). If required, remove any remaining liquid from the tip of the syringe and attach a Minisart® to the luer lock or luer slip connector on the filled syringe. Sterile blister units are opened by peeling off the protective backing. Hold the opened blister part on the outside and connect the inside sterile syringe filter to the syringe (c). Apply consistent pressure to press in the plunger of the syringe in order to filter the liquid through the Minisart® filter into a suitable vial for collection (d). Afterwards, press the plunger all the way in so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter. As a result, this will reduce the hold-up volume.

Caution!

1. If sterile Minisart® blister units are used please check that the blister packing is not damaged!
2. Take care when using syringes with a volume of less than 10 ml, as they can easily generate a pressure greater than the maximum recommended pressure resistance of 6 bar (87 psi) for Minisart®. Therefore, slowly press in the plunger on a syringe with a volume of less than 10 ml. As soon as you detect considerable resistance (= max. filter pressure capacity), do not use force while continuing to press the plunger exceeding 6 bar. Otherwise, you can damage the filter, and filtration will no longer guarantee an ultrapure filtrate!
3. Minisart® is designed for bidirectional use. However, once you have selected one direction of filtration, be sure to maintain this direction. Never use the same syringe filter for both directions!

Note

1. If you need to spike a closure first before injecting the filtrate, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile needle (for a male luer slip) to the outlet of the Minisart® and remove the protective cap before injection.
2. If you would like to achieve maximum recovery of your sample to be filtered, detach the syringe from the syringe filter and draw air once more into the syringe. After reconnecting the syringe and filter, use the air cushion created in this way to press out the remaining liquid. To do so, you will have to apply pressure beyond the bubble point of the particular Minisart® membrane incorporated (see Table 1). This will reduce the hold-up volume to the lowest possible amount.
3. Our filters are leading in terms of purity, so prerinsing is usually not necessary. For very sensitive methods or samples, you may need to additionally prerinse the particular syringe filter you use with 1–3 ml ultrapure water or buffer. This way, you can be sure that you have removed even the slightest traces of extractables.
4. Use a Minisart® for only one sample to reliably prevent carrying over residues of one sample to the next.
5. You will find an overview of chemical compatibilities in the 'Minisart® Chemical Compatibility Guide' on our website.

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Table 1

Specifications for Minisart® High Flow NML with 28 mm filtration area Ø, Air with 15 mm Ø, HY Acticosart with 26 mm Ø										
Housing material	MBS (Methacrylate butadiene styrene)									
Membranes	High Flow: – PES = Polyethersulfone – NML: (SF)CA = (Surfactant-free) Cellulose Acetate – Air HY: PTFE = Polytetrafluoroethylene – Acticosart: PTFE and a dome reservoir filled with active carbon									
Application limits High Flow	Max. recommended operating pressure 6.0 bar 87 psi									
Application limits NML HY Air	Max. recommended operating pressure 4.5 bar 65 psi									
Application limits Acticosart	Max. recommended operating pressure 1 bar 14.5 psi									
Housing burst pressure	>7 bar 102 psi (not determined for Acticosart)									
Max. temperature	60°C not autoclavable!									
Sterilization	Non-sterile Minisart® High Flow NML can be sterilized by ethylene oxide (EO) or Gamma sterilization Non-sterile Minisart® Air HY Acticosart can be sterilized by ethylene oxide (EO) sterilization*									
Minisart ³ type	PES 0.1 µm	PES 0.2 µm	PES 0.45 µm	CA 0.65 µm	CA 0.8 µm	CA 1.2 µm	CA 5.0 µm	PTFE 0.2 µm	Actico-sart	PTFE 1.0 µm
Bubble point (≥)	with water 5.0 bar 73 psi	with water 3.2 bar 46 psi	with water 2.0 bar 29 psi	with water 1.3 bar 19 psi	with water 0.8 bar 12 psi	with water 0.7 bar 10 psi	with water 0.4 bar 6 psi	with ethanol 1.1 bar 16 psi	with ethanol 0.9 bar 13 psi	with ethanol 0.5 bar 7 psi
Flow rate: High Flow NML, 28 mm Ø = 6.2 cm ² filter area Air, 15 mm Ø = 1.7 cm ² filter area HY Acticosart, 26 mm Ø = 5.3 cm ² filter area										
Hold-up volume¹: High Flow NML: 100-150 µl Air HY Acticosart: n.a.										
with water at 1 bar:	40 ml/ min	140 ml/ min	220 ml/ min	250 ml/ min	400 ml/ min	500 ml/ min	600 ml/ min	– ³	– ³	– ³
with air at 0.1 bar:	– ²	– ²	– ²	– ²	– ²	– ²	– ²	2.0 l/min	2.3 l/min	4 l/min
Water penetration point	–	–	–	–	–	–	–	>3.0 bar 44 psi ³	n.a.	>1.5 bar 22 psi
Sterile filtration capability ⁴ acc. to BCT	yes	yes	no	no	no	no	no	yes	n.a.	no
Main applications	hydrophilic filters with PES membrane for fast sterilization or clarification of liquids			hydrophilic filters for particle removal or clarification of aqueous solutions			hydrophobic filters for cleaning or ultracleaning (Actico-sart) of air and gases and venting of containers			

¹ Hold-up volume after air purge. For minimum hold-up please see "Note | Hinweis | Remarques | Nota | Note 1"

² Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

³ Hydrophobic membranes cannot be wetted with aqueous solutions unless you overcome their water penetration point

⁴ According to bacterial challenge test (BCT) with 10⁷ *Brevundimonas diminuta*.

All non-sterile Minisart³ types can be sterilized according to the sterilization processes listed above.
n.a. = not applicable

* Minisart⁶ Air can be sterilized by Gamma sterilization according to the following parameters:
Range 25–40 kGy (validated with 50 kGy).

Table 2

Sartorius Syringe Filters for Liquid Filtration and Venting

Minisart® High Flow | NML | Air | HY | Acticosart

Ø mm	Pore Size	Outlet	Qty. Pk	Order Number	Sterility Status
Minisart® High Flow (PES)					
28 mm	0.1 µm	Male Luer Lock	50	16553-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Lock	50	16532-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Lock	50	16532-----GUK	single-blister-packed, Gamma sterile
28 mm	0.2 µm	Male Luer Lock	500	16532-----Q	non-sterile
28 mm	0.2 µm	Male Luer Slip	50	16541-----K	single-blister-packed, EO sterile
28 mm	0.2 µm	Male Luer Slip	500	16541-----Q	non-sterile
28 mm	0.45 µm	Male Luer Slip	50	16533-----K	single-blister-packed, EO sterile
28 mm	0.45 µm	Male Luer Slip	50	16533-----GUK	single-blister-packed, Gamma sterile
28 mm	0.45 µm	Male Luer Slip	500	16533-----Q	non-sterile
28 mm	0.45 µm	Male Luer Lock	50	16537-----K	single-blister-packed, EO sterile
28 mm	0.45 µm	Male Luer Lock	500	16537-----Q	non-sterile
Minisart® NML (SFCA)					
28 mm	0.65 µm	Male Luer Slip	50	16569-----K	single-blister-packed, EO sterile
28 mm	0.8 µm	Male Luer Lock	50	16592-----K	single-blister-packed, EO sterile
28 mm	0.8 µm	Male Luer Lock	50	16592-----GUK	single-blister-packed, Gamma sterile
28 mm	0.8 µm	Male Luer Lock	500	16592-----Q	non-sterile
28 mm	1.2 µm	Male Luer Lock	50	17593-----K	single-blister-packed, EO sterile
28 mm	1.2 µm	Male Luer Lock	500	17593-----Q	non-sterile
28 mm	5.0 µm	Male Luer Lock	500	17594-----Q	non-sterile
Minisart® Air (hydrophobic PTFE)					
15 mm	0.2 µm	Male Luer Slip	500	1751A-----Q	non-sterile
15 mm	0.2 µm	Needle	50	16596-----HMK	single-blister-packed, Gamma sterile
Minisart® HY (hydrophobic PTFE)					
26 mm	1.0 µm	Male Luer Lock	500	1659A-----HYQ	non-sterile
26 mm	1.0 µm	Male Luer Slip	500	1659B-----HYQ	non-sterile
Minisart® Acticosart with Dome Reservoir + hydrophobic PTFE					
26 mm	0.45 µm	Male Luer Slip	500	17840-----Q	non-sterile, with active carbon
26 mm	0.45 µm	Male Luer Slip	500	40080-----Q	Non-sterile, empty reservoir

- Would you like to filter solvents, acids or bases?
- Do you need to filter small volumes?
- Would you like to use PP housings and other membranes?
Please refer to Minisart® RC, Minisart® NY, or Minisart® SRP for highest chemical compatibility;
also available in 4 mm or 15 mm filtration area
- Do you need Minisart® with pre-filters for filtration of highly particle laden samples?
Please refer to Minisart® NML Plus and Minisart® NY Plus which incorporate a glass fiber pre-filter
- Are you looking for CE-certified Minisart®?
Please request CE-certified article numbers for Minisart® NML, Minisart® HY and Minisart® SRP

Quality Assurance Certificate

Minisart is manufactured by Sartorius in accordance with the Applicable Good Manufacturing Practice Standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Sterile, individually packed products were sterilized with ethylene oxide gas using a validated process following DIN/EN ISO 11135 regulations, or were sterilized with gamma irradiation using a validated process following DIN/EN ISO 11137 regulations.

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no. alone or lot no. with one of the following: pore size, ID no., material type, order no. are printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biosafe according to the USP Class VI Biological tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.72 for non-fiber-releasing filters.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

Non-Sterile and Sterile Products
 Burst Pressure Test
 Bubble Point Test
 Endotoxine Test
 Pressure Hold Test
 Flow Rate Performance

Additional Tests for 0.2 µm Pore Size
 Sterile Filtration Capability (Bacterial Challenge Test)

Additional Tests for Sterile Products
 Sterility

09.03.2015

Date

Hartmut Hennig

Dr. Hartmut Hennig,
 Site Manager QA Filtration Products

Manufactured by
 Sartorius Stedim Biotech GmbH
 August-Spindler-Strasse 11
 37079 Goettingen, Germany
 www.sartorius.com

Symbols

- | | | |
|--|---|--|
| <p> Nicht zur Wiederverwendung
 Single Use
 Ne pas réutiliser
 No reutilizable
 Monouso</p> | <p>LOT Chargenbezeichnung
 Lot Number Code du lot
 Numéro de lote
 Numero di lotto</p> | <p> pyrogenfrei
 non-pyrogenic
 non pyrogène
 Apirógeno
 Apirogeno</p> |
| <p>REF Bestellnummer
 Order Number
 Référence du catalogue
 Número de referencia
 Codice d'ordine</p> | <p>STERILE ETO Sterilisation mit Ethylenoxid
 sterilized, ETO
 Méthode de stérilisation
 utilisant de l'oxyde d'éthylène
 Esterilización por ETO
 Sterilizzazione con ETO</p> | <p> Porengröße
 Pore size
 Taille de pore
 Tamaño de poro
 Porosità</p> |
| <p> Achtung! Warning!
 Attention! Atención!
 Attenzione!</p> | <p>STERILE R Sterilisation durch Bestrahlung
 sterilized, irradiated
 Méthode de stérilisation
 utilisant l'irradiation
 Esterilización por radiación
 Sterilizzazione per irradiazione</p> | <p> Temperaturbegrenzung
 maximum and
 minimum temperature
 Limite de température
 Temperatura máxima
 y mínima
 Temperatura massima
 e minima</p> |
| <p> Hersteller Manufacturer
 Fabricant Fabricante
 Produttore</p> | <p>STERILE I Sterilisation mit Dampf
 sterilized, hot steam
 Méthode de stérilisation
 utilisant la vapeur ou
 la chaleur sèche
 Esterilización por vapor caliente
 Sterilizzazione con vapore caldo</p> | |
| <p>IVD In-Vitro-Diagnostika
 In-Vitro-Diagnostics
 Dispositif médical
 de diagnostic in vitro
 Para uso en diagnóstico
 in vitro
 Per uso diagnostico
 in vitro</p> | <p> Gebrauchsanweisung beachten
 Consult instructions for use
 Consulter le mode d'emploi
 Consulte las instrucciones de uso
 Consultare le istruzioni per l'uso</p> | |
| <p> verwendbar bis
 Use before Utiliser jusqu'au
 Fecha caducidad
 Data di scadenza</p> | | |

Sartorius Stedim Biotech GmbH
 August-Spindler-Strasse 11 • 37079 Goettingen, Germany
 Phone +49.551.308.0 • Fax +49.551.308.3289
 www.sartorius.com