Minisart® NML | HY | SRP

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Sartorius Minisart® filters with hydrophilic cellulose acetate membrane are suitable for particle removal, clarification, bioburden reduction and sterile filtration of aqueous liquids. Filters with hydrophobic PTFE membrane are suitable for sterilization of air and gases and for sterile filtration of oily and aggressive liquids.

Minisart® is free of toxic substances and has been tested for in-vivo usage according to applicable regulations and directives.

Applications

- Minisart® NML with hydrophilic cellulose acetate is virtually free of extractables and is optimized for clarification and sterilization of aqueous samples and drugs.
- Minisart® HY with hydrophobic PTFE is designed for the sterilization and cleaning of air and gases and is optimized for the venting of bags and containers and the filtration of oily solutions. The housing material MBS can be glued with tubing by using suitable glue or residue-free solvents.
- Minisart® SRP with hydrophobic PTFE is designed for the sterilization and cleaning of air and gases and could be used for the venting of bags, containers and fermenters. The PP housing makes it also suitable for the filtration of aggressive solutions.

General Function

- Minisart® NML: 16534—GUK, 16534—K, 16534—Q, 17597—K, 17597—Q

Custom-tailored order number: Minisart® NML: 17597—GIO

Filters designed as infusion accessories, for sterile filtration and for ultracleaning of liquids for eventual infusion, administration or introduction into the human body. Membrane material: CA (cellulose acetate), pore size: 0.2 µm.

- Minisart® HY: 16596—HYK, 16596—HYQ, 16599—HYQ, 40078—Q
- Minisart® SRP: 17575—ACK

Custom-tailored order numbers: Minisart® HY: 16596—ZSK, Minisart® SRP: 17575—ZSK

Medical-grade vent filter for sterile filtration and ultracleaning of air and gases for eventual infusion, administration or introduction into the human body.

Membrane material: PTFE (polytetrafluoroethylene), pore size: 0.2 µm

- Minisart® NML: 16555—GUK, 16555—K, 17598—K

Medical-grade injection filter and vent filter for removal of particles from sterile liquids and from sterile dry air and gases for eventual infusion, administration or introduction into the human body. Membrane material: CA (cellulose acetate), pore size: 0.45 µm


- Minisart® NML: 16555—Q, 17598—Q

Medical-grade injection filter and vent filter for removal of particles from liquids as well as reduction of microbes and/or removal of particles from dry air and gases for eventual infusion, administration or introduction into the human body. Membrane material: CA (cellulose acetate), pore size: 0.45 µm.

- Minisart® NML: 17594—K, 17594—GJR

Medical-grade injection filter and vent filter for removal of particles from sterile liquids and from sterile dry air and gases for eventual infusion, administration or introduction into the human body. Membrane material: CA (cellulose acetate), pore size: 5.0 µm

Instructions for Use

If possible, draw a slight amount of air (at least 1 ml) into a syringe (a). Then fill the syringe with the liquid to be filtered: pull the plunger upwards to draw this liquid from the tip of the syringe and attach a Minisart® to the luer lock or luer slip connector on the filled syringe. [Note: A tight connection is achieved with half a turn or by plugging into the inlet with normal force.] Sterile blister units are opened by peeling off the protective backing. Hold the opened blister part on the outside and connect the inside sterile syringe filter to the syringe (c). Apply consistent pressure to press in the plunger of the syringe in order to filter the liquid through the Minisart® filter into a suitable vial for collection (d). Afterwards, press the plunger all the way in so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter. As a result, this will reduce the hold-up volume.

Caution!

1. If sterile Minisart® blister units are used please check that the blister packing is not damaged. Do not use if you identify any external or internal damage!
2. Minisart® filters are designed for one-time usage only. Never use more than once or for more than one sample to avoid any cross contamination!
3. Do not use Minisart® outside of the technical specifications!
4. Typical filtration volumes of liquids range from 5 – 100 ml. An accurate determination of reachable filter volumes is not possible due to the different particle load of solutions. A decreasing flow through at normal effort is clear evidence for approaching maximum filter pressure capacity. Do not exceed 6 bar to guarantee particle retention according to given pore size!
5. Take care when using syringes with a volume of less than 10 ml, as they can easily generate a pressure greater than the maximum recommended pressure resistance of 6 bar (87 psi) for Minisart®. This could damage the syringe filter and clarification or sterile filtration is no longer guaranteed! We recommend using 10 ml syringes and bigger. For syringes smaller than 10 ml, press down plunger uniform and slowly! If a clear resistance is reached (= maximum filter pressure capacity), please do not continue filtration exceeding 6 bar!
6. The Minisart® unit may be connected to an active medical device in Class IIa or a higher class (e.g. to a drug delivery pump). Take care that the pressure in the system is not greater than the maximum operation pressure of the Minisart® unit.
7. Minisart® is designed for bidirectional use. However, once you have selected one direction of filtration, be sure to maintain this direction. Never use the same syringe filter for both directions!

Note

1. If you need to spike a closure first before injecting the filtrate, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile needle (for a male luer slip) to the outlet of the Minisart® and remove the protective cap before injection.
2. Hold the syringe with filter and needle pointing upwards. Slowly press the liquid upwards. This assures proper wetting of the membrane filter and prevents air bubbles.
3. If you are working with radioactive or infectious material, discard the Minisart® according to the regulations for disposal.
4. Minisart® NML, non-sterile, can be sterilized with ethylene oxide gas (acc. to ISO 11135) or by gamma irradiation (acc. to ISO 11137). Minisart® with a PTFE membrane can be sterilized with ethylene oxide gas (acc. to ISO 11135).
5. For integrity testing of a used Minisart® filter, a bubble point test could be performed. Please ask for our validation guide or download from www.sartorius.com.
6. You will find an overview of chemical compatibilities in the 'Minisart® Chemical Compatibility Guide' on our website.

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Table 1

Specifications for Minisart [®] NML HY SRP, NML with 28 mm membrane filtration area Ø, HY with 26 mm Ø, SRP with 25 mm Ø	
Housing material	- NML HY: MBS (Methacrylate butadiene styrene) - SRP: PP (Polypropylene)
Membranes	- NML: (SF)CA = (Surfactant-free) Cellulose Acetate - HY SRP: PTFE = Polytetrafluoroethylene
Application limits	Max. recommended operating pressure 4.5 bar 65 psi
Housing burst pressure	> 7 bar 102 psi
Max. temperature	- NML HY: 60°C not autoclavable! - SRP: 121°C, 30 min
Sterilization	- Non-sterile Minisart [®] NML can be sterilized by ethylene oxide (EO) or Gamma sterilization - Non-sterile Minisart [®] HY can be sterilized by ethylene oxide (EO) - Non-sterile Minisart [®] SRP can be autoclaved or sterilized by ethylene oxide (EO)

Minisart [®] type	NML, 0,2 µm	NML, 0,45 µm	NML, 5,0 µm	HY SRP, 0,2 µm
Bubble point (z)	with water 3.2 bar 46 psi	with water 2.0 bar 29 psi	with water 0.4 bar 6 psi	with ethanol ¹ 1.4 bar 20 psi
Hold-up volume	100–150 µl	100–150 µl	100–150 µl	100–200 µl
Flow rate				
with water at 1 bar	60 ml/min	180 ml/min	600 ml/min	- ²
with methanol at 1 bar	- ¹	- ¹	- ¹	HY: - ¹ SRP: 160 ml/min
with air at 0.1 bar	- ³	- ³	- ³	HY: 2.0 l/min SRP: 1.8 l/min
Water penetration point ²	-	-	-	>4.0 bar 58 psi
Sterile filtration capability ⁴ acc. to BCI	yes	no	no	yes
Pyrogen-free according to USP	yes	yes	yes	yes
Cytotoxicity	No inhibition with MRC-5 (human lung cells) and L929			

¹ Minisart[®] with MBS housing have a limited stability with solvents! After 1 min exposure liquid-stress-cracking is observed. Bubble point testing of Minisart[®] HY should be performed by a visual bubble point test in less than 60 sec, after pre-wetting of the PTFE membrane with ethanol.

² Hydrophobic membranes cannot be wetted with aqueous solutions unless you overcome their water penetration point.

³ Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

⁴ According to bacterial challenge test (BCI) with 10⁷ *Brevundimonas diminuta*.

Non-sterile Minisart[®] types need to be sterilized before use for sterile filtration.

Ordering Information | Bestellinformation | Données de Commande |
 Datos de Pedido | Dati d'Ordinazione

Table 2
 Standard Sartorius Syringe Filters for Liquid Filtration and Venting
 Minisart® NML | HY | SRP

Ø mm	Pore Size	Outlet	Qty. Pk	Order Number	Sterility Status	CE
Minisart® NML Plus (GF 0.7* µm + CA)						
28 mm	0.2 µm	Male Luer Lock	50	16534-----K	single-blister-packed, EO sterile	Yes
28 mm	0.2 µm	Male Luer Lock	50	16534-----GUK	single-blister-packed, Gamma sterile	Yes
28 mm	0.2 µm	Male Luer Slip	500	16534-----Q	non-sterile	Yes
28 mm	0.2 µm	Male Luer Slip	50	17597-----K	single-blister-packed, EO sterile	Yes
28 mm	0.2 µm	Male Luer Slip	500	17597-----Q	non-sterile	Yes
28 mm	0.2 µm	Male Luer Slip	50	17528-----K	single-blister-packed, EO sterile	Yes
28 mm	0.45 µm	Male Luer Slip	50	16555-----K	single-blister-packed, EO sterile	Yes
28 mm	0.45 µm	Male Luer Slip	50	16555-----GUK	single-blister-packed, Gamma sterile	Yes
28 mm	0.45 µm	Male Luer Slip	500	16555-----Q	non-sterile	Yes
28 mm	0.45 µm	Male Luer Lock	50	17598-----K	single-blister-packed, EO sterile	Yes
28 mm	0.45 µm	Male Luer Lock	500	17598-----Q	non-sterile	Yes
28 mm	5.0 µm	Male Luer Lock	50	17594-----K	single-blister-packed, EO sterile	Yes
28 mm	5.0 µm	Male Luer Lock	1000	17594-----GJR	single-blister-packed, Gamma sterile	Yes
Minisart® HY (PTFE)						
26 mm	0.2 µm	Male Luer Lock	50	16596-----HYK	single-blister-packed, EO sterile	Yes
26 mm	0.2 µm	Male Luer Lock	500	16596-----HYQ	non-sterile	Yes
26 mm	0.2 µm	Male Luer Lock	50	16598-----HYQ	non-sterile, with luer slip inlet	Yes
26 mm	0.2 µm	Tube connectors	500	40078-----Q	non-sterile	Yes
Minisart® SRP (PTFE)						
25 mm	0.2 µm	Male Luer Slip	50	17575-----ACK	single-blister-packed, EO sterile	Yes

- Would you like to filter solvents, acids or bases?
- Do you need to filter small volumes?
- Would you like to use PP housings and other membranes?
 Please refer to Minisart® RC, Minisart® NY, or Minisart® SRP for highest chemical compatibility;
 also available in 4 mm or 15 mm filtration area
- Do you need Minisart® with pre-filters for filtration of highly particle laden samples?
 Please refer to Minisart® NML Plus and Minisart® NY Plus which incorporate a glass fiber pre-filter
- Are you looking for other pore sizes (0.8 µm, 1.0 µm, 1.2 µm)?
 Please request corresponding article numbers for Minisart® NML

DOMINIQUE DUTER S.A.S

Quality Assurance Certificate

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no., alone or lot no. with one of the following: pore size, ID no., material type, order no., etc. printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biocompatible according to the USP Class VI Biological Tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.22 for non-fiber-releasing filters.

In addition Minisart is classified in class IIa of the Medical Devices Directive 93/42/EEC.

12.03.2014
Date

Minisart is manufactured by Sartorius in accordance with the Applicable Good Manufacturing Practice Standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Sterile, individually packed products were sterilized with ethylene oxide gas using a validated process following DIN/EN ISO 11135 regulations, or were sterilized with gamma irradiation using a validated process following DIN/EN ISO 11137 regulations.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

Non-Sterile and Sterile Products
Burst Pressure Test
Bubble Point Test
Endotoxine Test
Pressure Hold Test
Flow Rate Performance

Additional Tests for 0.2 µm Pore Size
Sterile Filtration Capability (Bacterial Challenge Test)

Additional Tests for Sterile Products
Sterility

Hartmut Hennig
 Dr. Hartmut Hennig,
 Site Manager QA Filtration Products

Manufactured by
 Sartorius Stedim Biotech GmbH
 August-Spindler-Strasse 11
 37079 Goettingen, Germany
 www.sartorius.com

Symbols

<p> Nicht zur Wiederverwendung Single Use Ne pas réutiliser No reutilizable Monouso</p> <p> Bestellnummer Order Number Référence du catalogue Número de referencia Codice d'ordine</p> <p> Achtung! Warning! Attention! Atención! Attenzione!</p> <p> Hersteller Manufacturer Fabricant Fabricante Produttore</p> <p> In-Vitro-Diagnostika In-Vitro-Diagnostics Dispositif médical de diagnostic in vitro Para uso en diagnóstico in vitro Per uso diagnostico in vitro</p> <p> verwendbar bis Use before Utiliser jusqu'au Fecha caducidad Data di scadenza</p>	<p> LOT Chargenbezeichnung Lot Number Code du lot Número de lote Numero di lotto</p> <p> STERILE EO Sterilisation mit Ethylenoxid sterilized, ETO Méthode de stérilisation utilisant de l'oxyde d'éthylène Esterilización por ETO Sterilizzazione con ETO</p> <p> STERILE R Sterilisation durch Bestrahlung sterilized, irradiated Méthode de stérilisation utilisant l'irradiation Esterilización por radiación Sterilizzazione per irradiazione</p> <p> STERILE I Sterilisation mit Dampf sterilized, hot steam Méthode de stérilisation utilisant la vapeur ou la chaleur sèche Esterilización por vapor caliente Sterilizzazione con vapore caldo</p> <p> Gebrauchsanweisung beachten Consult instructions for use Consulter le mode d'emploi Consultare le istruzioni per l'uso</p>	<p> pyrogenfrei non-pyrogenic non pyrogène Apirógeno Apirogeno</p> <p> Porengröße Pore size Taille de pore Tamaño de poro Porosità</p> <p> Temperaturbegrenzung maximum and minimum temperature Limite de température Temperatura máxima y mínima Temperatura massima e minima</p>
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------