## Minisart® NML Plus | NY Plus

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## Minisart® NML Plus | NY Plus

High-capacity Sartorius Syringe Filters

Sartorius Minisart® filters with glass fiber (GF) pre-filters are the optimal solution for highly-particle-laden samples. Minisart® Plus filters contain a GF filter only or a GF pre-filter on top of a filter membrane. The GF depth-filter on top of a membrane protects the syringe filter against early blockage.


Minisart® NML Plus is optimized for clarification and sterile filtration of aqueous solutions. The binder-free GF depth-filter has a high loading capacity and highly increases reachable filter volumes.

Minisart® NY Plus with ultrapure quartz pre-filters are virtually free of extractables for a pure filtrate and are optimal for sample preparation prior to analytics.

### Applications

- **Minisart® NML Plus** for cleaning of aqueous difficult-to-filter liquids
- **Minisart® NY Plus** for ultracleaning of difficult-to-filter aqueous and organic samples

### Instructions for Use

If possible, draw a slight amount of air (at least 1 ml) into a syringe (a). Then fill the syringe with the liquid to be filtered: pull the plunger upwards to draw this liquid from a suitable container into the syringe (b). If required, remove any remaining liquid from the tip of the syringe and attach a Minisart® to the luer lock or luer slip connector on the filled syringe. Sterile blister units are opened by peeling off the protective backing. Hold the opened blister part on the outside and connect the inside sterile syringe filter to the syringe (c). Apply consistent pressure to press in the plunger of the syringe in order to filter the liquid through the Minisart® filter into a suitable vial for collection (d). Afterwards, press the plunger all the way in so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter. As a result, this will reduce the hold-up volume.

### Caution!

1. Minisart-Plus units are for laboratory use only, not for in-vivo applications.
2. If sterile Minisart® blister units are used please check that the blister packing is not damaged!
3. Take care when using syringes with a volume of less than 10 ml, as they can easily generate a pressure greater than the maximum recommended pressure resistance of 6 bar (87 psi) for Minisart®. Therefore, slowly press in the plunger on a syringe with a volume of less than 10 ml. As soon as you detect considerable resistance (= max. filter pressure capacity), do not use force while continuing to press the plunger. Otherwise, you can damage the filter, and filtration will no longer guarantee an ultrapure filtrate!
4. Filtration must always be from the female luer lock inlet side (printed housing part) to the outlet side (NML Plus: colored housing part, NY Plus: non-printed housing part)!

### Note

1. If you need to spike a closure first before injecting the filtrate, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile needle (for a male luer slip) to the outlet of the Minisart® and remove the protective cap before injection.
2. If you would like to achieve maximum recovery of your sample to be filtered, detach the syringe from the syringe filter and draw air once more into the syringe. After reconnecting the syringe and filter, use the air cushion created in this way to press out the remaining liquid. To do so, you will have to apply pressure beyond the bubble point of the particular Minisart® membrane incorporated (see Table 1). This will reduce the hold-up volume to the lowest possible amount.
3. Our filters are leading in terms of purity, so prerinsing is usually not necessary. In case of very sensitive methods or critical solutions, filters can be pre-washed with 1-3 ml ultrapure water or buffer. This way, you can be sure that you have removed even the slightest traces of extractables.
4. Use a Minisart® for only one sample to reliably prevent carrying over residues of one sample to the next.
5. You will find an overview of chemical compatibilities in the 'Minisart® Chemical Compatibility Guide' on our website.

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Table 1

Specifications for Minisart® NML Plus   NY Plus, NML Plus with 28 mm, NY Plus with 25 mm membrane filtration area Ø						
Housing material	– NML Plus: MBS (Methacrylate butadiene styrene) – NY Plus: Polypropylene (PP)					
Glass fiber pre-filter	– NML Plus: Binder-free GF, 0.7 <sup>1</sup> µm pore size – NY Plus: Ultrapure quartz, 0.7 <sup>1</sup> µm pore size					
Membranes	– NML Plus: (SF)CA = (Surfactant-free) Cellulose Acetate – NY Plus: NY = Polyamide					
Application limits	Max. recommended operating pressure 4,5 bar   65 psi					
Housing burst pressure	>7 bar   102 psi					
Max. temperature	– NML Plus: 60°C not autoclavable! – NY Plus: 121°C, 30 min (autoclavable)					
Sterilization	– Non-sterile Minisart <sup>1</sup> NML Plus can be sterilized by ethylene oxide (EO) or Gamma sterilization – Non-sterile Minisart <sup>5</sup> NY Plus can be autoclaved or sterilized by ethylene oxide (EO) sterilization					
Minisart <sup>3</sup> type	GF+CA 0,2 µm	GF+CA 0,45 µm	GF+CA 1,2 µm	GF 0,7 <sup>1</sup> µm	GF+NY 0,2 µm	GF+NY 0,45 µm
Bubble point (≥)	with water 3.2 bar 46 psi	with water 2.0 bar 29 psi	with water 0.7 bar 10 psi	with water 0.5 bar 7 psi	with water 3.0 bar 44 psi	with water 2.0 bar 29 psi
<b>Flow rate, NML Plus, 28 mm Ø = 6,2 cm<sup>2</sup> filter area</b>						
<b>NY Plus, 25 mm Ø = 4,8 cm<sup>2</sup> filter area</b>						
Hold-up volume <sup>2</sup> : 200–250 µl						
with water at 1 bar	60 ml/min	160 ml/min	350 ml/min	450 ml/min	50 ml/min	100 ml/min
with methanol at 1 bar	–	–	–	–	70 ml/min	200 ml/min
with air at 0,1 bar	– <sup>3</sup>	– <sup>3</sup>	– <sup>3</sup>	– <sup>3</sup>	– <sup>3</sup>	– <sup>3</sup>
Sterile filtration capability <sup>4</sup> acc. to BCF	yes	no	no	no	yes	no
Main applications	hydrophilic filters with GF pre-filter for difficult-to-filter aqueous samples with high particle load				hydrophilic filters with ultra-pure quartz pre-filter for sample preparation of aqueous and organic solutions	

<sup>1</sup> 0.7 µm is the nominal particle retention ≠ pore size of the used GF pre-filter in Minisart® NML Plus and NY Plus

<sup>2</sup> Hold-up volume after air purge. For minimum hold-up please see “Note|Hinweis|Précisions|Nota|Note 1”

<sup>3</sup> Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

<sup>4</sup> According to bacterial challenge test with 10<sup>7</sup> *Brevundimonas diminuta*. All non-sterile Minisart® types listed above can be sterilized according to the sterilization recommendation in this table.

**Table 2**

**High Capacity Syringe Filters**

Minisart® NML Plus | NY Plus

Ø mm	Pore Size	Outlet	Qty.   Pk	Order Number	Sterility Status
<b>Minisart® NML Plus (GF 0.7* µm + CA)</b>					
28 mm	0.2 µm	Male Luer Lock	50	17823-----K	single-blister-packed, sterile
28 mm	0.2 µm	Male Luer Lock	500	17823-----Q	non-sterile
28 mm	0.45 µm	Male Luer Lock	50	17829-----K	single-blister-packed, sterile
28 mm	0.2 µm	Male Luer Lock	500	17829-----Q	non-sterile
28 mm	1.2 µm	Male Luer Lock	500	17825-----Q	non-sterile
<b>Minisart® NML Plus (with GF filter only)</b>					
28 mm	0.7* µm	Male Luer Lock	50	17824-----K	non-sterile
28 mm	0.7* µm	Male Luer Lock	500	17824-----Q	non-sterile
<b>Minisart® NY Plus (GF 0.7* µm + NY)</b>					
25 mm	0.2 µm	Male Luer Lock	50	1784B-----K	ultrapure for HPLC sample prep, non-sterile
25 mm	0.2 µm	Male Luer Lock	500	1784B-----Q	ultrapure for HPLC sample prep, non-sterile
25 mm	0.45 µm	Male Luer Lock	50	1784C-----K	ultrapure for HPLC sample prep, non-sterile
25 mm	0.45 µm	Male Luer Lock	500	1784C-----Q	ultrapure for HPLC sample prep, non-sterile

\* Particle retention of GF depth filter

- You don't need a glass fiber pre-filter?  
Please refer to Minisart® NML and Minisart® NY
- Would you like to use other membranes?  
Please refer to Minisart® High Flow
- You need a high chemical compatibility?
- You need to filter small volumes?  
Please refer to Minisart® RC, Minisart® NY, or Minisart® SRP for highest chemical compatibility;  
also available in 4 mm or 15 mm filtration area Ø
- Are you looking for CE-certified Minisart®?  
Please request CE-certified article numbers for Minisart® NML, Minisart® HY and Minisart® SRP

## Quality Assurance Certificate

Minisart is manufactured by Sartorius in accordance with the Applicable Good Manufacturing Practice Standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Sterile, individually packed products were sterilized with ethylene oxide gas using a validated process following DIN/EN ISO 11135 regulations.

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no. alone or lot no. with one of the following: pore size, ID no., material type, order no. are printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biosafe according to the USP Class VI Biological tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.72 for non-fiber-releasing filters.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

Non-Sterile and Sterile Products  
 Burst Pressure Test  
 Bubble Point Test  
 Pressure Hold Test  
 Flow Rate Performance

Additional Tests for 0.2 µm Pore Size  
 Sterile Filtration Capability (Bacterial Challenge Test)

Additional Tests for Sterile Products  
 Sterility

09.03.2015

Date

*Hartmut Hennig*

Dr. Hartmut Hennig,  
 Site Manager QA Filtration Products

Manufactured by  
 Sartorius Stedim Biotech GmbH  
 August-Spindler-Strasse 11  
 37079 Goettingen, Germany  
 www.sartorius.com

## Symbols

- |  |   |  |
|--|---|--|
| <p> Nicht zur Wiederverwendung<br/>                 Single Use<br/>                 Ne pas réutiliser<br/>                 No reutilizable<br/>                 Monouso</p>  | <p><b>LOT</b> Chargenbezeichnung<br/>                 Lot Number   Code du lot<br/>                 Numéro de lote<br/>                 Numero di lotto</p>   | <p> pyrogenfrei<br/>                 non-pyrogenic<br/>                 non pyrogène<br/>                 Apirógeno<br/>                 Apirogeno</p>   |
| <p><b>REF</b> Bestellnummer<br/>                 Order Number<br/>                 Référence du catalogue<br/>                 Número de referencia<br/>                 Codice d'ordine</p>   | <p><b>STERILE ETO</b> Sterilisation mit Ethylenoxid<br/>                 sterilized, ETO<br/>                 Méthode de stérilisation<br/>                 utilisant de l'oxyde d'éthylène<br/>                 Esterilización por ETO<br/>                 Sterilizzazione con ETO</p>  | <p> Porengröße<br/>                 Pore size<br/>                 Taille de pore<br/>                 Tamaño de poro<br/>                 Porosità</p>  |
| <p> Achtung!   Warning!<br/>                 Attention!   Atención!<br/>                 Attenzione!</p>   | <p><b>STERILE R</b> Sterilisation durch Bestrahlung<br/>                 sterilized, irradiated<br/>                 Méthode de stérilisation<br/>                 utilisant l'irradiation<br/>                 Esterilización por radiación<br/>                 Sterilizzazione per irradiazione</p>                                  | <p> Temperaturbegrenzung<br/>                 maximum and<br/>                 minimum temperature<br/>                 Limite de température<br/>                 Temperatura máxima<br/>                 y mínima<br/>                 Temperatura massima<br/>                 e minima</p> |
| <p> Hersteller   Manufacturer<br/>                 Fabricant   Fabricante<br/>                 Produttore</p>  | <p><b>STERILE I</b> Sterilisation mit Dampf<br/>                 sterilized, hot steam<br/>                 Méthode de stérilisation<br/>                 utilisant la vapeur ou<br/>                 la chaleur sèche<br/>                 Esterilización por vapor caliente<br/>                 Sterilizzazione con vapore caldo</p> |  |
| <p><b>IVD</b> In-Vitro-Diagnostika<br/>                 In-Vitro-Diagnostics<br/>                 Dispositif médical<br/>                 de diagnostic in vitro<br/>                 Para uso en diagnóstico<br/>                 in vitro<br/>                 Per uso diagnostico<br/>                 in vitro</p> | <p> Gebrauchsanweisung beachten<br/>                 Consult instructions for use<br/>                 Consulter le mode d'emploi<br/>                 Consulte las instrucciones de uso<br/>                 Consultare le istruzioni per l'uso</p>  |  |
| <p> verwendbar bis<br/>                 Use before   Utiliser jusqu'au<br/>                 Fecha caducidad<br/>                 Data di scadenza</p>  |   |  |