Minisart® NML Plus | NY Plus

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Minisart® NML Plus | NY Plus

High-capacity Sartorius Syringe Filters

Sartorius Minisart® filters with glass fiber (GF) pre-filters are the optimal solution for highly-particle-laden samples. Minisart® Plus filters contain a GF filter only or a GF pre-filter on top of a filter membrane. The GF depth-filter on top of a membrane protects the syringe filter against early blockage.


Minisart® NML Plus is optimized for clarification and sterile filtration of aqueous solutions. The binder-free GF depth-filter has a high loading capacity and highly increases reachable filter volumes.

Minisart® NY Plus with ultrapure quartz pre-filters are virtually free of extractables for a pure filtrate and are optimal for sample preparation prior to analytics.

Applications

- **Minisart® NML Plus** for cleaning of aqueous difficult-to-filter liquids
- **Minisart® NY Plus** for ultracleaning of difficult-to-filter aqueous and organic samples

Instructions for Use

If possible, draw a slight amount of air (at least 1 ml) into a syringe (a). Then fill the syringe with the liquid to be filtered: pull the plunger upwards to draw this liquid from a suitable container into the syringe (b). If required, remove any remaining liquid from the tip of the syringe and attach a Minisart® to the luer lock or luer slip connector on the filled syringe. Sterile blister units are opened by peeling off the protective backing. Hold the opened blister part on the outside and connect the inside sterile syringe filter to the syringe (c). Apply consistent pressure to press in the plunger of the syringe in order to filter the liquid through the Minisart® filter into a suitable vial for collection (d). Afterwards, press the plunger all the way in so that the air cushion initially created will discharge any liquid remaining in the inlet and outlet of the filter. As a result, this will reduce the hold-up volume.

Caution!

1. Minisart-Plus units are for laboratory use only, not for in-vivo applications.
2. If sterile Minisart® blister units are used please check that the blister packing is not damaged!
3. Take care when using syringes with a volume of less than 10 ml, as they can easily generate a pressure greater than the maximum recommended pressure resistance of 6 bar (87 psi) for Minisart®. Therefore, slowly press in the plunger on a syringe with a volume of less than 10 ml. As soon as you detect considerable resistance (= max. filter pressure capacity), do not use force while continuing to press the plunger. Otherwise, you can damage the filter, and filtration will no longer guarantee an ultrapure filtrate!
4. Filtration must always be from the female luer lock inlet side (printed housing part) to the outlet side (NML Plus: colored housing part, NY Plus: non-printed housing part)!

Note

1. If you need to spike a closure first before injecting the filtrate, use a suitable disinfectant, such as 70% alcohol, to disinfect the site to be pierced. If necessary, attach a sterile needle (for a male luer slip) to the outlet of the Minisart® and remove the protective cap before injection.
2. If you would like to achieve maximum recovery of your sample to be filtered, detach the syringe from the syringe filter and draw air once more into the syringe. After reconnecting the syringe and filter, use the air cushion created in this way to press out the remaining liquid. To do so, you will have to apply pressure beyond the bubble point of the particular Minisart® membrane incorporated (see Table 1). This will reduce the hold-up volume to the lowest possible amount.
3. Our filters are leading in terms of purity, so prerinsing is usually not necessary. In case of very sensitive methods or critical solutions, filters can be pre-washed with 1-3 ml ultrapure water or buffer. This way, you can be sure that you have removed even the slightest traces of extractables.
4. Use a Minisart® for only one sample to reliably prevent carrying over residues of one sample to the next.
5. You will find an overview of chemical compatibilities in the 'Minisart® Chemical Compatibility Guide' on our website.

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Table 1

| Specifications for Minisart® NML Plus NY Plus, NML Plus with 28 mm, NY Plus with 25 mm membrane filtration area Ø | | | | | | |
|--|--|------------------------------|------------------------------|-----------------------------|---|------------------------------|
| Housing material | – NML Plus: MBS (Methacrylate butadiene styrene) – NY Plus: Polypropylene (PP) | | | | | |
| Glass fiber pre-filter | – NML Plus: Binder-free GF, 0.7 ¹ µm pore size – NY Plus: Ultrapure quartz, 0.7 ¹ µm pore size | | | | | |
| Membranes | – NML Plus: (SF)CA = (Surfactant-free) Cellulose Acetate – NY Plus: NY = Polyamide | | | | | |
| Application limits | Max. recommended operating pressure 4,5 bar 65 psi | | | | | |
| Housing burst pressure | >7 bar 102 psi | | | | | |
| Max. temperature | – NML Plus: 60°C not autoclavable! – NY Plus: 121°C, 30 min (autoclavable) | | | | | |
| Sterilization | – Non-sterile Minisart ¹ NML Plus can be sterilized by ethylene oxide (EO) or Gamma sterilization – Non-sterile Minisart ⁵ NY Plus can be autoclaved or sterilized by ethylene oxide (EO) sterilization | | | | | |
| Minisart ³ type | GF+CA 0,2 µm | GF+CA 0,45 µm | GF+CA 1,2 µm | GF 0,7 ¹ µm | GF+NY 0,2 µm | GF+NY 0,45 µm |
| Bubble point (≥) | with water 3.2 bar 46 psi | with water 2.0 bar 29 psi | with water 0.7 bar 10 psi | with water 0.5 bar 7 psi | with water 3.0 bar 44 psi | with water 2.0 bar 29 psi |
| Flow rate, NML Plus, 28 mm Ø = 6,2 cm² filter area | | | | | | |
| NY Plus, 25 mm Ø = 4,8 cm² filter area | | | | | | |
| Hold-up volume ² : 200–250 µl | | | | | | |
| with water at 1 bar | 60 ml/min | 160 ml/min | 350 ml/min | 450 ml/min | 50 ml/min | 100 ml/min |
| with methanol at 1 bar | – | – | – | – | 70 ml/min | 200 ml/min |
| with air at 0,1 bar | – ³ | – ³ | – ³ | – ³ | – ³ | – ³ |
| Sterile filtration capability ⁴ acc. to BCF | yes | no | no | no | yes | no |
| Main applications | hydrophilic filters with GF pre-filter for difficult-to-filter aqueous samples with high particle load | | | | hydrophilic filters with ultra-pure quartz pre-filter for sample preparation of aqueous and organic solutions | |

¹ 0.7 µm is the nominal particle retention ≠ pore size of the used GF pre-filter in Minisart® NML Plus and NY Plus

² Hold-up volume after air purge. For minimum hold-up please see "Note|Hinweis|Précisions|Nota|Note 1"

³ Hydrophilic membranes can filter dry air or gas but become impermeable to air or gas when wetted!

⁴ According to bacterial challenge test with 10⁷ *Brevundimonas diminuta*. All non-sterile Minisart® types listed above can be sterilized according to the sterilization recommendation in this table.

Table 2
High Capacity Syringe Filters
 Minisart® NML Plus | NY Plus

| Ø mm | Pore Size | Outlet | Qty. Pk | Order Number | Sterility Status |
|---|-----------|----------------|-----------|--------------|---|
| Minisart® NML Plus (GF 0.7* µm + CA) | | | | | |
| 28 mm | 0.2 µm | Male Luer Lock | 50 | 17823-----K | single-blister-packed, sterile |
| 28 mm | 0.2 µm | Male Luer Lock | 500 | 17823-----Q | non-sterile |
| 28 mm | 0.45 µm | Male Luer Lock | 50 | 17829-----K | single-blister-packed, sterile |
| 28 mm | 0.2 µm | Male Luer Lock | 500 | 17829-----Q | non-sterile |
| 28 mm | 1.2 µm | Male Luer Lock | 500 | 17825-----Q | non-sterile |
| Minisart® NML Plus (with GF filter only) | | | | | |
| 28 mm | 0.7* µm | Male Luer Lock | 50 | 17824-----K | non-sterile |
| 28 mm | 0.7* µm | Male Luer Lock | 500 | 17824-----Q | non-sterile |
| Minisart® NY Plus (GF 0.7* µm + NY) | | | | | |
| 25 mm | 0.2 µm | Male Luer Lock | 50 | 1784B-----K | ultrapure for HPLC sample prep, non-sterile |
| 25 mm | 0.2 µm | Male Luer Lock | 500 | 1784B-----Q | ultrapure for HPLC sample prep, non-sterile |
| 25 mm | 0.45 µm | Male Luer Lock | 50 | 1784C-----K | ultrapure for HPLC sample prep, non-sterile |
| 25 mm | 0.45 µm | Male Luer Lock | 500 | 1784C-----Q | ultrapure for HPLC sample prep, non-sterile |

* Particle retention of GF depth filter

- You don't need a glass fiber pre-filter?
Please refer to Minisart® NML and Minisart® NY
- Would you like to use other membranes?
Please refer to Minisart® High Flow
- You need a high chemical compatibility?
- You need to filter small volumes?
Please refer to Minisart® RC, Minisart® NY, or Minisart® SRP for highest chemical compatibility;
also available in 4 mm or 15 mm filtration area Ø
- Are you looking for CE-certified Minisart®?
Please request CE-certified article numbers for Minisart® NML, Minisart® HY and Minisart® SRP

Quality Assurance Certificate

Minisart is manufactured by Sartorius in accordance with the Applicable Good Manufacturing Practice Standards and in a facility whose Quality Management System is certified for compliance with the DIN EN ISO 9001 Quality System Standard and DIN EN ISO 13485.

Sterile, individually packed products were sterilized with ethylene oxide gas using a validated process following DIN/EN ISO 11135 regulations.

Each unit is tested 100% during manufacture for membrane and housing integrity by a leakage test under automatic conditions. The lot no. alone or lot no. with one of the following: pore size, ID no., material type, order no. are printed on top of the housing and/or Tyvek. Before packing each unit has been checked by visual inspection.

Minisart is biosafe according to the USP Class VI Biological tests for Plastics. It also complies with the Title 21 Code of Federal Regulations, Section 210.3(b)(6) and 211.72 for non-fiber-releasing filters.

Following 100% control, each lot has been sampled, tested and released by QA Department for the following characteristics:

Non-Sterile and Sterile Products
 Burst Pressure Test
 Bubble Point Test
 Pressure Hold Test
 Flow Rate Performance

Additional Tests for 0.2 µm Pore Size
 Sterile Filtration Capability (Bacterial Challenge Test)

Additional Tests for Sterile Products
 Sterility

09.03.2015

Date

Hartmut Hennig

Dr. Hartmut Hennig,
 Site Manager QA Filtration Products

Manufactured by
 Sartorius Stedim Biotech GmbH
 August-Spindler-Strasse 11
 37079 Goettingen, Germany
 www.sartorius.com

Symbols

- | | | |
|--|---|--|
| <p> Nicht zur Wiederverwendung Single Use Ne pas réutiliser No reutilizable Monouso</p> | <p>LOT Chargenbezeichnung Lot Number Code du lot Numéro de lote Numero di lotto</p> | <p> pyrogenfrei non-pyrogenic non pyrogène Apirógeno Apirogeno</p> |
| <p>REF Bestellnummer Order Number Référence du catalogue Número de referencia Codice d'ordine</p> | <p>STERILE ETO Sterilisation mit Ethylenoxid sterilized, ETO Méthode de stérilisation utilisant de l'oxyde d'éthylène Esterilización por ETO Sterilizzazione con ETO</p> | <p> Porengröße Pore size Taille de pore Tamaño de poro Porosità</p> |
| <p> Achtung! Warning! Attention! Atención! Attenzione!</p> | <p>STERILE R Sterilisation durch Bestrahlung sterilized, irradiated Méthode de stérilisation utilisant l'irradiation Esterilización por radiación Sterilizzazione per irradiazione</p> | <p> Temperaturbegrenzung maximum and minimum temperature Limite de température Temperatura máxima y mínima Temperatura massima e minima</p> |
| <p> Hersteller Manufacturer Fabricant Fabricante Produttore</p> | <p>STERILE I Sterilisation mit Dampf sterilized, hot steam Méthode de stérilisation utilisant la vapeur ou la chaleur sèche Esterilización por vapor caliente Sterilizzazione con vapore caldo</p> | |
| <p>IVD In-Vitro-Diagnostika In-Vitro-Diagnostics Dispositif médical de diagnostic in vitro Para uso en diagnóstico in vitro Per uso diagnostico in vitro</p> | <p> Gebrauchsanweisung beachten Consult instructions for use Consulter le mode d'emploi Consulte las instrucciones de uso Consultare le istruzioni per l'uso</p> | |
| <p> verwendbar bis Use before Utiliser jusqu'au Fecha caducidad Data di scadenza</p> | | |