
sartorius stedim  
biotech

Directions for Use | Gebrauchsanleitung | Mode d'emploi  
Instrucciones para el uso | Istruzioni per l'uso

## Midisart® 2000

17804 | 17805 | 17877

82521-060-50

# Directions for Use of Midisart 2000

## 17804|17805|17877

### 1. Product Description and Applications

Midisart 2000 ready-to-connect filtration units consist of a hydrophobic (PTFE) membrane filter in a polypropylene housing, and are supplied in sterile or unsterile condition. The filtration unit can be repeatedly autoclaved. They are available with 0.2  $\mu\text{m}$  and 0.45  $\mu\text{m}$  pore size membranes, and are ideally suitable for all small-scale air/gas sterilizing purposes, such as:

- Sterile venting of filling vessels and fermentation carboys including culture vessels
- Venting of holding tanks for sterile, distilled water and liquid culture media
- Autoclave venting
- In-line sterilization of and particulate removal from air and gases

The excellent chemical stability of these practical units also makes them useful for removal of particles from chemicals and organic solvents.

The female luer connection in the hose nipple allows the unit to be firmly fitted onto an appropriate syringe.

### 2. Order Numbers

Fore size in $\mu\text{m}$	Inlet and outlet connectors	Box of 12, sterile, individually packaged units	Box of 25, sterile, individually packaged units	Box of 100, bulk-packaged, non-sterile units
0.2 $\mu\text{m}$	Hose nipple	17805-E	17805-G	17805-UPN
0.2 $\mu\text{m}$	Modified hose nipple	-	-	17877-UPN
0.45 $\mu\text{m}$	Hose nipple	17805-E	17804-G	
0.2 $\mu\text{m}$	1/8" NPT thread	17805-NPE	17805-NPG	
0.45 $\mu\text{m}$	1/8" NPT thread	17804-NPE	17804-NPG	

### 3. Directions for Use

3.1 Remove a Midisart from the box.

3.2 Before using "sterile" Midisarts, check the bag to make sure that it is completely intact.

3.3 Cut open the packaging on one side.

3.4 Hold the packaging in one hand and fit tubing, or a syringe, to the Midisart hose nipple without touching it.

Midisart units can also be used for the filtration of aqueous liquids, but the hydrophobic PTFE membrane filter must first be wetted with a water-miscible solvent, such as ethanol.

#### 4. Specifications

Filter material: PTFE reinforced with polypropylene gauze  
Housing material: Polypropylene  
Filtration area: 20 cm<sup>2</sup>  
Housing diameter: 64 mm  
Priming volume: approx. 3 ml  
Hold-up volume: approx. 1 ml before, 0.5 ml after bubble point  
Max. operating pressure: 300 kPa (3 bar)  
Water penetration point: approx. 300 kPa (0.45 µm), 400 kPa (0.2 µm)  
Max. temperature (autoclave): 134°C  
Presterilized with: Ethylene oxide  
Biosafety: Materials pass USPXXII Plastics class VI  
Bubble point (isopropanol, 60%):  
0.45 µm: ≥0.9 bar  
0.2 µm: ≥1.1 bar  
Connections: Inlet and outlet: Conical hose nipple for hoses with an inner diameter of 6–12 mm, modified hose nipple for hoses with an inner diameter of 4–12 mm, as well as a female luer inlet connector, for use on syringes;  
or  
Inlet and outlet: 1/8" NPT thread

#### Typical Air Flow Rates

Flow rate for air at	Pore size 0.2 µm	Pore size 0.45 µm
Δp = 0.02 bar	1.1 l/min	1.8 l/min
Δp = 0.05 bar	2.9 l/min	4.6 l/min
Δp = 0.10 bar	5.0 l/min	8.5 l/min

1 bar = 100 kPa = 14.5 psi

#### 5. Chemical Compatibility

The materials (PP and PTFE) used in Midisart give it excellent compatibility with the products listed below. However, its compatibility can be affected by various factors, such as temperature, concentration, composition, etc. We therefore recommend that you test whether Midisart is compatible with the particular medium you wish to filter.

#### Media

Acetic acid (concentrated)  
Acetone  
Acetonitrile  
n-Butanol  
Cellosolve (ethyl)  
Chloroform  
Diethylacetamide  
Dimethyl formamide  
Dimethylsulfoxide  
Dioxane  
Ethyl acetate  
Ethanol  
Ethylene glycol  
Freon TF  
Gasoline  
1 N hydrochloric acid  
Hexane  
Isobutanol  
Isopropanol  
Methanol  
Methylene chloride  
Methyl ethyl ether  
Methyl ethyl ketone  
Pentane  
Sodium hydroxide (5%)  
Tetrahydrofuran  
Toluene  
Trichloroacetic acid  
Trichloroethane  
Water  
Xylene

In the interest of the further development of Sartorius Stedim Biotech products, we reserve the right to make changes without notice.

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Phone +49.551.308.0  
Fax +49.551.308.3289  
www.sartorius-stedim.com

Copyright by  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany.

All rights reserved. No part  
of this publication may  
be reprinted or translated in  
any form or by any means  
without the prior written  
permission of Sartorius Stedim  
Biotech GmbH.

The status of the information,  
specifications and illustrations  
in this manual is indicated  
by the date given below.  
Sartorius Stedim Biotech GmbH  
reserves the right to make  
changes to the technology,  
features, specifications and  
design of the equipment  
without notice.

Status:  
February 2012,  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany

Printed in Germany on paper that has  
been bleached without any use of chlorine  
W1\_Midisæt 2000 - KT  
Publication No.: SLD6017-p12021  
Ver. 02 | 2012