Operating Instructions | Bedienungsanleitung | Mode d'emploi
Istruzioni per l'uso | Instrucciones de manejo

Laboratory Vacuum Pumps | Labor Vakuumpumpen | Pompes à vide de laboratoire | Pompe a vuoto da laboratorio | Bombas de vacío para laboratorio

16694-2-50-22
16694-1-60-22

85032-535-32

English Page 3
Deutsch Seite 11
Français page 19
Italiano pagina 27
Español página 35

DOMINIQUE DUTSCHER SAS



List of Contents

1.	Description, Operating Conditions	4	You have selected a high-quality Sartorius product; the following tips will help you operate it safely, and reliably over a long period of time. Carefully study operating instructions before using the pumps and observe at all times the relevant instructions to avoid dangerous situations.
2.	Safety	4	
3.	Installation	4	
4.	Operation	5	
5.	Shutdown	5	
6.	Servicing	5	
7.	Trouble Shooting	7	
8.	Replacement Parts	7	
9.	Specifications	8	

Symbols

- ① Position in the illustration
- Important point
- ❶ Task
- Advise to users
- ⚠ Warning

Sartorius reserve the right to make changes

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1. Description, Operating Conditions

The pump transfer, and evacuate
100 % oil-free.

1.1 Electrical Equipment

- See section 9 for full electrical data.
- The protection class of standard versions is IP 44.
- The pumps are fitted with a thermal – switch to protect against overloading.

1.2 Operating Conditions

- Handling air, gases, and vapours at temperatures between +5°C ... +40°C.
- For maximum permissible operating pressure, ultimate vacuum, and flow capacity: see chapter 9.

 The pumps must not be used in areas where there is a danger of explosion.

- Before pumping a medium, the compatibility of materials of pump head, diaphragm and valves with the medium must be checked (for pump materials: see chapter 9).
- The pumps must not be used for liquids. You will find suitable liquid pumps in our product programme.

If your potential application lies outside the above limits discuss it with our technical adviser (see last page for contact telephone number).

1.3 Ambient Condition

When the pumps are operating the following ambient conditions must be maintained:

- Ambient temperature during operation: between +5°C ... +40°C.
- During operation an adequate supply of air for cooling must be provided.
- The pumps must not be used in areas where there is a danger of explosion.

1.4 Pump materials

See chapter 9.

2. Safety

- Note that the pumps may only be used for their intended purpose (see section 1).

The pumps must not be used in areas where there is a danger of explosion.

- The discharge of gas on the pressure side of the pump must be safely disposed.

Components connected to the pumps must be designed to withstand the pneumatic performance of the pumps (see chapter 9).

- Plug the pump only into properly installed grounded outlets.

When the operation of the pump is interrupted by the thermal switch, it will re-start automatically after cooling down. Take all care necessary to prevent this leading to a dangerous situation.

Specific safety instructions and measures for the media being handled must be observed.

- Parts of the casing marked with the sign below are only allowed to be opened after pulling out the plug (disconnecting the power source).

- Use only original Sartorius replacement parts.
- The pumps correspond to IEC 664:
 - the overvoltage category II
 - the pollution degree 2.

3. Installation

- Choose a safe location (flat surface) for the pump.

- Install the pump so as to ensure adequate flow of air cooling.

- Fit the pump at the highest point in the system, so that condensate cannot collect in the head of the pump – that prolongs working life.

- Connect the suction and pressure lines (tube ID 9). For flow direction see marking on the pump head.

- Arrange the suction and pressure lines so that condensate cannot run into the pump (sloping lines).

- Plug the pump only into properly installed grounded outlets.

- Compare the supply data with the electrical data of the pump. The voltage must not vary by more than +6%, and -10% from that shown on the type-plate.

4. Operation

- The pump can be switched on and off using the power switch.
-  Specific safety instructions for the media being handled must be observed.
- Before pumping a medium, the compatibility of materials of pump head, diaphragm, and valves with the medium must be checked (for pump materials: see chapter 9).
- The pumps must not start against pressure or vacuum. This also applies when the pump restarts after the power has been cut off for a short period. If a pump starts against pressure or vacuum, then the thermal switch will be activated and switch the pump off.

 When the operation of the pump is interrupted by the thermal switch, the pump will re-start automatically after cooling down. Take all care necessary to prevent this leading to a dangerous situation.

 The maximum permissible operating pressure (1 bar g) must not be exceeded.

- To prevent the maximum permissible operating pressure being exceeded, restriction or control of the air or gas flow should only be carried out in the suction line.
- If restriction or control of the air or gas flow is made on the pressure side ensure that the maximum permissible operating pressure is not exceeded.
- Diaphragm and valve plates | sealings are the only parts subject to wear. Wear is usually indicated by a drastic reduction in the pneumatic performance (vacuum, pressure, delivery). When replacing parts proceed as described in section 6.
- Ambient conditions: see chapter 1.3.

5. Shutdown

- ① Clean the containers.
- ② Purge the complete pump system including the connecting hoses with air for about 5 minutes.
- ③ Shutdown the pump power.
- ④ Disconnect the power source.

6. Servicing

 Before working on the pumps, disconnect it from the electricity supply by pulling the plug out.

When structured diaphragms are changed, valve plates | sealings should also be replaced. If the structured diaphragms are not changed in both heads at the same time or structured diaphragms and valve plates | sealings are not changed at the same time the nominal performance of the pump is not guaranteed after the service.

-  If a pump has been used for toxic substances or other types of substances which are hazardous, hazardous to health, or injurious, the following points must be observed:
- 1.) Clean the pump and its components before servicing.
 - 2.) Ensure that the service personnel is not subject to a health hazard. Apply the safety and protection measures that are necessary for the medium that has been handled by the pump (example: the use of protective gloves).
 - 3.) Ensure that discarded parts and materials are safely and correctly disposed of.
 - Use only original Sartorius replacement parts.

Required tools and material

- Service Set (see section 8)
- Philips-head screwdriver No. 2

Preparatory Steps

- ① Shut down the pump (see section 5) including disconnecting from the power source.
- ② Remove tubing from the inlet and outlet connectors of the pump.
- ③ If the pump is integrated in a vacuum system: Remove pump from the Baseplate.

Removing the Pump Head

- ① Loosen the 6 outer screws ① on each pump head.
- ② Carefully remove the pump head.

Changing the Diaphragm

- ① Carefully unscrew the structured diaphragm ⑥ anti-clockwise using both hands.
- ② Place the spacers ⑦ and ⑧ onto the screw thread of the new structured diaphragm (same number and order).
- ③ Screw in the new structured diaphragm ⑥ using both hands (do not overtighten).
- ④ Press the lip on the edge of the diaphragm into the groove in the housing.

Changing the Valve Plates | Sealings

- ① Carefully remove the top plate ② exposing the valve plates | sealings ③.
- ② Carefully remove old valve plates | sealings ③.
- ③ Insert new valve plates | sealings ③.

Refitting the Pump Head

- ① Position the top plate ② onto intermediate plate ④ according to guide pin ⑤; securely tighten screws ① in a crosswise pattern.
- ② Remount the pump to the Baseplate (if applicable).
- ③ Reconnect tubing.
- If the pump does not reach the desired vacuum after changing diaphragms and | or valve plates | sealings:
 - ① Check whether the spacers ⑦ and ⑧ have been replaced onto the diaphragm screw thread.
 - ② Check the interconnecting pipework connection between both pump heads as well as the tubing for leaks.
 - ③ Possibly the screws on one of the pump heads (or both heads) are insufficiently tightened (carefully tighten them again clockwise).
 - If you have any questions about servicing call our technical adviser (see last page for contact telephone number).


Fig. 3: Pump head (exploded drawing)

Specification

- ① Screw
- ② Top plate
- ③ Valve plate | sealing
- ④ Intermediate plate
- ⑤ Pin guide
- ⑥ Structured diaphragm
- ⑦ Spacer (thick)
- ⑧ Spacer (thin)
- ⑨ Dampening ring
- ⑩ Adapter
- ⑪ Screw

7. Troubleshooting

Before working on the pump, disconnect it from the electricity supply by pulling the plug out.

The following tips for fault-finding are best employed in the sequence shown.

- Flow, pressure, or vacuum too low
 - Compare the actual performance with the figures in chapter 9.
 - There is pressure on the pressure side, and at the same time vacuum, or a pressure above atmospheric, on the suction side.
- The pump is not designed for this condition.
- The cross-section of pneumatic lines, or connected components is too small, or they are restricted.
- To measure the performance, disconnect the pump from the system (small diameter tubing or a valve can significantly affect performance).
 - There is a leak at a connector, in a line, or in the pump head.
 - Diaphragm or valve plates|sealings are worn, or dirt is in the head.
- Section 6 Servicing.

8. Replacement Parts

A Service Set contains all replacement parts needed for one complete service.

- Service Set contains:
1 diaphragm and 2 valve plates| sealings.

No.: 1ED----0055

- Pump is switched on, but does not run, the on | off-switch on the pump is not lit.
- Pump is not connected with the power source.
- No voltage in the power source.
- Fuse in the pump is defective (only authorized | qualified personnel should investigate this problem; disconnect the pump from the power source, the fuse is accessible after opening the marked cover at the bottom of the pump. Dimension of the fuse: see section 9).
- Pump is switched on, but does not run, the on | off-switch on the pump is lit
- The heat sensitive circuit breaker has cut in because of the pump overheating (disconnect the pump from the power source, let the pump cool down, investigate the reason for the overheating and irradiate it).
- If the pump or the system does not operate properly and you cannot find any of the above faults, send them to the Sartorius Service Department.
- In order for Sartorius to repair the pump, the customer must provide a statement on the media which were pumped and on pump cleaning. Please fill out the corresponding Sartorius Stedim Biotech form, and submit it together with the pump. A sample statement for copying can be found in the Appendix of these operating instructions.

9. Specifications

Maximum capacity:	22 l/min
Ultimate vacuum:	≤100 mbar abs
Maximum permissible operating pressure:	1 bar ⁰

Electrical Type		
Voltage [V]	115	230
Frequency [Hz]	60	50
Power Consumption [W]	145	130
Fuse (2x) [A]	3.5	2.0

Pump materials		
Pump head:	Aluminium alloy	
Structured diaphragm:	CR	
Valves:	NBR	

Appendix: Customer statement for repair order (sample statement for copying)

- In order for Sartorius to repair the pump, the customer must provide a statement on the media which were pumped and on pump cleaning. Please fill out the corresponding Sartorius form, and submit it together with the pump.

Statement | Certificate

We confirm that the pump model listed below (please specify)

Serial-No. (please specify)

was used to pump the following media:

and that the pump listed above was cleaned. There are no poisonous, aggressive, biological, radioactive or other dangerous media in the pump.

Company _____

Date | Signature _____

9 Information and Instructions on Disposal

Attention must be paid to the applicable laws, rules and regulations when disposing of contaminated components.

If the packaging is no longer needed, it can be disposed of by local waste disposal authorities. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

arium® advance – including accessories – should not be disposed of as regular household waste.

EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and several other countries, Sartorius itself assumes responsibility for the return and conformant disposal of its electronic and

electrical products. Such equipment may not be thrown out with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Lab Instruments GmbH & Co. KG
Service Center
Weender Landstrasse 94–108
37075 Goettingen

WEEE reg. no. DE 89907997

In countries that are not members of the European Economic Area (EEA) or where no Sartorius subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and / or scrapping of the components, any batteries – if present – should be removed and disposed of in local collection boxes.

Non-hazard Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment | component arrives. An additional copy must accompany the equipment | component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilisable and | or repairable used electrical and electronic equipment and no longer utilisable and | or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by the supervisory authorities (used electrical and electronic equipment).

Model:

Serial No.:

I | we guarantee that the following measures have been taken:

- The equipment | component has been freed of hazardous materials to ensure that persons involved with the handling | repair are not exposed to any hazard or danger whatsoever.
- The equipment has been packaged for safety and the complete labeling has been affixed to it.
- The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.

The person sending in the equipment | component is hereby notified that said person shall be liable to Sartorius, its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius, its affiliates, its subsidiaries, dealers or distributors involved with the handling | repair of the equipment | component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.

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in the form as delivered fulfils all the relevant provisions of the following European Directives – including any amendments valid at the time this declaration was signed – and meets the applicable requirements of the harmonized European Standards listed below:

2004/108/CE 2004/108/CE	Compatibilità elettromagnetica Electromagnetic compatibility EN 61326-1:2013
2006/42/CE 2006/42/CE	Macchine Machines EN 60304-1: 2007
2011/65/EU 2011/65/EU	Restrizione dell'uso di determinate sostanze pericolose nelle apparecchiature elettriche ed elettroniche (RoHS) <i>Restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS)</i> EN 50581:2012

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This declaration certifies conformity with the above mentioned EC and EU Directives, but does not guarantee product attributes. Unauthorized product modifications make this declaration invalid. The safety information in the associated product documentation must be observed.

Sartorius Lab Instruments GmbH & Co. KG
International Certification Management
D-37079 Goettingen, Germania

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Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

Copyright by Sartorius, Goettingen,
Germany. All rights reserved.
No part of this publication may
be reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius. The status of the
information, specifications and
illustrations in this manual is
indicated by the date given below.
Sartorius reserves the right to
make changes to the technology,
features, specifications and design
of the equipment without notice.

Status:
Oktober 2015,
Sartorius Lab Instruments
GmbH & Co. KG,
Goettingen, Germany

Specifications subject to change without notice.
Copyright Sartorius Lab Instruments GmbH & Co. KG.
Printed in the EU on paper bleached without chlorine.
Publication No.: SL-6155-p151002
Ver. 10 | 2015