



sartorius

arium[®] mini – Ultrapure Water System
More Time for Your Research Thanks to
Unique Bagtank Technology



turning science into solutions

Bagtank Technology – Innovation for Efficient Research

Ultrapure water is the basis for every successful analysis. Yet analytical methods are becoming more and more sensitive and thus more susceptible to interference. Consistently high water quality is thus decisive for the reproducibility of results and to avoid time-consuming repeat analyses.

You can always rely on the arium® mini and its unique Bagtank system to deliver consistently high water quality that even exceeds ASTM Type 1 requirements. Because the results you achieve are more reliable, you save valuable time while increasing your efficiency.

As a pioneering and market leading manufacturer of disposable products for biopharmaceutical production, our success story continues: we have now integrated our innovative Bagtank technology into microdevices for daily laboratory use.

The closed Bagtank system reliably protects the pure water from secondary contamination by atmospheric components. Unlike conventional tank systems, the closed Bagtank system also prevents the formation of a permanent biofilm.

Depending on your needs, you can easily replace the bag (in just 5 minutes) and start always with a new basis every 6 months. This not only saves you immense effort and time spent cleaning and rinsing the system, you also avoid the use of chemicals.



Your arium® mini Benefits

- **Reliable**

Delivers consistently high water quality for reliable and reproducible results

- **Efficient**

Avoids tank cleaning thanks to its unique Bagtank technology

- **Intuitive**

Can be operated even with gloves via a multi-colored touch display

- **Compact**

With a footprint of only 28 cm in width, it fits into any laboratory environment



Your Requirements in the Focus – Comfortable Operation and Highest Quality Standards

In addition to the Bagtank, as our core feature we put a special focus on user-friendliness. With this in mind, we've equipped the arium® mini with a color touch screen featuring clear icons so that you can comfortably control the system – even when wearing lab gloves. This not only speeds up your workflows, it also makes them more secure thanks to active fault prevention – all in all, your daily lab work becomes easier and more efficient.

To ensure the best ultrapure water quality and the reliability of your process, we have also integrated the highest quality standards into all functional areas. For example, you can set a limit for conductivity. If this limit is exceeded, further dispensing of water will either be stopped or accompanied by a warning, as configured.

Integrated Illumination for Easier Visual Support

- Pulsating light during water dispensing
- Impulse signals accompany error messages
- Dimmed light in standby mode

arium® Scientific Pack for Highest Ultrapure Water Quality

- Highest cartridge capacity thanks to efficient filler materials
- Optimal flow over the resin mixed bed
- Designed for quality standard ASTM Type 1

arium® Sterile Filter for Effective Sterility

- Validated according to HIMA and ASTM F-838-05
- Meets quality standards in accordance with USP
- Safely removes microorganisms using 0.2-µm membrane




Current bag filling volume displayed

Color touch display with easy-to-clean surface

Product water quality and temperature displayed

Menu | Home function for settings and care programs



Manual dispensing

Volume-controlled dispensing in 50-ml increments

Last volume dispensed is saved by "Favorites" function



Intuitive Menu Navigation for Optimal User Comfort

- Icon-supported interface with touch display
- Display of current measured values and messages
- Favorites function to dispense recurring volumes

arium® UV Lamp for Critical Analyses (185 | 254 nm*)

- Optimal temperature gradient thanks to horizontal installation
- Degradation of organic components (TOC) by irradiation (< 5 ppb)
- Inactivates microorganisms by damaging their DNA

arium® Pre-treatment Cartridge for Effective Purification

- Effective adsorption using high-quality activated carbon
- Special catalyst removes oxidants
- Containment retained by reverse osmosis

Integrated Bagtank technology – A revolution in Laboratory Water Purification

arium® mini is the only ultrapure water system with integrated Bagtank technology. The 5-liter bag developed for the arium® mini is made of the proprietary material S71. Originally developed for the pharmaceutical industry, this multilayer, specialized film conforms to the biocompatibility standards of ISO and USP. Its quality and composition allows you to optimally store your pre-treated pure water for the production of Type 1 ultrapure water.

The bag is laterally integrated into the system. This enables a footprint-reducing and compact device width of only 28 cm. The system can thus be integrated into any laboratory environment almost independently of the available space. The place of installation can also be chosen flexibly.

The integrated Aqua Stop system lets you enjoy active leakage protection while ensuring maximum safety in your laboratory.




Intelligent Technology – Systematic Ultrapure Water Purification

With intelligent technology the arium® mini permanently protects against secondary contamination and consistently ensures the highest water quality.

The closed Bagtank system automatically regulates pressure equalization so that no CO₂ can penetrate. This additionally prevents ion contamination and maintains consistently low conductivity. As an added benefit, you profit from particularly high cartridge capacity because there is less demand placed on the ion exchangers.

Moreover, you can prevent the formation of a permanent biofilm and thereby additional TOC loading. Depending on your requirements, you replace the bag every 6 months and simply start with a new basis.



5-liter Bag for Minimal Maintenance

- Easy bag replacement in less than 5 minutes
- No use of cleaning agents needed
- Practical and safe connection using quick connector

Innovative Technology Prevents Secondary Contamination

- Optimum water quality by preventing permanent biofilm
- Closed system guarantees purity
- Stable quality thanks to the special high-gas barrier multilayer S71 film

Sophisticated Power Backup Functions

- Integrated pump for increasing water pressure
- Filling arium® mini with pre-treated water: 2 minutes
- Filling arium® mini plus with RO permeate: 8 liters per hour

Intuitive Menu Navigation for Optimal User Comfort

- Icon-supported interface with touch display
- Display of current measured values and messages
- Favorites function to dispense recurring volumes

arium® UV Lamp for Critical Analyses (185 | 254 nm*)

- Optimal temperature gradient thanks to horizontal installation
- Degradation of organic components (TOC) by irradiation (< 5 ppb)
- Inactivates microorganisms by damaging their DNA

arium® Pre-treatment Cartridge for Effective Purification

- Effective adsorption using high-quality activated carbon
- Special catalyst removes oxidants
- Containment retained by reverse osmosis



The Right arium® mini Version for Your Everyday Applications

arium® mini

Standard Version – Without Direct Water Connection

The device is fed via the integrated bag, which is filled by means of a built-in pump. Pre-treated water is automatically transferred from a vessel in the bag and thus serves as feed water for the production of Type 1 ultrapure water.

arium® mini plus

Advanced Version – with Direct Connection to the Drinking Water Supply Line

Combined with a pre-treatment stage, the device is connected directly to the potable water supply. The processed water is then produced directly into the bag. The final production of Type 1 ultrapure water is then carried out with the specially developed arium® Scientific purification cartridge.

Water Quality	arium® mini	arium® mini plus
Type 1 Water	■	■
Type 3 Water		■

Applications	arium® mini mini plus	arium® mini UV mini plus UV
Buffer and media preparation	■	■
Histology	■	■
Preparation of reagents, blank samples etc.	■	■
ELISA (enzyme-linked immunosorbent assay)	■	■
AAS GF-AAS (atomic absorption spectrometry)	■	■
Photometry	■	■
SPE (solid phase extraction)		■
IC (ion chromatography)		■
LC-MS (liquid chromatography–mass spectrometry)		■
GS-MS (gas chromatography coupled with mass spectrometry)		■
HPLC (high-performance liquid chromatography)		■

Services available for your arium® mini

Optimize the performance of your ultrapure water system and look forward to many years of trouble-free operation by availing yourself of the professional installation, extended warranty and preventive maintenance services we offer.

Our service specialists undergo continuous and thorough training in order to guarantee you'll always receive the best possible on-site service.

- Comprehensive feedwater analysis and corresponding device configuration
- Installation and startup for quick and safe usability
- Maximum safety by extended warranty
- Extended service life through regular checks tailored to your needs



Sales and Service Contacts

For further contacts, visit www.sartorius.com

Europe

Germany
Sartorius Lab Instruments
GmbH & Co. KG
Weender Landstrasse 94–108
37075 Goettingen
Phone +49.551.308.0
Fax +49.551.308.3289

France Et Suisse Romande
Sartorius France
2, rue Antoine Laurent de Lavoisier
ZA de la Gaudriée
91410 Dourdan
Phone +33.1.70.62.50.00
Fax +33.1.64.59.76.39

Austria
Sartorius Austria GmbH
Modecenterstrasse 22
1030 Vienna
Phone +43.1.7965760.0
Fax +43.1.7965760.24

Belgium
Sartorius Belgium N.V.
Rue Colonel Bourg 105
1030 Bruxelles
Phone +32.2.756.06.90
Fax +32.2.481.84.11

Finland Et Baltics
Sartorius Biohit Liquid Handling Oy
Laippatie 1
00880 Helsinki
Phone +358.9.755.951
Fax +358.9.755.95.200

Hungary
Sartorius Hungária Kft.
Kagyló u. 5.
2092 Budakeszi
Phone +3623.457.227
Fax +3623.457.147

Ireland
Sartorius Ireland Ltd.
Unit 41, The Business Centre
Stadium Business Park
Ballycoolin Road
Dublin 11
Phone +353.1.8089050
Fax +353.1.8089388

Italy
Sartorius Italy S.r.l.
Viale A. Casati, 4
20835 Muggiò (MB)
Phone +39.039.4659.1
Fax +39.039.4659.88

Netherlands
Sartorius Netherlands B.V.
Phone +31.30.60.53.001
Fax +31.30.60.52.917
info.netherlands@sartorius.com

Poland
Sartorius Poland sp.z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn
Phone +48.61.6473830
Fax +48.61.6473839

Russian Federation
LLC "Sartorius RUS"
Uralskaya str. 4, Lit. B
199155 St. Petersburg
Phone +7.812.327.53.27
Fax +7.812.327.53.23

Spain Et Portugal
Sartorius Spain, S.A.
Avda. de la Industria, 32
Edificio PAYMA
28108 Alcobendas (Madrid)
Phone Spain +34.902.123.367
Phone Portugal +351.800.855.800
Fax Spain +34.91.358.96.23
Fax Portugal +351.800.855.799

Switzerland
Sartorius Mechatronics Switzerland AG
Ringstrasse 24a
8317 Tagelswangen (ZH)
Phone +41.44.746.50.00
Fax +41.44.746.50.50

U.K.
Sartorius UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9QQ
Phone +44.1372.737159
Fax +44.1372.726171

Ukraine
LLC "Biohit"
Post Box 440 "B"
01001 Kiev, Ukraine
Phone +380.44.411.4918
Fax +380.50.623.3162

Americas

USA
Sartorius Corporation
5 Orville Drive, Suite 200
Bohemia, NY 11716
Phone +1.631.254.4249
Toll-free +1.800.635.2906
Fax +1.631.254.4253

Argentina
Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires
Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil
Sartorius do Brasil Ltda
Avenida Senador Vergueiro 2962
São Bernardo do Campo
CEP 09600-000 – SP– Brasil
Phone +55.11.4362.8900
Fax +55.11.4362.8901

Canada
Sartorius Canada Inc.
2179 Dunwin Drive #4
Mississauga, ON L5L 1X2
Phone +1.905.569.7977
Toll-Free +1.800.668.4234
Fax +1.905.569.7021

Mexico
Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
Ciudad Satélite
53100, Estado de México
México
Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia
Sartorius Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175
Phone +61.3.8762.1800
Fax +61.3.8762.1828

China
Sartorius (Shanghai) Trading Co., Ltd.
3rd Floor, North Wing, Tower 1
No. 4560 Jinke Road
Zhangjiang Hi-Tech Park
Pudong District
Shanghai 201210, P.R. China
Phone +86.21.6878.2300
Fax +86.21.6878.2882

Hong Kong
Sartorius Hong Kong Ltd.
Unit 1012, Lu Plaza
2 Wing Yip Street
Kwun Tong
Kowloon, Hong Kong
Phone +852.2774.2678
Fax +852.2766.3526

India
Sartorius Weighing India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra,
Nelamangala Tq
562 123 Bangalore, India
Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan
Sartorius Japan K.K.
4th Fl., Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku, Tokyo, 140-0001 Japan
Phone +81.3.3740.5408
Fax +81.3.3740.5406

Malaysia
Sartorius Malaysia Sdn. Bhd
Lot L3-F-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia
Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore
Sartorius Singapore Pte. Ltd
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528
Phone +65.6872.3966
Fax +65.6778.2494

South Korea
Sartorius Korea Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, Bundang-Gu
SeongNam-Si, GyeongGi-Do, 463-400
Phone +82.31.622.5700
Fax +82.31.622.5799

Thailand
Sartorius (Thailand) Co. Ltd.
129 Rama 9 Road,
Huaykwang
Bangkok 10310
Phone +66.2643.8361-6
Fax +66.2643.8367



◀ www.sartorius.com