sartorius

Installations- und Bedienungsanleitung | Installation and Operating Instructions |  
Notice d'installation et d'utilisation | Instrucciones de instalación y manejo |  
Istruzioni per l'installazione e l'uso | 取扱説明書 | Инструкция по установке и эксплуатации |  
安裝和操作說明 | Instrukcja montażu i obsługi | Instruções de instalação e operação | 설치 및 사용 설명서

## Secura<sup>®</sup>, Quintix<sup>®</sup>

26, 35, 65, 125, 125D, 225D, 324

Laborwaagen | Laboratory Balances | Balances de laboratoire | Balanzas de laboratorio | Bilance da laboratorio |  
研究室用電子天びん | Лабораторные весы | 实验室电子天平 | Wagi laboratoryjne | Balanças de laboratório |  
실험실용 저울

1000038784

**English – page 3**

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DOMINIQUE DUTSCHER SAS

# Contents

<b>1 User Information</b> .....	3
1.1 Warning / Danger Symbols .....	3
1.2 Explanation of Symbols .....	3
1.3 About These Instructions .....	3
1.4 Advice on Application / Technical Support .....	3
<b>2 Safety Precautions</b> .....	3
2.1 Guidelines and General Information .....	3
2.2 Installation Instructions .....	4
2.3 Intended Use .....	4
<b>3 Getting Started</b> .....	4
3.1 Unpacking and Equipment Supplied .....	4
3.2 Installing the Balance .....	4
3.3 Choosing a Location .....	5
3.4 AC Power Supply .....	5
3.5 Leveling .....	6
3.6 Anti-theft Locking Device .....	6
<b>4 Transporting the Balance</b> .....	7
<b>5 Shipping the Balance</b> .....	7
<b>6 Care and Maintenance</b> .....	8
<b>7 Disposal</b> .....	8
<b>8 Technical Data</b> .....	9
Model-specific Data .....	81
CE Declaration for Electromagnetic Compatibility .....	83

## Other Documents

Apart from these installation instructions, an up-to-date user manual containing additional information about the device is also available as a PDF file, which can be downloaded from the Sartorius website: [www.sartorius.com](http://www.sartorius.com)

You need Acrobat Reader to read and print the PDF file, which you can download for free from the Adobe website ([www.adobe.com](http://www.adobe.com)).

## 1 User Information

### 1.1 Warning / Danger Symbols

These notes identify hazards which have a high probability of resulting in death or serious physical injury if not avoided.

These notes identify hazards that can result in moderate or mild injuries if not avoided.

These notes identify hazards associated with the risk of material damage.

## 1.2 Explanation of Symbols

This symbol identifies useful information and tips.

Note on use in legal metrology. For instruments the conformity assessment procedure has been performed (verified balances)

The following symbols are used in these instructions:

- ▶ Indicates required steps
- ▷ Describes what happens after you have performed a particular step

Perform steps in the specified order:

1. First action
  2. Second action
  3. ...
- Indicates an item in a list.

### Conventions for these Instructions:

The figures in these instructions are based on "standard" balances. On balances verified for use in legal metrology, some displays and reports may deviate slightly from the figures. Where this is significant for operation, the differences will be explained in the text.

### 1.3 About These Instructions

- ▶ Please read these instructions carefully and completely before putting the equipment into operation for the first time. Observe the safety instructions.
- ▶ These instructions are an important part of the product. Keep these instructions in a safe place. If you give the device to others to use, give them these instructions too.
- ▶ If these instructions are lost, please contact Sartorius for a replacement or download the latest manual from our website: [www.sartorius.com](http://www.sartorius.com)

### 1.4 Advice on Application / Technical Support

Contact addresses for application advice and our technical support can be found online at: <http://www.sartorius.com>

## 2 Safety Precautions

### 2.1 Guidelines and General Information

- The balance complies with EU Directives and standards for electrical safety and electromagnetic compatibility\*. Improper use or handling can, however, result in damage and / or injury. Any improper use or operation of the balance, i.e., that is not consistent with the instructions, will result in forfeiture of all claims under the manufacturer's warranty.
  - Personnel need to have read and understood these instructions, including the safety instructions.
  - In the event of use in systems and ambient conditions which have greater safety requirements, you must observe the requirements and provisions applicable in your country.
  - Always keep the equipment and balance freely accessible.
- \*See Chapter "Technical Specifications", 9.

#### Danger of explosion

Do not use this equipment in hazardous areas in which explosive materials are present.

Make sure that the voltage rating printed on the AC adaptor is identical to your local mains voltage.

## 2.2 Installation Instructions

Do not operate the balance if its housing or AC adaptor is damaged. Immediately disconnect the damaged device from the power by pulling the plug.

Do not expose the balance, its AC adaptor or the accessories supplied by Sartorius to extreme temperatures, aggressive chemical vapors, moisture, shocks, vibrations or strong electromagnetic fields. Observe the conditions of operation described in the "Technical Data"!

The operator shall be solely responsible for any modifications to the equipment and for connecting any cables or equipment not supplied by Sartorius! Information on operational quality is available upon request from Sartorius. You should only use peripherals supplied by Sartorius!

Note the IP protection class of the balance and its AC adaptor! Do not allow liquid penetration. The protection class specifies the suitability of equipment for various environmental conditions (moisture, foreign bodies).

Before cleaning the AC adaptor or the balance: Unplug from the mains power.

The balance may only be opened by specialized personnel trained by Sartorius. Do not open the AC adaptor.

### Warnings at the Place of Installation and During Operation:

If glass breaks, there is a risk of injury posed by cuts on glass edges.

Lay the cables where they pose no risk of causing someone to trip.

Observe the additional safety and danger information in the following chapters.

## 2.3 Intended Use

This high-resolution balance is exclusively for use in laboratories and indoors under normal atmospheric conditions. It was developed specifically for the exact determination of the mass of materials in liquid, paste, powder or solid form.

Appropriate containers must be used for each type of sample material.

## 3 Getting Started

The balance must be disconnected from the mains power for all assembly work.

### 3.1 Unpacking and Equipment Supplied

- ▶ Open the packaging, making sure to remove of all parts carefully.
- ▶ After unpacking the device, check it immediately for any external damage.
- ▷ If you detect any damage, proceed as directed in the "Care and Maintenance" chapter.
- ▷ Save all parts of the original packaging for any future transportation. During shipment, please do not leave cables plugged in!

The following parts are included in the equipment supplied:

Components Model	Secura®		Quintix®		
	26	225D 125	324	125D 65	35
Draft shield (4 panels)	x	x	x	x	x
Collection container	x	x	x	x	x
Weighing pan: Ø 80 mm		x	–	x	x
Weighing pan: Ø 90 mm	–	*	x	*	*
Weighing pan: Ø 50 mm	x	–	–	–	–
Shield disk (normal/high)	x	–	–	–	–
AC adaptor with country-specific power plug	x	x	x	x	x
In-use dust cover	x	x	x	x	x
Dust cover	x	x	x	x	x
Installation and Operating Instructions	x	x	x	x	x
Application guide	x	x	x	x	x

\* Optional

### 3.2 Installing the Balance


#### Collection Container Installation

- ▶ Insert the collection container into the balance.


#### Assembling the Panels


- ▶ Press against the front of the cover from both sides and fold it upwards.


- ▶ Inserting the front panel of the draft shield:
  - Insert the panel on the top into the guide rails of the draft shield.
  - Press the panel gently upwards, while simultaneously pressing downwards into the slot of the draft shield.
  - The panel must latch into the slot.

- ▶ Inserting the side panels:
  - Mount the guide carriages (1) of the side panel by inserting them into the guide rails on the top.
  - Insert the side panel at the bottom behind the guide rod (2) of the draft shield.


- ▶ Make sure that:
  - the guide carriage (1) is properly and fully positioned on the guide rail.
  - the side panels are positioned behind the guide rod (2).
  - the panels can be moved smoothly and freely.

- ▶ Inserting the upper panel:
  - Insert the top panel on the top of the draft shield guide.
  - Push the panel carefully backwards.
  - Gently lift the panel by its grip in order to push it below the guide rod (3) of the guide.
- ▶ Close the cover by folding it downwards completely until it locks.

- ▶ Insert the weighing pan into the opening of the collection container.
- ▶ Push the upper and side panels forward to close them.

- Secura® 26**
- ▶ Position one of the two shield disks (1 or 2) on the collection container.
  - ▶ Insert the weighing pan (3) into the opening of the collection container.
  - ▶ Push the upper and side panels forward to close them.

### 3.3 Choosing a Location

Select the right setup location:

- Set up the device on a stable, even surface that is not exposed to vibrations (e.g., a weighing stone).
- Maintain free access to the device at all times.

Choose a location that is not subject to the following negative influences:

- Heat (heater or direct sunlight)
- Drafts from open windows, AC systems and doors
- Vibrations during weighing
- Heavy traffic areas (personnel)
- Excessive moisture
- Electromagnetic fields
- Extremely dry air

### Acclimatization

Condensation from humidity can form on the surfaces of a cold device when it is brought into a substantially warmer area. You should therefore let a device that has been disconnected from its power source acclimatize for approximately 2 hours before reconnecting it to the power.

### 3.4 AC Power Supply


#### AC Adaptor Assembly


**Fatal electric shocks can be caused by use of the incorrect power plug adapter or improper use of the power plug adapter.**

- ▶ Attach the country-specific power plug adapter to the AC adapter. The power plug adapter must be suitable for the wall outlet at the installation location.
- ▶ Do not insert the power plug adapter into the socket without an AC adapter.

Item number on packaging	Power supply/country-specific power plug adapter (packed in PE bag with printed country code, e.g. EU)	Illustration (from left to right)
YEPS01-PS4	Power supply with connection cable	
	China (CN) South Africa (ZA) Australia (AU)	
	India (IN) USA and Japan (US+JP) Europe (EU)	
	United Kingdom (UK)	
YEPS01-PS5	Power supply with connection cable	
	Argentina (AR) Brazil (BR) Korea (KR)	

- ▶ Select the country-specific power plug adapter. The power plug adapter must be suitable for the wall outlet at the installation location.
- ▶ Push the power plug adapter into the AC adapter's holder. The ribbed button must be facing forward.
- ▶ Push the power plug adapter all the way in until it audibly engages.
- ▶ Check that the power plug adapter is firmly locked in place. To do this, gently pull on the power plug adapter.
- ▶ If the power plug adapter cannot be moved then it is locked in place.


#### Removing the power plug adapter

- ▶ Press on the ribbed button from above while sliding the power plug adapter backward.
- ▶ Push the power plug adapter out of the AC adapter and remove it.

#### Connecting the AC Adapter

- ▶ Check the voltage rating on the AC adapter's type plate. Make sure that the voltage rating printed on this unit matches the local supply voltage at the place of installation.
- ▶ If the stated supply voltage does not comply with the local supply voltage or there is no suitable AC adapter available: Do not use the AC adapter. Contact Sartorius Service.
- ▶ Only use original Sartorius AC adapters.

- ▶ Connect the angle plug (2) to the device's power supply socket (operating voltage connection).
- ▶ Connect the power cable (1) to the wall outlet (supply voltage) at the installation location.

#### Warm-up Time

To ensure accurate results are delivered, the balance must warm up for at least 30 minutes after initial connection to the power supply. Only after this time will the device have reached the required operating temperature.

**M**

When a verified balance of the accuracy class  $\text{I}$  for use in legal metrology is connected to the mains power, it must warm up for at least one hour before operation.

### 3.5 Leveling

**i**

Leveling the balance compensates for slant or unevenness at the place of installation. The balance must be perfectly horizontal to ensure consistent, reproducible weighing results. The balance needs to be re-leveled and then adjusted each time its setup location is changed.


#### Secura Balances:

When "Level" appears on the display, you must level the balance.

- ▶ Make sure that the weighing pan is unloaded.
- ▶ Press "LEVEL".
- ▶ Follow the instructions on the interactive display.

#### Quintix® Balances:

- ▶ Turn the two feet as shown in the illustration until the air bubble is centered within the circle of the level indicator.
  - Air bubble at "12 o'clock": turn the two feet clockwise.
  - Air bubble at "3 o'clock": turn the left foot clockwise and the right foot counterclockwise.
  - Air bubble at "6 o'clock": turn the two feet counterclockwise.
  - Air bubble at "9 o'clock": turn the left foot counterclockwise and the right foot clockwise.

#### Instructions for Use in Legal Metrology

Additional Display Elements:

Metrology line with additional information about minimum load (Min) and verification scale interval (e).

Warning symbol to indicate calculated values, negative gross values or values that are smaller than the minimum sample quantity (SQ<sub>min</sub>).

**M**


**Setup Menu "Settings":**

On balances verified for use in legal metrology, not all functions and settings are available.


#### Interface Specification (Data Output)

Settings "PC – SBI" and "USB RS232 SBI" for the identification of differentiated scale divisions for use in legal metrology: In the setting "SBI", there is no automatic identification of display digits not verified for use in legal metrology. Corresponding measures or settings must be carried out on the peripheral device.

### 3.6 Anti-theft Locking Device

- ▶ A "Kensington" key lock can be installed at the fastening point on the back of the balance if required.

- ▶ The balance can also be secured at the installation location with a chain and lock, for example.

## 4 Transporting the Balance

Before moving the balance, switch it to the standby mode, then unplug the power connector and remove the weighing pan from the balance.

### Moving the Balance over Short Distances in the Lab

Avoid glass breakage, shocks and vibrations: Never lift or carry the balance by its draft shield!

- ▶ Hold the balance under the housing, lift it up carefully and carry it to its new location.

The balance needs to be re-leveled and then adjusted each time its setup location is changed.

## 5 Shipping the Balance


Prior to shipping, switch the balance to standby mode (see 8) and then pull the plug.

Please use the original packaging for shipping. Packaging can be ordered through the Sartorius Service Center if required.


### Detaching the Draft Shield / Weighing Pan

**IMPORTANT**

#### Danger of damage to the side panels!


The guide carriages of the side panels are provided with fine steel filaments. Removal of these steel filaments or carrying out any other modifications to the guide carriages is forbidden.

- ▶ Press against the front of the cover from both sides and fold it upwards.
- ▶ Push the upper panel forwards and pull it out of the draft shield.
- ▶ Push the side panels slightly backwards and remove them from the draft shield.
- ▶ Pull the front panel over the groove of the draft shield and remove the panel.
- ▶ Pull the front panel over the groove on the bottom of the draft shield and remove the panel.


- ▶ Pull the weighing pan out of the collection container.
- ▶ Remove the collection container.

#### Secura® 26

- ▶ Pull the weighing pan out of the collection container.
- ▶ Take the shield disk off the collection container.
- ▶ Lift out the collection container.


### Packing the Balance

- ▶ Place the balance in the bottom part of the packaging.

- ▶ Place the remaining components of the balance into the recesses provided in the packaging part:
  - Collection container
  - Weighing pan
  - Right and left side panels
  - Top panel
  - Front panel
  - Shield disks (Secura® 26 only)

- ▶ Seal the packaging part with the components and place in the recess provided.

- ▶ Put the top part of the packaging over the device parts.
- ▶ Next, put the packaged balance into the cardboard box and seal the box.

## 6 Care and Maintenance

### Service

To ensure the continued accuracy of your balance, we recommend scheduling regular servicing at least once a year. Sartorius Service offers different service contracts with maintenance intervals that are tailored to your needs. A calibration certificate should always be issued as part of every maintenance session. Safety inspections of the AC adaptor and its connections must be performed at appropriate intervals by a qualified electrician (e.g. every 2 years).

### Repairs

Repair work must only be carried out by trained service technicians. The device must be unplugged during repair work. Unplug from the mains power. In addition, the measurement accuracy of your balance may be affected and serious risks to the user may occur. Contact Sartorius Service or a Sartorius dealer for proper repairs.

### Cleaning the Balance

#### Cleaning the Control Panel

Switch the display to standby mode to avoid modifying the settings for operation during cleaning.

- ▶ Touch the Menu key to switch to the application selection.

- ▶ If you then touch , the display goes off.

- ▶ To switch the display back on again: Touch  at the bottom left of the display.

#### Cleaning the Balance Housing

#### Contaminated equipment:

Health risk posed by product contamination due to product deposits and collection of residue with microbial contamination.

Health risk posed by biological or microbiological substances. Observe cleaning specifications. Examine the cleaning results closely.

#### Disconnect from the power supply:

Disconnect the AC adaptor from the power supply. If necessary, disconnect the data cable connected to the balance.

Never open the balance or the AC adaptor. The parts contained in these cannot be cleaned, repaired or replaced by the operator.

- Make sure that no liquid or dust gets into the balance or the AC adaptor.
- Never use cleaning agents that contain solvents or abrasive ingredients, which can ultimately damage the equipment.


Do not clean the following parts with acetone or aggressive cleaning agents:  
Control panel, mains plug, data interface, labels, verified seal and all other plastic parts.

#### Detaching the Draft Shield / Weighing Pan

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If required, remove the rear wall for cleaning.

- ▶ Press the locking mechanism above the rear wall upwards and remove the rear wall.

- ▶ Clean the housing / the panels and then dry the balance with a soft cloth.
- ▶ After cleaning put together the components in reverse order:
  - Reinsert the rear wall. Ensure that this engages in the locking mechanism above the rear wall.
  - Observe the notes in the Chapter “Getting Started – Installing the Balance” 4.

## 7 Disposal

The packaging is made from environmentally-friendly materials that can be used as secondary raw materials. Packaging that is no longer needed can be disposed of in Germany free of charge via the VfW dual system (contract number D-59101-2009-1129).

Otherwise, you dispose of the material at the local waste management facility according to the applicable regulations. The device, including accessories and batteries, must not be disposed of in general household waste and must be recycled similar to electrical and electronic devices. For disposal or recycling, please contact our service staff on-site. In addition, the partners listed on the following website can be used for disposals within the EU:

1. Select <http://www.sartorius.com>.
2. Click “Service”.
3. Select “Information on Disposal”.
4. Addresses for local Sartorius disposal contacts can be found in the PDF files given on this webpage.

Sartorius will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal.

#### Service address for disposal:

Please refer to our website ([www.sartorius.com](http://www.sartorius.com)) or contact the Sartorius Service Center for more detailed information regarding repair service addresses or the disposal of your device.

## 8 Technical Data

### General Data

(Model-specific data see page 81)

Information	Unit	Value
<b>AC adapter</b>		
Sartorius AC Adapter		YEPS01-PS4 or YEPS01-PS5 with replaceable country-specific power plug adapters
Primary		100 – 240 V~, ± 10 %, 50 – 60 Hz ± 5 %, 0.2 A
Secondary		15 V <sub>DC</sub> , ± 5%, 530 mA (max.) / 8 Watt (max.): 0°C to +40°C 15 V <sub>DC</sub> , ± 5%, 330 mA (max.) / 5 Watt (max.): +40°C to +50°C
Further details		Protection class <b>II</b> in accordance with EN/IEC 60950-1   up to 3000 m above sea level IP40 in accordance with EN 60529/IEC 60529
<b>Balance</b>		
Supply voltage		Only via Sartorius AC adapter YEPS01-PS4 or YEPS01-PS5 with replaceable country-specific power plug adapters
Power consumption	W	4.5 (typically)
Further details		IP43 in accordance with EN 60529/IEC 60529
<b>Ambient conditions</b>		
The technical data apply under the following ambient conditions:		
Environment		For indoor use only
Ambient temperature*	°C	+10 to +30
Operational capability	°C	Guaranteed between +5 and +45.
Storage and shipping	°C	-10 to +60
Height above sea level (NN)	m	Up to 3000
Relative humidity**		15 – 80% for temperatures up to 31°C non-condensing, decreasing linearly up to 50% relative humidity at 40°C and 20% at 50°C.
<b>Safety of Electrical Equipment</b>		
As per EN 61010-1 / IEC61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements		
<b>Electromagnetic Compatibility</b>		
As per EN 61326-1 / IEC61326-1 Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements		
Defined immunity to interference		Suitable for use in industrial areas
Interference emission		Class B Suitable for use in residential areas and areas that are connected to a low voltage network that also supplies residential buildings.
<b>Standard Equipment</b>		
Selectable weight units <sup>1</sup>		Gram, kilogram, carat, pound, ounce, troy ounce, Hong Kong teal, Singapore teal, Taiwan teal, grain, pennyweight, milligram, parts per pound, China teal, mommas, Austrian carat, tola, baht, mesghal and Newton
Available application programs		<b>Secura® Models</b> Counting, weighing in percent, density determination, peak hold, unstable condition, checkweighing, mixing, components (totalizing), statistics, conversion, S0min and identifiers <b>Quintix® Models</b> Counting, weighing in percent, density determination, peak hold, unstable condition, checkweighing, mixing, components (totalizing), statistics, conversion

Verified balances in accordance with EU requirements comply with the requirements of Council Directive 2014/31/EU with EN 45501:2015 and OIML R76:2006.

\* For verified balances in accordance with EU requirements, refer to the information on the balance.

\*\* For verified balances in accordance with EU requirements, the legal regulations apply.

<sup>1</sup> Depending on the country-specific model version, not all weight units listed may be available.

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### FCC Rules

#### Safety Precautions:

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by

turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver
- Consult the dealer or an experienced radio/TV technician for help.

#### Shielded Cables:

Connections between the devices and peripherals must be made using shielded cables in order to maintain compliance with FCC radio frequency emission limits.

#### Modifications:

Any modifications made to this devices that are not approved by Sartorius Weighing Technology GmbH may void the authority granted to the user by the FCC to operate this equipment.

<sup>1)</sup> Possible terms for country-specific models:  
 x = S: Standard balances without country specific additions  
 x = SAR: Standard balances with country specific additions for Argentina only  
 x = SJP: Standard balances with country specific additions for Japan only  
 x = SKR: Standard balances with country specific additions for South Korea only

\* Upper edge of the weighing pan to the lower edge of the upper draft shield panel

## Technical Specifications

### Model-specific Data

Verified Models with Country-specific Type Approval Certificate

		Secura <sup>®</sup>			
Model		26-1x <sup>2)</sup>	225D-1x <sup>2)</sup>	125-1x <sup>2)</sup>	324-1x <sup>2)</sup>
Accuracy class		Ⓡ	Ⓡ	Ⓡ	Ⓡ
Type <sup>3)</sup>		SQP-H	SQP-F	SQP-F	SQP-G
Max	g	21	120 / 220	120	320
Scale interval d	g	0.000002	0.00001 / 0.0001	0.00001	0.0001
Verification scale interval e	g	0.001	0.001	0.001	0.001
Min	g	0.001	0.001	0.001	0.01
Min (only for Models ...-10IN)	g	–	0.1	0.1	0.1
Tare (subtractive)		< 100% of max. weighing capacity			
Minimum sample weight as per USP, typical	mg	4	40	40	160
Typical stabilization time	s	8	6 / 2	6	2
isoCAL:					
– Temperature change	K	1.5	1.5	1.5	1.5
– Time interval	h	4	4	4	4
Display result (depending on the set filter level)	s	0.2 / 0.4	0.2 / 0.4	0.2 / 0.4	0.2 / 0.2
Weighing pan size	mm	∅ 50	∅ 80 (optional ∅ 90)	∅ 80 (optional ∅ 90)	∅ 90
Weighing chamber height*	mm	218	218	218	218
Net weight, approx.	kg	8.0	7.8	7.8	7.9

		Quintix <sup>®</sup>		
Model		125D-1x <sup>2)</sup>	65-1x <sup>2)</sup>	35-1x <sup>2)</sup>
Accuracy class		Ⓡ	Ⓡ	Ⓡ
Type <sup>3)</sup>		SQP-F	SQP-F	SQP-F
Max	g	60 / 120	60	30
Scale interval d	g	0.00001 / 0.0001	0.00001	0.00001
Verification scale interval e	g	0.001	0.001	0.001
Min	g	0.001	0.001	0.001
Min (only for Models ...-10IN)	g	0.1	0.1	0.1
Tare (subtractive)		< 100% of max. weighing capacity		
Typical stabilization time	s	6 / 2	6	6
isoCAL:				
– Temperature change	K	1.5	1.5	1.5
– Time interval	h	4	4	4
Display result (depending on the set filter level)	s	0.2 / 0.4	0.2 / 0.4	0.2 / 0.4
Weighing pan size	mm	∅ 80 (optional ∅ 90)	∅ 80 (optional ∅ 90)	∅ 80 (optional ∅ 90)
Weighing chamber height*	mm	218	218	218
Net weight, approx.	kg	7.8	7.8	7.8

<sup>2)</sup> Possible terms for country-specific models:  
**x = CEU:** Verified balances with EC Type Approval Certificate D12-09-014 without country specific additions  
**x = CFR:** Verified balances with EC Type Approval Certificate D12-09-014 for France only  
**x = CIT:** Verified balances with EC Type Approval Certificate D12-09-014 for Italy only  
**x = CCH:** Verified balances with EC Type Approval Certificate D12-09-014 for Switzerland only  
**x = CN:** CMC Type Approval Certificate for China

**x = OJP:** Balance with Type Approval Certificate for Japan  
**x = OBR:** Balance with Type Approval Certificate for Brazil  
**x = ORU:** Balance with Type Approval Certificate for Russia  
**x = OIN:** Balance with Type Approval Certificate for India

<sup>3)</sup> All models with "...CN": type "SQP"

\* Upper edge of the weighing pan to the lower edge of the upper draft shield panel

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*The person authorised to compile the technical file:* Sartorius Lab Instruments GmbH & Co. KG  
International Certification Management  
37070 Goettingen, Germany

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Sartorius Lab Instruments GmbH & Co. KG  
Goettingen, 2016-04-20

Dr. Reinhard Baumfalk  
Vice President R&D

Dr. Dieter Klausgrete  
Head of International Certification Management

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*This declaration certifies conformity with the above mentioned EC and EU Directives, but does not guarantee product attributes. Unauthorised product modifications make this declaration invalid. The safety information in the associated product documentation must be observed.*

Doc: 2003687-02 SLI13CE001-02.de,en 1 / 1 PMF: 2018167 OP-113\_fo1\_2015.10.12

**M**

For verified balances, the conformity declaration supplied with the balance is valid for use in the EEA. Please retain this declaration.

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# Certificate of Compliance

**Certificate:** 2497246 **Master Contract:** 167555 (056628\_0\_000)  
**Project:** 70130821 **Date Issued:** 2017-04-28  
**Issued to:** Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
Goettingen, 37079  
GERMANY

*The products listed below are eligible to bear the CSA Mark shown with adjacent indicators 'C' and 'US' for Canada and US or with adjacent indicator 'US' for US only or without either indicator for Canada only.*

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## **PRODUCTS**

CLASS - C872105 - LABORATORY EQUIPMENT-Electrical

CLASS - C872185 - ELECTRICAL EQUIPMENT FOR LABORATORY USE-Certified to US Standards

Laboratory Scales, Practum, Quintix and Secura, Models AAAAAAABBBB-CDDD, rated: 15 Vdc, 3 W

(max.).

High Resolution Laboratory Scale Series models Quintix.....-..... and Secura.....-....., rated: 15 Vdc, 4.5 W (max.).

## **Notes:**

1. The external AC Adapter for use with the above models is considered Equipment Class II, Pollution Degree 2.
2. To be used with external AC Adapter (Sartorius part 6971790), UL approved (120 Vac) and manufactured by FRIWO GERAETEBAU GMBH (type 153779). Rated output 15 Vdc, 530 mA.  
The equipment can also be used with any other power supply approved for Canada and USA under the 60950-1 standard and which provides Reinforced insulation between mains and secondary circuits. It must be rated for 3000m and up to 50°C and must provide an output of 15Vdc with a minimum of 300mA.
3. The above Laboratory scales were evaluated for 50 C ambient, 3000 m elevation.

**Certificate:** 2497246  
**Project:** 70130821

**Master Contract:** 167555  
**Date Issued:** 2017-04-28

**APPLICABLE REQUIREMENTS**

CAN/CSA-C22.2 No. 61010-1-04  
UL Std. No. 61010-1 (2nd Edition)

Safety Requirements for Electrical Equipment for Measurement,  
Control, and Laboratory Use, Part 1: General Requirements  
Safety Requirements for Electrical Equipment for Measurement,  
Control, and Laboratory Use - Part 1: General Requirements

### Supplement to Certificate of Compliance

Certificate: 2497246

Master Contract: 167555 (056628\_0\_000)

*The products listed, including the latest revision described below, are eligible to be marked in accordance with the referenced Certificate.*

#### Product Certification History

Project	Date	Description
70130821	2017-04-28	CSA c/us report update (2497246) for modifying the external AC Adapter detail of Laboratory Scales, Practum, Quintix and Secura, Models AAAAAAABBBB-CDDD and High Resolution Laboratory Scale Series models Quintix..... and Secura.....
2757579	2014-10-15	Update to Report to cover new models - high resolution scales Secura and Quintix, addition of Peltier element, new PCB and TFT display.
2604694	2013-02-22	Report correction - model name from Quintex to Quintix.
2497246	2012-04-13	CSA Certification (c us Mark) for Scales Secura, Quintex and Practum.

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen, Germany

Phone: +49.551.308.0  
Fax: +49.551.308.3289  
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:  
08 | 2017