Operating Instructions | Bedienungsanleitung | Mode d'emploi |  
Manuale d'uso | Instrucciones de manejo

## Microsart® e.jet Transfer Pump

Laboratory Vacuum Pump  
166MP-4

# Contents | Inhalt | Table des matières | Indice | Índice

English.....	Page 2
Deutsch.....	Seite 11
Français.....	page 20
Italiano.....	pagina 29
Español.....	página 38

## Table of contents

<b>1 Delivery contents.....</b>	<b>3</b>
1.1. Microsart® e.jet Transfer Pump.....	3
<b>2 Description.....</b>	<b>4</b>
2.1. General notes.....	4
2.2. Areas of application.....	4
2.3. Features of the Microsart® e.jet.....	4
2.4. Operating principles.....	4
<b>3 General safety instructions.....</b>	<b>4</b>
<b>4 Transport and storage.....</b>	<b>5</b>
<b>5 Setting up.....</b>	<b>5</b>
5.1. General advice.....	5
5.2. Operating conditions.....	5
<b>6 Installation.....</b>	<b>6</b>
6.1. Installation.....	6
6.2. Installing the Microsart® e.jet.....	6
6.3. Operation.....	6
6.4. Removing the pump.....	6
6.5. Troubleshooting.....	6
<b>7 Technical Data.....</b>	<b>7</b>
7.1. Microsart® e.jet.....	7
<b>8 Spare parts.....</b>	<b>7</b>
8.1. Spare parts list Microsart® e.jet.....	7

You have selected a Sartorius Fluid Transfer pump. The following tips will help you to operate it safely and reliably over a long period. It is essential that you read these Operating Instructions before putting the pump into service. Follow them for all applications as they help prevent dangerous situations and damage to the pump and system.

## Symbols

① Position in the illustration

– Important point

● Task

• Advise to users

Warning

Sartorius reserve the right to make changes

## 1 Delivery contents

### 1.1 Microsart® e.jet Pump with Quick Connections

Power adapter with country specific power plugs are included.

## 2 Description

### 2.1 General notes

- It is essential that you read these Operating Instructions before putting the transfer pump into service. Follow them for all applications as they help prevent dangerous situations and damage.
- The Microsart® e.jet transfer pumps have been developed by Sartorius especially for laboratory and industrial applications and distinguish in reliability and simple operation.

### 2.2 Areas of application

Transfer of neutral and the most aggressive liquids in laboratories. Especially for vacuum filtration in Microbiology sample preparation.

### 2.3 Features of the Microsart® e.jet

- They are self-priming and can operate dry.
- By using chemically resistant materials such as PTFE, ETFE, EPDM and PP for parts that come into contact with the liquid, almost all neutral and aggressive liquids can be handled. Chemical Compatibility is given for example for the following disinfectants: Ethanol, Isopropanol, Peracetic Acid and Hydrogen Peroxide.
- Due to the compact construction it very little working surface required.
- Maintenance free

### 2.4 Operating principles

The Microsart® e.jet transfer pump is an oscillating positive displacement pump. An eccentric converts the rotary motion of the drive shaft into an oscillating movement of a connecting rod, which in turn transmits its motion to the diaphragm. In combination with inlet and outlet valves, this diaphragm motion provides the pumping action.

## 3 General safety instructions

- Ensure that the pump is only used for its intended purpose.
- The unit may not be used in areas in which there is a danger of explosion.
- Only connect the pump into fused sockets.
- Observe all relevant safety and accident prevention regulations.

### Safety Notice

„Before working on the pump, disconnect it from the mains supply!“

- The parts of the equipment that contain live parts may only be opened by qualified persons.
- The Operating Instructions should always be kept readily available and near to the equipment.
- When cleaning the unit, make sure that no liquid enters the housing.
- Use only genuine spare parts from Sartorius.

**If you return your Sartorius pump for repair, please state for what medium it has been used for, especially if this is very aggressive | dangerous.**

**If the pump has been used for dangerous or highly aggressive media, please clean it before dispatch.**

## 4 Transport and storage

- When packing the pump, make sure that it is properly secured inside the package.
- The package must be stable enough to withstand any transport damages.

**Use original packaging**

## 5 Setting up

### 5.1 General advice

As you set up the pump please pay attention to the following points.

- Ensure that the pump is on a stable surface and is standing upright on its 4 rubber feet.
- The pump can be mounted from underneath with four screws.
- All special safety instructions for the liquid | gas being handled must be observed.
- The head materials, hose connectors and hoses must be resistant against the liquids being pumped.

① Check that the connections between hoses and hose barbs are connected securely.

**It is sensible to secure the hoses on to the hose barbs with a clip | clamp. This is even more important when pumping aggressive liquids or when operating at high pressures.**

### 5.2 Operating conditions

- Permissible ambient temperature +5 to +40°C
- The presence of solids in the medium being handled may lead to faulty operation of the pump.
- Permissible temperature of the medium being handled is +5 to +80°C.
- The pump should not operate against a closed system.
- The pump is limited against which pressure | vacuum it can start running. The ideal conditions are when the pump can start at atmospheric conditions as then the priming is guaranteed.

Typical set-up: Microsart® e.jet (166MP-4) connected to a 3-branch manifold (16824-CS)

## 6 Installation

### 6.1 Installation

Refer to the points in section 5.

### 6.2 Installing the Microsart® e.jet

- The pump is delivered ready to use.
- ❶ Either push use the installed quick connections and push the hose over the thinner part of the quick connection or push the hose directly over the hose bars (Microsart® e.jet hose diameter 10 mm)

### 6.3 Operation

- ❷ Plug the adaptor in.
- ❸ Start the pumping by pressing the ON | OFF button
- ❹ Note: The pump is adjusted to maximum flow. The generated vacuum is 0.4 bar abs..

### 6.4 Removing the pump

- ❶ After using the pump rinse the pump and system with a neutral liquid and then pump it empty.

### 6.5 Troubleshooting

Before working on the pump ensure the mains supply is isolated.

#### Motor fails to start

- Pump is not connected to power supply.
- Power supply is switched off.
- The pressure in the delivery line from the pump is too high.
- Pump is on "STANDBY" ON- | OFF button not pressed.
- External controlling is connected without input signal.

#### Pump not sucking liquid in

- Suction side of pump is not connected.
- Hose connections are leaking.
- A system valve is closed, or blocked.
- The pump head is filled with gas. The gas pressure after the pump can not be overcome because of its compressibility in the working chamber.
- Diaphragm and | or valves are worn out or dirty.

#### Flow rate, suction height, or pressure height is insufficient.

- The periphery installation contains components (hoses, valves, filters etc) which create too much resistance.
- Hose connectors are not tight.
- Diaphragm or valves are worn out or are covered in dirt.

## 7 Technical Data

### 7.1 Microsart® e.jet

		Microsart® e-jet
Dimensions	W×L×H	12×17×19 cm
Flow rate <sup>1)</sup>	l/min	> 4.0
Suction height	mWg	6.0
Max Pressure	bar	1.0
Weight		
– of pump	g	1,425.3
– of power supply	g	242.6
Max. ambient temperature	°C	+5...+40
Max. temperature of the liquid	°C	+5...+80
Max. viscosity of the liquid	cSt	max. 150
Protection type		IP64
Protection class		Ⓜ
Power supply		100–240 V AC, 47–63 Hz 800 mA UL Listed, USA & Canada
Pump		24 V DC   1.96 A
Materials in direct contact with liquid		PTFE, ETFE, Polypropylen, EPDM, PSU, POM

<sup>1)</sup> Measured with pressure 0 mWg

### 7.2 CE-Marking

the Microsart® e.jet transfer pumps are in accordance with the following standards.

- EN 61010-1
- EN 60529
- EN 61326-1
- EN 55022
- EN 61000-4-6
- EN 61000-4-3
- EN 61000-4-4
- EN 61000-4-2
- EN 61000-4-5
- EN 61000-4-11
- EN 61000-3-2
- EN 61000-3-3

## 8 Spare parts & Accessories

### 8.1 Spare parts list Microsart® e.jet

Description	Order No.
Service Kit Microsart® e.jet Complete Kit for Pump head exchange	1EP---0003
Foot switch	1ZE---0053
Power supply adapter	1EE---0012
Threaded fitting DN 10	1EAF--0020
Microsart® e.jet Quick Connection Coupling (PSU) and Nipple (POM)	1EAS--0027

## Information and Instructions on Disposal

The equipment, including consumables and empty non-rechargeable and rechargeable batteries, does not belong in your regular household waste; this equipment is manufactured from high-grade materials which can be recycled and reused. European Directive 2012/19/EU (WEEE) on Waste Electrical and Electronic Equipment (WEEE) requires that electrical and electronic equipment be collected and disposed of separately from other unsorted municipal waste, with the aim of recycling it.

The crossed-out waste bin symbol indicates that separate collection is required. In Germany and several other countries, Sartorius itself assumes responsibility for the return and legally compliant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators. Please contact Sartorius Service for this purpose.

### Procedure

- ▷ Detailed information with service addresses for the disposal of your device can be found on our website: [www.sartorius.com](http://www.sartorius.com).
- ▷ In countries that are not members of the European Economic Area or in which there are no Sartorius subsidiaries: Contact the local authorities or the local waste management company.

## Non-hazardous Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment | component arrives. An additional copy must accompany the equipment/component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilizable and | or repairable used electrical and electronic equipment and no longer utilizable and | or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by supervisory authorities (used electrical and electronic equipment).

Model: \_\_\_\_\_ Serial No.: \_\_\_\_\_

I | we guarantee that the following measures have been taken:

- The equipment | component has been freed of hazardous materials to ensure that persons involved with the handling/repair are not exposed to any hazard or danger whatsoever.
- The equipment has been packaged for safety and the complete labeling has been affixed to it.
- The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.

The person sending in the equipment | component is hereby notified that said person shall be liable to Sartorius, its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius, its affiliates, its subsidiaries, dealers or distributors involved with the handling/repair of the equipment | component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.

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*The person authorised to compile the technical file:*

Sartorius Lab Instruments GmbH & Co. KG  
International Certification Management  
37070 Goettingen, Germany

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Dr. Reinhard Baumfalk  
Vice President R&D

*i.v. Klausgrete*

Dr. Dieter Klausgrete  
Head of International Certification Management

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*This declaration certifies conformity with the above mentioned EC and EU Directives, but does not guarantee product attributes. Unauthorised product modifications make this declaration invalid. The safety information in the associated product documentation must be observed.*

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Sartorius Lab Instruments GmbH & Co. KG  
Goettingen, 20/04/2016

Dr Reinhard Baumfalk  
Vice President R&D

Dr Dieter Klausgrete  
Head of International Certification Management

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Phone: +49.551.308.0  
Fax: +49.551.308.3289  
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:  
04 | 2016