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Instructions for Use
Polycarbonate In-Line Filter Holder 50 mm

1. Specifications

Material:
Polycarbonate, except for
Polypropylene filter supports

Sealing:
Silicone O-ring, 40×5 mm
(6980110)

Filter diameter: 50 mm
Prefilter diameter: 40 mm
Filtration area: 12.5 cm²

Connections inlet and outlet:
Female thread with connectors
6985004, M 12×1 male
thread/nipple for ~10 mm hose.
Hose nipple at inlet may be replaced
with connector 16880 for syringe.
Max. operating pressure:
700 kPa (7 bar)
Sterilization: Autoclave at 121°C

2. Cleaning and Care of Filter Holder

- 2.1 Completely disassemble filter holder (Fig. 1).
- 2.2 Using a soft brush, wash all parts with warm water and a laboratory detergent for plastics.
- 2.3 Rinse all parts with warm water first, then with distilled water.
- 2.4 Dry all parts with compressed air or drip dry. Do not dry using towels or cloth since they may leave lint particles on filter holder.

3. How to Assemble Filter Holder and Insert Membrane Filter

- 3.1 Screw hose nipple onto bottom portion of filter holder (Fig. 2).
- 3.2 With the corrugated side face-up, place frit (Fig. 3) in bottom portion and then place membrane filter on top of frit.
- 3.3 Place sealing ring on top of membrane filter. If a prefilter is required for a particular filtration procedure, place it in the center of ring (fits exactly inside sealing ring) (Fig. 4).
- 3.4 Screw second hose nipple onto upper portion of filter holder and press second frit with corrugated side face-down in upper portion (Fig. 5).

3.5 Place upper portion on top of bottom portion so that they fit exactly together and screw both parts together (Fig. 6).

3.6 The filter holder is ready to use.

4. Clarifying Filtration or Particle Removal

- 4.1 Directions for assembling filter holder and inserting membrane filter and prefilter, if desired, are given above in steps 3–3.5 and illustrated in Figs. 2–6.
- 4.2 For filtration of liquids or gases, you may use filter holder either integrated into a piping system or connected to a peristaltic pump, or you may supply the filter holder with the medium to be filtered through a hose connected to a pressure tank.
- 4.3 The 16508 in-line filter holder may be converted into a syringe filter holder by unscrewing the hose nipples and replacing them with connectors 16881 and 16880 (Fig. 7).

5. Sterilizing Filtration

- 5.1 Directions for assembling filter holder, inserting membrane filter and prefilter, if desired, are given above in steps 3–3.5 and are illustrated in Figs. 2–6.
- 5.2 Wrap filter holder containing filter in aluminum foil or waterproof paper and autoclave at 121°C for 30 minutes.
- 5.3 After sterilizing, connect to the system as under 4.2, 4.3, observing usual precautions in order to maintain sterility.

6. Microbiological Testing

6.1 For microbiological testing, follow the same procedure given for sterilizing filtration. After filtering, remove membrane filter from filter holder while observing all precautions to avoid contamination. Incubate filter on a culture medium, i.e. on agar or on a Sartorius Nutrient Pad (Fig. 8).

7. Recommended Accessories and Spare Parts, Accessories

16680	Stainless steel connector, M 12×1 male thread/ male Luer slip
16685-----2	Dosing syringe, 0,5 to 5,0 ml in 0,5 ml graduations
16647-E	Disposable plastic syringe (12 per pkg.)
16639	3-way valve for syringe filter holders
16624	Stainless steel forceps
1324	Stainless steel needle for syringes with male Luer slip

Spare Parts

6980110	Silicone O-ring, 40×5 mm
6980232	Polypropylene filter support (2 per pkg.)

Instructions for Use
Polycarbonate In-Line Filter Holder 50 mm

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Printed and copyrighted by Sartorius
Sartorius AG