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Lire attentivement ce mode d'emploi et conserver proche de l'appareil.
Please read this operating instructions carefully and keep it near the equipment.

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Table of contents	Page
Structure of the device	5
Extent of delivery	7
Warranty	7
schuett solaris Loop Sterilizer	7
Safety instructions	9
Instructions for set-up / Line/Mains connection / Putting into operation	9
Operation	11
Cleaning	11
Trouble shooting	13
Replacing the heater element	15
Technical Data	17
Accessories	17
Spare-parts	17
Laboratory report - Expertise	19

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
Structure of the device

Fig. 1: Front view


- (a) Main switch
- (b) Sensor
- (c) Frontal opening – Quartz tube
- (d) Head section of heat sink

Fig. 2: Side view

- (e) Knurled head screws
- (f) Heat sink
- (m) Cylinder pin

Fig. 3: Back view

- (g) IEC socket with fuse carrier
- (h) Bore for extraction of the Quartz tube
- (i) Covering plate
- (j) Back part of the sterilizer
- (k) Supply circuit

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Extent of delivery

- schuett solaris Loop Sterilizer
- Quartz tube with bottom (**before putting into operation see page 7**)
- Cylinder pin, stainless steel (for removal of Quartz tube, see page 9)
- Power cord
- Operating instructions
- Guarantee card with reply card

Warranty

The manufacturer grants a one year guarantee on the schuett solaris Loop Sterilizer. To establish your right for warranty coverage, send in the reply card accompanying the Warranty Card to the manufacturer within 14 days. Device faults and malfunctions resulting from improper use are not covered by the warranty. Only original accessories and spare-parts supplied by the manufacturer may be used.

schuett solaris Loop Sterilizer

The schuett solaris Loop Sterilizer is ideal for sterilizing inoculating loops, forceps and instruments made from temperature resistant material, e. g. stainless steel or platinum-iridium alloy in a laboratory surrounding.

This application with the electrical driven sterilizer replaces a gas driven Bunsen Burner. Providing maximal safety because the sterilizer works without any open flame or flammable gases. The schuett solaris Loop Sterilizer is operable directly after switched on the device, there is no need of preheating.

The schuett solaris Loop Sterilizer offers an integrated timer. The sterilizer detects the approaching inoculating loop by means of a contact-free sensor and automatically starts the reproducible sterilization cycle for 5 or 7 seconds. After the sterilization time has elapsed, the heating is automatically stopped which is readily identifiable by observation of the red illumination inside the sterilizer. Following a 3 second holding period after the inoculation loop is pulled out the sterilizer is ready for use again.

Only reusable inoculating loops or instruments made of heat resistant wire must be used. No responsibility is taken for using the device for any other purpose.

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Safety instructions

The schuett solaris Loop Sterilizer is only to be used for the intended operation as sterilizer and only in a faultless condition. The manufacturer assumes no liability for inappropriate applications.

Never use the schuett solaris Loop Sterilizer close to flammable and combustible substances. Do not use the sterilizer in an explosion hazardous surrounding.

Instructions for set-up / Line/Mains connection / Putting into operation

Take the schuett solaris Loop Sterilizer out of its shipping package and place it in its intended location. **Attention: The Quartz tube is pre-installed. Do not operate the schuett solaris without quartz tube as this will cause severe damage to the device!** Plug the plug of the power cord into the IEC socket and connect the mains plug to a suitable, grounded socket.

The schuett solaris Loop Sterilizer will be delivered with the maximum working angle. If you wish to set the sterilizer to a smaller angle, unscrew the knurled head screws (e) for one turn, take the two front knurled screws and pull the heat sink (f) to the front. Lower the central part, as desired. At the lowest position the sterilization feature is blocked. Use the four working angle positions for adjustment. Retighten the front knurled screws and the back knurled screws.

Switch the Loop Sterilizer on by pressing the main switch (a). The device starts a sterilization cycle for sterilization of the protecting tube for one time and is ready for use.

!! Caution !!

Putting the schuett solaris into operation without the quartz tube in place will cause severe damage to the high gloss reflector and adversely affect its sterilization characteristics.

Use cleaned quartz tubes only – otherwise the sterilization temperature may decrease.

Use only original quartz tubes supplied by the manufacturer.

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Operation

Insert the inoculating loop together with the nut of the loop holder to the frontal opening (c) of the sterilizer. The sensor (b) recognizes the approach of the instrument to the sterilizer and starts the sterilization cycle. The internal space of the sterilizer lights up in red color. Leave the inoculating loop inside the sterilizer until the red light is turned off by the timer.

Pull the sterilized inoculating loop out of the sterilizer. Following a 3 second holding period the sterilizer is ready for the next sterilization cycle.

Setting of sterilization time:

The schuett solaris is set to a sterilization time of 5 seconds upon delivery. This time is optimized for inoculating loops. If you should like to sterilize instruments of higher mass the sterilization time can be increased to 7 seconds.

Unplug the power plug. Put the device upside down with the back part (i) facing you. Please look for the opening in the bottom plate and use a small screwdriver etc. to adjust the slide switch downward.

Caution:

During the sterilization cycle there will be a temperature of approx. 1300 °C inside the sterilizer.

The red cooling body and the covering plate at the back of the sterilizer could be heated up to 200°C during the intermittent operation. Please avoid to touch these parts of the sterilizer when using the sterilizer in continuous operation.

Cleaning

All surfaces of the schuett solaris Loop Sterilizer are made of aluminium or stainless steel. Use a wet cloth and standard detergents for stainless steel for cleaning the device. Please switch off the device! Take care that no liquid enters the inside of the sterilizer.

It is possible to disinfect the device by spraying suitable disinfectants onto the device or wiping it with them. Do not sterilize the schuett solaris Loop Sterilizer! If you have questions, please contact the manufacturer.

Please clean the Quartz tube (c) of the loop sterilizer regularly. Unplug the sterilizer and let it cool down approx. 10 minutes. Insert the cylinder pin (m) supplied with the schuett solaris or a suitable object (e.g. the shaft of a loop holder) into the bore at the backside (h) for the extraction of the tube and push the Quartz tube through the front opening of the sterilizer.

Fill the tube with an appropriate cleaning solution and let it soaking, for example over night, to loosen the burned in residues. After the cleaning procedure rinse the tube carefully with water. **Prior to reinserting the Quartz tube please make sure it is clean and dry, both on the inside and on the outside!**

The quartz tube is held on the bottom side by a spring. Use slight force to push the tube towards its final position until it is completely surrounded by the heat sink. (see page 3, Fig. 1).

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Trouble shooting

A simple fault may often be the reason for a malfunction and as such can easily be corrected by the user himself.

Malfunction	Measures to eliminate the fault
After switching the sterilizer on, the red LED ring light round the main switch (a) is not lit	<ul style="list-style-type: none">▪ No power supply. Check the correct plugging of the power cord and if the power socket is supplied with power.▪ Check the fuse (2AmT). In order to check the fuse, unplug the IEC plug and pull the fuse holder out of the IEC socket (g) with an appropriate tool (screwdriver).
After insertion of the inoculating loop into the sterilizer, the internal space of the sterilizer does not light up in red color, the sterilizer stays cold.	<ul style="list-style-type: none">▪ Clean both windows of the sensor (b).▪ Defect of the heater element. Please exchange the heater element.
Inoculating loop is not sterile.	<ul style="list-style-type: none">▪ Quartz tube is dirty. Please clean and <u>dry</u> the Quartz tube.

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Replacing the heater element

Let the heat sink cool down for a minimum of 10 min. Disconnect the power cord. Turn the schuett solaris upside down and while the red heat sink faces you, the edge of the solaris top plate is in line with the front of your work table. Remove 2 screws using 2 mm Allen wrench screwdriver to loosen the head section of the heat sink (d). Remove both knurled head screws (e) while blocking the axle from turning with your fingers. Lower the heat sink onto the table. Power is connected to the upper socket of the heater element by means of a wire covered with ceramic tubing. While progressing, please pay attention to not damage this ceramic tubing. Tilt the head section of the heat sink while turning it slightly counterclockwise until the heater element becomes visible.

Fig. 4 und 5: Tilting and turning of the head section for disassembly of the heater element:

Remove and replace the heater element. When replacing, make sure the back portion of the heater element slides into the lower socket without using force. Tilt and turn the head section of the heat sink back into its former position. Please make sure the upper contact portion of the heater element clicks nicely into its socket. The head section of the heat sink should show now a spring-loaded feeling. Return the Allen head screws and readjust the heat sink to its proper working angle. Replace and tighten the knurled head screws and return the schuett solaris to its working position. Now, reconnect the device to power and test for proper function.

If the measures listed above do not restore the device to its operational state or if the fault is not listed above, please contact the manufacturer.

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Accessories

	Cat.-No.
Inoculating loop E 1 made of stainless steel wire, loop-Ø 1 mm (qty.=10)	3.686 112
Inoculating loop E 2 made of stainless steel wire, loop-Ø 2 mm (qty.=10)	3.686 122
Inoculating loop E 3 made of stainless steel wire, loop-Ø 3 mm (qty.=10)	3.686 132
Inoculating loop E 5 made of stainless steel wire, loop-Ø 5 mm (qty.=10)	3.686 152
Inoculating loop P-I 1 made of platin-iridium wire, loop-Ø 1 mm (qty.=1)	3.686 212
Inoculating loop P-I 2 made of platin-iridium wire, loop-Ø 2 mm (qty.=1)	3.686 222
Inoculating loop P-I 3 made of platin-iridium wire, loop-Ø 3 mm (qty.=1)	3.686 232
Inoculating loop P-I 5 made of platin-iridium wire, loop-Ø 5 mm (qty.=1)	3.686 252
Inoculating loop-holder 1 based on Kolle, 160 mm length	3.686 362
Inoculating loop-holder 2 based on Kolle, 240 mm length	3.686 442

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Department of Medical Microbiology
Drinking Water and Hygiene laboratory

Dr.med. Dipl.-Chem. Dipl.-Ing.(FH) Ulrich Schmelz – Head of the Lab
www.hygiene-goettingen.de

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Datum: 29. Juli 2013

Laboratory report – Expertise

Microbiological validation of the Inoculation loop sterilizer "schuett solaris" by means of aerobic bacterial endospores

1. Question:

In the context of this report the physical microbiozide of the device "schuett solaris" with the claim of sterilization by test germs will be checked.

Intended purpose of the device is the thermal sterilization of inoculation loops in the laboratory. The device operates with an infrared heat-source with the power of 500 watts, which focuses the infrared light onto a quartz cuvette via a mirror system.

For sterilization the inoculation loop is inserted into the quartz cuvette. In the moment of insertion infrared light is activated for 5 seconds. Within this time the inoculation loop, which is located in the cuvette in the focus of the infrared light, is heated up to ca. 750°C. Sterilization occurs through devitalization of the microorganisms on the inoculation loop.

Because soiling at the inner walls of the quartz cuvette are possible during operation, an examination of the sterilizing power also takes place in a partially darkened cell.

The test germ *Bacillus subtilis* is used in form of an endospores suspension. The test germ has an extreme resistance against inactivation and is usually used as test germ for hot air sterilization.

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The spore suspension is applied to Caso broth, so that besides the test germ, an organic load of the liquid of 22.5 g/L in form of proteins (casein), and carbohydrates, or fats (soy lecithin) is charged, which represents additional load of the thermal sterilization process.

2. Methodology:

2.1 Application of the spore suspension:

Bacillus subtilis as a test organism is purchased from the DSMZ, Braunschweig. DSM-No. 1088.

The lyophilized pellet is suspended in 5 mL of sterile physiological saline. Afterwards a three-eyelet smear on Caso-Agar is applied ("Caso" means caseine and soya-lecithine-nutrient media, also known as Tryptone-Soya-Agar (TSA)).

After incubation for 48h at a temperature of 36°C the colonies are examined macroscopically. Rough, cockade-shaped colonies are the sign of bacterial endospores formers.

Colony material is removed, gram stained and a biochemical differentiation (BBL-Crystal-test) is conducted in parallel to confirm the identity of the test germ.

A subculture confirmed as Bacillus subtilis colonies on Caso agar is created by homogeneous smears. After incubation for 48 hours at 36°C the subculture is stored for 14 days at room temperature. In this way the vegetative bacteria pass over in the form of endospores.

After this time, the culture is washed off with 80% v/v ethanol. 80% ethanol eliminates vegetative bacteria and fungi. Furthermore it conserves bacterial endospores.

After sedimentation of the endospores the alcohol is decanted. The endospores are rinsed with sterile physiological saline. After the initial flushing two more rinsing processes are following.

The obtained flushed endospores of Bacillus subtilis are used for the microbiozide test.

For examination the endospores enrichment is filled with Caso broth up to 10mL and homogenized with a vortex shaker.

The homogenate is the test suspension, which includes endospores as well as a load of protein, carbohydrate and fat.

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2.2 Determination of density of germ number in the test suspension

To determine the initial bacterial count, the test suspension is diluted in powers of ten up to a factor 10^{-5} in sterile physiological saline (serial dilution 1mL in 9mL). From each dilution step the creation of a culture on Caso agar is made by dripping 0.1 mL of the suspension of the respective dilution and by homogeneous smearing with the Drigalski spatula.

The cultures were evaluated after 48 h incubation at 36 °C. The nutrient medium of that dilution which shows 10 to 30 colonies is used for the calculation of the initial bacterial load.

The initial bacterial load is calculated as follows:

$$\left(\frac{N_{\text{Culture medium}} * 10}{0.1} \right) * VS = N_{\text{total}}$$

Explanation:

$N_{\text{culture medium}}$ = Germ number on the culture medium used for count

VS = Dilution level at 10^{-2} is the value for VS = 100

N_{total} = Total germ number in the initial approach of 10mL

2.3 Implementation of the inactivation experiment:

The test suspension described in 2.1 is homogenized by shaking. Following it is transferred to the test suspension with an annular curved inoculating loop. The inoculating loop has an approximate volume of 10µL liquid.

Then the inoculating loop is transferred to the sterilizer and the sterilization process is performed for 5 sec.

After this, the inoculating loop is dipped into sterile Caso broth. Microorganisms still adhering to the inoculation loop are transferred into the broth and can be cultivated by this way.

The inoculated broths were incubated at 56 ° C for maximum 10 days in the incubator. By this method still viable, but damaged endospores are allowed to reactivate.

Non-sterility is indicated by cloudiness of the suspension.

The experiment is performed a total of 10 times repeatedly.

Parallel a second try of inactivation is conducted in the same way with the difference that about 15% of the quartz cuvette's surface is darkened previously with a non-translucent, metal-loaded paint pen. This is to simulate routine soiling of the cuvette. For this purpose, 4 lines are arranged axially (90 °) on the cuvette that has a 17 mm outer diameter and a circumference of 5.34 cm = 53.4 mm with the metal-loaded paint pen. The line thickness is 2 mm, i.e. the overall darkening is about 8 mm of 53.4 mm overall size.

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3. Results:

3.1 Results of the determination of the initial bacterial count

Diluti.level	CFU/Plate	Diluti.factor	V _{total}	V _{Plate}	CFU total
10 ⁰	> 100	0	10mL	0.1mL	n.a.
10 ⁻¹	83	10	10mL	0.1mL	n.a.
10 ⁻²	12	100	10mL	0.1mL	12*10 ⁶ CFU/10mL
10 ⁻³	2	1000	10mL	0.1mL	n.a.
10 ⁻⁴	0	10000	10mL	0.1mL	n.a.

The 10⁻² dilution was evaluable (12 CFU/nutrient medium), this results in a germ number of 12 * 10⁶ CFU/10mL.

At an average volume of 10µL = 0.010mL of an inoculation loop the germ number per inoculation loop is approximately 12 CFU in the form of bacterial endospores.

The organic load is at 22.5 g/L (compare point 1) with 10µL theoretically 0.000225 g, corresponding to 0.225 mg per loop.

3.2 Result of the inactivation experiment

Nr.:	solaris	solaris, darkened
01	-	-
02	-	-
03	-	-
04	-	-
05	-	-
06	-	-
07	-	-
08	-	-
09	-	-
10	-	-
Pos.contr.	+	+

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4. Interpretation:

It was proven that the device "schuett solaris" leads to safe sterilization of inoculating loops (and diluting loops) in 5 sec. The test germ *Bacillus subtilis* was safely inactivated in all 10 tested samples. A reactivation of the spore-activation (56°C 10 days) at special conditions of stocks failed to appear.

Besides, it was proven that dirt appearing with daily usage has no negative impact on the success of sterilization. Nevertheless, it is recommended to clean the quartz-cuvette daily, so that the optical lucency and the effectiveness of the process of sterilization are given at any time.

Dr. Ulrich Schmelz, M.D.
Head of the Drinking Water and Hygiene Lab of the University Medical Center Goettingen
Medical doctor (specialist in Medical Microbiology); Food Chemist and Graduate Engineer

DOMINIQUE DUTSCHER SAS

schuett-biotec GmbH
Rudolf-Wissell-Straße 13
D-37079 Göttingen
Germany

Fon +49 (0) 551/5 04 10-0
Fax +49 (0) 551/5 04 10-99

info@schuett-biotec.de
www.schuett-biotec.de