Gebrauchsanleitung TPP - Zellkulturflaschen / *Instructions for Use TPP Tissue Culture Flasks*

Anleitung

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Instructions

TPP tissue culture flasks are intended for cell culture in manual use. The growth area is opto-mechanically treated for an optimal cell adhesion and proliferation. The chanted neck decreases the risk of media contact at the inside of the screw cap. It also enables a better access for pipettes and scrapers to reach the corners of the flask.

- *Open product packaging and inspect for intact bag and product.*
- *Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the flask.*
- *Close flask:*
 1. *VENT screw cap for open culture systems, by opening the cap the gas exchange takes place*
 - *Introduce screw cap until it clicks in the VENT position.*
 - ✓ *Position for defined aeration of medium: rectangle on 12 o'clock.*
 - *Further turn screw cap over the click about 90° to close flask gas tight.*
 - ✓ *Position for gas tight: rectangle on 3 o'clock*
 2. *Filter screw cap with integrated hydrophobic membrane of pore size 0.22 µm for sterile gas exchange after the cap is screwed down tightly.*
 - *Screw cap down tightly.*
- *Incubate: Stapling ring and rim at the angled neck result in a stable and secure stand in the incubator.*
- *Harvest cells:*
 - *Harvest adherent cells by trypsinization or scrape with cell scraper.*
 - *Harvest suspension cells by exchange of medium.*

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Advice and Recommendation

For homogenous cultivation of the cells it is most important to level the incubator and the continuous supply of medium. Sources of vibrations must be eliminated to avoid circular growth pattern. The necessary gas exchange takes place through the screw cap.

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Materialien

Flasche PS
 Kappe PE
 Membrane PTFE,
 Porengrösse 0.22 µm

Technical Data

Materials

Flask PS
 Screw cap PE
 Membrane PTFE,
 pore size 0.22 µm

| Abmessungen | Measurement | 90025 | 90026 | 90075 | 90076 | 90150 | 90151 | 90300 | 90301 |
|---------------------------------|-----------------------------|-------|--------|--------|--------|---------|--------|---------|--------|
| Kappe | Cap | VENT | | VENT | | VENT | | VENT | |
| Kappe | Cap | | Filter | | Filter | | Filter | | Filter |
| Membrane µm | Membrane µm | | 0.22 | | 0.22 | | 0.22 | | 0.22 |
| Höhe mm | Heights mm | 29 | | 40 | | 50 | | 50 | |
| Breite mm | Width mm | 51 | | 87 | | 122 | | 140 | |
| Länge mm | Length mm | 92 | | 155 | | 210 | | 275 | |
| Wachstumsfläche cm ² | Growth area cm ² | 25 | | 75 | | 150 | | 300 | |
| Empf Volumen ml | Rec. volume ml | 3 – 8 | | 8 – 22 | | 15 – 45 | | 30 – 85 | |
| Max Volumen ml | Max. volume ml | 15 | | 65 | | 165 | | 410 | |