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### Instructions

TPP tissue culture flasks are intended for cell culture with necessary sterile handling techniques. The growth area is opto-mechanically treated for optimal cell adhesion and proliferation. The angled neck decreases the risk of media contact with the inside of the screw cap.

- Open the product packaging and inspect for intact bag and product.
- Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the flask.
- Tightly screw closed the filter screw cap with the integrated hydrophobic membrane for sterile gas exchange
- Incubate: Lay the flask on the incubator shelves for short or long term storage.
- **Open flask with peel-off foil:** Peel-off foil is peeled off before access to the cells. Grasp the flab and peel off the foil. Re-closure is not possible.
- Harvest cells with large access from above to flask chamber
  - Harvest adherent cells by trypsinization or scrape with the cell scraper.
  - Harvest suspension cells by exchange of medium.

### Advice and Recommendation

- Stable ring enables slip-free stacking of multiple flasks also in combination with the standard TPP tissue culture flasks.
- Barrier creates a growth area of 115 cm<sup>2</sup>
- Both flasks are available with Filter screw cap with a hydrophobic PTFE membrane with a pore size of 0.22 µm only. If a «VENT» screw cap is used, the build-up of inside pressure of more than 0.03 bar could press open the flask.
- A controlled level base surface and regular use of feeding medium are essential to uniform tissue growth.
- Sources of vibrations must be eliminated to avoid circular growth pattern.

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