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Instruction

The centrifuge tubes are intended for centrifugation of suspensions (cells, bacteria, precipitate etc). Tubes are not intended to be used for freezing use.

- *Introduction of medium: Place tube into appropriate rack, use a pipette allow the medium to flow down the sides of the tube.*
- *Tube closure: Screw down screw cap tightly to ensure gas- and liquid tightness.*
- *Centrifugation: Use applicable centrifuge rotor or adapter, carefully balance the load in the centrifuge and check rotor and lid locking mechanism.*

Advice and Recommendation

- *Alcohol <90 % plus rubbing can dissolve and/or blur the white color print on the tube.*
- *Maximum centrifugal acceleration force can be influenced by following criteria: use of organic solutions, temperature below 0 °C, centrifuging without centrifuge adapter, etc. Conduct a trial run to determine proper conditions before any critical use begins.*
- *The RCF or gx (relative centrifugal force) ratings for the disposable centrifuge tubes have been established at room temperature using tubes filled to nominal capacity with water and spun in a horizontal rotor (swing-out) centrifuge for 5 minutes. The centrifuge must be equipped with the recommended carriers, adapters, and cushions that fully support the tubes. If an angle head rotor is used or proper support is not provided, gx values will be lower.*
- *Do not confuse speed or revolutions per minute (RPM) with relative centrifugal force (RCF). Instructions for centrifuging a sample at a given RPM and time are incomplete unless the rotor or radius is specified. Protocols should always state the time and RCF value for centrifuging a sample.*
- *When spun in a high speed, deformation of tubes might occur, thus the tubes might no more fit into the standard rack.*

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Technical Data

Materialien

Schraubkappe PE
 Röhrchen PP / PS

Materials

Screw Cap PE
 Tubes PP / PS

Abmessungen	Measurements	91015	91016	91017	91019
Länge mm	<i>Length mm</i>	120	105	120	100
Durchmesser mm	<i>Diameter mm</i>	17.1	17.1	17.1	17.1
Volumen grad. ml	<i>Volume grad. ml</i>	15	13	15	12
Volumen max. ml	<i>Volume max. ml</i>	15	14	15	13
(g) max	<i>(gx) max</i>	15'500	15'500	15'500	15'500
Form	<i>Form</i>	konisch / <i>conical</i>	rund + kurz / <i>round + short</i>	rund + lang / <i>round + long</i>	flach / <i>flat</i>

Abmessungen	Measurements	91050	91051	91056	91115 (PS)
Länge mm	<i>Length mm</i>	115	115	115	120
Durchmesser mm	<i>Diameter mm</i>	30	30	30	17.1
Volumen grad. ml	<i>Volume grad. ml</i>	50	50	50	14
Volumen max. ml	<i>Volume max. ml</i>	50	50	50	15
(g) max	<i>(gx) max</i>	15'500	15'500	15'500	1'700
Form	<i>Form</i>	konisch / <i>conical</i>	konisch mit Stehrand / <i>conical with rim</i>	rund / <i>round</i>	konisch / <i>conical</i>