



## T5 Exonuclease

(recombinant protein)

Cat. No.	size
E1130-01	1 000 units
E1130-02	5 000 units

### Unit Definition:

One unit of T5 Exonuclease catalyzes the release of 1 nmol of acid-soluble nucleotides from double-stranded calf thymus DNA in 30 minutes at 37°C.

### Storage Conditions:

Store at -20°C.

### Inactivation:

No heat inactivation.

### Concentration:

10 U/μl

### References:

1. Kaliman, A.V. et al. (1986) *FEBS*. 195, 61-64.
2. Sayers, J.R. et al. (1991) *Nucleic Acids Res.* 19, 4127-4132.
3. Sayers, J.R. et al. (1990) *J. Biol. Chem.* 265, 18311.
4. Sayers, J.R. et al. (1996) *Analytical Biochem.* 241, 186.
5. Sayers, J.R. et al. (1997) *Nucleic Acids Res.* 25 (19): 3801-3807

T5 Exonuclease is a product of T5 phage D15 gene derived from *E. coli* (1). T5 Exonuclease degrades DNA in a 5'→3' direction and also exhibits endonucleolytic activity on flap structures (branched duplex DNA containing a free single-stranded 5'-end) (2, 5). T5 Exonuclease has been shown to have no effect on double stranded circular DNA (2).

### Application:

- Plasmid mutagenesis methods.
- Removal of denatured DNA from alkaline-based plasmid purification procedures for improved cloning (4).
- Generation of plasmid sequencing templates.
- Degradation of linear ssDNA, dsDNA or nicked plasmid DNA while maintaining supercoiled plasmid DNA.

### Storage Buffer:

50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 0.1 mM EDTA, 0.1% Triton X-100, 1 mM DTT and 50% glycerol.

### Quality Control:

T5 Exonuclease is free of detectable RNases or contaminating DNA endonuclease activities. Preparations are greater than 95% pure, as judged by SDS polyacrylamide gel electrophoresis.

### Example Reaction:

#### 1. Mix:

Component:	Amount:
DNA	x μl (up to 1 μg)
10 x Buffer ONE	5 μl
T5 Exonuclease	1 μl
Water	to 50 μl

2. Incubate for 15–30 min at 37°C.
3. Add 10 mM EDTA to stop the reaction.
4. Clean-up treated samples by column purification or phenol/chloroform extraction followed by ethanol precipitation.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039  
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23