
20000005402  
MT 40\_100\_072016

**Disposable Milling Tube MT 40**  
**Disposable Milling Tube MT 100**

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4424900C

## Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Refer to the operating instructions for the drive.
- Addition of the milling unit or the related accessories can get into the medium you are working on.
- Do not use any damaged tubes e.g. with hairline cracks or splits.
- The tubes are appropriate only for a single-use. Multiple-use can destroy the tube and the milling unit.
- Excessive size and excessive hardness can damage the tubes (see "Technical data").
- Twist the mill tube on securely until it clicks home audibly.
- Take care not to fill the vessel above the mark.
- Ensure that the tube is firmly attached to the bayonet lock connector on the drive unit prior to operating the device.
- After the milling operation, the tubes may be hot; please either allow it to cool down or wear appropriate personal protective equipment when detaching it from the drive.
- The processing time is limited to three minutes. The tubes must be allowed to cool down for a few minutes before processing the next sample.
- The temperature of the material must not exceed 100 °C. **Caution!** A sample of dry maize can reach this temperature just by milling for three minutes at maximum speed.
- The only access medium for directly cooling the sample in the sample vessel is not permitted.
- Only process material that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do not use the tubes in explosive atmospheres, with hazardous substances or under vacuum.
- Do not process materials that are radioactive, flammable or explosive, or materials that react chemically with each other (releasing large amounts of energy).

## Production information

- The milling tubes are for use in conjunction with the **Tube Mill control** for milling dry samples up to 5 Mohs hardness. The sample pieces must be larger than 1 cm<sup>2</sup>, if the sample pieces are larger than this, they must be pre-milled before placing them in the tubes. The harder the sample pieces, the smaller they must be. If for instance the hardness is 5 Mohs the edge length must be no longer than 0.5 cm (corresponding to a piece size of approx. 0.125 cm<sup>2</sup>). The tubes can be used for preparing samples for foodstuffs quality control. Commodities suitable for processing include grains, sugar, coffee and tea etc..
- Dry ice can be used to extend the application range to include moist and resilient samples such as fruits, leaves, jelly babies etc.
- The tubes are ideal for applications, with which no cross contamination is permissible.
- There is no need for laborious cleaning.
- The sample can be stored in the tubes.
- The tubes are stackable with the caps (MTC 40-100), available as an accessory.
- **Materials:**
- **Sample vessel and coupling cover:** PP
- **Sealing ring:** PUR and viscose
- **Coupling and blade attachment:** PA
- **Blade:** Stainless steel (AISI 301)
- **Delivery scope:**
  - 10 Disposable Milling Tube MT 40
  - 100 Disposable Milling Tube MT 40
  - Operating instructions
- **MT 40-100:**
  - 10 Disposable Milling Tube MT 40 (gamma sterilized)
  - Operating instructions
- **MT 40-100 sterile:**
  - 10 Disposable Milling Tube MT 40 (sterilisation gamma)
  - Operating instructions
- **MT 100-10:**
  - 10 Disposable Milling Tube MT 100
  - Operating instructions
- **MT 100-50:**
  - 50 Disposable Milling Tube MT 100
  - Operating instructions
- **Accessories:**
  - MTC 40-100: 100 caps to close the vessel; of the MT 40 and MT 100 (make the vessel stackable).

## Application instruction

- Store the tubes in a dry place, in their packaging and at room temperature.
- Please note the limit on the usable life of sterile variants.
- To open the tubes, press the locking cam and twist the coupling cover anticlockwise.
- Do not fill the tubes above the maximum mark on the sample vessel. Milling less material may give better results (finer eventual particle size).
- If the samples are resilient or contain residual moisture, add a teaspoonful of dry ice (dry snow, not compacted dry ice) directly into the sample vessel.
- To close the tubes, twist the coupling cover clockwise until it clicks home audibly, see Fig. 1.
- The tube is attached by twisting it clockwise on to the bayonet coupling of the drive. The tube is correctly attached when it has clicked home audibly, see Fig. 2.
- When the coupling cover has been correctly screwed on to the tube and the tube is correctly attached to the drive, one of the grip cams on the sample vessel will point directly to the rear, see Fig. 3.
- Longer processing times will cause high temperatures in the milled material. For instance a 40 ml sample of maize milled for three minutes at maximum speed can reach a temperature of approx. 100 °C.
- If the sample is temperature-sensitive, run the mill for only a short time. Cooling the sample with dry ice snow can also prevent it overheating.
- Likewise the optimal process duration and rotating frequency must be determined by attempt.
- Before you open the tubes it is advisable to tap lightly on the mounting so that all the milled material falls into the sample vessel.
- To open the tubes, release cam and twist the coupling cover anticlockwise. Open the tube always upside down.

**Attention:** Don't open the tube in the operation position.

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