


Perfect™ 100-1000 bp DNA Ladder

| Cat. No. | size |
|----------|--------|
| E3141-01 | 50 µg |
| E3141-02 | 250 µg |

Storage Conditions: Store at +4°C.
For long term storage, store at -20°C.



DNA ladder with 100 bp increments for sizing small-to-medium DNA fragments.

Description:

- Ideal for sizing linear double-stranded DNA fragments from 100 to 1000 bp.
- Contains 10 bands with fragments of the following sizes: 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000 bp.
- For easy reference, band at 500 bp is three times brighter than other bands in the ladder.
- Can be 5'-end labeled with radioisotopes and T4 Polynucleotide Kinase for visualization by autoradiography after a dephosphorylation step.

Storage Buffer:

10 mM Tris-HCl (pH 8.0 at 22°C), 1 mM EDTA, dye.

Loading:

The recommended amount of size marker to load on a gel is 0.5 µg per lane. Mix well after thawing.

Concentration:

The Perfect™ 100-1000 bp DNA Ladder is supplied at 50 µg/ml in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

5 µl load would load 240 ng DNA.*

| Band size | DNA conc. |
|-----------|-----------|
| 1000 bp | 20 ng |
| 900 bp | 20 ng |
| 800 bp | 20 ng |
| 700 bp | 20 ng |
| 600 bp | 20 ng |
| 500 bp | 60 ng |
| 400 bp | 20 ng |
| 300 bp | 20 ng |
| 200 bp | 20 ng |
| 100 bp | 20 ng |

*not for use in order to quantify the amount of DNA in bands

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

EURx Ltd. 80-297 Gdańsk Poland ul. Przyrodników 3, NIP 957-07-05-191, KRS 0000202039
www.eurx.com.pl, orders@eurx.com.pl, tel. +48 58 524 06 97, fax +48 58 341 74 23