

thermo scientific

# UV-Vis Spectrophotometer Product Guide

Powerful spectroscopy for accurate analysis

**ThermoFisher**  
SCIENTIFIC



# Understanding UV-Vis

## Why UV-Vis?

The power of UV-Vis lies in its simplicity. Ultraviolet-Visible (UV-Vis) spectrophotometry is literally as simple as “abc” since most analyses will follow Beer’s Law: Absorbance = abc. “a” is unique to each molecule at a particular wavelength; “b” is the pathlength; and “c” is the concentration of the analyte.

## The significance of spectral bandwidth

Bandwidth is like a pupil in the human eye. The wider it is, the more light the instrument can see. But sometimes, when the bandwidth is too wide, the instrument may not see enough detail in your sample. Therefore, if you know that your samples have wide absorption peaks—large bandwidths are okay (and less expensive!), but if you want to look at narrow peaks—small bandwidths may be required. Check with our product specialist when in doubt.

## Cell material

Even though glass looks perfectly transparent to our vision, it is not transparent to ultraviolet radiation (that is why you cannot get a suntan behind a closed window). For this reason “quartz” (aka “silica”) cells need to be used to measure samples in the ultraviolet (UV) region; if you are always going to be measuring above ~325 nm, then “glass” (aka “optical glass”) or plastic cells are okay and less costly to use.

## Beam configurations

Single beam spectrophotometers are simple and inexpensive. However, they are somewhat more prone to drift and should be “zeroed” once every hour. Dual-beam instruments offer better stability and can be used to monitor and subtract instrument and solvent variations, as well as provide better precision and stability for longer term measurements.

## Fixed wavelength measurement vs. scanning

UV-Vis measurements for quantitative analysis are typically performed at a single wavelength. A blank solution is measured to establish “100 %T” and then the sample is measured and the absorbance is calculated. If a calibrated method has been set up by measuring one or more standards the instrument may report concentration directly. UV-Visible spectroscopy is ideally suited for this kind of measurement. Some quantitative analyses require measurement at more than one wavelength, however, and these require a more sophisticated instrument equipped to move across the spectrum during the experiment or to measure the entire spectrum at once. These “scanning” instruments can also be used to find the location and size of absorbance peaks that can be used to identify analytes and as the best wavelength at which to establish a quantitative measurement method for the analyte. All Thermo Scientific™ spectrophotometers offer scanning capability.

## Visible vs. UV-Visible spectrophotometers

A spectrophotometer equipped with a tungsten-halogen lamp and conventional glass optics can measure wavelengths no lower than 325 nm because the lamp does not emit and the optics do not transmit below this wavelength. These instruments are generally referred to as visible spectrophotometers even though they can measure slightly into the UV and near-infrared ranges. If all the measurements made in your laboratory are done at wavelengths between ~350 and 1100 nm then a visible spectrophotometer will be a good choice.

In practice, glass optics begin absorbing light at all wavelengths less than ~350 nm, so for UV measurements below 350 nm an instrument equipped with quartz optics and a lamp that emits in the UV region is recommended. Below 325 nm it is required. Even if your sample appears transparent it may absorb ultraviolet light. Many pharmaceuticals, organic chemicals, DNA, RNA and proteins absorb in the UV range. Established methods for quantifying their concentration in solution rely on UV-Visible spectrophotometers to perform the analysis.

## The significance of pathlength


Pathlength is the distance that the light travels through the sample. Doubling the pathlength will double the measured absorbance. A 10 mm pathlength is most commonly used for UV-Vis measurements. Special long-path cells with 20, 40, 50 or 100 mm pathlengths are available for applications where the absorbance of the sample is very low, or the sample solution is very dilute. For highly absorbing or more concentrated solutions, short-path cells of 5, 2 or 1 mm pathlengths can help to keep measured absorbance within the measurement range of the instrument.

## Instruments and software featured in this guide

- Thermo Scientific™ BioMate™ 160 UV-Vis spectrophotometer
- Thermo Scientific™ Evolution™ 201 and 220 UV-Vis spectrophotometer
- Thermo Scientific Evolution 260 Bio UV-Vis spectrophotometer
- Thermo Scientific Evolution 350 UV-Vis spectrophotometer
- Thermo Scientific™ GENESYS™ 30 Visible spectrophotometer
- Thermo Scientific GENESYS 40 Visible/GENESYS 50 UV-Vis spectrophotometer
- Thermo Scientific GENESYS 140 Visible/GENESYS 150 UV-Vis spectrophotometer
- Thermo Scientific GENESYS 180 UV-Vis spectrophotometer
- Thermo Scientific™ SPECTRONIC™ 200 Visible spectrophotometer
- Thermo Scientific™ CUE™ software
- Thermo Scientific™ INSIGHT™ software
- Thermo Scientific™ VISION*lite*™ software

## Education, Basic Quality Control

Instrument	SPECTRONIC 200	GENESYS 30	GENESYS 40/140
Catalog Numbers	222-265800 (UK) 222-265900 (EU)	840-277400 (UK) 840-277300 (EU)	GENESYS 40: 840-297400 (UK) GENESYS 40: 840-297300 (EU) GENESYS 140: 840-308400 (UK) GENESYS 140: 840-308300 (EU)



<b>Spectrum</b>	Visible		
<b>Spectral Range</b>	340–1100 nm	325–1100 nm	
<b>Type</b>	Single Beam		Dual Beam
<b>Bandwidth</b>	≤4.0 nm	5 nm	
<b>Lamp(s)</b>	Tungsten-Halogen		
<b>Detector</b>	CCD	Photodiode	
<b>Standard Functions</b>	A/%T/Concentration/Scan/multi-wavelength/Quant (Up to 4 STDs)	A/%T/Concentration/Scan/Fixed with calcs/Quant (up to 6 STDs)/OD600	A/%T/Concentration/Scan/Fixed with calcs/Quant/Rate/OD600
<b>User Interface</b>	On-board or Computer Control		On-board Touchscreen Tablet or Computer Control
<b>Method/Data Storage</b>	Quant Methods to USB/All to computer	Methods On-board, Data to USB	On-board and Computer
<b>Display</b>	Color graphical LCD, 320 × 240 pixels, 7 × 5 cm, tiltable	5" diagonal, 32-bit color display, 800 × 480 pixels	7" color touchscreen, HD, 800 × 1280 pixels, fixed/tiltable
<b>Printing</b>	External Printer	Snap-on Printer Accessory	Snap-on Printer Accessory, USB, Network or Wi-Fi to page printers
<b>Optional Specialized Functionality</b>	Computer control with VISION/ite. Adds kinetics, more sophisticated scanning and quant, plus PV capabilities		Solution color with VISION/ite ColorCalc software
<b>Validation IQ/OQ</b>	No		Available soon
<b>Capacity</b>	Single Cell		Standard: Single Cell Optional: 8 or 4 cell (140 only)
<b>Liquid Samples – Cuvettes</b>	Rectangular or Cylindrical to 100 mm		
<b>Liquid Samples – Test Tubes</b>	13 mm × 100 mm	12–25 mm dia, up to 150 mm tall	12–25 mm dia, up to 150 mm tall (40 only) N/A (140 only)
<b>Temperature Control</b>	N/A	Water thermostatted single cell holder	Water thermostatted single-cell holder (40 only) Peltier thermostatted cell holder (140 only)
<b>Solid Samples</b>	Thin solid samples such as films and filters		
<b>Warranty</b>	1 year	2 year (plus 1 additional year with registration)	
<b>Benefits and Competitive Advantage</b>	<ul style="list-style-type: none"> <li>• Very simple operation</li> <li>• Auto-zero at all wavelengths will save time, reduce error</li> <li>• Live measurements every 2 seconds</li> <li>• Scan 400–900 nm in 10 seconds</li> <li>• Quant with 4 standards</li> </ul>	<ul style="list-style-type: none"> <li>• Larger, removable sample compartment</li> <li>• Large, color screen</li> <li>• Familiar keypad controls</li> <li>• Export data files to USB</li> <li>• Print data to snap-on thermal printer</li> </ul>	<ul style="list-style-type: none"> <li>• 7", HD color touchscreen</li> <li>• Easy-to-use user interface</li> <li>• Large sample compartment with front access simplifies sample handling</li> <li>• Network equipped – print or store data to a USB stick or a network</li> <li>• Go wireless – optional Wi-Fi adaptor</li> </ul>

### Featured Accessories for SPECTRONIC 200



**100 mm Rectangular Cell Holder**  
840-250500



**Test Tube Holder**  
840-314800\*  
\*Available soon


### Featured Accessories for GENESYS 30/40




**Test-Tube Holder**  
840-277500




**Tall Test-Tube Cap**  
840-277600



**Long Path Rectangular Cell Holder**  
840-277700



**5–100 mm Path Cylindrical Cell Holder**  
840-277800



**Film/Filter Holder**  
840-278100



**Thermostat Rectangular Cell Holder (no baseplate)**  
335079-000

### Featured Accessories for GENESYS 140




**8 Position Carousel for 10 mm Pathlength Cells**  
840-303400



**Peltier Thermostatted Single Cell Holder**  
840-306800 (UK)  
840-306700 (EU)



**Integrated Sipper System**  
840-307100 (UK)  
840-307000 (EU)




**4 × 10–50 mm Rectangular Cell Holder**  
840-314700



**Long Path Rectangular Cell Holder**  
840-303800



**5–100 mm Path Cylindrical Cell Holder**  
840-303100



**Film/Filter Holder**  
840-303200

## Routine QC and Research

Instrument	GENESYS 40/50	GENESYS 140/150	GENESYS 180
Catalog Numbers	GENESYS 40: 840-297400 (UK) GENESYS 40: 840-297300 (EU) GENESYS 50: 840-298400 (UK) GENESYS 50: 840-298300 (EU)	GENESYS 140: 840-308400 (UK) GENESYS 140: 840-308300 (EU) GENESYS 150: 840-300400 (UK) GENESYS 150: 840-300300 (EU)	840-309400 (UK) 840-309300 (EU)



<b>Spectrum</b>	Visible/UV-Visible		UV-Visible
<b>Spectral Range</b>	325–110 nm/190–1100 nm		190–1100 nm
<b>Type</b>	Dual Beam		Double Beam
<b>Bandwidth</b>	2.0 nm		
<b>Lamp(s)</b>	Xenon		
<b>Detector</b>	Dual Photodiodes		
<b>Standard Functions</b>	A/%T/Concentration/Scan/Rate/Quant, Fixed with Equations		
<b>User Interface</b>	On-board Touchscreen or Computer Control		
<b>Method/Data Storage</b>	On-board and Computer		
<b>Display</b>	7" color touchscreen, HD, 800 × 1280 pixels, fixed	7" color touchscreen, HD, 800 × 1280 pixels, tiltable	
<b>Printing</b>	USB, Network and Wireless printing to page printers. Snap-on thermal printer option.		
<b>Optional Specialized Functionality</b>	Solution color with VISION/ite ColorCalc software		
<b>Validation IQ/OQ</b>	No		
<b>Capacity</b>	Single Cell	Single Cell (standard), Multi-cell up to 8 (optional)	Single Cell and Multi-cell up to 8 (standard)
<b>Liquid Samples</b>	Up to 100 mm cuvettes; Up to 25 mm test tubes	Up to 100 mm cuvettes	
<b>Temperature Control</b>	Water thermostatted single cell holder	Single cell Peltier thermostatted cell holder	
<b>Solid Samples</b>	Thin solid samples such as films and filters		
<b>Warranty</b>	2 years (plus 1 additional year with registration)		
<b>Benefits and Competitive Advantage</b>	<ul style="list-style-type: none"> <li>7", HD color touchscreen</li> <li>Easy-to-use user interface increases speed and reliability for routine analysis</li> <li>Large sample compartment with front access simplifies sample handling</li> <li>Network equipped – print or store data to a USB stick or a network</li> <li>Go wireless – optional Wi-Fi adaptor</li> </ul>	<ul style="list-style-type: none"> <li>Tiltable 7", HD color touchscreen tablet to avoid glare</li> <li>Compatible with automated cell changers, Peltier, sipper and fiber optic probe accessory options</li> <li>Network equipped – print or store data to a USB stick or a network</li> <li>Go wireless – optional Wi-Fi adaptor</li> </ul>	<ul style="list-style-type: none"> <li>Double-beam optics</li> <li>Ideal solution when a reference beam is required, such as for kinetics</li> <li>8-cell carousel standard</li> <li>Make measurements outside the sample compartment and without cuvettes using fiber optic coupler and probe</li> </ul>

### Featured Accessories for GENESYS 40/50



**Test Tube Holder**  
840-277500




**Tall Test Tube Cap**  
840-277600




**Long Path Rectangular Cell Holder**  
840-277700



**Printer**  
840-278000



**Long Path Cylindrical Cell Holder**  
840-277800



**Film Filter Holder**  
840-278100

### Featured Accessories for GENESYS 140/150/180




**8 Position Carousel for 10 mm Pathlength Cells**  
840-303400



**Wi-Fi Adaptor for GENESYS 40-180**  
840-309900



**Long Path Rectangular Cell Holder**  
840-303800



**Film/Filter Holder**  
840-278100




**5–100 mm Path Cylindrical Cell Holder**  
840-277800



**Printer**  
840-278000



**Integrated Sipper System**  
840-307100 (UK)  
840-307000 (EU)



**4 × 10–50 mm Rectangular Cell Holder**  
840-314700

## Advanced QC, Research and Pharma

Instrument	Evolution 201	Evolution 220	Evolution 350
Catalog Numbers	Computer Control: 840-210800 Tablet Control: 912A0883	Computer Control: 840-210600 Tablet Control: 912A0884	840-310800



<b>Spectrum</b>	UV-Visible		
<b>Spectral Range</b>	190–1100 nm		
<b>Type</b>	Double Beam		
<b>Bandwidth</b>	1.0 nm	1.0 nm, 2.0 nm & AFBG	0.5, 1.0, 1.5, 2.0 & 4.0 nm
<b>Lamp(s)</b>	Xenon		
<b>Detector</b>	Dual Photodiodes		
<b>Standard Functions</b>	A/%T/Quant/Fixed Wavelength with Advanced Calculations/Scan/Rate/Performance Verification		
<b>User Interface</b>	Computer or Tablet Control		Computer Control
<b>Method/Data Storage</b>	On-board Tablet and Computer		Computer
<b>Display</b>	Optional Tablet Control Accessory, tiltable		N/A
<b>Printing</b>	Windows® Compatible Printers		
<b>Optional Specialized Functionality</b>	21 CFR Part 11, Customized Analyzer, Color, Biochemical Methods including DNA Melting, Autosampler		
<b>Validation IQ/OQ</b>	Yes		
<b>Capacity</b>	Multi-cell (up to 8)		
<b>Liquid Samples</b>	Up to 100 mm cuvettes; Up to 25 mm tubes and vials		
<b>Temperature Control</b>	Peltier accessories and a range of thermostatted cell holders with liquid recirculators		
<b>Solid Samples</b>	Thin solid samples	Thin and thick solids, powders, reflecting surfaces, films, DRA, SRA	Thin and thick solids, specular surfaces and films, SRA
<b>Warranty</b>	1 year		2 years
<b>Benefits and Competitive Advantage</b>	<ul style="list-style-type: none"> <li>• Available with Tablet Control Module</li> <li>• Double-beam configuration</li> <li>• Xenon lamp requires no warm-up for instant measurements and has 3-year warranty</li> <li>• Application Focused Beam Geometry (AFBG) matches hardware to your application</li> <li>• INSIGHT software is straightforward and powerful</li> <li>• CUE software scripting capability to simplify complex assays with the touch of a button</li> <li>• Specialized, easy-to-use accessories for increased productivity</li> </ul>		<ul style="list-style-type: none"> <li>• Powered by INSIGHT software</li> <li>• Double-beam configuration with large sample compartment and 21 cm beam separation</li> <li>• Xenon lamp requires no warm-up for instant measurements and has 3-year warranty</li> <li>• 0.5 nm bandwidth for highest resolution data</li> <li>• Smart Accessories™ with unique serial number recorded for comprehensive event logs and audit trails</li> <li>• Automated performance verification with CVC accessory</li> </ul>

### Featured Accessories for Evolution Spectrophotometers



**8-Cell Peltier System**

699-131200



**Single Cell Peltier System**

699-131100



**Combination Filter, Vial & Square 10 mm Cell Holder**

840-211800



**100 mm Rectangular Cell Holder**

10011101



**Smart Thermostatted Rotary 7-Cell Changer**

840-214900



**Smart Thermostatted Linear 8-Cell Changer**

840-219000



**Standard Fiber Optic Dip Probe with 10 mm Tip**

840-305000



**Temperature Probe Hub and Temperature Probes**

840-214600



**ISA-220 Integrating Sphere**

222-269400



**Powder Cell Holder**

840-289300

## Life Science

Instrument	GENESYS 30	BioMate 160	Evolution 260 Bio
Catalog Numbers	840-277400 (UK) 840-277300 (EU)	840-301400 (UK) 840-301300 (EU)	840-211000



<b>Spectrum</b>	Visible	UV-Visible	
<b>Spectral Range</b>	325–1100 nm	190–1100 nm	
<b>Type</b>	Single Beam	Dual Beam	Double Beam
<b>Bandwidth</b>	5.0 nm	2.0 nm	1.0 nm, 2.0 nm & AFBG
<b>Lamp(s)</b>	Tungsten-Halogen	Xenon	
<b>Detector</b>	Silicon Photodiode	Dual Silicon Photodiodes	
<b>Standard Functions</b>	A/%T/Concentration/Scan/OD600/ Fixed with eqns/Quant	A/%T/Concentration/Scan/OD600/ Fixed with eqns/Quant/Rate/DNA & Protein methods	A/%T/Concentration/Scan/Rate/ Temperature kinetics/DNA Melting
<b>User Interface</b>	On-board keyboard and screen or Computer	On-board Touchscreen Tablet or Computer Control	Computer Control with Desktop PC or optional tablet
<b>Method/Data Storage</b>	On-board: Methods USB: Method and data Computer: Method and data	On-board: Method and Data Computer: Method and Data	Computer or tablet: method and data
<b>Display</b>	5" diagonal, 32-bit color display, 800 × 480 pixels	7" color touchscreen, HD, 800 × 1280 pixels, tiltable	Optional, tiltable
<b>Printing</b>	Optional GENESYS printer only	Optional GENESYS printer, Page printers by USB, Network or Wi-Fi (with optional transmitter)	Windows printers
<b>Optional Specialized Functionality</b>	Optional accessory for measuring test tubes and nephelo flasks	Optional: Beckman cell holder, 50 µL disposable cell holder, 8 position cell holder (not thermo- statted), Peltier thermostatted single cell holder (20–60 °C)	Temperature ramping with Peltier thermostatted cell holders Advanced measurement automation with thermostatted multiple cell holders
<b>Validation IQ/OQ</b>		No	Yes
<b>Capacity</b>	Single Cell	Single Cell (standard), Multi-cell up to 8 (optional)	
<b>Liquid Samples</b>	Up to 100 mm cuvettes; Up to 25 mm test tubes	Up to 100 mm cuvettes	
<b>Temperature Control</b>	Optional water thermostatted single cell holder	Optional Peltier thermostatted single cell holder (20–60 °C)	Peltier accessories and a range of thermostatted cell holders with liquid recirculators
<b>Warranty</b>	2 year (plus 1 additional year with registration)		1 year
<b>Benefits and Competitive Advantage</b>	<ul style="list-style-type: none"> <li>• Easy to learn, easy to use interface with 5" high resolution color screen and tactile rubber keypad</li> <li>• Large sample compartment with front access simplifies sample handling</li> <li>• Best choice for OD600 measurements</li> </ul>	<ul style="list-style-type: none"> <li>• Tiltable 7", high-resolution touchscreen tablet to avoid glare</li> <li>• Compatible with automated cell changers, Peltier, sipper and fiber optic probe accessory options</li> <li>• Network equipped – print or store data to a USB stick or a network</li> <li>• Go wireless – optional Wi-Fi adaptor</li> </ul>	<ul style="list-style-type: none"> <li>• Double-beam optics</li> <li>• Ideal solution when a reference beam is required, such as for kinetics</li> <li>• Widest range of accessory options</li> <li>• Best option for temperature control</li> <li>• Temperature kinetics and DNA melting with single or 8-cell Peltier systems</li> </ul>

### Featured Accessories for GENESYS 30



**Test Tube Holder**  
840-277500



**Tall Test Tube Cap**  
840-277600



**GENESYS Printer  
Accessory**  
840-278000

### Featured Accessories for BioMate 160



**Beckman Style  
Cell Holder**  
840-299000



**8 Position Carousel  
for 10 mm  
Pathlength Cells**  
840-303400



**Microcell Holder**  
840-303300



**Peltier Thermostatted  
Single-Cell Holder**  
840-306600



**Printer Accessory**  
840-278000

### Featured Accessories for Evolution 260 Bio



**8-Cell Peltier System**  
699-131200



**Single Cell Peltier  
System**  
699-131100



**Temperature Probe  
Hub and  
Temperature Probes**  
840-214600

# Replacement Guide

Original Unit	Replace with	Original Unit	Replace with	Original Unit	Replace with	Original Unit	Replace with	Original Unit	Replace with
<b>Agilent/HP/Varian</b>		<b>GE Healthcare<sup>1</sup>, Pharmacia/LKB<sup>1</sup></b>		<b>Jenway</b>		<b>Shimadzu</b>		<b>Thermo Scientific/SPECTRONIC/Milton Roy</b> – continued	
HP8450A/HP8452A	Evolution 201	Novaspec i/ii/iii	SPECTRONIC 200	6051/6300/6320D/7200	SPECTRONIC 200	UV 100-01/110-01	Evolution 201	Helios Alpha	GENESYS 180 or Evolution 201
Agilent 8453/8454	Evolution 201	Novaspec Plus	GENESYS 30	6100	GENESYS 30	UV 120-01/120-02	Evolution 201	Helios Beta	GENESYS 150
DMS 70/80/90/100	Evolution 201	Ultrospec/Ultrospec III	BioMate 160	6105/6305/7205	GENESYS 50	UV-140/150	Evolution 201	Helios Gamma	GENESYS 50 or 150
DMS 200/300	Evolution 350	Ultrospec 10 Cell Density	GENESYS 30	6310	GENESYS 30 or 40	UV-160A/UV-160U	Evolution 201	Helios Delta	GENESYS 30 or 40
Cary 50/60	Evolution 201	Ultrospec 500pro	BioMate 160	6400	GENESYS 140	UV-1201V	GENESYS 30	Helios Epsilon	GENESYS 30
Cary 1, 1E, 100 / 3, 3E, 300	Evolution 350	Ultrospec 1000/2000	BioMate 160	6405/6505/6705	GENESYS 150	UV-1201UV	GENESYS 50	Helios Zeta	Evolution 201
<b>Beckman<sup>1</sup></b>		Ultrospec 3000	Evolution 201	6500/6700	GENESYS 40 or 140	UV mini-1240	GENESYS 150	Helios Omega	GENESYS 150
DU-62	GENESYS 30	Ultrospec 1100pro	BioMate 160	6715	GENESYS 150 or Evolution 201	UV mini-1240V	GENESYS 40	UV1/UV2/UV-10	GENESYS 180 or Evolution 201
DU-63/64/65	BioMate 160	Ultrospec 2100pro/2300pro	Evolution 201	6800	Evolution 201	UV-1250PC	Evolution 201	UV3	Evolution 220
DU 520/720	GENESYS 50	Ultrospec 3100pro/3300pro	Evolution 201	6850	Evolution 220	UV-1280	GENESYS 150 or BioMate 160	UV4/UV5	Evolution 350
DU 530	BioMate 160	Ultrospec 6300	Evolution 201	7300	GENESYS 30	UV-1601/UV-1650PC	Evolution 201	UV 300 Series (LC/PC)	Evolution 350
DU 640/650	BioMate 160 or Evolution 260 Bio	Ultrospec 7000	GENESYS 150	7305/7315	GENESYS 150	UV-1700/UV-1800	Evolution 201	UV 500 Series (LC/PC)	Evolution 350
DU 730	BioMate 160	Ultrospec 8000	Evolution 201	7310	GENESYS 30 or 40	UV-2100U	Evolution 201	Evolution 500/600	Evolution 350
DU 800	Evolution 260 Bio	Ultrospec 9000	Evolution 220	Genova	BioMate 160	UV 2400, 2450	Evolution 350	<b>Turner</b>	
<b>Biochrom/WPA</b>		GeneQuant/GeneQuant II	BioMate 160	Aquanova/Aquagem	GENESYS 40	UV 2500, 2550, 2600	Evolution 350	320/340	GENESYS 30
C07500	GENESYS 30	GeneQuant Calculator	BioMate 160	<b>Merck-Millipore</b>		BioSpec-1601	BioMate 160	390	GENESYS 50
S800	SPECTRONIC 200	GeneQuant pro	BioMate 160	Pharo100	GENESYS 30 or 40	BioSpec-Mini	BioMate 160	690	GENESYS 30
S1200	SPEC 200 or GENESYS 30	Biochrom 4060	BioMate 160	Pharo300	GENESYS 150	<b>SI Analytics/WTW/Secomam</b>		SP-830/830+	GENESYS 30
Lightwave/Lightwave II	GENESYS 50	<b>Hach</b>		Prove 100	GENESYS 30 or 40	UviLine SI 5000	SPECTRONIC 200	SP-850	GENESYS 30
Biowave II	BioMate 160	DR 2/3	GENESYS 30	Prove 300/600	GENESYS 150	UviLine SI 6000	GENESYS 30	SP-870/870+	GENESYS 40
Libra S50	GENESYS 150	DR 2000/2010/3000	GENESYS 30	<b>Mettler-Toledo</b>		UviLine SI 6100/ Photolab 7100/ UviLine 9300	GENESYS 40	SP-890	GENESYS 50
Libra S60	GENESYS 180	DR 3800/3900	GENESYS 30 or 40	UV5	GENESYS 150	UviLine SI 7000/ Photolab 9600	GENESYS 50	<b>Unico</b>	
Libra S70	Evolution 201	DR 4000 V	GENESYS 30	UV5 Bio	BioMate 160	UviLine SI 7100	GENESYS 50	1100/1100RS	SPECTRONIC 200
Libra S80	Evolution 220	DR 4000 U	GENESYS 150	UV7	Evolution 201	<b>Thermo Scientific/SPECTRONIC/Milton Roy</b>		1200/1205	GENESYS 30
<b>Bio Rad<sup>1</sup></b>		DR 5000/6000	GENESYS 150 or Evolution 201	<b>PerkinElmer</b>		SPECTRONIC 20/20+/20D	SPEC 200 or GENESYS 30	2100-S/2150	GENESYS 140
SmartSpec/SmartSpec 3000	BioMate 160	<b>Hitachi</b>		Lambda 11	GENESYS 50	SPECTRONIC 21 MV/DV	GENESYS 30	S-2150UV	GENESYS 150
<b>Biotek</b>		U-1000	GENESYS 30	Lambda 12/20/25	Evolution 201	SPECTRONIC 21DUV	GENESYS 50	SQ2802/SQ3802	GENESYS 150
DS-C	GENESYS 50	U-1100/U-1800	GENESYS 150	Lambda 14/40/45	Evolution 220 or 350	SPECTRONIC 301	SPECTRONIC 200	SQ2802S	GENESYS 150 or Evolution 220
<b>Eppendorf</b>		U-2000/U-0080D	Evolution 201	Lambda 30/35	Evolution 220	SPECTRONIC 401/501/601	GENESYS 30	SQ4802	GENESYS 180
BioPhotometer	BioMate 160	U-2800/U-2900	Evolution 220	Lambda 265	GENESYS 150	SPECTRONIC 1001	GENESYS 50	<b>VWR House Brand (all models OEM by Mapada)</b>	
BioPhotometer D30	BioMate 160	U-3010/U-3310	Evolution 350	Lambda 365	Evolution 220 or 350	SPECTRONIC 1201	GENESYS 150	V1200	GENESYS 30
BioSpectrometer	BioMate 160 or Evo 260 Bio	U-5100	GENESYS 150	Lambda 465	Evolution 201 or 350	SPECTRONIC 3000	GENESYS 180	M4/P4	GENESYS 50
BioSpectrometer Kinetic	BioMate 160 with Peltier Accy	UH-5300	Evolution 201	EZ150/EZ201	GENESYS 150	GENESYS 2/5	GENESYS 150	UV1600PC	GENESYS 50 or 150
<b>Implen</b>		<b>Jasco</b>		Lambda XLS/XLS+	GENESYS 150	GENESYS 10 UV	GENESYS 50 or 150	UV-3000-PC	GENESYS 30 or 40
NanoPhotometer C40	BioMate 160	V-530/V-630/V-730	Evolution 201	<b>Pye-Unicam<sup>2</sup></b>		GENESYS 10 Vis	GENESYS 30 or 40	UV-3100PC	GENESYS 50 or 150
<b>Jasco</b>		V-550/V-650/V-660	Evolution 350	PU-8600/PU-8660	GENESYS 150	GENESYS 20	GENESYS 30	UV-6300PC	Evolution 201
V-630 Bio/V-730 Bio	Evolution 260 Bio	V-630 Bio/V-730 Bio	Evolution 260 Bio	SP6-300	GENESYS 150	Evolution 100	GENESYS 150 or 180	<b>Footnote 1:</b> These brands/vendors have exited the UV-Vis market. Although the brands still exist selling other instruments, there is no formal successor brand to these models.	
V-750/V-760	Evolution 350	V-750/V-760	Evolution 350	PU-8800	Evolution 350	Evolution 160	BioMate 160	<b>Footnote 2:</b> The Turner and Pye-Unicam brands were acquired by Thermo Scientific. The current Thermo Scientific line represents the successor products for these brands.	

continued

Find out more at [thermofisher.com/uv-vis](http://thermofisher.com/uv-vis)

**ThermoFisher**  
SCIENTIFIC