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Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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Assembling of the dispersion tool:

- Open the knurled screw (pos. 10) so that the thread does not protrude into the hole.
- Insert the dispersing element (pos. 9) from the bottom into the drive unit as far as the stop. Ensure that the dispersion tool is aligned correctly when inserted (refer to picture): the tabs must project sideways when viewed from the front of the drive unit/dispersion tool. The shaft will audibly click into place when a small amount of force is applied. The dispersion tool is correctly mounted when the edge of the projecting part is aligned with the lower edge of the drive unit (see fig. 2).

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Safety instructions

- The dispersion tool may heat up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
- When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.

- No liquid may get into the drive unit of the dispenser.
- While working with the dispenser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- The work has to be interrupted immediately if you notice unusual noise and/or emission of liquid from the upper rinsing drilling. The seals must be controlled, and replaced if necessary.
- Note the operating instructions of the dispenser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- The gasket and bearings are made of PTFE and non-ferrous steel, the following points should therefore be noted: **Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorofluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.** (Source: Rompps Chemie-Lexikon and "Ulmann", Volume 19)

Product information

Rotor-Bearing: the shaft bearing of the dispersion tool is a sleeve bearing bushing (material PTFE). Output - laterally: the shaft of the rotor is inserted in the coupling of the drive unit. The plastic material used are approved by FDA (Food and Drug Administration Approval).

Application instructions

- Attention!** Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed. Always ensure that the dispersing element does not run dry in some cases should spout develop.
- The fill of the medium must be at least 10 mm above the lower overflow hole of the shaft.
- While working with the dispenser, the user must wear his personal protective equipment.
- Longer application times bring hardly any improvements, increase however the sample temperature substantially.
- The function of the dispersion tool depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the dispense effect diminishes.

Spare parts list

Item	Designation	Item	Designation
1	Start tube	1	Start tube
2	Shaft	2	Shaft
3	Rotor	3	Rotor
4	Grooved piston ring bearing	4	Stator
5	Shaft	5	Rotor
6	Shaft key	6	Shaft key

Some parts dispocan see front page. For orders of spare parts, please specify the type of dispersion tool and the designation of the spare part.

Maintenance and cleaning

Allowed procedures	Sterilization techniques	Chemical procedures
Damp heat	Autoclave to 121 °C at 2 bar positive pressure	By germ-killing disinfectants alcohol, phenol formalin...
Hot air	Germ-killing by hot air with 160 to 190 °C (approx. 30 min)	Disinfectant: remanents are to be removed with germ-free water

For cleaning the dispersion tool, operate it in a disinfectant for dissolving the residues of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residues of the substance don't cultivate unwanted bacterial cultures. For cleaning purposes the dispersion tool can be disassembled by the flat key and the shaft key, which are parts of the tool sets of the drive, how described below.

- 5 18 N-10 G:** With the flat key (locking of the rotor (3)) and the shaft key (turning of the shaft (2) clutch-laterally), pull the shaft (2) out upwards of the shaft tube (1). Now press the grooved piston ring bearing (4) out upwards of the shaft tube (1).

- 5 18 N-19 G:** Unscrew the shaft (1) from the shaft tube (1) by the flat key (Attention! Left-handed thread) and pull it out of the shaft tube (1). Operate the rotor (5) with the shaft key (3). Unscrew the rotor (5) from the shaft (3). Pull the shaft key on the end of the shaft and turn the shaft key. The assembly of the dispersion tool takes place in reverse order. **Attention!** The rotor must not be tightened too strongly, because otherwise the small teeth at the rotor become bent and so the dispersion tool gets destroyed.

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