Kunststoff - Dispergierwerkzeuge Plastic dispersion tools Outils dispersants

S18D – 14G – KS







S25D – 14G – KS

S18D – 10G – KS

S25D – 10G – KS

IKA®-Werke GmbH & Co. KG

Janke & Kunkel-Str. 10

D-79219 Staufen


Tel. +49 7633 831-0

Fax +49 7633 831-98

sales@ika.de

Technische Information
Technical information
Informations techniques

DE
EN
FR

3451600b

www.ika.com

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- Safety instructions**
- The dispersion tool may heat up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
 - When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.
 - Don't touch rotating parts during operation.
 - No liquid may get into the drive unit of the disperser.
 - Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
 - While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
 - Follow the health and safety regulations accident prevention regulations.
 - Do not use any damaged dispersion tool e.g. with hairline cracks or splits.
 - Note the operating instructions of the disperser unit.
 - When using disposable tube **518/25-ET50**, ensure that the dispersing element has fully engaged in the cover.
 - The maximum capacity is 30 ml.
 - The drive is only allowed in a vertical working position.
 - The distance between the stand and the dispersing instrument must not exceed 80 mm.

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Product information

- The dispersion tool comprises shaft tube, rotor shaft and two bearing shafts.
- Dispersion tools made of plastics are not sterile and not free of DNA or RNA.
- The tools are ideal for applications, with which no cross contamination is permissible.
- The tools are ideally appropriate for a single-use. However they can be used, depending upon the application, several times too.
- In this case please consider the cleaning references.
- Use disposable tube **518/25-ET50** for working with the plastic dispensing element in a closed system.

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Application instruction

- Attention!** Never run the dispersion tool dry. The bearing will be destroyed without cooling by the medium.
- Always ensure that the dispensing element does not run dry in some cases should sprouts develop.
- The fill level of the medium must be between the minimum and maximum marks.
- Medium may leak out of the overflow hole around the maximum mark.
- The dispersion tool are not suitable for continuous operation.
- Dispensing causes the medium to heat up.
- The dispersing element must not be operated at the maximum permitted speed of 25000 rpm for more than one minute.
- Store the dispensing elements in a dry place, in their packaging and at room temperature.
- The solid bodies in samples must not be frozen or too hard.
- Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient; maximum duration is one minute.
- Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispose of the tool after use according to the applicable regulations.

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Materials and chemical resistance

The shaft tubes are made of polycarbonate (PC), the rotor shafts and bearing shafts are made of polyetheretherketone (PEEK).
 The plastic materials used conform to the FDA (Food and Drug Administration Approval) - rules.
 The parts have a good stability against weak acids, chlorides, hypochlorites and many other chemicals.

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Cleaning, disinfection and sterilization

The dispersion tools are not sterile, but they can be sterilised before use.

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	Not allowed procedures	Not allowed procedures
Method of sterilisation	Autoclave up to 125 °C with 2 bar UV-irradiation Hot-air sterilisation 160-190 °C	UV-irradiation Hot-air sterilisation 160-190 °C
Biological disinfection	Alcohol (70%) Hypochlorite (5%) Formaldehyde (40%)	Phenol/Deviate Expanded use of quindin
Radioactive decontamination	Radiowash Compt - Off Water / ethyl alcohol / SDS	

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